



ICE is Easy A.H. Grange NERL, ESD, ECB



ICE is Easy

Written and Prepared by:

Andrew H. Grange, Ph.D.

Environmental Chemistry Branch

P.O. Box 93478, Las Vegas, NV 89193

grange.andrew@epa.gov

* Ion Composition Elucidation

Mt. Hood, Oregon

Click

Acronyms

MPPSIRD – Mass Peak Profiling from Selected Ion Recording Data

PGM - Profile Generation Model

ICE is Nice, which is available from grange.andrew@epa.gov, describes **MPPSIRD** and the **PGM**, and illustrates applications of **ICE**. It should be viewed before **ICE is Easy**.

During this presentation:

**right click to pause, then
left click to resume.**

Numerous column labels contain “M+1” or “M+2”. Fragment ions are also investigated using ICE. Ignore the “M” in these headers.

Determining an Ion Composition

Select an ion in the low resolution mass spectrum





Acquire Survey Data

Acquire Full Profile with 10,000 Resolution

**Enter Exact Mass into the PGM –
Select Hypothetical Composition**

Acquire Full or Partial Profile Data

**Enter Exact Masses and
Relative Abundances into the PGM**



Ion chromatograms do not track m/z 198

Click

Ion compositions will be determined for the apparent molecular ions from background subtracted mass spectra for 3 chromatographic peaks in the total ion chromatogram.

Prior Preparation ...

If MPPSIRD has not been performed recently, it is prudent to perform:

- 1. A magnetic calibration using PFK over a mass range that includes at least 50 Da below and 100 Da above all masses to be studied.**


.pfkcal 50,500,10,3

- 2. An electric scan calibration**

.aec

**If lock-on fails, perform both calibrations.
This generally solves the problem.**

These procedures are described in screen 37.



ENTER the Resolution, Center Masses & Times desired.

Allow at least 17 sec between each end time and the next start time.

ENTER <Cntl>G after last entry.

ENTER <Cntl>B after data has been acquired.

STEP Resolution: 1000

GC/Probe (G or P): G

SLIT Resolution: 3000

Scan Speed: 1.00

CI: EI

Precede the times with a ' to indicate the times are characters, not values.

file: demo1

LK: PFK

Center Mass

Start Time

End Time MID Grp

186.97

'17:35

'17:46

133.07

'20:50

'21:02

198.07

'29:01

'29:13

Click

▲
▼
▶
◀
+100
-100
+10
-10
+1
-1
+0.1
-0.1

-11.0 V **SCAN** **AUTOTUNE** **TUNE** **R** = 906 2 728 811 uV

p

1.00 mA
70.0 V

250 C
257 C


34
34

-62.0 V
-62.0 V

**Within the *tune* window
click **tune** to start
scanning across a single
mass peak profile.**

EI POS magnetic
U(a) 4772 V **u** U(esa) 727.4 V

Entr. Slit	(ENS)	340.0	
Exit Slit	(EXS)	390.0	
Sweep Mass	(MASS)	180.73	amu
Virt. Mass	(VMASS)	180.73	amu
Link Mass	(LMASS)	180.73	amu
Sweep Width	(SW)	1.00	%
Sweep Speed	(SS)	1.0	/s
U(a) Offset	(UOFS)	0.0	V
U(a)/U(esa)	(URATIO)	2.025	%
Multiplier	(EMULT)	1.60	kV
Dyn. Volt.	(DYNODE)	pos	
Ref. Inlet	(TREF)	150	C
Beam Rot.	(UROT)	-0.5	V
Focus Quad	(FQUAD)	-7.4	V



VG off VVR off
VV off VS2 off
 VS3 off

TUNE DESCRIPTOR
eiposag
 RESTORE

179.83 180.73 181.63

<< >> <<< >>> ^ v C V LOCK UNLK

TUNE : .midty166 **Enter .midty166**


Click

Navigation buttons:

- ▲
- ▼
- ▶
- ◀
- +100
- 100
- +10
- 10
- +1
- 1
- +0.1
- 0.1

```


GC-LC Type: HP 6890 GC          GC not ready          IC < GC
GC Descriptor
On Display: gc37                vmake 1;#gc
In GC      : gc37                vsend;#gc;#sc150
Current GC Temperature          doze4
GC elapsed time 0.0 Min.
GC run time 10.0 Min.
Stabilize time 0.0 Min.
  
```



The column temperature is set to 150°C for autotuning. The ion source conditions are then more similar to those when analytes elute.

```

Column: Constant Flow 0.1 ml/min
Inlet: Splitless Mode
Purge time 1.0 Purge flow 70.0 Saver time 3.0 Saver flow 15.0
  
```




START GC
 RESTORE
 LOAD
 CYCLE
 New Program

GC: _

Click

▲
▼
▶
◀
+100
-100
+10
-10
+1
-1
+0.1
-0.1

-11.0 V **SCAN** **AUTOTUNE** **TUNE** **R** = **2 532** 599 869 uV




EI POS magnetic
U(a) **4532** U(esa) 727.4 V

Enter
E:
S:
V:
L:
S:
S:
U:
U:
M:
D:
Re...
Beam Rot. (UROT) -0.5 V
Focus Quad (FQUA

VG off VVR off
VV off VS2 off
 VS3 off

TUNE :



```
vmake 1;#tune
ratio 1.05:doze 3
tune:doze .5
vsend;#tune;#m192.98882
doze 1
vsend;#tune;#w1
lock:doze 9
sr 3000
vsend;#tune;#w0.100
lock:doze 9
vsend;#tune;#r
doze 2
```

If the tune display “locks up”, simply click on TUNE to restart *midtune* from the lock-up point.

The PFK lock mass for the last SIR group becomes the center mass of the tune display.

] **LOCK** UNLK

Click

NOT ACQUIRING

ELAPSED :

12 Jan 2001 10:31


END :

```

FILE : demol          STUDY :          Time : MID run time
SAMP : demol          AMT : 0.00      VOL : 0.00
OPER :                CLIENT :         INJ.VOL : 1.00
COMM : Clark Co. effluent

```

START



METHOD LIST :

init(1)

SAVE

FILE	AUTOSAMP	VIAL	D:INLET	ICL PROC.	DS PROC.	S/T/M	DLY
demol	pm	1	GC:gc37	simchr	mpp.prn	\$	mid37

```

vmake 1;#analysis
doze 2
vsend;#analysis;#f1^demol^a^pm^v^1
vsend;#analysis;#g^gc37^i^simchr
vsend;#analysis;#d^mpp.icl $^m^mid37
doze 5

```

To run the correct GC temperature program, save the GC program desired as GC37.

ANALYSIS :

***The puzzle for running this program automatically has not been solved. Click**

midtune.icl

```

:
:
:
vsend;#analysis;#d^mpp.icl^m^mid37
doze 5
vmake 1;#mid
doze 2
mpp1
vmake 1;#tune
doze 2
autotune

```

The *mid* view will appear and subroutine *mpp1.icl* will be run to prepare a SIR descriptor for each eluting GC peak.

```

reset 3000
vsend;#tune;#r
vmake 4;#tune;#gc;#mid;#analysis
vsend;#tune;#w0.5
ratio 1
lock:doze 9
vsend;#tune;#w0.050
vsend;#tune;#r
doze 1
vsend;#gc;#sc90
vmake 1;#mid:doze 3
repeat 5
.#click >& enter mpp2, if > 9 groups
doze 1:##
doze 1:end
vmake 4;#tune;#gc;#mid;#analysis

```

Note that subroutines such as *mpp1.icl* can be run independently of *midtune.icl*. This will be illustrated in the next screen.

MID Set Up Parameters

```

MID File          mid37
Measure/lock ratio (X)      1
Set Damping relay (T)      TRUE
Width first lock (A)       1.00 amu
Electric jump time (E)     10 ms
Magnetic jump time (D)    60 ms
Offset (O)                1 cts
Electric range (R)        300 %
Sweep peak width (W)      5.00
Acq mode (C|P)           Cent. mode
MID mode (J|M|L|N)       Lock mode
    
```

MID Time Windows

#	Start	Measure	End	Cycle	Time
1	29:07	0.31	29:38 min	1:00	sec
2	0:00				
3					
4					
5					
6					
7					
8					
9					

Similarly, the m/z ratios and group times would appear for Groups 2 and 3.

- Clear Menu
- Clear Times
- Clear Masses
- Start MID
- SAVE
- Main

MID Masses for Time Window 1

#	mass	F	int	gr	time (ms)
8	186.8204		1	1	20.48
9	186.8391		1	1	20.48
10	186.8578		1	1	20.48
11	186.8765		1	1	20.48
12	186.8952		1	1	20.48
13	186.9139		1	1	20.48
14	186.9326		1	1	20.48
15	186.9513		1	1	20.48
16	186.9700		1	1	20.48
17	186.9887		1	1	20.48
18	187.0074		1	1	20.48
19	187.0261		1	1	20.48
20	187.0448		1	1	20.48
21	187.0635		1	1	20.48
22	187.0822		1	1	20.48
23	187.1009		1	1	20.48
24	187.1196		1	1	20.48
25	187.1383		1	1	20.48
26	187.1570		1	1	20.48
27	192.9760		1	1	20.48
28	192.9824		1	1	20.48
29	192.9888	c	1	1	20.48
30	192.9953		1	1	20.48
31	193.0017		1	1	20.48

- Lock Mass
- Cali Mass

MID: .mpp13 **Enter .mpp1**

Click


Navigation buttons:

- ▲
- ▼
- ▶
- ◀
- +100
- 100
- +10
- 10
- +1
- 1
- +0.1
- 0.1

-11.0 V **SCAN** **AUTOTUNE** **TUNE** **R** = **3 077** 362 264 uV

p -62.0 V

The ion source is autotuned and subroutine *reset 3000* is run to adjust the resolution to 3000 x (1 to 1.05). The resolution is then shown.



**vmake 1:.#tune
doze 2
autotune
reset 3000
vsend;.#tune;.#r**

Entr. Slit	(ENS)	183.5	
Exit Slit	(EXS)	189.8	
Sweep Mass	(MASS)	192.989	amu
Virt. Mass	(VMASS)	192.989	amu
Link Mass	(LMASS)	192.989	amu
Sweep Width	(SW)	0.10	%
Sweep Speed	(SS)	1.0	/s
U(a) Offset	(UOFS)	0.0	V
U(a)/U(esa)	(URATIO)	2.025	%
Multiplier	(EMULT)	1.60	kV
Dyn. Volt.	(DYNODE)	pos	
Ref. Inlet	(TREF)	150	C
Beam Rot.	(UROT)	-0.5	V
Focus Quad	(FQUAD)	-7.4	V

VG off VVR off
VV off VS2 off
 VS3 off

TUNE DESCRIPTOR
eiposag
 RESTORE

192.892 192.989 193.085

Control buttons: << >> <<< >>> ^ v C V LOCK UNLK

TUNE :

Click

GC-LC Type: HP 6890 GC GC ready IC < GC
 GC Descriptor: GC elapsed time 0.0 Min
 On Display: gc37 GC run time 32.0 Min
 In GC : gc37 Stabilize time 0.0 Min



MID Set Up Parameters

MID File mid37
 Measure/lock ratio (X) 1
 Set Damping relay (T) TRUE
 Width first lock (A) 1.00 amu
 Electric jump time (E) 10 ms
 Magnetic jump time (D) 60 ms
 Offset (O) 1 cts
 Electric range (R) 300 %
 Sweep peak width (W) 5.00
 Cent mode 9
 Lock mode 10


MID Masses for Time Window 4

#	mass	F int	gr time(ms)
1		0.0000	
2			
3			
4			
5			
6			
7			
8			
9			
10			
11			
12			
13			
14			
15			
16			
17			
18			
19			
20			
21			
22			
23			
24			

You are now ready to collect:

MPPSIRD

Double left click **START** in the *analysis* view.



NOT ACQUIRING 5 Jun 2001 11:32 ELAPSED: END:

FILE:demo1 STUDY: Time: MID run time
 SAMP:demo1 AMT: 0.00 VOL: 0.00
 OPER: CLIENT: INJ.VOL: 1.00
 COMM:Clark Co. effluent

START

pm gc37 mpp mpp.icl mid37


METHOD LIST:
 init(1)
 SAVE

FILE AUTOSAMP VIAL D;INLET ICL PROC. DS PROC. S/T/H DLY
 1>demo1 pm 1 GC:gc37 mpp mpp.icl mid37

ANALYSIS:

Click

GC-LC Type: HP 6890 GC **GC method is running** IC < GC
 GC Descriptor GC elapsed time 0.0 Min
 On Display: gc37 GC run time 23.0 Min
 In GC * 8077 GC stabilize time ^ ^ Min



Simchr.icl sets the filament current to 0 until shortly before the first SIR group start time.

awaits the first SIR descriptor start time.

0.00 mA
70.0 V
250 C
251 C

EI POS U(a) 4769
 Entr. Slit
 Exit Slit
 Sweep Mas:
 Virt. Mas:
 Link Mass
 Sweep Widl
 Sweep Speed (SS) 1.0 /s
 U(a) Offset (UOFS) 0.0 V
 U(a)/U(esa) (URATIO) 2.025 %
 Multiplier (EMULT) 1.60 kV
 Dyn. Volt. (DYNODE) pos
 Ref. Inlet (TREF) 150 C
 Beam Rot. (UROT) -0.5 V
 Focus Quad (FQUAD) -7.4 V

VC off WVR off TUNE DESCRIPTOR
 VS2 off eiposag
 VV off VS3 off RESTORE

192,941 192,989 193,03

MULTIPLE ION DETECTION

Lock Mass

Int 5
Mass 181.0

STATUS

MID File mid37
 # Windows 3
 Measuring # 0
 Mode Lock
 MID is **not running**
 Elapsed time 23.0 min
 ECORR 1.000219
 EDACC 0.999820
 EDACZ 1435

Cali. Mass

Int 11
Mass 193.0

Remarks:

Start MID SAVE Setup

NOT ACQUIRING

After the last SIR group data is recorded, the MID vic

Signal is recorded for the interval between the start and end times specified by the user.

Done reading defaults file

Done reading defaults file

FILE QuanF write 1_tbl 2_tbls Rep_1 P


↑ ↓ N ← → DA All s&rt rt&a 1000 mmu 50 mmu 2 mmu d_all

QUAN: andreas.q0
Data file(s) used
Fill: demol 13-Dec-00
Samp: Vial 1 demol
Comm: Clark Co. effluent
Mode: EI +VE +LMR ESCAN LR NRM
Oper: Inlet

CHRO: demol 13-Dec-00 Elapse: 21:42.4 215

Disp: CHRO
Amount = (area * ref amnt) / (ref area * response fa
Response factors from quantitation file

Quan Masses	Height	Area	Actual Time
192.9727	45270.	509144.	28:55
192.9791	62917.	722844.	28:55
192.9855	67157.	798735.	28:57
192.9920	57144.	672876.	29:00
192.9984	38237.	449699.	29:00
197.8686	2.	6.	29:07
197.9082	3.	7.	29:04
197.9479	4.	14.	29:01
197.9677	802.	7095.	28:54
197.9875	2593.	30654.	28:54
198.0073	539.	3886.	29:01
198.0271	9.	31.	28:58
198.0271	8289.	18875.	29:05
198.0469	66315.	155028.	29:05
198.0667	38507.	91424.	29:05
198.0865	35.	90.	28:54
198.0865	161.	561.	29:05
198.1261	70.	337.	28:58
198.1261	203.	600.	29:05
198.1459	172.	869.	29:05
198.1855	8.	20.	29:05
198.2450	3.	5.	28:59
198.2648	1.	2.	28:55
204.9718	27088.	310320.	28:54
204.9787	39473.	453776.	28:54
204.9857	53908.	617290.	28:54
204.9926	40418.	473261.	28:54
204.9996	26507.	332648.	28:54



```










DECTerm
File Edit Commands Options Print Help
finnigan@icis5>
finnigan@icis5>
finnigan@icis5>
finnigan@icis6> cd ~/icl
finnigan@icis7> alias m csh mpp.icl
finnigan@icis8> m
finnigan@icis9>

```

Enter cd ~/icl
The cursor line reappears
after mpp.csh is mpp.icl.

Click


The adjacent plots are printed for archiving.



Click

ENTER the Resolution, Center Masses & Times desired.
Allow at least 17 sec between each end time and the next start time.
ENTER <Cntrl>G after last entry. ENTER <Cntrl>B after data has been acquired.
For all 3 SIR groups, the most abundant ion was an
analyte ion. Had $\frac{10000}{10000}$ FK ion had the greatest
abundance the user would estimate the mass of the
and <Cntrl>Z enters 10,000 for both enter it manually.
STEP and SLIT resolutions.

Center Mass	Start Time	End Time
186.07040	17:35	17:47
133.06490	20:52	21:04
198.05522	29:00	29:12



Enter <Cntrl>I to increment the Data Filename.
wide window 5 sec later, because the chromatographic peak maximum was missed.

Click

Survey Data Mass Increment

#	mass
1	180.9767
2	180.9828
3	180.9888
4	180.9948
5	181.0009
6	186.7830
7	186.8017
8	186.8204
9	186.8391
10	186.8578
11	186.8765
12	186.8952
13	186.9139
14	186.9326
15	186.9513
16	186.9700
17	186.9887
18	187.0074
19	187.0261
20	187.0448
21	187.0635
22	187.0822
23	187.1009
24	187.1196

**Calibrant ion of known mass
33 ppm**

**Analyte ion of unknown mass
100 ppm**

10,000 Resolution Mass Increment

#	mass
1	180.9852
2	180.9870
3	180.9888
4	180.9906
5	180.9924
6	186.9452
7	186.9471
8	186.9490
9	186.9508
10	186.9527
11	186.9546
12	186.9564
13	186.9583
14	186.9602
15	186.9621
16	186.9639
17	186.9658
18	186.9677
19	186.9695
20	186.9714
21	186.9733
22	186.9751
23	186.9770
24	186.9789

**Calibrant ion of known mass
10 ppm**

**Analyte ion of unknown mass
10 ppm**

1000 resolution

Click

**Review screens 7 through 17 to perform
MPPSIRD at 10,000 resolution.**

***Midtune.icl* contains two
additional instructions:**

***fq*
*ur***

**These are subroutines that autotune the
focus quadrupoles and the rotational
quadrupoles.**

**These adjustments are only important at high
mass resolution and are not performed by
*Autotune.***

Sheet Window Help
 EPA ICE is Easy A.H. Grange
 NERL, ESD, ECB
 Enter D to return to
 Enter 4P for the hypothetical composition.

m/z 186.95955 ± 6 ppm

C10 H22 N10 O10 P10 S10 F10 CL2

Enter the Hypothetical Composition #: P

#	Err(mmu & ppm)	RDB	Range	Composition	%M+1	%M+2
1	-0.6	-3.1	2.5	C2 H2 N4 O F CL2	3.74	64.21
2	-0.1	-0.4	1.5	C3 H6 N2 O P CL2	4.17	64.22
3	+0.8	+4.1	2.0 3.0	C4 H4 N O2 F CL2	4.91	64.45
4	-0.4	-2.0	6.0	C7 H3 N O CL2	8.16	64.45

The "P" automatically executes a modified C:\QB\1.BAS file to provide calculated mass peak profiles, exact masses, and relative abundances for the hypothetical composition.

Click

ENTER the Resolution, Center Masses & Times desired.

Allow at least 17 sec between each end time and the next start time.

ENTER <Cntrl>G after last entry.

ENTER <Cntrl>B after data has been acquired.

STEP Resolution: 10000
SLIT Resolution: 10000


GC/Probe (G or P): G
Scan Speed: 1.00
El or Cl: EI
Datafile: demo2
PFK or ALK: PFK



Center Mass	Start Time	End Time	MID Grp
186.95917	17:35	17:47	1
187.96234	20:52	21:04	2
188.95626	29:00	29:12	3

The 3 exact masses were entered into the mass list. Enter <Cntrl>Y, if an error was made (e.g. the wrong SIR group was considered) to remove the last masses entered.

Click



Enter 3P for the hypothetical composition.

m/z 133.06411 \pm 6 ppm

C11 N9 O8 F7 P4 S4 Si4 H132

Enter the Hypothetical Composition #: 3P

#	Err(mmu & ppm)	RDB	Range	Composition	%M+1	%M+2
1	-0.3 -2.3	2.5		C2 H6 N6 F	4.50	0.08
2	+0.2 +1.5	1.5		C3 H10 N4 P	4.93	0.10
3	-0.1 -0.8	6.0		C7 H7 N3	8.92	0.35

Again, the apparent molecular ion was prominent, suggesting aromatic character. Only the last composition has 4 or more rings and double bonds.

Click

ENTER the Resolution, Center Masses & Times desired.

Allow at least 17 sec between each end time and the next start time.

ENTER <Cntl>G after last entry. ENTER <Cntl>B after data has been acquired.

STEP Resolution: 10000
SLIT Resolution: 10000

GC/Probe (G or P): G
Scan Speed: 1.00
EI or CI: EI
Datafile: demo2
PFK or ALK: PFK



Center Mass	Start Time	End Time	MID Grp
186.95917	17:35	17:47	1
187.96234	20:52	21:04	2
188.95626	29:00	29:12	3
133.06400			
134.06661			

Only 2 exact masses were entered into the mass list.

Click

Sheet Window Help
 ICE is Easy A.H. Grange NERL, ESD, ECB
 ^demo2

Enter 16P for the hypothetical composition.

m/z 198.05576 ± 6 ppm

C16 N14 O12 F10 P6 S6 SI7 H196

Enter the Hypothetical Composition #: 16 P

#	Err(mmu & ppm)	RDB	Range	Composition	%M+1	%M+2
1	+0.8 +4.3	-1.0	0.0	C5 H18 O2 S SI2	16.77	12.67
2	+0.2 +1.0	0.0		C5 H14 O6 SI	11.01	5.01
3	-0.4 -2.0	1.0		C5 H13 N2 F2 P SI	11.50	3.86
4	+0.4 +2.2	5.5		C5 H8 N7 S	8.99	4.80
5	+0.1 +0.6	0.0		C6 H17 F P2 SI	11.93	3.91
6	+1.1 +5.4	-1.0	0.0	C6 H18 O S2 SI	13.56	13.17
7	-1.2 -5.8	0.0	1.0	C6 H15 O2 F S SI	12.77	8.84
8	+0.4 +2.2	0.0	1.0	C6 H14 O5 S	7.80	5.70
9	+0.4 +2.2	1.5		C6 H11 N O F3 SI	12.25	4.15
10	-0.9 -4.7	0.0	1.0	C7 H15 O F S2	9.56	9.47
11	+0.7 +3.4	1.5		C7 H11 N F3 S	9.04	4.79
12	-0.5 -2.3	6.0	7.0	C7 H7 N4 O2 F	9.36	0.79
13	+0.1 +0.3	5.0	6.0	C8 H11 N2 O2 P	9.79	0.83
14	-0.7 -3.5	5.5		C9 H10 N F2 SI	15.50	4.38
15	+0.9 +4.5	5.5	6.5	C9 H9 N O3 F	10.53	1.10
16	-0.3 -1.3	9.5	10.5	C12 H8 N O2	13.79	1.28

Of these, the apparent composition (m/z 198) was the least and is cheap, as the hypothetical composition. Only the composition with the fewest atoms is more correct in the majority of cases.

Click

A B C D E F G H I J K L

ENTER the Resolution, Center Masses & Times desired.

Allow at least 17 sec between each end time and the next start time.

ENTER <Cntl>G after last entry.

ENTER <Cntl>B after data has been acquired.

STEP Resolution: 10000

SLIT Resolution: 10000

GC/Probe (G or P): G

Scan Speed: 1.00

El or Cl: EI

Datafile: demo2

PFK or ALK: PFK




10.0 pts/prof

5 5 ST ET 5

Center Mass	Start Time	End Time	MID Grp
186.96392	17:35	17:46	1
133.0639	20:50	21:02	2
198.05522	29:01	29:13	3
133.064			3
134.06661			3
198.0555			3
199.05877			3
200.06129			3

The 3 exact masses for the third analyte were added to the mass list.

Click



ENTER the Resolution, Center Masses & Times desired.

Allow at least 17 sec between each end time and the next start time.

ENTER <Cntrl>G after last entry.

ENTER <Cntrl>B after data has been acquired.

STEP Resolution: 10000

GC/Probe (G or P): G

SLIT Resolution: 10000

Scan Speed: 1.00

Use 1.9 kV on the SEM to obtain adequate signals for the less abundant +1 and +2 profiles.

Review screens 7 through 17 to perform MPPSIRD at 10,000 resolution.

**MO end times into to assign ORK group numbers.
an to m/z 200. Release the button.**

Click

ICE is Easy A.H. Grange
NERL, ESD, ECB

Lock Mass

The first two
mpp1.lcl as p

Multiple analy


m/z 133

+1

Cal. Mass

MID Masses for Time


#	mass	F	int	Time (ms)
1	130.9894	1		21.84
2	130.9907	1		21.84




ors written by
e.icl are shown.

ll be monitored.

10	133.0633			21.84
11	133.0647			21.84
12	133.0660			21.84
13	133.0673			21.84
14	133.0687			21.84
15	133.0700			21.84
16	134.0606			21.84
17	134.0619			21.84
18	134.0633			21.84
19	134.0646			21.84
20	134.0659			21.84
21	134.0673			21.84
22	134.0686			21.84
23	134.0700			21.84
24	134.0713			21.84






25	134.0726			21.84
26	142.9891			21.84
27	142.9906			21.84
28	142.9920 c			21.84
29	142.9934			21.84
30	142.9948			21.84
31				



MID Masses for Time

#	mass	F	int	Time (ms)
1	180.9852	1		20.48
2	180.9870	1		20.48
3	180.9888			
4	180.9906			
5	180.9924			
6	186.9536			
7	186.9554			
8	186.9573			
9	186.9592			
10	186.9610	1		20.48
11	186.9629	1		20.48
12	186.9648	1		20.48
13	187.9567	1		20.48
14	187.9586	1		20.48
15	187.9605	1		20.48
16	187.9623	1		20.48
17	187.9642	1		20.48
18	187.9661	1		20.48
19	187.9680	1		20.48
20	188.9506	1		20.48
21	188.9525	1		20.48
22	188.9544	1		20.48
23	188.9563	1		20.48
24	188.9581	1		20.48
25	188.9600	1		20.48
26	188.9619	1		20.48
27	192.9850	1		20.48
28	192.9869	1		20.48
29	192.9888 c	1		20.48
30	192.9908	1		20.48
31	192.9927	1		20.48



Click

#	RDB	Range	Composition	M	M+1	M+2	%M+1	(%M+1 Range)	%M+2	(%M+2 Range)
44										
45										
46										
47	1	-2.0 0.0	H4 N O3 F3 S2	.95847	.95780 X	.95441 X	0.62	(0.15-1.22) X	8.27	(6.39-10.09) X
48	2	-0.5 0.5	N2 O2 F5 S	.96006	.95854 X	.95611	0.65	(0.17-1.33) X	4.80	(3.76-5.96) X
49	3	2.5 3.5	H2 N4 O3 F P2	.95862	.95591 X	.96279 X	0.16	(0.01-0.40) X	0.03	(0.00-0.13) X
50	4	1.5 4.5	H3 N4 O4 S2	.95957	.95812 X	.95561	1.03	(0.24-2.14) X	9.19	(7.39-11.10) X
51	5	-0.5 0.5	C H O3 F5 P	.95835	.96181	.96260 X	1.26	(1.05-1.47) X	0.02	(0.00-0.14) X
52	6	2.0 3.0	C H N O10	.96005	.96283	.96429 X	1.79	(1.37-2.26) X	0.01	(0.00-0.18) X
53	7	1.5 2.5	C H6 N2 O3 P3	.95913	.96082 X	.96333 X	1.51	(1.04-2.04) X	0.01	(0.00-0.09) X
54	8	3.0 5.0	C N3 O4 F2 P	.95945	.95996 X	.96362 X	1.54	(0.91-2.28) X	0.01	(0.00-0.11) X
55	9	1.5	C H4 N4 F S3	.95821	.95778 X	.95403 X	1.44	(0.46-2.68) X	10.77	(7.41-13.89) X
56	10	2.5 3.5	C H2 N4 O2 F P CL	.95879	.95867 X	.95590	1.23	(0.66-1.92) X	32.36	(25.37-39.99) X
57	11	2.0	C H4 N5 P CL2	.95814	.95727 X	.95519	0.88	(0.39-1.50) X	63.99	(51.37-77.73)
58	12	2.0 4.0	C2 H4 N O4 F P2	.95996	.96299	.96423 X	2.70	(2.18-3.26) X	0.00	(0.00-7.00) X
59	13	2.0 4.0	C2 H2 N O6 F S	.95869	.96074 X	.95537	2.93	(2.07-3.81) X	4.89	(4.09-5.81) X
60	14	0.5	C2 H8 N2 P S3	.95872	.95955 X	.95455 X	3.21	(1.76-4.82) X	11.71	(8.35-14.92) X


Again, a blank row between the last and preceding compositions indicate that **only the last composition was consistent with all 5 measured values**. All other possible compositions based on the Exact Mass of the m/z 187 ion were rejected by comparing the measured and calculated values of Exact Masses and Relative Abundances for the +1 and +2 partial profiles. At least 1 X is seen in each other composition's row.

77	31	5.0	C6 H5 N S3	.95841	.96064 X	.95435 X	7.92	(5.90-9.94)	11.48	(8.27-14.51) X
78	32	6.0 7.0	C6 H3 N O2 P CL	.95899	.96214	.95613	6.99	(5.98-8.02)	32.27	(25.36-39.79) X
79	33	6.5							4.10	(3.21-4.97) X
80	34	6.5							31.91	(25.08-39.30) X
81	36	10.5 1							89)	0.01(0.00-0.12) X
82	37	10.5							50)	32.74(25.20-41.04) X
83										
84	35	6.0	C7 H3 N O CL2	.95917	.96234	.95626	8.05	(6.90-9.24)	64.01	(50.11-79.05)
85										
86	Experimental Values:			.95915	.96232	.95616	7.39		64.3	
87										



The **hypothetical composition** was correct: **C₇H₃NOCl₂**

CM = 134.06661 PM = 134.06718
M+1 PROFILE

CM = 135.06910 PM = 135.07007
M+2 PROFILE



With a step resolution of 10,000, the 5% levels were not monitored for these broadened profiles. A step resolution of 6700 would usually ensure observation of 5% of profile maximum levels with a slit resolution of 10,000. Because both +1 and +2 partial profiles were monitored for two analytes, 10,000 resolution was preferred.



A	EI	EJ	EK	EL	EM	EN	EO	EP	EQ	ER	ES
43	m/z = 198.0	5538 ± 6 ppm			C16 N14 O12 F10 P6 S6 SI7 H196				Resolution: 10000		
44											
45	#	RDB Range	Composition	M	M+1	M+2	%M+1 (%M+1 Range)	%M+2 (%M+2 Range)			
46											
47	1	1.0	H11 N8 O P SI	.05627	.05504 X	.05337 X	3.92 (1.51-7.18) X	0.11 (0.00-0.40) X			
48	2	2.0	H7 N10 F S	.05599	.05355 X	.05177 X	0.99 (0.13-2.40) X	0.01 (0.00-0.16) X			
49	3	1.0	C H11 N8 P S	.05650	.05518 X	.05233 X	3.03 (1.36-5.26) X	0.02 (0.00-0.24) X			
50	4	2.0	C H8 N8 O F P	.05427	.05282 X	.05762 X	0.98 (0.38-1.78) X	0.10 (0.03-0.18) X			
51	5	-0.5 0.5	C2 H16 N3 O2 SI3	.05503	.05495 X	.05237 X	8.77 (4.09-14.27)	0.13 (0.00-0.56) X			
52	6	1.0	C2 H12 N6 O P2	.05478	.05527 X	.05833 X	2.46 (1.46-3.63) X	0.16 (0.08-0.25) X			
53	7	2.5 3.5	C2 H6 N7 O2 F2	.05510	.05521 X	.05888 X	2.56 (1.42-3.94) X	0.33 (0.17-0.50) X			
54	8	6.0	C2 H6 N10 SI	.05462	.05409 X	.05164 X	3.75 (1.60-6.49) X	0.01 (0.00-0.09) X			
55	9										
56	10										
57	11										
58	12										
59	13										
60	14										
61	15										
62	16										
63	17										
64	18										
65	19										
66	20										
67	21										
68	22										
69	23										
70	24										
71	25										
72	26										
73	27	1.0									
74	28	2.5									
75	29	6.0 7.0									
76	30	5.0 6.0									
77	31	5.5									
78											
79	32	9.5 10.5	C12 H8 N O2	.05550	.05877	.06130	13.73 (11.79-15.68)	1.25 (1.02-1.49)			
80											
81	Experimental Values:			.05538	.05860	.06125	13.27	1.26			

Alternatively, one could repeat data acquisition over a narrower time window to discriminate against the interfering ion.

The hypothetical composition was correct: $C_{12}H_8NO_2^+$

Click






**Lock
Mass**

m/z 198

+1

+2

**Cal.
Mass**



mass

mass

151	192.98496	543265		198.04956	1510919		199.05280	207608		200.05529	23122		204.98472	323446
152	192.98689	633411		198.05154	2453982		199.05479	329295		200.05729	39002		204.98677	441194
153	192.98882	760773		198.05352	3256842		199.05678	427492		200.05929	54867		204.98882	495096
154	192.99075	643767		198.05550	3535016		199.05877	467209		200.06129	61469		204.99087	432224
155	192.99268	452116		198.05748	2952665		199.06076	387838		200.06329	57528		204.99292	313807
156				198.05946	2132532		199.06275	282175		200.06529	43771			
157	192.98882	606444.5	192.98871	198.06144	1244582		199.06474	161460		200.06729	28612		204.98882	404451.5 204.98879
158														
159						198.05538			199.05860				200.06154	
160														
161						100 %			13.27 %				1.72 %	
162														


Click

C:\WINNT\System32\CMD.exe



c:\viva\HL123

	EI	EJ	EK	EL	EM	EN	EO	EP	EQ	
3	m/z = 152.99798 ± 6 ppm	C12 N10 O9 F8 P4 S4 CL4 H151						Resolution:		
5	#	RDB Range	Composition	M	M+2			%M+2	(%M+2 Range)	
7	1	2.5 3.5	C H3 N4 O2 F P	.99777	.00168	X		0.43	(0.36-0.50) X	
8	2	2.0	C H5 N5 P CL	.99711	.99416	X		32.01	(25.70-38.89)	
9	3	1.5 2.5	C2 H7 N2 O2 P2	.99828	.00245	X		0.44	(0.36-0.51) X	
0	4	3.0 4.0	C2 H N3 O3 F2	.99860	.00269	X		0.65	(0.54-0.76) X	
1	5	2.5	C2 H3 N4 O F CL	.99794	.99504			32.24	(25.89-39.15)	
2	6	1.5	C3 H7 N2 O P CL	.99845	.99556			32.25	(25.90-39.17)	
3	7	6.5	C3						(4.00-5.21) X	
4	8	1.5	C4					5	(51.44-77.84)	
5	9	2.5	C5						(4.00-5.21) X	
6	10	7.0 8.0	C5						(0.49-0.70) X	
7	11	6.0 7.0	C6 H4 N O2 P	.99797	.00286	X		0.62	(0.52-0.73) X	
8	12	6.5	C7 H2 O3 F	.99880	.00383	X		0.87	(0.73-1.02) X	
9	14	10.5	C10 H O2	.99765	.00336	X		0.97	(0.81-1.13) X	
1	13	6.0	C7 H4 N O CL	.99814	.99532			32.48	(26.09-39.43)	
3	Experimental Values:			.99798	.99505			31.14		

A single composition is now consistent with the 3 measured values.



+2 Profiles



However, for the second composition, a good estimate of the Exact Mass would still be provided. Monitoring the +2 profile would distinguish between these compositions.

UCLR USER LIST FILE (4x500): init

1:	2802.00000	999.64839	0.00000	0.00000
2:	2787.00000	999.65017	0.00000	0.00000
3:	2787.00000	999.65023	0.00000	0.00000
4:	0.00000	0.00000	0.00000	0.00000
5:	0.00000	0.00000	0.00000	0.00000
6:	0.00000	0.00000	0.00000	0.00000
..	^ ^ ^ ^ ^	^ ^ ^ ^ ^	^ ^ ^ ^ ^	^ ^ ^ ^ ^

Scan Parameter Display and Editor

Acquisition Mode: HCENT
 Scan Type: Electric
 Ionization Mode: EI POS
 Scan Function: TUNE sweep


SCAN File:
 SCAN TUNE BREAK SRANGE
 ASCAN ASTAT SAVE LRANGE

First Mass (FH): 50.0
 500.0
 10.00
 0.20 s
 10.0000 s
 181.0
 1000
 14
 10
 200
 4
 1.00 /s
 850

Channel Plate Voltage (UCHAN): 850

Scan Cycle Time: 10.200 s
 Sampling Frequency: 0.200 kHz
 Sweep Samples: 200
 100 Ions per Peak: 7368 uV

Gas Control:
 VG off WVR off VC1 off
 VW off VS2 off VAUX off
 VS3 off VC2 off



Three electric calibrations are made between the m/z 181 and 231 PFK ions and averaged, and the results are saved.

17:	v.vvvvvv	v.vvvvvv	v.vvvvvv	v.vvvvvv
18:	0.00000	0.00000	0.00000	0.00000
19:	0.00000	0.00000	0.00000	0.00000
20:	0.00000	0.00000	0.00000	0.00000
LINK	NONE	NONE	NONE	NONE
SIZE	3.00000	3.00000	0.00000	0.00000
MEAN	2792.00000	999.64959	0.00000	0.00000
SUM	8376.00000	2998.94878	0.00000	0.00000
S.D.	8.66025	0.00105	0.00000	0.00000
MAX	2802.00000	999.65023	0.00000	0.00000
MIN	2787.00000	999.64839	0.00000	0.00000
ULIST:				

-13.0 V SCAN AUTOTUNE TUNE R = 1000 1064551 uV

1.00 mA 3300 2817 2393 1379
 70.0 V 250 C 259 C 3300 2983 2407 1621

EI POS electric
 U(a) 4767 V U(esa) 763.8 V

Entr. Slit (ENS) 340.0
 Exit Slit (EXS) 390.0
 Sweep Mass (MASS) 180.99 amu
 Virt. Mass (VMASS) 180.99 amu
 Link Mass (LMASS) 180.99 amu
 Sweep Width (SW) 0.58 %
 Sweep Speed (SS) 1.0 /s
 U(a) Offset (UOFS) 0.0 V
 U(a)/U(esa) (URATIO) 1.925 %
 Multiplier (EMULT) 1.60 kV
 Dyn. Volt. (DYNODE) pos
 Ref. Inlet (TREF) 150 C
 Beam Rot. (UROT) -0.5 V
 Focus Quad (FQUAD) -7.4 V

VC off WVR off TUNE DESCRIPTOR
 VS2 off eiposag
 VW off VS3 off RESTORE

180.47 180.99 181.5

Peak Register Display Copy result to ULIST

REG 1/2: 230.963314 REG 1/3: REG 2/3:

MASS : 180.91101 230.86402
 REF MASS : 180.98882 230.98563
 PEAK WIDTH [ppm]: 699.99804 624.99114

Active Register: 2 0.30 1.00

N REG RCOPY1 RCOPY2 RCOPY3 LIMIT CENTER SPREAD ERASE RESUME PA

PEAK:

Click

A		B	C	D	E	F	G	H	I	JK	L
---	--	---	---	---	---	---	---	---	---	----	---

ENTER the Resolution, Center Masses & Times desired.

Allow at least 17 sec between each end time and the next start time.

ENTER <Cntl>G after last entry. ENTER <Cntl>B after data has been acquired.

STEP Resolution: 10000
SLIT Resolution: 10000

GC/Probe (G or P): G
Scan Speed: 1.00
El or Cl: El
Datafile: demo2
PFK or ALK: PFK





Center Mass	Start Time	End Time	MID Grp
186.96392	17:35	17:46	1
133.0639	20:50	21:02	2
198.05522	29:01	29:13	3
133.064			3
134.06661			3
198.0555			3
199.05877			3
			3
			3
			3
			3

The program exits back to the spreadsheet.

Click

These manuals discuss all of the MPPSIRD and the Profile Generation Model code.



However, the MPPSIRDware manual describes Lotus 123 v2.2 and WordPerfect 5.1 code no longer used. The manual has not been updated to describe the new Lotus 123 v9.x code currently used. The PGM manual in Quick Basic fails to describe only a few features that have been added to permit its use directly from the spreadsheet.

Articles about ICE

- Grange, A.H.; Osemwengie, L.; Brilis, G.; Sovocool, G.W. "Ion Composition Elucidation (ICE): An Investigative Tool for Characterization and Identification of Compounds of Regulatory Importance" *International J. of Environmental Forensics*, **2001**, 2, 61-74.
- Snyder, S.A.; Kelly, K.L.; Grange, A.H.; Sovocool, G.W.; Snyder, E.M. Giesy, J.P. "Pharmaceuticals and Personal Care Products in the Waters of Lake Mead, Nevada" Book Chapter, ACS Publications, August, **2001**.
- Grange, A.H.; Sovocool, G.W. "Mass Peak Profiling from Selected Ion Recording Data (MPPSIRD) as a Tool for Regulatory Analyses" *J. AOAC Int'l*, **1999**, 82, 1443-1457.
- Grange, A.H.; Sovocool, G.W. "Identifying Endocrine Disruptors by High Resolution Mass Spectrometry" *Advances in Analysis of Environmental Endocrine Disruptors*, **1999**, Ch. 9, 133-145.
- Grange, A.H.; Sovocool, G.W. "Determination of Elemental Compositions by High Resolution Mass Spectrometry without Mass Calibrants" *Rapid Commun. Mass Spectrom.*, **1999**, 13, 673-686.
- Grange, A.H.; Sovocool, G.W.; Donnelly, J.R.; Genicola, F.A.; Gurka, D.F. "Identification of Pollutants in a Municipal Well Using High Resolution Mass Spectrometry" *Rapid Commun. Mass Spectrom.*, **1998**, 12, 1161-1169.
- Grange, A.H.; Brumley, W.C.; Sovocool, G.W. "Powerful New Tools for Analyzing Environmental Contaminants: Mass Peak Profiling from Selected-Ion-Recording Data and a Profile Generation Model" *American Environmental Laboratory*, **1998**, 10, 1, 6-7.
- Grange, A.H.; Brumley, W.C. "A Mass Peak Profile Generation Model to Facilitate Determination of Elemental Compositions of Ions Based on Exact Masses and Isotopic Abundances", *J. Amer. Soc. for Mass Spectrom.* **1997**, 8, 170-182.
- Grange, A.H.; Donnelly, J.R.; Brumley, W.C.; Sovocool, G.W. "Determination of Elemental Compositions from Mass Peak Profiles of the Molecular (M), M+1 and M+2 Ions" *Anal. Chem.* **1996**, 68, 553-560.
- Grange, A.H.; Brumley, W.C. "Mass Spectral Determination" *Environmental Testing & Analysis* **1996**, March/April, 22-26.
- Grange, A.H.; Brumley, W.C. "Determining Elemental Compositions from Exact Masses and Relative Abundances of Ions" *Trends in Analytical Chemistry*, **1996**, 15(1), 12-17.
- Grange, A.H.; Brumley, W.C. "Identification of Ions Produced from Components in a Complex Mixture by Determination of Exact Masses and Relative Abundances Using Mass Peak Profiling" *LC•GC*, **1996**, 14, 478-486.
- Grange, A.H.; Donnelly, J.R.; Brumley, W.C.; Billets, S.; Sovocool, G.W. "Mass Measurements by an Accurate and Sensitive Selected-Ion-Recording Technique" *Anal. Chem.* **1994**, 66, 4416-4421.
- Grange, A.H.; Brumley, W.C. "Plotting Mass Peak Profiles from Selected Ion Recording Data" *Rapid Comm. in Mass Spectrom.* **1992**, 6, 68-70.

ICE Posters at <http://www.epa.gov/nerlesd1/chemistry/ecb-posters.htm>

Well Pollutants Identified With A New Mass Spectrometric Technique

Ion Composition Elucidation (ICE) of Ions from Trace Levels of Pharmaceuticals and Disinfection Byproducts in Water Supplies

Characterizing Hazardous Waste Constituents: A New Tool

Identification of Analytically Problematic Pollutants with a New Mass Spectrometric Technique

Mass Determination of Intact α -Chain Hemoglobin Adducts to within 0.2 Da Using MPPSIRD with Electrospray Ionization

Deconvoluting Overlapping Isotopic Patterns Using Mass Peak Profiling from Selected Ion Recording Data (MPPSIRD)

Determination of Elemental Compositions by High Resolution Mass Spectrometry without Mass Calibrants





Mention of trade or commercial products in this presentation does not constitute endorsement or recommendation for use by the Office of Research and Development (ORD) or the Environmental Protection Agency (EPA).

Note: although the ORD has peer reviewed this PowerPoint presentation, it does not necessarily reflect the views or opinions of ORD or EPA.

The End