

## METHOD 9078

### SCREENING TEST METHOD FOR POLYCHLORINATED BIPHENYLS IN SOIL

#### 1.0 SCOPE AND APPLICATION

1.1 The method may be used to determine the amount of PCB (polychlorinated biphenyl) contamination in soils such as sand, gravel, loam, sediment, and clay, assuming that PCBs are the sole source of organic halogens in the sample.

1.2 This electrochemical method is designed to provide quantitative field results over a range of 2 to 2000 µg/g PCBs, significantly cutting down on the number of samples requiring laboratory testing.

1.3 Chlorines are removed from the PCB molecule using an organo-sodium reagent. The resulting chloride ions are measured using a chloride specific electrode. Analysts must identify the type of Aroclor contamination in order to use this as a quantitative method.

1.4 This method is restricted to use by or under the supervision of trained analysts. Each analyst must demonstrate the ability to generate acceptable results with this method.

#### 2.0 SUMMARY OF METHOD

A sample of the soil to be tested is extracted with a hydrocarbon based solvent. The resulting extract is filtered to remove moisture and inorganic salts. The dried extract is reacted with metallic sodium and a catalyst to strip chloride from any PCB that may be present. The resulting chloride ions are extracted into an aqueous buffer solution where they are detected using a chloride ion specific electrode.

**CAUTION:** Some of the reagents used with this testing procedure contain flammable solvents, dilute acids, and metallic sodium. Wear gloves and safety glasses while performing tests. Read all MSDS and warnings included with the instrument before starting testing procedure.

#### 3.0 INTERFERENCES

3.1 This procedure is sensitive to any chlorinated compound that is preferentially soluble in a non-polar solvent. When analyzing for PCBs, the presence of other chlorinated organics will result in a high bias. Iodine and bromine containing compounds will affect results if present in significant quantities. Wet or dry samples may be run, but results for all samples are calculated on a wet-weight basis. In one evaluation study (Table 1), 1.4% of the measurements were false negatives.

3.2 Inorganic chlorides should not interfere using this method if the sample is extracted with organic solvent.

## 4.0 APPARATUS AND MATERIALS

Electrochemical PCB test kit: L2000® PCB/Chloride Analyzer, (Dexsil Corporation, One Hamden Park Drive, Hamden, CT), or equivalent. Each commercially available test kit will supply or specify the apparatus and materials necessary for successful completion of the test.

## 5.0 REAGENTS

Each commercially available test kit will supply or specify the reagents necessary for successful completion of the test. Reagents should be labeled with appropriate expiration dates.

## 6.0 SAMPLE COLLECTION AND HANDLING

6.1 See the introductory material to this chapter, Organic Analytes, Sec. 4.1.

6.2 Soil samples may be contaminated, and should therefore be considered hazardous and handled accordingly. All samples should be collected using a sampling plan that addresses the considerations discussed in Chapter Nine.

6.3 To achieve accurate analyses, soil samples should be well homogenized prior to testing. PCBs are generally not evenly distributed in a soil sample and extensive mixing must be done to assure consistency.

## 7.0 PROCEDURE

Follow the manufacturer's instructions for the test kit being used. Those test kits used must meet or exceed the performance specifications indicated in Tables 1 and 2.

## 8.0 QUALITY CONTROL

8.1 Follow the manufacturer's instructions for quality control procedures specific to the test kit used. Additionally, guidance provided in Chapter One should be followed.

8.2 Use of replicate analyses, particularly when results indicate concentrations near the action level, is recommended to refine information gathered with the kit.

8.3 Method 9078 is intended for field or laboratory use. The appropriate level of quality assurance should accompany the application of this method to document data quality.

## 9.0 METHOD PERFORMANCE

9.1 146 soil samples from a PCB contaminated site were analyzed. There were 114 individual samples and 32 field duplicates. Each sample was analyzed using both the L2000 and GC/MS. The L2000 analyses were performed on-site in a mobile lab and the PCBs were analyzed as Aroclor 1242. Laboratory analyses were performed on splits of the same samples. The results from the analyses are presented in Table 1.

9.2 After applying accepted statistical methods to account for the detection limit difference between the two methods the data were evaluated to determine the acceptability of the L2000 method. A matched-pair students t-test performed on the L2000 and CLP GC/MS data results in a t value of 0.2141. This is well below the critical value (1.645 @ 0.05) for rejecting the null hypothesis indicating that there is no statistical difference between the data pairs. An analysis of the data for outliers identified only 2 data points whose residuals were greater than 3 standard deviations (10 and 5 respectively). Both points were determined to be in error using other evidence and were eliminated from the data set. A linear regression analysis of the remaining data results in a correlation coefficient of 0.95 and a positive intercept of 10.98  $\mu\text{g/g}$ . The slope of 0.985 was not statistically different from 1 and the intercept was not statistically different from 0.

9.3 The relative percent difference (RPD) calculated from all valid duplicates greater than the L2000 detection limit of 2  $\mu\text{g/g}$  for each method resulted in a mean RPD of 19% for the L2000 data and a mean RPD of 43% for the CLP GC-MS method. A Dunnett's test shows that this is statistically significant.

9.4 In a second study, soil samples contaminated with Aroclor 1260 were taken during a site cleanup. The samples were split and sent for lab analysis by Method 8082 as well as analysis by the L2000 in the field. The results are reported in Table 2. A linear regression analysis of the data resulted in a correlation coefficient of 0.995, a slope of 1.048 and an intercept of -1.48  $\mu\text{g/g}$  indicating that the L2000 is accurate compared to the lab method. A calculation of the relative percent difference for data, where duplicates were run within a method, results in a lower RPD for the L2000 indicating a tighter data spread and better repeatability.

## 10.0 REFERENCE

1. Griffin, Roger D. Application of a New PCB Field Analysis Technique for Site Assessment. Proceedings of Hazmacon '92 March - April 1992.

TABLE 1

COMPARISON OF L2000 AND GC/MS RESULTS FROM SPLIT SAMPLES  
Summary of Results

| Sample Number | L2000 (µg/g) | GC/MS (µg/g) | Results Agree? |
|---------------|--------------|--------------|----------------|
| 1             | ND           | 0.593        | Yes            |
| 3             | ND           | 0.114        | Yes            |
| 4             | 23.6         | 6.71         | Yes            |
| 6             | ND           | 0.679        | Yes            |
| 7             | ND           | 0.552        | Yes            |
| 8             | 3.9          | 2            | Yes            |
| 9             | 6.9          | 1.3          | Yes            |
| 10            | 5.1          | 0.172        | Yes            |
| 11            | 2.7          | 1.15         | Yes            |
| 15            | 9.4          | 9.13         | Yes            |
| 15D           | 12.5         | 9.84         | Yes            |
| 16            | 484          | 2110         | Yes            |
| 17            | 6.5          | 2.55         | Yes            |
| 18            | 382          | 45.4         | Yes            |
| 19            | 71.1         | 6.7          | False Pos.     |
| 23            | 48.8         | 20.8         | Yes            |
| 25            | 3.5          | 11.7         | Yes            |
| 32            | 36           | 47.6         | Yes            |
| 33            | ND           | 6            | Yes            |
| 34            | 14.4         | 34           | Yes            |
| 36            | >2000        | 816          | Yes            |
| 38            | 778          | 1030         | Yes            |
| 40            | 5.7          | 4.25         | Yes            |
| 43            | 4.1          | 1.69         | Yes            |
| 43D           | 3.6          | 1.74         | Yes            |
| 50            | ND           | 3.6          | Yes            |
| 50D           | ND           | 4.4          | Yes            |
| 52            | 9.3          | 4.21         | Yes            |
| 53            | 25.7         | 0.958        | False Pos.     |
| 54            | 5.1          | 0.516        | Yes            |
| 55            | 4.4          | 2.4          | Yes            |
| 59            | ND           | 7.9          | Yes            |
| 60            | 2.3          | 0.624        | Yes            |
| 60D           | 4.4          | 0.577        | Yes            |
| 61            | 549          | 580          | Yes            |

TABLE 1 (cont.)

| Sample Number | L2000 (µg/g) | GC/MS (µg/g) | Results Agree? |
|---------------|--------------|--------------|----------------|
| 62            | 111          | 2.35         | False Pos.     |
| 64            | 172          | 19           | Yes            |
| 65            | ND           | 3.1          | Yes            |
| 66            | 2.1          | 1.98         | Yes            |
| 67            | 7.5          | 0.081        | Yes            |
| 68            | 8            | 0.504        | Yes            |
| 69            | 5.8          | ND           | Yes            |
| 69D           | 4.4          | ND           | Yes            |
| 73            | 37           | 15.8         | Yes            |
| 74            | 22           | 13.3         | Yes            |
| 75            | 61           | 23           | Yes            |
| 76            | 82           | 46.7         | Yes            |
| 78            | 21           | 2.27         | Yes            |
| 79            | 148          | 42.8         | Yes            |
| 80            | ND           | 3.8          | Yes            |
| 84            | 7.6          | 1.16         | Yes            |
| 84D           | 10.9         | 1.08         | False Pos.     |
| 85            | 593          | 428          | Yes            |
| 85D           | 596          | 465          | Yes            |
| 88            | ND           | 2.7          | Yes            |
| 88D           | ND           | 1.77         | Yes            |
| 89            | ND           | 45           | False Neg.     |
| 90            | 2            | 1.01         | Yes            |
| 90D           | ND           | 1.4          | Yes            |
| 91            | 1650         | 1630         | Yes            |
| 91D           | 1608         | 1704         | Yes            |
| 92            | 3.14         | 1.21         | Yes            |
| 92D           | 3.4          | ND           | Yes            |
| 95            | 20.6         | 17.5         | Yes            |
| 95D           | 20.1         | 31.2         | Yes            |
| 100           | 384          | 177          | Yes            |
| 100D          | 363          | 167          | Yes            |
| 101           | 8.3          | 1.21         | Yes            |
| 102           | 6.3          | 293          | False Neg.     |
| 102D          | 5            | 1.77         | Yes            |
| 103           | 75.2         | 40.3         | Yes            |
| 104           | 4.1          | 7.66         | Yes            |

TABLE 1 (cont.)

| Sample Number | L2000 (µg/g) | GC/MS (µg/g) | Results Agree? |
|---------------|--------------|--------------|----------------|
| 107           | 161          | 14.1         | Yes            |
| 108           | 6.1          | 3.84         | Yes            |
| 109           | P            | ND           | Yes            |
| 109D          | 10.3         | ND           | False Pos.     |
| 111           | 20           | ND           | False Pos.     |
| 112           | 240          | 315          | Yes            |
| 113           | 21.8         | 14.9         | Yes            |
| 114           | 107          | 66.3         | Yes            |

NOTE: 75 out of 146 samples are reported in Table 1. Samples that were found to be ND for both the L2000 kit and the GC/MS determination were not reported. The determination of a "false negative" result for the L2000 technique is based on an action level of 10 µg/g. If another action limit is chosen, the rate of false negative results may differ. Similarly, a "false positive" result for the L2000 technique is indicated when the L2000 results are above 10 µg/g and the GC/MS results are "ND" or below 10 µg/g, or when the results of the L2000 techniques are higher than the GC/MS results by more than two orders of magnitude.

ND = Not detected

6 False positives: ND - 14.1 ppm by GC/MS

2 False negatives: 2.7 - 293 ppm by GC/MS

71 Non-detects: ND - 2.5 ppm by GC/MS

TABLE 2

COMPARISON OF L2000 AND GC/EC RESULTS FROM SPLIT SAMPLES  
Summary of Results

| Sample Number | Method 8082 (µg/g) | L2000 Results (µg/g) |
|---------------|--------------------|----------------------|
| 1             | 83                 | 79/76                |
| 2             | 21                 | 22                   |
| 3             | 12                 | 14                   |
| 4             | 300/375            | 357/326/327          |
| 5             | 29                 | 27                   |
| 6             | 106/134            | 116/117              |
| 7             | 3                  | 7.6                  |
| 8             | 9.3                | 7.2                  |
| 9             | 1.5                | 5.2                  |
| 10            | 99                 | 93                   |
| 11            | 7/9                | 13                   |
| 12            | 3.6                | 12                   |
| 13            | 4.2/6.2            | 2.9                  |
| 14            | 290                | 254/265              |