

## 2 Objective

The objective of this study was to demonstrate that method FN-002-S05-02 ("BYI08330: Analytical method for the determination of BYI08330 and its metabolites BYI08330-enol, BYI08330-keto-hydroxy and BYI08330-MA-amide in soil and sediment by LC/MS/MS") can be performed with acceptable recoveries for determination of the compounds BYI08330, BYI08330-enol, BYI08330-keto-hydroxy and BYI08330-MA-amide at an independent laboratory having no prior experience with the method. The method was developed by Bayer CropScience LP, Stilwell, USA, and validated in that laboratory with results reported as Method FN-002-S05-02, by D.J. Netzband & J.M. Wade, in report dated January 25, 2006. Soil Höfchen (Germany) and soil Laacher Hof (Germany) were chosen as representative matrices for validation within the present study.

This study was performed in accordance with EC Guidance Document on Residue Analytical Methods, SANCO/825/00 rev. 7 of March 17, 2004, Commission Directive 96/46/EC amending Council Directive 91/414/EEC of July 16, 1996, and BBA Guideline on Residue Analytical Methods for Post-Registration Control Purposes of July 21, 1998 US EPA Ecological Effects Test Guidelines, OPPTS 850.7100 Data Reporting for Environmental Chemistry Methods, EPA 712-C-96-348, April 1996 and US EPA Residue Chemistry Test Guidelines, OPPTS 860.1340 Residue Analytical Method, EPA 712-C-96-174, August 1996.


## 3 Materials

### 3.1 Test and Reference Items

The structures for BYI08330 and its metabolites BYI08330-enol, BYI08330-keto-hydroxy and BYI08330-MA-amide are presented below:

#### **BYI08330**

Structural formula:



Chemical code:	BYI08330, AE 1302943
Chemical name (CAS):	cis-3-(2,5-Dimethylphenyl)-8-methoxy-2-oxo-1-azaspiro[4.5]dec-3-en-4-yl ethyl carbonate
CAS No.:	203313-25-1
Empirical formula:	C <sub>21</sub> H <sub>27</sub> N O <sub>5</sub>
Molecular weight:	373.45 g/mol


## Reference standard:

Certificate of analysis: AZ 12545  
 Purity: 99.2 %  
 Expiry date: May 2008  
 Origin: Bayer CropScience GmbH, PT – Analytics Frankfurt,  
 D-65926 Frankfurt am Main, Germany

**<sup>13</sup>C<sub>3</sub>-BYI08330**

(used as internal standard for the active ingredient BYI08330):

## Structural formula:




Code name: [azaspirodecenyl-2,3,4-<sup>13</sup>C<sub>3</sub>] cis BYI08330  
 Chemical name (CAS): 3-(2,5-Dimethylphenyl)-8-methoxy-2-oxo-1-azaspiro[4.5]dec-3-en-4-yl-2,3,4-<sup>13</sup>C<sub>3</sub> ethyl carbonate  
 Empirical formula: C<sub>18</sub> <sup>13</sup>C<sub>3</sub> H<sub>27</sub> N O<sub>5</sub>  
 Molecular weight: 376.41 g/mol

## Reference standard:

Batch no.: KML3386-1-2  
 Purity: 95.1%  
 Origin: Bayer CropScience AG, R-PT Isotope Chemistry  
 D-42096 Wuppertal, Germany

**BYI08330-enol**

## Structural formula:



Chemical code: AE 1302944  
 Chemical name (CAS): cis-3-(2,5-Dimethylphenyl)-4-hydroxy-8-methoxy-1-azaspiro[4.5]dec-3-en-2-one  
 CAS No.: 203312-38-3

Empirical formula:  $C_{18}H_{23}NO_3$   
 Molecular weight: 301.39 g/mol


## Reference standard:

Certificate of analysis: AZ 12333  
 Purity: 99.4 %  
 Expiry date: February 2008  
 Origin: Bayer CropScience GmbH, PT – Analytics Frankfurt,  
 D-65926 Frankfurt am Main, Germany

 **$^{13}C_3$ -BYI08330-enol**

(used as internal standard for the test item BYI08330-enol):

## Structural formula:




Code name: [azaspirodecenyl-2,3,4- $^{13}C_3$ ] FHN13777  
 Chemical name (CAS): cis-3-(2,5-Dimethylphenyl)-4-hydroxy-8-methoxy-1-azaspiro[4.5]dec-3-en-2-one-2,3,4- $^{13}C_3$   
 Empirical formula:  $^{13}C_3C_{15}H_{23}NO_3$   
 Molecular weight: 304.35 g/mol

## Reference standard:

Batch no.: KML 3384-1-1  
 Purity: 96.3%  
 Origin: Bayer CropScience AG, R-PT Isotope Chemistry  
 D-42096 Wuppertal, Germany

**BYI08330-keto-hydroxy**

## Structural formula:



Chemical code: AE 1422479

Chemical name (CAS): cis-3-(2,5-Dimethylphenyl)-3-hydroxy-8-methoxy-1-azaspiro[4.5]decane-2,4-dione  
 Empirical formula: C<sub>18</sub> H<sub>23</sub> N O<sub>4</sub>  
 Molecular weight: 317.38 g/mol


## Reference standard:

Certificate of analysis: AZ 11970  
 Purity: 92 %  
 Expiry date: September 2006  
 Origin: Bayer CropScience GmbH, PT – Analytics Frankfurt, D-65926 Frankfurt am Main, Germany

**<sup>13</sup>C<sub>3</sub>-BYI08330-keto-hydroxy**

(used as internal standard for the test item BYI08330-keto-hydroxy):

## Structural formula:




Code name: [azaspirodecane-2,3,4-<sup>13</sup>C<sub>3</sub>] FHN14066  
 Chemical name (CAS): 3-(2,5-Dimethylphenyl)-3-hydroxy-8-methoxy-1-azaspiro[4.5]decane-2,4-dione-2,3,4-<sup>13</sup>C<sub>3</sub>  
 Empirical formula: <sup>13</sup>C<sub>3</sub> C<sub>15</sub> H<sub>23</sub> N O<sub>4</sub>  
 Molecular weight: 320.35 g/mol

## Reference standard:

Batch no.: KML 3387-1-8  
 Purity: 99.4%  
 Origin: Bayer CropScience AG, R-PT Isotope Chemistry D-42096 Wuppertal, Germany

**BYI08330-MA-amide**

## Structural formula:



Chemical code: AE 1786350


Chemical name (CAS): cis-1-[[[(2,5-Dimethylphenyl)hydroxyacetyl]amino]-4-methoxycyclohexanecarboxylic acid  
 Empirical formula: C<sub>18</sub> H<sub>25</sub> N O<sub>5</sub>  
 Molecular weight: 335.40 g/mol

Reference standard:  
 Certificate of analysis: AZ 11525  
 Purity: 96.1 %  
 Expiry date: April 2006  
 Origin: Bayer CropScience GmbH, PT – Analytics Frankfurt, D-65926 Frankfurt am Main, Germany

**<sup>13</sup>C<sub>3</sub> -BYI08330-MA-amide**

(used as internal standard for the test item BYI08330-MA-amide):

Structural formula:



Code name: [acetyl-<sup>13</sup>C<sub>2</sub>, carboxylic acid-<sup>13</sup>C] FHN14065  
 Chemical name (CAS): 1-[[[(2,5-Dimethylphenyl)hydroxyacetyl-<sup>13</sup>C<sub>2</sub>]amino]-4-methoxycyclohexanecarboxylic-<sup>13</sup>C acid  
 Empirical formula: <sup>13</sup>C<sub>3</sub> C<sub>15</sub> H<sub>25</sub> N O<sub>5</sub>  
 Molecular weight: 338.37 g/mol

Reference standard:  
 Batch no.: KML 3389-1-4  
 Purity: 88.1%  
 Origin: Bayer CropScience AG, R-PT Isotope Chemistry D-42096 Wuppertal, Germany

**3.2 Test System**

The method was validated using two German soils *Höfchen* and *Laacher Hof*. Two different soils were used in order to assess a possible influence of different soil characteristics. The soil samples were classified according to DIN and/or USDA specifications. Soil characteristics of the used soils are summarised in Table 1. Complete soil parameterisation is reported in Table 9 and Table 10.

Table 1: Soil Types

Soil	Texture of Soil	Organic Matter [%]
Höfchen	silt loam (USDA)	1.58
Laacher Hof	sandy loam (USDA)	2.06

## 4 Experimental

### 4.1 Analytical Method

The recovery data for the study were generated using the following method, which gives full details of preparing the analytical sample extracts and the conditions for high performance liquid chromatography (HPLC):

Number of the method: FN-002-S05-02  
Title of the method: BY108330: Analytical Method for the Determination of BY108330 and its Metabolites BY108330-enol, BY108330-keto-hydroxy and BY108330-MA-amide in Soil and Sediment by LC/MS/MS  
Author of the method: Derek J. Netzband  
Jami M. Wade  
Bayer CropScience LP  
17745 S Metcalf Avenue  
Stilwell, Kansas 66085  
Reference: Method FN-002-S05-02  
Limit of quantitation: 0.005 mg/kg

The following sample sets were analysed:

Table 2: Level and Number of Recoveries per Fortification Level

Soil	Control sample	Level 0.005 mg/kg	Level 0.05 mg/kg
Höfchen	2	5	5
Laacher Hof	2	5	5

Additionally, two solvent blanks with internal standards added were analysed.

#### 4.1.1 Outline of the Method

Residues of BY108330 and its metabolites BY108330-enol, BY108330-keto-hydroxy and BY108330-MA-amide are extracted from soil using an acidic extraction solution in the presence of cysteine hydrochloride and utilizing microwave extraction. The extraction solvent consists of a mixture of water (containing 8 g/L cysteine hydrochloride), acetonitrile, ethyl acetate and formic acid. An aliquot of the final extract is analyzed by LC/MS/MS. Quantification of residues is based on the use of isotopically labeled internal standards and comparison of peak areas with those of known standards.

The final quantitative detection of BY108330 and its metabolites BY108330-enol, BY108330-keto-hydroxy and BY108330-MA-amide is accomplished by LC/MS/MS.

---

#### 4.1.2 Instruments

- Microwave Extractor:      MLS-Ethos  
                                  MWS Vertriebs GmbH  
                                  88299 Leutkirch, Germany
- Balances:                    PC 4400, PM 4800 and AT 261  
                                  Mettler Instruments GmbH  
                                  35387 Giessen, Germany
- Ultrasonic Bath:            Transsonic 890/H  
                                  Heinrich Faust  
                                  51145 Cologne, Germany
- Liquid Chromatograph:    HP 1100 Column Compartment G1316A  
                                  HP 1100 Binary Pump G1312A  
                                  HP 1100 Isocratic Pump G1310A  
                                  HP 1100 Degasser G1322A  
                                  Agilent  
                                  40880 Ratingen, Germany
- Autosampler:                HTC PAL System  
                                  CTC Analytics AG  
                                  4222 Zwingen, Switzerland
- Mass Spectrometer:         API 4000 with turbo-ionspray interface  
                                  mass selective detector (MS/MS)  
                                  Perkin Elmer Sciex Instruments  
                                  64331 Weiterstadt, Germany
- Note:                        Some mass spectrometric conditions are instrument specific.  
                                  The spectrometric conditions were optimised by a competent  
                                  operator prior to analysis.

#### 4.1.3 Reagents and Equipment

- Column (HPLC):             Luna 3u C8(2) 100A, size 50 mm x 2 mm  
                                  Part No. 00B4248-B0  
                                  Phenomenex  
                                  63741 Aschaffenburg, Germany
- Magnetic stirring bar:     plain (large, e.g. 40 x 8 mm [length x i.d.])
- Syringe filter:              Acrodisc® CR 13 mm syringe filter with 0.45 µm PTFE  
                                  membrane, parts no. 514-4008  
                                  VWR International GmbH  
                                  64301 Darmstadt, Germany

---

Acetonitrile:	for HPLC, super gradient grade Riedel de Haen, No. 34998 30926 Seelze, Germany
Ethyl acetate:	for residue analysis Promochem 46485 Wesel, Germany
Formic acid:	p.a. (100%) Merck, No.1.00264.0100 64271 Darmstadt, Germany
L-Cysteine hydrochlorid:	for biochemistry Merck, No.1.02839.1000 64271 Darmstadt, Germany
Water:	purified in a Milli-Q unit Milli-Pore GmbH 65731 Eschborn, Germany

Volumetric flasks, pipettes and other equipment commonly used in the laboratory.

#### **4.1.4 Chromatographic Conditions and Mass Spectrometric Parameters**

Liquid chromatographic conditions were identical to those described in Appendix 1 of the original method report FN-002-S05-02.

MS/MS parameter settings were optimized for the instrument being used and therefore not identical with those reported in method FN-002-S05-02.

#### **4.1.5 Calculation**

For calculation of the concentrations, calibration curves were used. These curves were calculated automatically after each sequence run with the Applied Biosystem quantitation software Analyst (vers. 1.4) using linear regression. Further calculations were performed using the software EXCEL 2002 (Office 2002®).

Matrix effects for BY108330 and its metabolites are eliminated by using an internal standard solution of the isotopically labelled reference items. Generally, the concentration of the internal standards should be in the range of the concentration of the analytes in the sample solutions.

The linear equation is expressed as:

$$y = \text{Intercept} + \text{Slope} \cdot x$$
$$y = \text{Area}, x = \text{Concentration}$$

If an internal standard is used:

$$y = \frac{\text{Area}_{\text{Standard}}}{\text{Area}_{\text{Internal Standard}}} = \text{Int. Ratio} \quad \text{and} \quad x = \frac{\text{Conc}_{\text{Standard}}}{\text{Conc}_{\text{IS}}} = \text{Conc}_{\text{Ratio}}$$

<i>Int. Ratio:</i>	<i>intensity ratio</i>
<i>Conc<sub>Standard</sub>:</i>	<i>concentration of standard solution [µg/L]</i>
<i>Conc<sub>IS</sub>:</i>	<i>concentration of internal standard solution [µg/L]</i>
<i>Conc<sub>Ratio</sub>:</i>	<i>concentration ratio</i>

Since the concentrations of the isotopically labelled internal standards were the same in all sample and standard solutions that were injected into the HPLC instrument, their concentrations can be neglected for calculations. However, the unit (µg/L or ng/mL) of the internal standard concentration has always to be considered in the calculation while the value is set to 1 (see example calculation). In cases where the concentration is taken into account in the formula for  $\text{Conc}_{\text{Ratio}}$ , it has to appear in the formula for  $\text{Conc}_{\text{Soil}}$ , too. By means of the linear equation, the compounds concentration in soil can be calculated as follows:

$$\text{Dilution}_{\text{Factor}} = \frac{\text{Volume}_{\text{Extraction}}}{\text{Weight}}$$

$$\text{Conc}_{\text{Analyte}} = \frac{\text{Int. Ratio} - \text{Intercept}}{\text{Slope}}, \quad \text{Int. Ratio} = \frac{\text{Area}_{\text{Analyte}}}{\text{Area}_{\text{Internal Standard}}}$$

$$\text{Conc}_{\text{Soil}} = \text{Conc}_{\text{Analyte}} \times \text{Dilution}_{\text{Factor}} \times \text{Conc}_{\text{IS}}$$

<i>Volume<sub>Extraction</sub>:</i>	<i>volume of the extraction solvent [L]</i>
<i>Weight:</i>	<i>weight of the soil sample [kg]</i>
<i>Intercept:</i>	<i>intercept of the linear regression curve</i>
<i>Slope:</i>	<i>slope of the linear regression curve</i>
<i>Area<sub>Analyte</sub>:</i>	<i>area of the analyte in the sample solution</i>
<i>Conc<sub>Soil</sub>:</i>	<i>concentration of the analyte in soil (sediment) [µg/kg]</i>

The recovery is calculated according to the following equation:

$$\text{Recovery} = \frac{\text{Conc}_{\text{Soil (Sediment)}} \times 100\%}{\text{Conc}_{\text{Soil (Sediment) Spiked}}}$$

<i>Conc<sub>Soil (Sediment) Spiked</sub>:</i>	<i>concentration of the reference item spiked [µg/kg]</i>
---	---