Evaluation of High-throughput Genotoxicity Screening Assays Used in Profiling the 320 US EPA ToxCast™ Chemicals.












A.W. Knight, L. Birrell and R.M. Walmsley. Gentronix Limited, CTF Building, 46 Grafton Street, Manchester, M13 9NT UK, andrew.knight@gentronix.co.uk

Introduction

A key aim of the U.S. Environmental Protection Agency's ToxCast project is to investigate state of the art high-throughput screening (HTS) assays to provide a biologically informed basis for predicting toxicity and prioritizing chemicals for further testing1. The assessment of genotoxicity is a key element of the ToxCast project and an important component in toxicity risk assessment. Three HTS genotoxicity assays were used to analyze the collection of 320 predominantly pesticide active compounds in Phase I of the ToxCast research project: GreenScreen HC GADD45a-GFP (Gentronix Ltd.)², CellCiphr p53 activation (Cellumen Inc.)³ and CellSensor p53RE-bla activation (Invitrogen Corp.)4, the latter performed by the National Institutes of Health Chemical Genomics Center (NCGC). These assays each invoke components of cellular responses to a genotoxic challenge, which can lead to DNA damage, mis-repair, mis-segregation, mutations and, ultimately, to tumorigenesis and carcinogenesis. The assays measure 2 gene expression endpoints, the induction of p53 and the induction of GADD45a in a p53 competent cell line. p53 is known to act as a 'gate keeper' to ensure genetic and cellular integrity during the cell cycle. GADD45a (growth arrest and DNA damage) mediates the cell's response to genotoxic stress.

The HTS results were compared with comprehensive tumorogenicity data from the US EPA ToxRefDB database and, where available, historical mutagenicity (Ames) data. The protocols and performance characteristics of the genotoxicity assays have also been assessed as to their suitability for deployment in HTS campaigns. The aim of the present study was to evaluate the utility of HTS assays to identify potential genotoxicity hazard in the larger context of the ToxCast project, and to aid prioritization of environmentally relevant chemicals for further testing. The purpose of these HTS assays is not to replace the use of Ames and other established genotoxicity tests, but to increase testing efficiency to enable earlier and more efficient screening of larger sets of chemicals of interest.

p53 and GADD45α Cellular Pathways



Complex regulation of GADD45a and p53.

GreenScreen HC

Positive Results - and overlap between HTS assavs

Assay	Genotoxicity		Cytotoxicity	
	Number	%	Number	%
GreenScreen HC	32	10.4	231	74.8
CellCiphr p53	27	8.7	171	55.3
CellSensor p53	36	11.7	-	ı

CellCiphr p53

CellSensor p53

The HTS - Genotoxicity Assays

gentronix

GreenScreen HC (GADD45a-GFP)

- TK6 human derived cell line with a GADD45a-GFP (green fluorescent protein) gene reporter. A second non-fluorescent strain is used to control for compound autofluorescence.
- Three serial dilutions of 200, 100 and 50 μM of each of 12 test compounds per 96-well microplate.
- Measurement at 24 and 48 hours by microplate spectrophotometer.
- Induction of cellular fluorescence is indicative of genotoxicity. Cytotoxicity is quantified by reduction in cellular proliferation, measured by optical absorbance.




CELLUMEN CellCiphr p53

- The p53 assay is one of 10 cellular response endpoints in HepG2 cells.
- Nine serial dilutions from 200 μM for each of 16 compounds per 384-well microplate.
- Measurement of differential staining by a High Content Imaging platform.
- Employs an anti-p53 antibody.



invitrogen CellSensor p53RE-bla

- ➤ HCT-116 cells with stably integrated beta-lactamase reporter gene and p53 response element.
- > 15 inter-plate dilutions from 92 μM for each of 1408 compounds over 15 1536-well microplates.
- Fluorescence resonance energy transfer (FRET) endpoint using a microplate spectrometer.



Replicate Samples – Randomly distributed in the compound library

Number of compounds giving positive

genotoxicity results in each HTS assay.



References and Links

- J.S. EPA lational Center for Computational Toxicology ToxCast Project. www.epagevircotto. [20 Kinglit AM, Breat LL, Wahnley BM, 2006 Development and validation to higher throughout screening approach to genotoxicity testing using the GADD45a-GIFP GreenScreen HC assay. 20 Journal of Biomisecular Screening 14, 16-30.
 Celloffort Cytotoxicity Profiling Assays. www.millipore.com/drugdscovery/dd3/cellofphr (4) Cell Sersor Cell Lines. www.infortgen.com

Kev Results and Observations

- > 9 12% of compounds were positive for genotoxicity in the HTS assays. The results of the three approaches only partially overlapped thus giving support to the proposed strategy of combining data from a battery of assays.
- The HTS assays were highly reproducible giving consistent results for replicate compounds.
- Overall the HTS assays demonstrated low sensitivity for rodent tumorigens likely due to; screening from a low, fixed concentration, coverage of selected genotoxic mechanisms, lack of metabolic activation in the assay protocols employed and difficulty in detecting non-genotoxic carcinogens.
- Conversely, HTS results demonstrated high specificity for rodent tumorigens (> 88%).
- The highest concordance with tumorogenicity data (>74% for HTS assays, c.f. 60% for Ames) was for compounds producing tumors in rodents at multiple sites and thus more likely genotoxic carcinogens.
- The highest 'positive predictivity' of multi-site rodent tumorigens was demonstrated by the GreenScreen HC assay (3 fold, vs. 1.41 - 1.65 for Ames and the HTS p53 assays).

Concordance with Multiple Site Rodent Tumorogenicity Data

	GreenScreen +	GreenScreen -		CellSensor +	CellSensor -	
Rodent +	13	45	Rodent +	10	48	
Rodent -	16	199	Rodent -	23	192	
	CellCiphr +	CellCiphr -		Ames +	Ames -	
Rodent +	8	50	Rodent +	14	14	
Rodent -	18	197	Rodent -	22	40	
		GreenScreen	CellSensor	CellClphr	Ames	
Number of com	parisons	273	273	273	90	
Sensitivity (% c	orrect positives)	22.4	17.2	13.8	50.0	
Specificity (% c	orrect negatives)	92.6	89.3	91.6	64.5	
Concordance		77.7	74.0	75.1	60.0	
Balanced Accur	acy	57.5	53.3	52.7	57.3	
Relative Predici	tivity (Positives)	3.01	1.61	1.65	1.41	
Relative Predict	ivity (Negatives)	1.19	1.08	1.06	1.29	

Comparison of Performance Characteristics

	GreenScreen HC	CellCiphr p53	CellSensor p53
Endpoints	Genotoxicity / Cytotoxicity	Genotoxicity / Cytotoxicity	Genotoxicity
Microplate	96 well	384 well	1536 well
Number of compounds / plate	12	16	1408
Number of dilutions / compound	3	10	15 (in 15 microplates)
Incubation time / hr	24 and 48	24 to 72	24
Compounds per week (typical)	720	100 - 500	30,000 - 100,000
Equipment required	Conventional microplate reader (Tecan Infinite F500)	Image analysis equipment (ArrayScan)	Microlitre liquid handling automation on ultra-high throughput robotic platforms
Data interpretation	Straightforward and automated - using proprietary Excel software template	Specialist image analysis software - 1 of 10 endpoints therefore data extraction required	Straightforward and automated - using proprietary software. Some user input to identify artifacts
Reagents	Proprietary cell lines and media from Gentronix	Proprietary cell lines and media from Millipore	Proprietary cell lines and media from Invitrogen

Concluding Remarks

Data from these HTS genotoxicity assays alone are not proposed as surrogates for the comprehensive regulatory genotoxicity assays or as accurate predictors of animal tumorogenicity. Rather data from the HTS assays can be applied in the following areas within the wider context of the aims of the ToxCast project:

- High-throughput screening for compounds capable of causing genetic damage in vitro and thus for 'hazard identification', i.e. highlighting chemicals with the potential to induce carcinogenicity.
- As part of a weight of evidence assessment of the likelihood of a compound's adverse effect for humans.
- To help determine the mode of action for carcinogenicity.
- To aid prioritisation of a compound for follow up in vitro and in vivo testing.