

**Validating Chlorinated Herbicides  
GC, SW-846, Method 8151A**



Prepared by: George Karras Date: 12/05/06  
George Karras, Chemist HWSS

Peer Reviewed by: Russell Arnone Date: 12/05/06  
Russell Arnone, Chemist HWSS

Concurred by: Linda Manuel Date: 12/5/06  
Linda Manuel Chief HWSS

Approved by: Robert Runyon Date: 12/11/06  
Robert Runyon, Chief HWSS

Annual Review

Reviewed by: \_\_\_\_\_ Date: \_\_\_\_\_

Reviewed by: \_\_\_\_\_ Date: \_\_\_\_\_

[illegible]

YES NO N/A

## 1.0 Traffic Reports and Laboratory Narrative

1.1 Are Traffic Report Forms present for all samples?         

ACTION: If no, contact lab for replacement of missing or illegible copies.

1.2 Do the Traffic Reports or SDG Narrative indicate any problems with sample receipt, condition of the samples, analytical problems or special circumstances affecting the quality of the data? \_\_\_\_ [ ] \_\_\_\_

ACTION: If any sample analyzed as a soil, other than TCLP, contains 50%-90% water, all data should be qualified as estimated (J). If a soil sample, other than TCLP, contains more than 90% water, all data should be qualified as unusable (R).

ACTION: If samples were not iced (4°C) upon receipt at the laboratory, flag all positive results "J" and all non-detects "UJ".

## 2.0 Holding Times

2.1 Has the technical holding times, determined from date of sample receipt to date of extraction, been exceeded? [ ]

Note: Samples may be analyzed for herbicide ester and acid. Check Laboratory SDG Narrative.

Note: Aqueous samples must be extracted within 7 days. Extracts must be analyzed within 40 days following extraction. Soil/Concentrated Waste samples must be extracted within 14 days and extracts analyzed within 40 days following extraction.

ACTION: If technical holding times are exceeded, flag all positive results and non-detects(U) as estimated ("J") and document in the narrative that holding times were exceeded. Samples extracted more than 28 days from sample receipt, either on the first analysis or upon re-analysis, flag all positive results as

YES NO N/A

### 3.0 Surrogate Recovery (Form II/Equivalent)

b. Soil [ ]

b. Soil/Concentrated Waste [ ]

Note: recommend surrogate is 2,4-Dichlorophenylacetic acid (DCAA)

YES NO N/A

[ ] \_\_\_\_\_

\_\_\_\_\_ [ ] \_\_\_\_\_

[ ]      \_\_\_\_\_

                          

[ ]

                          

Herbicides- 4 -

YES NO N/A

5.0                      Blanks (Form IV/Equivalent)

5.2 Frequency of Analysis: has a reagent/method blank been analyzed for each SDG or every 20 samples of similar matrix or concentration or each extraction batch, whichever is more frequent? [ ]

5.3 Has a Herbicide instrument blank been analyzed at the beginning of every analytical sequence of 10 samples ? [ ]

5.4 Chromatography: review the blank raw data - chromatograms, quant reports or data system printouts.

ACTION: Use professional judgement to determine the effect on the data.

## 6.0 Contamination

Herbicides- 5 -

YES NO N/A

- \_\_\_\_\_ [ ] \_\_\_\_\_

- \_\_\_\_\_ [ ] \_\_\_\_\_

ACTION: Follow the directions in the table below to qualify TCL results due to contamination. Use the largest value from all the associated blanks.

NOTE: If gross blank contamination exists, all data in the associated samples should be qualified as unusable (R).

- [ ]

## 7.1 Are the Gas Chromatograms and Data Systems

YES NO N/A

[ ]

[ ]

[ ]

[ ]

[ ]

[ ]

                          

[ ]

Herbicides- 7 -

YES NO N/A

7.7 Are the linearity criteria for the Initial Calibration analyses within limits for both columns? (% RSD must be < 20.0% for all analytes).

[ ]

7.8 Are there any transcription/calculation errors between raw data and Form VII - Herbicides-2?

\_\_\_\_\_ [ ] \_\_\_\_\_

7.9 Is the resolution between any two adjacent peaks in the QC Reference Check Mixture > 60.0% for both columns? (Form VI-Herbicides- 4)

                                

7.10 Is Form VII -Continuing Calibration present and complete for each analytical sequence for both columns?

                          

Herbicides- 8 -



YES NO N/A

7.11 Have all samples been injected within a 24 hr. period beginning with the injection of the first standard?

[ ]      \_\_\_\_\_      \_\_\_\_\_

7.12 Do all analyte retention times for the Mid-concentration Check standard (Form VII Herb-2) fall within the windows established by the initial calibration sequence? [

                                                                        

7.13 Are RPD values for all verification calibration standard compounds < 25.0%

                          

```
If %D is 25 -50% qualify as "J"
If %D is 51-100% qualify as "NJ"
If %D is >100% qualify as "R"
If %D is >100% with visible interferences/qualify as "JN"
```

8.1 Is Form VIII present and complete for each column

YES NO N/A

[ ] \_\_\_\_\_

[ ]

[ ]

\_\_\_\_\_ [ ] \_\_\_\_\_

                          

[ ]

Herbicides- 10 -

YES NO N/A

[ ]

| % Difference | Qualifier |
|--------------|-----------|
|--------------|-----------|

NOTE: The lower of the two values is reported on Form I. If using professional judgement, the reviewer determines that the higher result was more acceptable, the reviewer should replace the value and indicate the reason for the change in the data assessment.

[ 1 ]

ACTION: Use professional judgement to decide if the compound should be reported.

## 10.0 Compound Quantitation and Reported Detection Limits

[ ]

NOTE: The reviewer should use professional judgement to decide whether a much larger concentration obtained on one column versus the other indicates the presence of an interfering compound. If an interfering compound is indicated, the lower of the two values should be reported and qualified as presumptively present at an approximated quantity (NJ). This necessitates a determination of an estimated concentration on the confirmation column. The narrative should indicate the presence of interferences during the evaluation of the second column confirmation.

YES NO N/A

10.2 Are the CRQLs adjusted to reflect sample dilutions and, for soils, % moisture?

[ ]      \_\_\_\_\_

ACTION: If errors are large, call lab for explanation/resubmittal, make any necessary corrections and document effect in data assessments.

ACTION: When a sample is analyzed at more than one dilution, the lowest CRQLs are used (unless a QC exceedance dictates the use of the higher CRQL data from the diluted sample analysis). Replace concentrations that exceed the calibration range in the original analysis by crossing out the "E" value on the original Form I and substituting it with data from the analysis of diluted sample. Specify which Form I is to be used, then draw a red "X" across the entire page of all Form I's that should not be used, including any in the summary package.

**ACTION:** Quantitation limits affected by large, off-scale peaks should be qualified as unusable (R). If the interference is on-scale, the reviewer can provide an approximated quantitation limit (UJ) for each affected compound.

10.3 Have all data (Forms and associated chromatograms and quantitation reports) been submitted for original, diluted or re-extraction/re-analysis samples? [

                                

## 11.0 Chromatogram Quality

### 11.1 Were baselines stable?

\_\_\_\_\_

11.2 Were any electropositive displacement (negative peaks) or unusual peaks seen?

\_\_\_\_\_ [ ] \_\_\_\_\_

ACTION: Address comments under System Performance of data assessment. Explain use of professional judgement where used to qualify data.

12.0 Field Duplicates

12.1 Were any field duplicates submitted for  
Herbicides analysis? ☐ ☐ ☐

Note: Check whether SAS Client Request required  
field duplicates.

ACTION: Compare the reported results for  
field duplicates and calculate the  
relative percent difference.

ACTION: Any gross variation between field  
duplicate results must be addressed  
in the reviewer narrative. However, if  
large differences exist, identification  
of field duplicates should be confirmed  
by contacting the sampler.