

SOP HW-36
Revision 0
April 2006

SOP NO. HW-36/Pesticide Data Validation
USEPA Contract Laboratory Program
Statement of Work for Organic Analysis of Low/Medium
Concentration of Pesticide Organic Compounds SOM01.1



Prepared by: George Karas
George Karas, Chemist
Hazardous Waste Support Section

Date: 12/05/06

Peer Reviewed by: Russell Arnone
Russell Arnone, Chemist
Hazardous Waste Support Section

Date: 12/05/06

Concurred by: Linda Mauel
Linda Mauel, Chief
Hazardous Waste Support Section

Date: 12/5/06

Approved by: Robert Runyon
Robert Runyon, Chief
Hazardous Waste Support Branch

Date: 12/11/06

Annual Review

Reviewed by: _____
Name

Date: _____

Reviewed by: _____
Name

Date: _____

TABLE OF CONTENTS

INTRODUCTION	1
Scope and Applicability	1
Summary	1
Data Qualifiers	1
Lab Qualifiers	2
Reviewer Qualifications	2
PACKAGE COMPLETENESS AND DELIVERABLES	3
1. <u>Chain of Custody and Sampling Trip Reports</u>	3
2. <u>Data Completeness and Deliverables</u>	3
3. <u>Cover Letter SDG Narrative</u>	4
4. <u>Data Validation Checklist</u>	5
PART A: VOA ANALYSES	5
1. <u>Sample Conditions/Problems</u>	5
2. <u>Holding Times</u>	5
3. <u>Deuterated Monitoring Compound (DMC) Recovery (Form II)</u>	6
4. <u>Matrix Spike/Matrix Spike Duplicate Recovery (Form III)</u>	8
5. <u>Method Blanks (Form IV)</u>	9
6. <u>Contamination</u>	10
7. <u>GC/ECD Instrument Performance Check (Form V)</u>	11
8. <u>Analytical Sequence Check (Form VIII)</u>	16
9. <u>Florisil Cartridge and GPC Checks (Form IX)</u>	17
10. <u>Pesticide Identification (Form X)</u>	18
11. <u>Target Pesticide List (TCL)</u>	20
12. <u>Compound Quantitation and Reported Detection Limits</u>	20
13. <u>Field Duplicates</u>	21
Definitions	22
References	23

INTRODUCTION

Scope and Applicability

This SOP offers detailed guidance in evaluating laboratory data generated according to the method in the "USEPA Contract Laboratory Program Statement of Work for Organics Analysis Multi-Media, Multi-Concentration, SOM01.1, May 2005". The validation procedures and actions discussed in this document are based on the requirements set forth in the "USEPA Contract Laboratory Program National Functional Guidelines for Superfund Organic Methods Data Review, January 2005". This document attempts to cover technical problems specific to low/Medium concentration of Pesticide compounds. Situations may arise where data limitations must be assessed based on the reviewer's own professional judgement.

In addition to technical requirements, contractual requirements may also be covered in this document. While it is important that instances of contract non-compliance be addressed in the Data Assessment, the technical criteria are always used to qualify the analytical data.

Summary

To ensure a thorough evaluation of each result in a data case, the reviewer must complete the checklist within this SOP, answering specific questions while performing the prescribed "ACTIONS" in each section. Qualifiers (or flags) are applied to questionable or unusable results as instructed. The data qualifiers discussed in this document are as follows:

Data Qualifiers

- U - The analyte was analyzed for, but was not detected above the reported sample quantitation limit.
- J - The analyte was positively identified; the associated numerical value is the approximate concentration of the analyte in the sample.
- N - The analysis indicates the presence of an analyte for which there is presumptive evidence to make a "tentative identification."
- JN - The analysis indicates the presence of an analyte that has been "tentatively identified" and the associated numerical value represents its approximate concentration.

- UJ - The analyte was not detected above the reported sample quantitation limit. However, the reported quantitation limit is approximate and may or may not represent the actual limit of quantitation necessary to accurately and precisely measure the analyte in the sample.
- R - The sample results are rejected due to serious deficiencies in the ability to analyze the sample and meet quality control criteria. The presence or absence of the analyte cannot be verified.

Lab Qualifiers:

- D - The positive value is the result of an analysis at a secondary dilution factor.
- B - The analyte is present in the associated method blank as well as in the sample. This qualifier has a different meaning when validating inorganic data.
- E - The concentration of this analyte exceeds the calibration range of the instrument.
- P - Pesticide/Aroclor target analytes when the % Difference between the analyte concentrations obtained from the two dissimilar GC columns is greater than 25%.

The reviewer must prepare a detailed data assessment to be submitted along with the completed SOP checklist. The Data Assessment must list all data qualifications, reasons for qualifications, instances of missing data and contract non-compliance.

Reviewer Qualifications:

Data reviewers must possess a working knowledge of the USEPA Statement of Work SOM01.1 and National Functional Guidelines mentioned above.

```

USEPA Region II                                     Date: April 2006
Method: CLP/SOW, SOM01.1/Pesticide                SOP HW-36/Pesticide, Revision 0
S))))))))))))))))))))))))))))))))))))))))))))))))))))))))))))))))))))))))))))))))
      . . . . .                                     YES    NO    N/A

```

CASE NUMBER: _____ LAB: _____

SITE NAME: _____ **SDG No(s) . :** _____

1.1 Are the Traffic Reports/Chain-of-Custody Records present for all samples?

1.2 Is the Sampling Trip Report present for all samples? [] _____

2.1 Have any missing deliverables been received and added to the data package? _____ ☐ _____

2.2 Was SMO/CLASS CCS checklist included with the package? []

7

STANDARD OPERATING PROCEDURE

USEPA Region II

Date: April 2006

Method: CLP/SOW, SOM01.1/Pesticide

SOP HW-36/Pesticide, Revision 0

S))))))

YES NO N/A

• • • • •

- | | | | |
|--|-----------|-----|-----|
| a.) Form VI Pest-5/Pesticide Resolution Check Mix | <u> </u> | ___ | ___ |
| b.) Form VI Pest-6/Performance Evaluation Mixture | <u> </u> | ___ | ___ |
| c.) Form VI Pest-7/Individual Standard Mixture A | <u> </u> | ___ | ___ |
| d.) Form VI Pest-8/Individual Standard Mixture B | <u> </u> | ___ | ___ |
| e.) Form VI Pest-9/Individual Standard Mixture C | <u> </u> | ___ | ___ |
| f.) Form VI Pest-10/Individual Standard Mixture C | <u> </u> | ___ | ___ |
| g.) Form VII Pest-1/Calibration Verification | <u> </u> | ___ | ___ |
| h.) Were the appropriate GC columns used as specified on
page D-11/Pest, sections 6.26.1.3 to 6.26.1.3.2 in
SOM01.1? | <u> </u> | ___ | ___ |

7.2 The identification of a single component pesticide by GC method is based primarily on RT data. Were the following requirements met:

- a.) The chromatogram that results for PEM and Individual Standards Mixture analyses must display the analytes at > 10% full scale but < 100% full scale
- b.) The baseline of the chromatogram must return to below 50% of full scale before the elution of alpha-BHC, and return to below 25% of full scale after the elution time of alpha-BHC and before the elution time of decachlorobiphenyl

NOTE: If a chromatogram is replotted electronically to meet these requirements, the scaling factor used must be displayed on the chromatogram, and if standard, blank, etc chromatogram needs to be replotted electronically to meet these requirements, both the initial chromatogram and the replotted chromatogram(s) must be submitted in the data package.

ACTION: If all single component pesticides (SCP) are not clearly displayed on chromatograms for all Individual Standard Mixtures and PEM, notify the TOPO to obtain resubmittal of the necessary data.

7.3 Are there any transcription/calculation errors between raw data and the Forms? []

ACTION: If large errors exist, take action specified in section 3.1 above.

7.4 Resolution Check Mixture (Form VI Pest-5)

STANDARD OPERATING PROCEDURE

USEPA Region II

Date: April 2006

Method: CLP/SOW, SOM01.1/Pesticide

SOP HW-36/Pesticide, Revision 0

S))))))

YES NO N/A

Criteria [(Individual Standard Mixture (A and B))]		Criteria (Individual Standard Mixture C)	Action
Resolution Check Mixture % Resolution <60.0%	Resolution Check Mixture % Resolution <80.0% (primary column) % Resolution <50.0% (secondary column)		Detects: JN Non-detects: R
PEM % Resolution <90.0%			Detects: JN Non-detects : R
PEM: 4,4'-DDT % Breakdown >20.0% and 4,4'-DDT is detected			Detects for 4,4'-DDT: J Detects for 4,4'-DDD: J Detects for 4,4'-DDE: J
PEM: 4,4'-DDT % Breakdown >20.0% and 4,4'-DDT is not detected			Non-detects for 4,4'-DDT: R Detects for 4,4'-DDD: JN Detects for 4,4'-DDE: JN
PEM: Endrin % Breakdown >20.0% and Endrin is detected			Detects for Endrin: J Detects for Endrin aldehyde: J Detects for Endrin ketone: J
PEM: Endrin % Breakdown >20.0% and Endrin is not detected			Detects for Endrin: R Detects for Endrin aldehyde: JN Detects for Endrin ketone: JN
PEM: Combined % Breakdown > 30.0%			Apply qualifiers as described above considering degree of individual breakdown
Mid-point Individual Standard Mixtures (A and B) % Resolution <90.0%	Mid-point Individual Standard Mixture (C) % Resolution <80.0% (primary column) Mid-point Individual Standard Mixture (C) % Resolution <50.0% (secondary column)		Detects: JN Non-detects: R
PEM analysis not performed at the required frequency *			All results: R
Mid-point Individual Standard Mixtures analysis not performed at the required frequency **			All results: R

* The PEM is analyzed at the beginning (following the Resolution Check Mixture) and at the end of the initial calibration.

** Mid-point Individual Standard Mixture A and B: Analyzed as part of the initial calibration. The mid-point INDA and INDB must bracket one end of each 12-hour analytical period.

Mid-point Individual Standard Mixture C: Analyzed as part of the initial calibration. The mid-point INDC must bracket one end of each 12-hour analytical period.

