

S))q  
YES NO N/A

This Region II SOP document is based on SW846  
Method 8150A, Revision I, July,1992

1.0 Traffic Reports and Laboratory Narrative

- 1.1 Are Traffic Report Forms present for all ☐ ☐ ☐  
samples?

ACTION: If no, contact lab for replacement of  
missing or illegible copies.

- 1.2 Do the Traffic Reports or SDG Narrative indicate  
any problems with sample receipt, condition of  
the samples, analytical problems or special  
circumstances affecting the quality of the data? ☐ ☐ ☐

ACTION: If any sample analyzed as a soil, other  
than TCLP, contains 50%-90% water,  
all data should be qualified as estimated  
(J). If a soil sample, other than TCLP,  
contains more than 90% water, all data  
should be qualified as unusable (R).

ACTION: If samples were not iced upon receipt at  
the laboratory, flag all positive results  
"J" and all non-detects "UJ".

2.0 Holding Times

- 2.1 Have any technical holding times,determined from  
date of collection to date of extraction,  
been exceeded? ☐ ☐ ☐

Note: Water and soil samples for herbicide analysis  
must be extracted within 7 days of the date of  
collection. Extracts must be analyzed within 40  
days of the date extraction. However, the SAS Client  
Request takes precedence and the Holding Times  
specified in the SAS are the criteria used for

**S))Q**

Revision: 1.3

Herbicides- 2 -



Date: Nov., 1994  
Revision: 1.3

YES NO N/A

\_\_\_\_\_ [ ] \_\_\_\_\_

S))q  
YES NO N/A

ACTION: No qualification is done if the surrogate is diluted out. If recovery for the surrogate is below the contract limit, but above 10%, flag all results for that sample 'J'. If recovery is < 10%, qualify positive results 'J' and flag non-detects 'R'. If recovery is above the contract advisory limit qualify positive values 'J'.

3.5 Were surrogate retention times (RT) within the windows established during the initial 5-point calibration analysis? (see Form VI Herb-1) ☐ ☐ ☐

ACTION: If the RT limits are not met, the analysis may be qualified unusable (R) for that sample on the basis of professional judgement.

3.6 Are there any transcription/calculation errors between raw data and Form II? ☐ ☐ ☐

ACTION: If large errors exist, call lab for explanation/resubmittal. Make any necessary corrections and document effect in data assessments.

4.0 Matrix Spikes (Form III)

4.1 Is the Matrix Spike/Matrix Spike Duplicate Recovery Form (Form III) present? ☐ ☐ ☐

4.2 Were matrix spikes analyzed at the required frequency for each of the following matrices? (1 MS/MSD must be performed for every 20 samples of similar matrix or concentration level)

a. Low Water ☐ ☐ ☐

b. Soil ☐ ☐ ☐

Date: Nov., 1994  
Revision: 1.3

YES NO N/A

ACTION: If any matrix spike data are missing, take the action specified in 3.2 above.

**S)))))))))))))**

Revision: 1.3

YES NO N/A

- Soil

out of

- Soil

out of

5.0 Blanks (Form IV)

- 5.2 Frequency of Analysis: has a reagent/  
method blank been analyzed for each SDG or  
every 20 samples of similar matrix  
or concentration or each extraction batch,  
whichever is more frequent? [ ] [ ] [ ]

Herbicides- 7 -

Date: Nov., 1994  
Revision: 1.3

YES NO N/A

5.3 Has a Herbicide instrument blank been analyzed at the beginning of every analytical sequence of 10 samples ?



**S))Q**

Revision: 1.3

ACTION: If any blank data are missing, call lab for explanation/resubmittals. If missing deliverables are unavailable, document the effect in data assessments.

- 

\_\_\_\_\_ [ ] \_\_\_\_\_

\_\_\_\_\_ [ ] \_\_\_\_\_

Herbicides- 9 -

Date: Nov., 1994  
Revision: 1.3

YES NO N/A

day) may be used to qualify data. Blanks may not be qualified because of contamination in another blank. Field blanks must be qualified for surrogate, calibration, or any QC problems.

**S))Q**

Date: Nov., 1994

Revision: 1.3

YES NO N/A

ACTION: Follow the directions in the table below to qualify TCL results due to contamination. Use the largest value from all the associated blanks.

		Sample
conc > CRQL but < 5x blank	Sample conc < CRQL & is < 5x blank value	Sample conc > CRQL & > 5x blank value
Flag sample result with a "U";	Report CRQL & qualify "U"	No qualification is needed

NOTE: If gross blank contamination exists, all data in the associated samples should be qualified as unusable (R).

6.3 Are there field/rinse/equipment blanks associated with every sample?

ACTION: For low level samples, note in data assessment that there is no associated field/rinse/equipment blank. Exception: samples taken from a drinking water tap do not have associated field blanks.

## 7.0 Calibration and GC Performance

7.1 Are the Gas Chromatograms and Data Systems Printouts for both columns present for all samples, blanks, QC Check reference, and MS/MSD?

ACTION: If no, take action specified in 3.2 above.

7.2 Are Forms VI - Herbicides 1,2,4 present and complete for each column and each analytical sequence? [

ACTION: If no, take action specified in 3.2 above.

7.3 Are there any transcription/calculation errors

[illegible]

Revision: 1.3

YES NO N/A

\_\_\_\_\_ [ ] \_\_\_\_\_

**S)))))))))))))**

Revision: 1.3

YES NO N/A

7.4 Were the retention time windows calculated using the average absolute retention time (at least three measurements)  $\pm$  three times the standard deviation of the absolute retention time, for each standard? (Refer to Method 8000A, section 7.5). [ ]

7.5.2 If yes, was the surrogate recovery >50%?   1        

7.5.3 Was the QC check standard re-extracted/re-analyzed, if surrogate recovery was <50%, or any one analyte was < 40%, or two analytes < 70% ? [ ]

Action: If NO to any of the above, then qualify positive hits as estimated "J" and non-detects as rejected "R" in the original analysis of all samples in the associated analytical sequence.

7.6 Do all standard retention times, including each Herbicides in each level of Initial Calibration fall within the windows established during the initial calibration analytical sequence? (For Initial Calibration Standards, Form VI - Herbicides - 1).

ACTION: If no, all samples in the entire analytical sequence are potentially affected. Check to see if the chromatograms contain peaks within an expanded window surrounding the expected retention times. If no peaks are found and the surrogate is visible, non-detects are valid. If peaks are present

Date: Nov., 1994  
Revision: 1.3

YES NO N/A

Herbicides- 14 -

**S))Q**

Revision: 1.3

YES NO N/A

\_\_\_\_\_

[ ]

[ ]

\_\_\_\_\_

Herbicides- 15 -

Date: Nov., 1994  
Revision: 1.3

YES NO N/A



**S)))))))))))))**

Revision: 1.3

YES NO N/A

Herbicides- 17 -

Date: Nov., 1994  
Revision: 1.3

YES NO N/A

                          

ACTION: If no, take action specified in 3.2 above.



Date: Nov., 1994  
Revision: 1.3

YES NO N/A

RT shift. The reviewer should use professional judgement to assign an appropriate quantitation limit.

Date: Nov., 1994  
Revision: 1.3

YES NO N/A

ACTION: If the reviewer finds neither column shows interference for the positive hits, the data should be flagged as follows:

25-50 %	J
50-90 %	JN
> 90 %	R

9.5 Check chromatograms for false negatives.  
Were there any false negatives? [ ]

## 10.0 Compound Quantitation and Reported Detection Limits

10.1 Are there any transcription/calculation errors in Form I results? Check at least two positive values. Were any errors found? \_\_\_\_\_ [ ] \_\_\_\_\_

Herbicides- 21 -

**S**))**Q**

YES NO N/A

Revision: 1.3

YES NO N/A

S))q  
YES NO N/A

10.2 Are the CRQLs adjusted to reflect sample dilutions  
and, for soils, % moisture? ☐ ☐ ☐

ACTION: If errors are large, call lab for  
explanation/resubmittal, make any  
necessary corrections and document  
effect in data assessments.

ACTION: When a sample is analyzed at more than  
one dilution, the lowest CRQLs are used  
(unless a QC exceedance dictates the use  
of the higher CRQL data from the diluted  
sample analysis). Replace concentrations  
that exceed the calibration range in the  
original analysis by crossing out the "E"  
value on the original Form I and substituting  
it with data from the analysis of diluted  
sample. Specify which Form I is to be used,  
then draw a red "X" across the entire page  
of all Form I's that should not be used,  
including any in the summary package.

ACTION: Quantitation limits affected by large,  
off-scale peaks should be qualified as  
unusable (R). If the interference is  
on-scale, the reviewer can provide an  
approximated quantitation limit (UJ) for  
each affected compound.

10.3 Have all data (Forms and associated chromatograms and  
quantitation reports) been submitted for original,  
diluted or re-extraction/re-analysis samples? ☐ ☐ ☐

#### 11.0 Chromatogram Quality

11.1 Were baselines stable? ☐ ☐ ☐

11.2 Were any electropositive displacement  
(negative peaks) or unusual peaks seen? ☐ ☐ ☐

CHLORINATED HERBICIDES  
STANDARD OPERATING PROCEDURE  
U.S.EPA Region II

Date: Nov., 1994  
Revision: 1.3

**S))Q**  
YES NO N/A

ACTION: Address comments under System  
Performance of data assessment.  
Explain use of professional judgement  
where used to qualify data.



Date: Nov., 1994  
Revision: 1.3

**S**)))))))))))))  
YES NO N/A

## 12.0 Field Duplicates

12.1 Were any field duplicates submitted for Herbicides analysis?

                                

Note: Check whether SAS Client Request required field duplicates.

ACTION: Compare the reported results for field duplicates and calculate the relative percent difference.

ACTION: Any gross variation between field duplicate results must be addressed in the reviewer narrative. However, if large differences exist, identification of field duplicates should be confirmed by contacting the sampler.