









EXPLOSIVE RESIDUES  
STANDARD OPERATING PROCEDURE  
U.S.EPA Region II

Date: Sept., 1994  
Revision: 1.3

S))Q  
YES NO N/A

Holding Time Violations  
Table 1

(See Traffic Report)				
Sample ID	Date Sampled	Date Received	Date Lab Extracted	Date Analyzed
_____	_____	_____	_____	
_____	_____	_____	_____	
_____	_____	_____	_____	_____
_____	_____	_____	_____	_____
_____	_____	_____	_____	_____

ACTION: If technical holding times are exceeded, flag all positive results as estimated ("J") and sample quantitation limits as estimated ("UJ"), and document in the narrative that holding times were exceeded. The reviewer must use professional judgement to determine the reliability of the data and the effects of additional storage on the sample results. At a minimum, all results must be qualified "J", but the reviewer may determine that non-detect data are unusable (R).

## 8.0 Surrogate Recovery

8.1 Are the Organic Analysis Data Sheets (Form I) present and complete with surrogate recoveries



**S))Q**

Revision: 1.3

YES NO N/A

\_\_\_\_\_ [ ] \_\_\_\_\_

[ ]

**ACTION:** If the retention time shift (RTS) or % recovery of the surrogate in the field or QC samples is out of specification, then all associated sample data should be qualified. Use **Table 2** below as a guide.

## TABLE 2

SURROGATE	<u>NOT QUALIFIED</u>	<u>J</u>	<u>R</u>	<u>N</u>
% RECOVERY - FIELD SAMPLES				
Detects	50 - 125%	< 50%; > 125%		
Non-detects	\$ 10%	< 10%*	< 10%*	
% RECOVERY - BLANKS AND QC SAMPLES **				
Detects	50 - 125%	< 50%; > 125%		
Non-detects	\$ 20%	10 - 19%	< 10%	
RTS - field samples	± 1.0%		> 1.5%; < -1.5% *	± 1.1 - 1.5% *
RTS - QC samples	± 1.0%		> 1.5%; < -1.5% *	± 1.1 - 1.5% *

Explosive - 7 of 18 - Residues

Date: Sept., 1994  
Revision: 1.3

YES NO N/A

If the surrogate recovery in a QC sample or blank was less than 50 percent or greater than 125 percent, then the field sample data should be qualified if the surrogate in the field sample is outside 50 - 125 percent recovery.



Date: Sept., 1994  
Revision: 1.3

YES NO N/A

ACTION: If large errors exist, call lab for explanation/  
resubmittal, make necessary corrections and note  
errors in the data assessment.

9.1 Is the QC Check Reference Sample Recovery Form present? (Created by Lab) [ ]

NOTE: The QC Check Form has to be created by the Lab. Concentration of the spiking solution is 5-10 times the estimated quantitation limits for all analytes listed in section 12 of the SAS.

NOTE: QC Check Reference Sample information is important,  
data is used to judge extraction efficiency.

ACTION: If any QC Check data are missing, call the lab for explanation/resubmittal. If the lab cannot provide missing deliverables, document the effect on the validity of the data in the data assessment. Positive hits should be flagged "J" and non-detects should be flagged "R".

9.3 How many QC check recoveries are outside QC limits of 70-125%?

<u>Water</u>	<u>Soil</u>
out of	out of

**S))Q**

Revision: 1.3

ACTION: If any recovery is greater than 125%, positive results should be flagged "J" for the affected compound.

**S))Q**

Revision: 1.3

When two or more analytes have recoveries above 10% but below 70% all associated data should be flagged "J" (positive and non-detects).

Note in data assessment. If further assistance is required, contact SMO for instructions.

It should be noted for TPO action, if a laboratory fails to analyze a QC Check Reference Sample or if a laboratory consistently fails to generate acceptable recoveries.

## 10.1 Initial Calibration

The initial calibration standards must be analyzed daily prior to any sample analysis. The lab may have to create two sets of initial standard solutions if several analytes coelute.

The initial calibration curve must be injected in triplicate for each of the 5 levels.

10.2 Are the chromatograms and data system printouts (Quant reports) present for initial and triplicate calibrations?

\_\_\_\_\_

10.3 Are the modified Initial Calibration Forms (Pest VI-1, VI-2 Forms) present?

[ ]

10.4 Are modified Forms present and completed for each analytical sequence and each column?

\_\_\_\_\_

**S)))))))))))))**

Revision: 1.3

ACTION: If any Calibration Standard chromatograms and forms are missing, contact the lab for resubmittals.

\_\_\_\_\_ [ ] \_\_\_\_\_



**S))Q**

Revision: 1.3

**S))Q**

ACTION: If degradation is detected or suspected, qualify positive hits for tetryl and TNT as "J". If degradation is suspected, and tetryl is reported as non-detected, qualify the analyte as "UJ". Note in data assessment.

11.1 Are Forms (Modified Pest VII-2) present for each continuing calibration?

\_\_\_\_\_

11.2 Are there any transcription/calculation errors between raw data and forms?

\_\_\_\_\_ [ ] \_\_\_\_\_

ACTION: If large errors exist, call lab for explanation/resubmittals. Make any necessary corrections, and document effect in data assessment. Check at least 2-3 values from each calibration standard and more if errors are found.

11.3 Were continuing calibrations analyzed at the required frequency? (standard should bracket every 10 samples).

\_\_\_\_\_

ACTION: Criteria must be met. Determine effect on data. At a minimum, all data should be qualified "J", but the reviewer may determine that data are unusable "R". Document the data qualifications in data assessment.

Explosive - 14 of 18 - Residues

Date: Sept., 1994  
Revision: 1.3

YES NO N/A

                                

                                

[ ]      \_\_\_\_\_      \_\_\_\_\_

Explosive - 15 of 18 - Residues

**S))Q**

Revision: 1.3

for each initial calibration and subsequent analyses?

[ ]

ACTION: Flag "J" all data generated outside an acceptable twenty-four hour sequence starting after the calculation of the RT windows resulting from the initial calibration unless the daily calibration meets the calibration criteria for all target analytes. If the daily calibration standard does not meet the criteria for some analytes the reviewer might reject

## 13.0 Method Blanks

a. water

[ ]

b. soil

[ ]

Has a Method/Prep blank been analyzed at the same time for each batch of samples extracted.

[ ]



**S))Q**

Revision: 1.3

ACTION: If any method blank data are missing,  
call lab for explanation/resubmittal.  
If method blank data are not available,  
reject "R" all associated positive data.

[ ]

13.4 Chromatography: review the method blank raw data chromatograms, quant reports or data system printouts.

[ ] \_\_\_\_\_

Explosive - 17 of 18 - Residues

**S))Q**

Revision: 1.3

YES NO N/A

## 14.0 Contamination

NOTE: "Water blanks", "drill blanks", and distilled water blanks" are validated like any other sample, and are not used to qualify data. Do not confuse them with the other OC blanks discussed below.

[ ]

ACTION: Sample analysis results after the high concentration sample must be evaluated for carryover. Sample cross-contamination should be noted for TPO action if an effect on the data is suspected. An Instrument Blank is not required in the methodology.

[ ]

[ ] \_\_\_\_\_

ACTION: Note in data assessment that there is no associated field/rinse/equipment blank.

[ ]

ACTION: Prepare a list of the samples associated with each of the contaminated blanks.  
(Attach a separate sheet.)

NOTE: All field blank results associated with a

**S))Q**

Revision: 1.3

particular group of samples (may exceed one per case) must be used to qualify data. Blanks may not be qualified because of contamination in another blank. Field blanks must be qualified for surrogate recovery, instrument performance criteria, or calibration QC problems.

**S))Q**

Revision: 1.3

YES NO N/A

ACTION: Follow the directions in **Table 3** below to qualify TCL results due to contamination. Use the largest value from all the associated blanks. If any blanks are grossly contaminated, all associated data should be qualified as unusable (R) (see item # 12.3).

	should be qualified as unusable (R) (see Item # 12.3).	
Sample conc > CRQL but < 5x blank	Sample conc < CRQL & is < 5x blank value	Sample conc > CRQL & > 5x blank value
Flag sample result with a "U"	Report CRQL & qualify "U"	No qualification is needed

NOTE: Analytes qualified "U" for blank contamination are still considered as "hits" when qualifying for calibration or other QC criteria.

15.1 Are the Organic Analysis Data Sheets (Form I) present with required header information on each page, for each of the following:

- |    |   |                          |     |     |
|----|---|--------------------------|-----|-----|
| a. | Samples ?                                     | <input type="checkbox"/> | ___ | ___ |
| b. | Laboratory Control Samples?                   | <input type="checkbox"/> | ___ | ___ |
| c. | Blanks?                                       | <input type="checkbox"/> | ___ | ___ |
| d. | Matrix spikes and matrix spike<br>duplicates? | <input type="checkbox"/> | ___ | ___ |
| e. | Lab duplicate?                                | <input type="checkbox"/> | ___ | ___ |

Revision: 1.3

**S))Q**

YES NO N/A

**S))Q**

Revision: 1.3

YES NO N/A

a. Samples?

[ ]      [ ]      [ ]

b. Laboratory Control Samples?

\_\_\_\_\_

c.     Blanks?

[ ]      [ ]      [ ]

d. Matrix spikes and Matrix spike duplicates?

[ ]      [ ]      [ ]

e. Lab duplicate?

[ ]

Lab for resubmittals.

15.3 Are the response factors shown in the Quant Report?

\_\_\_\_\_

15.4 Is chromatographic performance acceptable with respect to:

Baseline stability?

\_\_\_\_\_

Resolution?

[ ]      [ ]      [ ]

Peak shape?

[ ]

Full-scale graph (attenuation)?

[ ]

Other:

[ ]

the acceptability of the data.

**S))Q**

Revision: 1.3

                                          

16.1 Is Form X completed for every sample in which an analyte was detected?

Explosive - 23 of 18 - Residues

**S))Q**

Revision: 1.3

**S))Q**

YES NO N/A

[ ]

[ ]

[ ]

<u>% Difference</u>	<u>Qualifier</u>
25-90%	J
>90%	R

Explosive - 24 of 18 - Residues



**S))Q**

Revision: 1.3

result was more acceptable, the reviewer should replace the value and indicate the reason for the change in the data assessment.

**S))Q**

Revision: 1.3

**S))Q**

YES NO N/A

ACTION: Use professional judgement to determine qualification of analytes.

17.1 Are there any transcription/calculation errors in Form I results? Check at least two positive values. Verify that the calculations were adjusted for percent moisture. Were any errors found?

\_\_\_\_\_ [ ] \_\_\_\_\_

Explosive - 26 of 18 - Residues



Date: Sept., 1994  
Revision: 1.3

YES NO N/A

ACTION: Flag all associated data with an \* for "out of control" duplicate.

## 19.0 Matrix Spikes



**S))Q**

Revision: 1.3

may use the matrix spike data results in conjunction with other QC criteria and determine the need for some qualification of the data.