

## METHOD 7091

### BERYLLIUM (ATOMIC ABSORPTION, FURNACE TECHNIQUE)

#### 1.0 SCOPE AND APPLICATION

1.1 See Section 1.0 of Method 7000.

#### 2.0 SUMMARY OF METHOD

2.1 See Section 2.0 of Method 7000.

#### 3.0 INTERFERENCES

3.1 See Section 3.0 of Method 7000 if interferences are suspected.

3.2 The long residence time and high concentrations of the atomized sample in the optical path of the graphite furnace can result in severe physical and chemical interferences. Furnace parameters must be optimized to minimize these effects.

3.3 In addition to the normal interferences experienced during graphite furnace analysis, beryllium analysis can suffer from severe nonspecific absorption and light scattering caused by matrix components during atomization. Simultaneous background correction is required to avoid erroneously high results.

#### 4.0 APPARATUS AND MATERIALS

4.1 For basic apparatus, see Section 4.0 of Method 7000.

4.2 Instrument parameters (general):

4.2.1 **Drying time and temp:** 30 sec at 125°C.

4.2.2 **Ashing time and temp:** 30 sec at 1000°C.

4.2.3 **Atomizing time and temp:** 10 sec at 2800°C.

4.2.4 **Purge gas:** Argon.

4.2.5 **Wavelength:** 234.9 nm.

4.2.6 **Background correction:** Required.

4.2.7 Other operating parameters should be set as specified by the particular instrument manufacturer.

NOTE: The above concentration values and instrument conditions are for a Perkin-Elmer HGA-2100, based on the use of a 20- $\mu$ L injection, continuous-flow purge gas, and nonpyrolytic graphite. Smaller sizes of furnace devices or those employing faster rates of atomization can be operated using lower atomization temperatures for shorter time periods than the above-recommended settings.

## 5.0 REAGENTS

5.1 See Section 5.0 of Method 7000.

### 5.2 Preparation of standards:

5.2.1 **Stock solution:** Dissolve 11.6586 g beryllium sulfate,  $\text{BeSO}_4$ , in Type II water containing 2 mL concentrated nitric acid and dilute to 1 liter. Beryllium metal can also be dissolved in acid. Alternatively, procure a certified standard from a supplier and verify by comparison with a second standard.

5.2.2 Prepare dilutions of the stock solution to be used as calibration standards at the time of analysis. The calibration standards should be prepared using the same type of acid and at the same concentrations as in the sample after processing (0.5% v/v  $\text{HNO}_3$ ).

## 6.0 SAMPLE COLLECTION, PRESERVATION, AND HANDLING

6.1 See Chapter Three, Section 3.1.3, Sample Handling and Preservation.

## 7.0 PROCEDURE

7.1 Sample Preparation: The procedures for preparation of the sample are given in Chapter Three, Section 3.2.

7.2 See Method 7000, Paragraph 7.3, Furnace Procedure. The calculation is given in Method 7000, Paragraph 7.4.

## 8.0 QUALITY CONTROL

8.1 See Section 8.0 of Method 7000.

## 9.0 METHOD PERFORMANCE

9.1 Precision and accuracy data are not available at this time.

9.2 The performance characteristics for an aqueous sample free of interferences are:

Optimum concentration range: 1-30 ug/L.

Detection limit: 0.2 ug/L.

## 10.0 REFERENCES

1. Methods for Chemical Analysis of Water and Wastes, EPA-600/4-82-055, December 1982, Method 210.2.

METHOD 7091  
BERYLLIUM (ATOMIC ABSORPTION, FURNACE TECHNIQUE)

