

## METHOD 7481

### MOLYBDENUM (ATOMIC ABSORPTION, FURNACE TECHNIQUE)

#### 1.0 SCOPE AND APPLICATION

1.1 See Section 1.0 of Method 7000.

#### 2.0 SUMMARY OF METHOD

2.1 See Section 2.0 of Method 7000.

#### 3.0 INTERFERENCES

3.1 See Section 3.0 of Method 7000 if interferences are suspected.

3.2 Molybdenum is prone to carbide formation. Use a pyrolytically coated graphite tube.

3.3 Memory effects are possible, and cleaning of the furnace may be required after analysis of more concentrated samples or standards.

#### 4.0 APPARATUS AND MATERIALS

4.1 For basic apparatus, see Section 4.0 of Method 7000.

4.2 Instrument parameters (general):

4.2.1 **Drying time and temp:** 30 sec at 125°C.

4.2.2 **Ashing time and temp:** 30 sec at 1400°C.

4.2.3 **Atomizing time and temp:** 5 sec at 2800°C.

4.2.4 **Purge gas:** Argon (nitrogen should not be used).

4.2.5 **Wavelength:** 313.3 nm.

4.2.6 **Background correction:** Required.

4.2.7 Other operating parameters should be set as specified by the particular instrument manufacturer.

4.2.8 **Pyrolytically coated graphite tube.**

NOTE: The above concentration values and instrument conditions are for a Perkin-Elmer HGA-2100, based on the use of a 20-uL injection, continuous-flow purge gas, and nonpyrolytic graphite. Smaller sizes of furnace devices or those employing faster rates of atomization can be operated using lower atomization temperatures for shorter time periods than the above-recommended settings.

## 5.0 REAGENTS

5.1 See Section 5.0 of Method 7000.

### 5.2 Preparation of standards:

5.2.1 **Stock solution:** Dissolve 1.840 g of ammonium molybdate,  $(\text{NH}_4)_6\text{Mo}_7\text{O}_{24} \cdot 4\text{H}_2\text{O}$  (analytical reagent grade), in Type II water and dilute to 1 liter; 1 mL = 1 mg Mo (1,000 mg/L). Alternatively, procure a certified standard from a supplier and verify by comparison with a second standard.

5.2.2 Prepare dilutions of the stock solution to be used as calibration standards at the time of analysis. The calibration standards should be prepared using the same type of acid and at the same concentrations as in the sample after processing (0.5% v/v  $\text{HNO}_3$ ).

## 6.0 SAMPLE COLLECTION, PRESERVATION, AND HANDLING

6.1 See Chapter Three, Section 3.1.3, Sample Handling and Preservation.

## 7.0 PROCEDURE

7.1 Sample preparation: The procedures for preparation of the sample are given in Chapter Three, Section 3.2.

7.2 See Method 7000, Paragraph 7.3, Furnace Procedure. The calculation is given in Method 7000, Paragraph 7.4.

## 8.0 QUALITY CONTROL

8.1 See Section 8.0 of Method 7000.

## 9.0 METHOD PERFORMANCE

9.1 Precision and accuracy data are not available at this time.

9.2 The performance characteristics for an aqueous sample free of interferences are:

Optimum concentration range: 3-60 ug/L.

Detection limit: 1 ug/L.

## 10.0 REFERENCES

1. Methods for Chemical Analysis of Water and Wastes, EPA-600/4-82-055, December 1982, Method 246.2.

METHOD 7481  
MOLYBDENUM (ATOMIC ABSORPTION, FURNACE TECHNIQUE)

