

VANADIUM (ATOMIC ABSORPTION, FURNACE TECHNIQUE)

1.0 SCOPE AND APPLICATION

1.1 See Section 1.0 of Method 7000.

2.0 SUMMARY OF METHOD

2.1 See Section 2.0 of Method 7000.

3.0 INTERFERENCES

3.1 See Section 3.0 of Method 7000 if interferences are suspected.

3.2 Background correction is required.

3.3 Vanadium is refractory and prone to form carbides. Consequently, memory effects are common, and care should be taken to clean the furnace before and after analysis.

3.4 Nitrogen should not be used as a purge gas.

4.0 APPARATUS AND MATERIALS

4.1 For basic apparatus, see Section 4.0 of Method 7000.

4.2 Instrument parameters (general):

4.2.1 **Drying time and temp:** 30 sec at 125°C.

4.2.2 **Ashing time and temp:** 30 sec at 1400°C.

4.2.3 **Atomizing time and temp:** 15 sec at 2800°C.

4.2.4 **Purge gas:** Argon (nitrogen should not be used).

4.2.5 **Wavelength:** 318.4 nm.

4.2.6 **Background correction:** Required.

4.2.7 Other operating parameters should be set as specified by the particular instrument manufacturer.

NOTE: The above concentration values and instrument conditions are for a Perkin-Elmer HGA-2100, based on the use of a 20- $\mu$ L injection, continuous-flow purge gas, and nonpyrolytic graphite. Smaller sizes of furnace devices or those employing faster rates of atomization can be operated using lower atomization temperatures for shorter time periods than the above-recommended settings.

## 5.0 REAGENTS

5.1 See Section 5.0 of Method 7000.

### 5.2 Preparation of standards:

5.2.1 **Stock solution:** Dissolve 1.7854 g of vanadium pentoxide,  $V_2O_5$  (analytical reagent grade), in 10 mL of concentrated nitric acid and dilute to 1 liter with Type II water. Alternatively, procure a certified standard from a supplier and verify by comparison with a second standard.

5.2.2 Prepare dilutions of the stock solution to be used as calibration standards at the time of analysis. The calibration standards should be prepared using the same type of acid and at the same concentrations as in the sample after processing (0.5% v/v  $HNO_3$ ).

## 6.0 SAMPLE COLLECTION, PRESERVATION, AND HANDLING

6.1 See Chapter Three, Section 3.1.3, Sample Handling and Preservation.

## 7.0 PROCEDURE

7.1 Sample preparation: The procedures for preparation of the sample are given in Chapter Three, Section 3.2.

7.2 See Method 7000, Paragraph 7.3, Furnace Procedure. The calculation is given in Method 7000, Paragraph 7.4.

## 8.0 QUALITY CONTROL

8.1 See Section 8.0 of Method 7000.

## 9.0 METHOD PERFORMANCE

9.1 Precision and accuracy data are available in Method 286.2 of Methods for Chemical Analysis of Water and Wastes.

9.2 The performance characteristics for an aqueous sample free of interferences are:

Optimum concentration range: 10-200 ug/L.

Detection limit: 4 ug/L.

## 10.0 REFERENCES

1. Methods for Chemical Analysis of Water and Wastes, EPA-600/4-82-055, December 1982, Method 286.2.

METHOD 7911  
VANADIUM (ATOMIC ABSORPTION, FURNACE TECHNIQUE)

