




# Information Collection Requirements Rule— Protozoa and Enteric Virus Sample Collection Procedures



## ABOUT THIS MANUAL

---

This manual is designed to be brought into the field by drinking water utility personnel when collecting source and finished water samples for protozoa and viruses. The sample collection steps in this manual are consistent with those demonstrated in the accompanying video. To further associate the steps in this manual with the sampling demonstration on the video, the photos for each step are taken directly from the video.

Several graphic conventions are used throughout the manual to differentiate steps or denote special actions:



A step icon is used at the beginning of each step. These steps are parallel to those in the accompanying video.



Actions denoted by this icon are critical to ensuring that the sample will be valid and uncontaminated, such as putting on fresh latex gloves before handling the filter.



Text denoted by this icon provides additional information to the samplers, but may not be part of the actual collection procedure.

Collecting protozoan and virus samples correctly under the Information Collection Requirements Rule can be challenging. Please watch the demonstration video before collecting the samples, and be sure to follow each step in this manual when in the field.



**CONTENTS**

---

The Information Collection Requirements Rule .....	1
Questions Commonly Asked by Drinking Water Utilities .....	3
Sample Collection Procedures for Detecting Protozoa in Water .....	7
Collecting Source Water Samples .....	11
Collecting Finished Water Samples .....	19
Sample Collection Procedures for Detecting Enteric Viruses in Water .....	27
Collecting Source Water Samples .....	33
Collecting Finished Water Samples .....	49
Credits and Acknowledgements .....	63

## PROTOZOAN AND ENTERIC VIRUS SAMPLE COLLECTION PROCEDURES AS DEFINED BY THE INFORMATION COLLECTION REQUIREMENTS RULE

---

This manual describes the procedures for collecting source water and finished water samples for protozoan and enteric virus monitoring under the Information Collection Requirements (ICR) rule. This manual and the accompanying video comprise a two-part set of instructional materials that provide public water supply systems with the information needed to properly collect samples for protozoan and virus monitoring. All water utility personnel involved with ICR monitoring should watch the video and review this manual before collecting any samples.

The protozoan collection procedures described in this manual and in the video are based on the procedures in the ICR Protozoan Method for Detecting *Giardia* cysts and *Cryptosporidium* Oocysts in Water by a Fluorescent Antibody Procedure. The total culturable virus collection procedures described in this manual and in the video are based on the procedures in the Virus Monitoring Protocol for the Information Collection Rule. Both of these methods can be requested by calling the Safe Drinking Water Hotline, at (800) 426-4791.



## QUESTIONS COMMONLY ASKED BY DRINKING WATER UTILITIES

---

### ***What is the purpose of the ICR rule?***

The ICR rule was developed by EPA to collect occurrence, exposure, and treatment data on drinking water pathogens and disinfectant by-products. The pathogen data are needed to determine whether current Surface Water Treatment Regulations should be revised to include new or more stringent treatment levels for some microbes. The disinfectant by-product data are needed to determine whether to regulate the chemical by-products that form when disinfectants react with organic chemicals in source water.

Although drinking water utilities will be involved in collecting both disinfectant by-product and waterborne pathogen data under the ICR rule, this manual describes the utility's role in collecting data on drinking water pathogen occurrence.

### ***What pathogens are monitored under the ICR rule?***

The ICR rule requires public water supply systems to monitor source water (and finished water in some cases) for the following pathogens:

- *Giardia* cysts
- *Cryptosporidium* oocysts
- Total culturable viruses
- Fecal coliform or *Escherichia coli* bacteria
- Total coliform bacteria

EPA is considering revising the current Surface Water Treatment Regulations because existing treatment levels for *Giardia* and viruses

may not be adequate to protect public health for systems supplied by poor-quality source water and because of the new threat posed by *Cryptosporidium*.

*Giardia* cysts in drinking water cause more reported waterborne disease outbreaks than any other single known pathogen. They also are more resistant to environmental stresses and disinfection than almost all other known waterborne pathogens.

*Cryptosporidium* oocysts in drinking water have caused major waterborne disease outbreaks in the U.S. and other countries and are even more resistant to disinfection than *Giardia*.

Several enteric viruses have caused waterborne disease and may be responsible for many, if not most, of the outbreaks where a causative agent was not identified (about half of all reported outbreaks). Adequate analytical methodology is not yet available for routine analysis for many enteric viruses, so EPA has required monitoring of total culturable viruses. Total culturable viruses are a group of enteric viruses commonly found in poor-quality waters and which EPA believes are at least somewhat representative of other pathogenic viruses. Monitoring for total culturable viruses is useful because this group contains pathogens and is a potential indicator of other viral pathogens.

Fecal coliforms, *E. coli*, and total coliforms have been used for decades to assess source water quality. Coliform bacteria are much more sus-

ceptible to environmental stress and disinfection than protozoa and viruses, and would be eliminated by any system that eliminated more resistant pathogens. However, the ICR rule requires drinking water utilities to submit coliform monitoring data as general indicators of water quality. Monitoring procedures for fecal coliform, *E. coli*, and total coliform densities have been established and are not addressed by this manual.

***Which drinking water utilities have to collect protozoan and virus samples?***

Public water supply systems that serve between 10,000 and 100,000 people and use surface water (or groundwater under the influence of surface water) are required to monitor their source water for *Giardia* cysts and *Cryptosporidium* oocysts.

Public water supply systems that serve more than 100,000 people and use surface water (or groundwater under the influence of surface water) are required to monitor their source water for *Giardia* cysts, *Cryptosporidium* oocysts, and total culturable viruses. If pathogen densities in the source water exceed 1 pathogen per liter during the first 12 months of monitoring, then public water supply systems also must sample finished water for the remaining months.

***How often must samples be taken?***

Public water supply systems that serve between 10,000 and 100,000 people must collect samples every two months for 12 months.

Systems that serve more than 100,000 people must take samples every month for 18 months.

However, these systems may discontinue monitoring if:

- Viruses are not detected in the source water during the first 12 months of monitoring, or
- Source water has been tested for either total coliforms or fecal coliforms at least five times per week for four months before and two months after the effective date of the ICR and the total coliform density is less than 100 colonies/100 mL or the fecal coliform density in 90 percent of all samples is less than 20 colonies/100 mL.

***Where should samples be collected?***

Samples must be taken at the intake of each treatment plant. If a plant has several sources of water, the system must sample the blended water from all sources. If this is not possible, the source with the highest expected pathogen concentration should be sampled.

***Who will analyze the samples?***

EPA has approved several laboratories to analyze the protozoan and virus samples. Before collecting samples, you must arrange to have them analyzed by an EPA-approved laboratory. If you have not already located an approved laboratory, notify:

ICR Laboratory Coordinator  
EPA Office of Ground Water & Drinking Water  
26 West Martin Luther King Drive  
Cincinnati, Ohio 45268.


EPA will provide you with a list of approved laboratories or other appropriate guidance.



**SAMPLE COLLECTION PROCEDURES  
FOR DETECTING PROTOZOA IN WATER**

---

# Sampling Train for Collecting Protozoa





**1** Each month, your laboratory will send you all of the equipment needed to collect samples for *Giardia* cyst and *Cryptosporidium* oocyst analyses. When you receive the sampling kit, check the contents of the carton. The sampling kit should contain the following items:

- Sampling train for collecting protozoa (left):
  - Inlet hose
  - Pressure regulator with pressure gauge
  - Fluid proportioning injector module, including an injector and pressure gauge\*
  - 1- $\mu$ m nominal porosity filter and holder made by Parker Hannifan or Filterite
  - Water meter
  - Effluent hose and flow control valve

\*Needed for finished water sample collection only

- ☐ Plastic sample bags
- ☐ Ice packs for shipping the collected samples
- ☐ Sample labels

If you are missing any items, contact your laboratory immediately. Do not attempt to collect the samples without a complete sampling kit.



Once you have verified the contents of the sampling kit, place the ice packs in the freezer and repack the box for later use.