

EPA/OPP MICROBIOLOGY LABORATORY
ESC, Ft. Meade, MD

Standard Operating Procedure
for
VITEK: Quality Control Procedures

SOP Number: QC-17-02

Date Revised: 08-22-02

Prepared By: _____ Date: ____/____/____

Print Name: _____

Reviewed By: _____ Date: ____/____/____

Print Name: _____

Technical Staff

_____ Date: ____/____/____

Print Name: _____

QA Officer

_____ Date: ____/____/____

Print Name: _____

Laboratory Director

Date Issued: ____/____/____

Withdrawn By: _____ Date: ____/____/____

Controlled Copy No.: _____

TABLE OF CONTENTS

<u>Contents</u>	<u>Page Number</u>
1.0 SCOPE AND APPLICATION.....	2
2.0 DEFINITIONS.....	2
3.0 HEALTH AND SAFETY.....	2
4.0 CAUTIONS.....	2
5.0 INTERFERENCES.....	2
6.0 PERSONNEL QUALIFICATIONS.....	2
7.0 SPECIAL APPARATUS AND MATERIALS.....	2
8.0 INSTRUMENT OR METHOD CALIBRATION.....	3
9.0 SAMPLE HANDLING AND STORAGE.....	4
10.0 PROCEDURE AND ANALYSIS.....	5
11.0 DATA ANALYSIS/CALCULATIONS.....	6
12.0 DATA MANAGEMENT/RECORDS MANAGEMENT.....	6
13.0 QUALITY CONTROL.....	7
14.0 NONCONFORMANCE AND CORRECTIVE ACTION.....	7
15.0 REFERENCES.....	7
16.0 FORMS AND DATA SHEETS.....	7

1.0 SCOPE AND APPLICATION:

- 1.1 This protocol describes the method for verifying the performance of the Vitek Gram-Positive Identification Card (GPI), *Bacillus* Biochemical Card (B), and the Gram-Negative Identification+ (GNI+) Card.

2.0 DEFINITIONS:

- 2.1 ATCC = American Type Culture Collection

3.0 HEALTH AND SAFETY:

- 3.1 Laboratory personnel should follow biosafety procedures appropriate for the organism being confirmed as outlined in SOP MB-01, Biosafety in the Laboratory.

4.0 CAUTIONS: None

5.0 INTERFERENCES:

- 5.1 Do not use Vitek cards beyond expiration date because it may result in erroneous readings.
- 5.2 Improper maintenance of stock cultures, subculturing, and filling of Vitek cards may result in inconsistent or erroneous biopatterns.

6.0 PERSONNEL QUALIFICATIONS:

- 6.1 Personnel are required to be knowledgeable of the procedures in this SOP. Documentation of training and familiarization with this SOP can be found in the training file for each employee.

7.0 SPECIAL APPARATUS AND MATERIALS:

- 7.1 Vitek 32 System for the automated identification of microorganisms.
- 7.2 Vitek 32 Identification cards (GPI, B, and GNI+).
- 7.3 Quality control organisms listed in Section 10.

8.0 INSTRUMENT OR METHOD CALIBRATION:

8.1 Factory Calibrations: Prior to shipment, the Vitek 32 Instrument met all acceptance test procedures stipulated by bioMerieux. The Field Service Engineer has performed a verification of the bioMerieux Vitek, Inc. factory calibration as a part of the installation procedure of this instrument. The instrument test analysis outlined on the installation checklist has been performed and has met the instrument engineering specifications as indicated on the Customer Calibration Verification certificate (see 16.2).

8.2 Internal monitoring of the Vitek reader/incubator module.

8.2.1 The Vitek reader/incubator module houses the card handling and scanning mechanism as well as the heater that maintains the cards at the required incubation temperature. The trays that hold the cards are mounted to a carousel that rotates once every 15 minutes to position the cards for data scanning and identification. A thermistor is located in the center of the carousel shaft and positioned to monitor any change of temperature in the carousel stack. A heater and fan on top of the carousel maintains the temperature at 35°C.

The incubation temperature is verified during the Vitek system test. Temperature deviations of $\pm 2^{\circ}\text{C}$ generate error messages at the data terminal module such as, "Reader Temperature High," or "Reader Temperature Low." The process cycle is aborted if the temperature varies $\pm 5^{\circ}\text{C}$ from the set temperature for more than one hour. The thermistor assembly in the Vitek Reader is used as an absolute reference in the Vitek System. Each thermistor assembly is tested at the Vitek manufacturing facility to Vitek's test specifications, using equipment calibrated to the National Bureau of Standards. In a final test, the Reader and the other Vitek modules are functionally tested as a system.

8.2.2 The optical system in the Vitek consists of Light Emitting Diodes (LED's) which generate light at a precise wavelength and photo transistors which detect that light. The LED array

is physically aligned in a plane 180 degrees with respect to the photo transistors. The LED/photo transistor array is formatted physically to the VITEK Card in a figure eight design for the identification segment of the card and in a row of five for the growth wells. The specifications of the LED's and the photo transistors, along with their corresponding tolerance, result in a transfer function which relates the LED excitation currents to the photo transistor detection current within a tolerance. That relationship is verified in a calibration algorithm which is coded into the firmware (software) of every Reader. If any part of that transfer function is out of tolerance, it is indicated with a "calibration failure" message; and processing will not resume until the calibration tests are completed successfully. Every LED/photo transistor pair is calibrated by this process, and this occurs during the production process as well as during actual operation of the Reader in the field. Calibration is completed before reading every tray.

In addition to the optical testing, a series of test protocols are performed on the Reader/Incubator system prior to shipment. The card handling system, carousel sensor, optical tray sensor, cars stop and horizontal in/out stops are verified and/or adjusted for optimal performance. Output voltages are verified on pins and connectors for the Intelligent Reader Board and Reader/Incubator Boards, and system is readied for a "burn-in" test. "Burn-in" consists of loading the instrument with non-inoculated but otherwise functional cards and then processing them for four hours. All messages are checked for correctness and the system is then secured for packaging and shipment.

- 8.2.3 A performance check of the VITEK unit will be performed on an annual basis by a trained bioMerieux technician. All service and recertification documents will be maintained in the VITEK Maintenance and Certification Records Book.

9.0 SAMPLE HANDLING AND STORAGE: Not applicable

10.0 PROCEDURE AND ANALYSIS:

- 10.1 Quality control procedures consist of performing standard VITEK analyses using a specified set of known microbes in conjunction with the appropriate card types. The QC will be performed on a semiannual basis. QC organisms and associated card types are listed in Table 1.
- 10.2 Lyophilized cultures will be purchased from the ATCC and will be stored according to the ATCC recommendations (in trypticase soy broth with 20% glycerol at -70°C). Cultures are propagated for storage according to the Product Information Sheets supplied by ATCC.
- 10.3 Stored cultures will be initiated for VITEK analysis according to the VITEK Pinset instructions outlined in the VITEK Pinset Folder.
- 10.4 For VITEK analysis, organisms are assigned a 7 digit identification number (format: 000000-0). The first six digits are the Reference Number for the organism as listed in Table 1. The seventh digit is assigned by VITEK automatically and is derived from the VITEK card lot being used. The six digit number is entered directly onto the VITEK card and will appear later on the VITEK QC Deviation Report. Follow the instructions for labeling cards outlined in the VITEK Pinset Folder. Identification numbers are recorded on the VITEK Quality Control Log in the VITEK QC Log Book. Quality Control organism tracking and Deviation Reports are maintained in the VITEK Quality Control Record Log Book.
- 10.5 Following analysis, the VITEK unit will automatically print a VITEK QC Deviations Report; a report which displays any deviations of "expected" reactions versus those "observed" during the analysis (see 16.3).

Table 1. Quality Control Organisms for VITEK 32 Automated Identification System

Organism	ATCC #	Card	Reference Number
<i>Bacillus licheniformis</i>	12759	B	600101
<i>Bacillus sphaericus</i>	4525	B	600102
<i>Proteus mirabilis</i>	7002	GNI+	900101
<i>Providencia alcalifaciens</i>	51902	GNI+	900102
<i>Klebsiella pneumoniae</i>	13883	GNI+	900103
<i>Plesiomonas shigelloides</i>	51903	GNI+	900104
<i>Bordetella bronchiseptica</i>	10580	GNI+	900105
<i>Serratia liquefaciens</i>	27592	GNI+	900106
<i>Leclercia adecarboxylata</i>	23216	GNI+	900107
<i>Burkholderia cepacia</i>	25608	GNI+	900108
<i>Enterococcus durans</i>	6056	GPI	800101
<i>Streptococcus equi</i>	9528	GPI	800102
<i>Streptococcus bovis</i>	9809	GPI	800103
<i>Erysipelothrix rhusiopathiae</i>	19414	GPI	800104
<i>Streptococcus pyogenes</i>	19615	GPI	800105
<i>Enterococcus faecalis</i>	29212	GPI	800106
<i>Staphylococcus xylosus</i>	29971	GPI	800107

11.0 DATA ANALYSIS/CALCULATIONS: None

12.0 DATA MANAGEMENT/RECORDS MANAGEMENT:

12.1 Data will be recorded promptly, legibly, and in indelible ink. Information on Quality Control will be recorded on the Quality Control Record Log Sheet (see 16.1). Quality Control organism tracking forms will be maintained in the VITEK Quality Control Record Book. Completed forms and reports are archived in notebooks kept in locked file cabinets in file

room D217. Only authorized personnel have access to the locked files. Archived data is subject to OPP's official retention schedule contained on SOP ADM-03, Records and Archives.

13.0 QUALITY CONTROL:

- 13.1 The OPP Microbiology Laboratory conforms to 40CFR Part 160, Good Laboratory Practices. Appropriate quality control measures are integrated into each SOP.

14.0 NONCONFORMANCE AND CORRECTIVE ACTION:

- 14.1 A VITEK QC Deviation Report is generated if an actual value differs from an expected value. If the QC Deviation Report is generated, the QC procedure should be repeated. Strict adherence to given parameters for that organism should be followed, specifically the instructions for maintaining the organism prior to inoculating the VITEK card (i.e., passing the organism the indicated number of times on the specified medium, the proper age organism, the proper optical density for that card, etc.). If the same results are returned on the QC Deviation Report, a bioMerieux technician should be notified and maintenance scheduled. The machine should not be used for confirmation until it is serviced and recertified.

15.0 REFERENCES:

- 15.1 bioMerieux VITEK, Inc. 1997. Industrial Automated Microbiology Systems Summary, Part Number 512312-2, REV 0597.
- 15.2 bioMerieux VITEK, Inc. Aug. 26, 1998. Pininsert, VITEK Software Programs and Data, VTK-R06.01

16.0 FORMS AND DATA SHEETS:

- 16.1 Quality Control Record Log Sheets
- 16.2 Customer Calibration Verification Certificate
- 16.3 Sample QC Deviation Report
- Vitek Quality Control Record Log Sheet

OPP Microbiology Laboratory

QUALITY CONTROL RECORD LOG					
Date/Init.	Reference No.	Card*	Card Lot	Organism	Miscellaneous Test Results ⁺

* Enter Vitek card type (B, GNI+, or GPI for Bacillus card, Gram Negative Identification, or Gram Positive Identification respectively)

+ Miscellaneous results include catalase test for Bacillus species and Staphylococcus, coagulase test for Staphylococcus species, and oxidase test for Gram negative rods.

Customer Calibration Verification Certificate US EPA/OPP/Microbiology Laboratory

Customer Calibration Verification	
Customer Name:	USEPA
Site #:	1874
Date:	08/18/98

<input checked="" type="checkbox"/>	This instrument has met all acceptance test procedures (ATP) before leaving bioMérieux Vitek's facility. The Field Service Engineer has performed a verification of the bioMérieux Vitek, Inc. factory calibration as a part of the installation procedure of this instrument. The instrument test analysis outlined on the installation checklist has been performed and has met the instrument engineering specifications as indicated on this document.
<input type="checkbox"/>	The bioMérieux Vitek, Inc. factory calibration of this instrument was verified by a bioMérieux Vitek, Inc. Field Service Engineer at installation and is verified at preventive maintenance intervals by the Field Service Engineer. This Field Service Report is your documentation that this instrument has met the preventive maintenance calibration specifications as required by bioMérieux Vitek, Inc. The Field Service Engineer has performed the service as indicated on the report.

FSE Signature:



Sample QC Deviation Report US EPA/OPP/Microbiology Laboratory

Date: 03/16/2000 11:09:58
KSTY-106.01

bioMerieux Titak
Titak QC Exception Report:

Test Type: bacII Card Lot: K32N Expiration: 09/07/2000
Test Date: 03/16/2000 Technologist:
Expected Organism ID: ATCC 12759 Bacillus licheniformis
Actual Organism ID: Bacillus licheniformis

Biochem
T T G I G I I W R S A I Z N T P J W A X W O Y I P S E Y
I U I A C H A S P A A I S W I C G O I W C C S A I A I S S
I C C G U I A I W P L S U E T E I R G Y N C Y S A I S E C W
Expected - + + + + -
Actual - + + + + -

Deviations: Nil