


# Performance Evaluation Program Data Management

Tim Hanley  
EPA - OAQPS



# PEP Data Management

- Background
- Data Validation
- Approval
- AIRS submittal
- Next Steps

# PEP Data Management Background

- Implemented Access database in each of two regional quality assurance laboratories in early '99
  - ▶ Athens, Georgia
  - ▶ Manchester, Washington
- First versions of database focused on acquisition of data
  - ▶ matching laboratory SOP's with need to track many types of data
  - ▶ acquiring sampler data through electronic data package


# Designing Current Version of PED to incorporate validation

- Need to distinguish Performance Evaluation (PE) data as "valid" or "rejected"
- Need to have appropriate manager approve data as "valid" or "rejected"
- Identify all criteria into a system that can be broken down by:
  - ▶ where data is produced
  - ▶ when data is produced
  - ▶ whether data is part of:
    - Critical criteria
    - Sample batch validation
    - Operational evaluation criteria

# PEP data validation

- Used PEP QAPP section on "Validation and verification of methods" as basis for setting up data validation process.
- Setup PED to identify for each performance evaluation the:
  - ▶ Critical criteria
  - ▶ Sample batch validation criteria
  - ▶ Operational evaluation criteria
- Compare associated PE data to acceptable parameters for each criteria

# PEP Validation Matrix



# **Data that feeds into the validation process**

- Laboratory data
- Field Data
  - ▶ Electronic data package
  - ▶ Information from the Field Data Sheets (FDS)
- Systematic Issues
  - ▶ Range and CV checks

# PEP Data Validation

- If all criteria "pass", PE identified as:
  - ▶ "Valid Data" and
  - ▶ "complete"
- If one or more criteria:
  - ▶ "fail", allow for correction or override:
    - Corrections may be made due to entry errors
    - Overrides may be made when non-critical criteria have minor flags associated with them
  - ▶ PE's that have a "fail" are not automatically marked as "complete"



## PEP Data Validation (*continued*)

- Re-run validation only for those samples that have not already been completed.
  - ▶ PE "Valid" and "Complete"
  - ▶ PE "Rejected"
  - ▶ PE can be marked "Complete" or more corrections and overrides can take place.

PM<sub>2.5</sub> Performance Evaluations

Record Last Changed: 5/15/00

PE Filter ID: T127630

Sampler ID: BG0217

AIRS Site Code: 220550005

POC: 1

Filter Lot: 99357

Sampler Serial No.: 0217

Model:

PE Cassette ID: 3543

Sampler Model: BG1PQ200A

Serial No.:

## PE Run Time

PE Start Date: 2/24/00 12:00:05 AM

Elapsed Time (days): 0.999

PE Stop Date: 2/24/00 11:59:05 PM

Filter Removed Date: 2/25/00 10:54:20 AM

PM<sub>2.5</sub> ConcentrationMeasured by PE Sampler ( $\mu\text{g}/\text{m}^3$ ): 8.537Measured by Site ( $\mu\text{g}/\text{m}^3$ ):

Percent Difference:

Field Flag:

Comment:

## Sample Flow/Volume

Flow Avg (L/min): 16.69

Flow CV: 0.23

Sample Vol. ( $\text{m}^3$ ): 24.013Ambient Sample Temperature ( $^{\circ}\text{C}$ )

Average: 21.3

Minimum: 17.9

Maximum: 26.5

## Barometric Pressure (mm Hg)

Average: 765

Minimum: 764

Maximum: 768

## Out of Spec.

Time: ☐Temperature: ☐Flow: ☐

## Field Verifications

BP: ☒Temperature: ☒Flow: ☒

## Field Blank

Cassette ID: 4604

Filter ID: T127631

Pre- Wt. (mg): 139.130

Post- Wt. (mg): 139.141

Wt. Change (mg): 0.011

## Trip Blank

Cassette ID:

Filter ID:

Pre- Wt. (mg):

Post- Wt. (mg):

Wt. Change (mg):

## Lab Blank

Cassette ID: N/A

Filter ID: T127625

Pre- Wt. (mg): 135.880

Post- Wt. (mg): 135.875

Wt. Change (mg): -0.005

## Collocated Samplers

Collocated: ☐

Count:

Std. Dev.:

Mean:

Max CV (%):

## Pre-sample (Unexposed)

## Post-sample (Exposed)

## Pre/Post Change

## Weigh Sessions

Equilibration Tray:

TR13

1/27/00 1:14:48 PM

TR02

3/3/00 3:47:36 PM

Equilibration Time (days):

20

4

Holding Time (days):	8	12	(Sample Mass)				
Filter Weight(mg):	139.286	2/16/00 9:38:02 AM	139.491	3/7/00 2:40:13 PM	0.205		
Weighing Room Conditions							
Avg. Temperature (°C):	21.8	21.8	(RH Difference) -2.100				
Temp. Std. Dev. (°C):	0.132	0.15					
Avg. Relative Humidity (%):	37.3	35.2					
RH Std. Dev. (%):	0.303	0.294					
Batch Duplicates							
Filter ID:	T127630	T127630					
Primary Weight (mg):	139.286	2/16/00 9:38:02 AM				139.491	3/7/00 2:40:13 PM
Duplicate Weight (mg):	139.285	2/16/00 9:49:11 AM				139.500	3/7/00 2:55:51 PM
Weight Change (mg):	-0.001	0.009					
Low Mass Balance Checks							
Verified Weight (mg):			(Max. Overall Wt. Change) 0.003				
Minimum Weight (mg):	99.999	2/16/00 9:49:54 AM				99.998	3/7/00 2:30:50 PM
Maximum Weight (mg):	100.001	2/16/00 9:35:05 AM				99.998	3/7/00 2:30:50 PM
Min/Max Wt. Change (mg):	0.002	0.000					
High Mass Balance Checks							
Verified Weight (mg):			(Max. Overall Wt. Change) 0				
Minimum Weight (mg):	199.998	2/16/00 9:35:50 AM				199.998	3/7/00 2:31:33 PM
Maximum Weight (mg):	199.998	2/16/00 9:35:50 AM				199.998	3/7/00 2:31:33 PM
Min/Max Wt. Change (mg):	0.000	0.000					
Batch Stability Tests (Lot Exposure Blanks)							
Filter Count:	3	3					
Average Wt. Diff(mg):	0.0033	0.003					
Maximum Wt. Diff. (mg):	0.005	0.004					

Lab Approval

Valid Data: ☒  
 PE Rejected: ☐  
 Complete: ☒

QA Manager Approval

Approval: ☐  
 Initials:  
 Date:

Transfer to AIRS

Approval: ☐  
 Initials:  
 Date:



PM<sub>2.5</sub> Performance Evaluations

Record Last Changed: 5/15/00

PE Filter ID: T127630

Sampler ID: BG0217

AIRS Site Code: 220550005

POC: 1

Quality Control Checks

Check Date: 5/15/00 4:10:06 PM

Check ID	Flag ID	Check Code	Description	Check Value	Auto QC		Override		Comment
					Pass	Fail	Pass	Fail	
Critical criteria									
1	2	DAM	no filter damage (visual defect)	True	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
2	3	EST	sample period 1380-1500 min	1439 min	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
3	8	FLR_1	flow rate $\leq \pm 5\%$ of 16.67 L/min	16.69 L/min	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
4	8	FLR_2	flow rate $\leq 2\%$ CV	0.23	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
5	13	FVL	no flow rate excursions $> \pm 5\%$ for $> 5$ min	True	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
6	33	HTE_1	sample recovery $\leq 96$ hours from sample end date	11 hours	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
7	33	HTE_2	post-sample weighing $\leq 10$ days at 25 deg C or $\leq 30$ days at 4 deg C	12 days, 0 deg C (cold pack: Frozen)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
8	34	ISP_1	pre-sample minimum 24 hour equilibration	476 hours	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
9	34	ISP_2	pre-sample mean temperature 20-23 degrees C	21.8 degrees C	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
10	34	ISP_3	pre-sample temperature control $\pm 2$ degrees C over 24 hours	0.132 degrees C	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
11	34	ISP_4	pre-sample mean RH 30-40%	37.3 percent	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
12	34	ISP_5	pre-sample RH SD control $\pm 5\%$	0.303 percent	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
13	34	ISP_6	post-sample minimum 24 hour equilibration	95 hours	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
14	34	ISP_7	post-sample mean temperature 20-23 degrees C	21.8 degrees C	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
15	34	ISP_8	post-sample temperature control $\pm 2$ degrees C over 24 hours	0.15 degrees C	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
16	34	ISP_9	post-sample mean RH 30-40%	35.2 percent	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
17	34	ISP_10	post-sample RH SD control $\pm 5\%$	0.294 percent	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
18	34	ISP_11	pre/postsample RH $\pm 5\%$	-2.1 percent	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Sample batch validation with major and minor flags

File Edit View View Tools Window Help						75% Close			
18	34	ISP_11	pre/postsample RH $\pm$ 5%	-2.1 percent	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Sample batch validation with major and minor flags									
19	26	FIS_1	pre-sample 100 mg balance check $\leq$ 3 ug	2 ug	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
20	26	FIS_2	pre-sample 200 mg balance check $\leq$ 3 ug	0 ug	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
21	28	FLD_1	pre-sample duplicate filter weight $\pm$ 15 ug	-1 ug	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
22	26	FIS_3	post-sample 100 mg balance check $\leq$ 3 ug	0 ug	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
23	26	FIS_4	post-sample 200 mg balance check $\leq$ 3 ug	0 ug	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
24	28	FLD_2	post-sample duplicate filter weight $\pm$ 15 ug	9 ug	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
25	27	FLB	lab blanks $\pm$ 15 ug	-5 ug	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
26	25	FFB	field blanks $\pm$ 30 ug	11 ug	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Operational evaluation criteria									
27	9	FLT	filter temperature, no excursions of $>$ 5 degrees C lasting longer than 30 min	True	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
28	6	FAT	temperature verification $\pm$ 2 degrees C	True	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
29	11	FPC	barometric pressure verification $\pm$ 10 mm Hg	True	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
30	12	FSC	flow rate verification $\pm$ 4%	True	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
31	33	HTE_3	(pre-sample) filter holding $\leq$ 30 days from pre-weigh	8 days	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
32	18	BDL	lower detection limit (PM2.5 conc. $\geq$ 2 ug/m3)	8.537 ug/m3	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
33			upper concentration limit (PM2.5 conc. $\leq$ 200 ug/m3)	8.537 ug/m3	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
34	24	FCS	collocated CV $\leq$ 10%	Not Collocated	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
35	32	FRW_1	pre-sample batch stability test $<$ 15 ug between weighings (each filter)	5 ug (max. diff.)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
36	32	FRW_2	pre-sample batch stability test $<$ 5 ug between weighings (average of at least 3 filters)	3.3 ug (3 filters)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
37	32	FRW_3	post-sample batch stability test $<$ 15 ug between weighings (each filter)	4 ug (max. diff.)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
38	32	FRW_4	post-sample batch stability test $<$ 5 ug between weighings (average of at least 3 filters)	3 ug (3 filters)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Page: 46

Ready

# QA Manager Approval

- Once Laboratory Analyst has PE's identified as either:
  - ▶ "Valid Data" and "Complete" or
  - ▶ "PE Rejected" and "Complete"
- PE's are submitted to the QA Manager for approval:
  - ▶ Can have "Approval" or
  - ▶ Returned to lab for further investigation

# AIRS Data Submittal

- Once data are approved by the QA manager AIRS transactions of those PE's can be generated.
- Example transactions from database output to text file:

942133000888101171059906111 118	
942043040188101171059906111 123	
942129000888101171059906171 125	
942101013688101171059906231 181	
942095002588101171059906081 266	
942101013688101171059906021 223	
942003029088101171059906021 220	

PM 2.5 Performance Evaluation Database - v 2.02 - [Output AIRS File]

File Edit View View Insert Format Records Tools Window Help

PE Primary Filter ID

PE Primary Cassette ID

AIRS Site ID

State

EPA Region

Evaluated after

Evaluated before

QA Approved

☐ Any value

☒ True

☐ False

AIRS Approved

☒ Any value

☐ True

☐ False

Write AIRS File

Reset Form

Exit Form

Form View

NUM



# Next Steps

- Full validation of '99 dataset in each lab
- Approval by each QA manager
- Decision on submittal to AIRS?
  - ▶ Current system does not provide the required statistics identified in part 58
  - ▶ Current system requires submittal of audit and State primary data in the same transaction; therefore need to have "Official" State value before submittal of PE values
  - ▶ Re-engineered system allows submission of audit data independent of State primary value