

BioSeek


Complex Human Cell Systems for Understanding Toxicity Mechanisms

Ellen L. Berg, PhD

May 23, 2007


Biological Complexity

Modular design / Networked architecture




Biological Complexity

Modular Design




- Great diversity from few components
 - A given component can contribute to “many” functions
- The ability to compartmentalize
 - Safety feature to limit damage




Biological Complexity

Networked Architecture




- **Rapid responses to environment**
 - Efficient information flow
- **Framework for control systems**
 - Feedback mechanisms, etc.
- **Many potential outcomes**
 - System “wiring” determines outcome




Biological Complexity

Control Systems



- Dominate in networks with demanding performance requirements
- Important for providing tolerance to perturbations (robustness)
- “Hidden nature”
 - Unexpected fragility




Modular, Networked Architecture

Consequences for Drug Discovery & Predictive Toxicology

- **Drug targets function in multiple biological processes**
 - Different pathways, cell types, different environments
 - Context, context, context
- **Hard to predict clinical trial efficacy / human exposure effects**
 - Redundant mechanisms are difficult to detect
 - Compensatory mechanisms often override drug effects
- **Drugs can have effects far downstream of the target**
 - Unexpected side effects are often only discovered in clinical testing
- **Problems become amplified since most drugs have multiple targets**


BioMAP® Technology

Bridging the Gap



BioMAP® Systems


Complex Cell-Based In Vitro Models of Human Disease



- ***BioMAP Systems integrate biological complexity***

Current BioMAP Systems

Inflammation / Autoimmune Disease / Respiratory / Cardiovascular

Primary Human Cell Types	Disease Relevance	BioMAP® System
Endothelial cells (EC)	Th1/ Th2 inflammation, allergy, asthma, dermatitis, angiogenesis, wound healing, restenosis, atherosclerosis (coronary art.)	
EC + Peripheral Blood Mononuclear Cells	Th1 inflammation, psoriasis, COPD fibrosis, monocyte and T cell responses	
EC + Macrophages	Macrophage responses, arthritis, COPD, fibrosis	
EC + Mast cells	Asthma, allergy, dermatitis, fibrosis	
EC + Smooth Muscle Cells	Vascular biology, restenosis, atherosclerosis	
EC + Th2 blasts	Allergy, asthma	
Fibroblasts	Arthritis, asthma, dermatitis, fibrosis, psoriasis, wound healing	
Myofibroblasts	Fibrosis, COPD, wound healing	
Keratinocytes	Psoriasis, dermatitis, wound healing	
Keratinocytes + Fibroblasts	Psoriasis, dermatitis, wound healing	
Bronchial Epithelial Cells	Th1 and Th2 inflammation Allergy, asthma, fibrosis, COPD	
Bronchial Epithelial Cells + Fibroblasts	Asthma, allergy, fibrosis, COPD	
Smooth Muscle Cells	Vascular inflammation, asthma, COPD, fibrosis (coronary artery SMC)	

Broad Coverage of Targets and Pathways

Detect and Distinguish Known and Novel Therapeutics

Inflammation / Autoimmune

- Calcineurin, TCR
- Glucocorticoid R
- Prostaglandin, leukotriene
- TNF- α
- IL-10
- NF κ B
- IL-1, IFN γ
- p38 kinase
- Jak/Stat
- Lck kinase
- p38 kinase
- JNK

Asthma/Allergy

- H1-Receptor
- β 2 Adrenergic
- cAMP/PDE
- PAF
- IL-4, IL-13

Cardiovascular

- ACE
- β 2 Adrenergic
- Ca++ Channel
- Cholesterol
- Hyperlipidemia
- Restenosis
- Antioxidant

Metabolism

- PPAR γ
- GR
- LXR
- FXR
- Estrogen receptor
- Androgen receptor
- HMG-CoA reductase
- AMPK
- GSK3 β



Cancer

- HDAC
- Hsp90
- Proteasome
- NF κ B
- PI-3K/AKT
- Mek
- CDK
- RAR/RXR
- Stat
- Ras/MAPK
- TGF β
- Microtubule
- Jak/Stat
- Tie2 R
- TGF β

....and many others (>2000 drugs/compounds in reference database)


What Makes BioMAP Systems Different

Why Primary Cells?




What Makes BioMAP Systems Different

Why Primary Cells?




**Primary cells retain regulatory networks
and feedback mechanisms**



What Makes BioMAP Systems Different

Why Activate Multiple Pathways?


One Pathway Active




What Makes BioMAP Systems Different

Why Activate Multiple Pathways?

One Pathway Active




Expression Level




What Makes BioMAP Systems Different

Why Activate Multiple Pathways?

Multiple Pathways Active




Expression Level




What Makes BioMAP Systems Different

Why Activate Multiple Pathways?



Multiple Pathways Active



Expression Level



Complex Assays Capture “System” Properties





- Constant level of expression over 2 orders of magnitude
- Complex assays demonstrate *robustness*

Is VCAM Special?


- VCAM is an adhesion molecule and ligand for the T cell integrin $\alpha 4\beta 1$

$\alpha 4\beta 1$ is the target of the multiple sclerosis drug, Tysabri



Example BioMAP Profile


Reference p38 Inhibitor



- BioMAP activity profiles are robust and reproducible
- Profiles retain shape over multiple concentrations

Example BioMAP Profile


Reference p38 Inhibitor



- No inhibition of VCAM under “3C” conditions
(IL-1 β + TNF- α + IFN γ)

Example BioMAP Profile


Reference p38 Inhibitor



- VCAM is differentially regulated
- Regulation of VCAM is “context dependent”


Classification of Compounds

Pairwise Correlation Analysis of BioMAP Profiles




BioMAP Classifies Drugs By Mechanism of Action

Multidimensional Scaling - Function Similarity Map




BioMAP Classifies Drugs By Mechanism of Action

Multidimensional Scaling - Function Similarity Map



BioMAP for Structure Activity Relationships

p38 kinase Inhibitors - Correlation with Biochemical Data




Kinase Profiling
227 Kinases
(Amit)

Amit Gene Symbol	SB202190	PD169316	SB220025	SB239063
ERK1	No hit	No hit	2.85	No hit
ERK2	No hit	No hit	2.1	No hit
JNK1	No hit	No hit	No hit	No hit
JNK2	0	0.15	1.65	0.1
JNK3	0.3	0.65	0.85	1.3
MAP4K4	No hit	11.2	7.05	5.1
MAP4K5	No hit	No hit	6.95	1.1
NLK	2.15	2.5	3.15	8.6
p38-alpha	0	0	0	0
p38-beta	0	0	0.4	0.3
p38-gamma	No hit	11.3	7.65	11.3
TGFRB1	15.35	5.8	10.1	No hit

- PD169316, 3.333 uM
- SB202190, 5.556 uM
- SB220025, 3.333 uM
- SB239063, 3.333 uM


BioMAP Classifies Toxic Agents By Mechanism

Pairwise Correlation Analysis




BioMAP Classifies Toxic Agents By Mechanism

Multidimensional Scaling - Function Similarity Map




Retinoic Acid



- **Cell-type selective cytotoxicity**
- **Unexpected fragility**

Terconazole



- Profile shows similarity to PI-3 kinase inhibitors
- Cytochrome P450 products have been shown to regulate fibroblast proliferation via PI-3 kinase (Nieves, 2006, J. Lipid Res., 47:2681)

Summary

- Addressing biological complexity in assay design is important for understanding compound mechanisms
 - Much of the regulatory architecture is hidden
 - “Complexity for simplicity”
 - We know much less than we think we know
- Complex human cell systems:
 - Practical format for evaluating drug / compound mechanisms
 - Broad range of target mechanisms
 - Useful platform for integrating molecular and clinical information

Acknowledgements

BioSeek

- Eric J. Kunkel
- Jonathan Bingham
- David Patterson
- Vangelis Hytopoulos
- Jian Yang

Stanford University

- Rob Tibshirani, Department of Statistics
- Trevor Hastie, Department of Statistics
- Paul Wender, Chemistry Department