

Microbial Source Tracking

Orin C. Shanks

Geneticist

US EPA Region 5

Genomics Training Workshop

April 28, 2005

EPA/ORD Source Tracking Group

**National Risk Management Research Laboratory
Water Supply & Water Resource Division
Microbial Contaminants Control Branch**

Cincinnati Team:

Jorge Santo Domingo
Orin Shanks
Jingrang Lu
Catherine Kelty
Mark Rodgers
Gina Lamendella

Collaborators:

Marirosa Molina (Athens)
Rich Haugland (Cincinnati)

Alumni:

Donald Reasoner
Joyce Simpson
Michelle Bonkoski
Jay Hua
Jean Tang

Microbial Source Tracking Presentation Outline

I. Overview

II. EPA Guide Document

III. Current Research



National Water Quality Inventory

2000 Report

Sample Area:

39% rivers/streams (269K miles)

45% lakes/ponds (7.7 million acres)


51% estuaries (15K square miles)

Rivers and Streams:


Most common biological contaminant

13% bacterial pollution

35% of reported problems




Protecting America's Public Health **PREVENTION** with **RISK ASSESSMENT**



Monitoring Fecal Pollution

Microbial “Fecal Indicators”

- Represents fecal pollution event
- Bacteria from animal intestine




Traditional Methods

- Presence/absence
- Count per unit volume

Feces Production in the U.S.

Human (0.7%)



1×10^{12} kg/year

Wildlife Contributions



Microbial Source Tracking

CONCEPT.... Match microbe from a polluted site and an animal source to suggest the origin of fecal pollution.

When are Microbial Source Tracking Methods Useful?

To supplement sanitary surveys:

- Identify sources of beach contaminants
- Identify sources of TMDL violations


For risk analyses:

- Human versus non-human
- Human versus domestic animal

Why Should Microbial Source Tracking Work?

**Intestinal microbes of animal groups are
expected to be different:**

- **Gut conditions**
 - Temperature
 - Diet
 - Digestive system
- **Natural selection**
 - Space
 - Nutrients



“Source Identifiers”

Definition . . . microbial populations that are particular to a specific animal host

Ideal Candidates:

- **Exhibit host-specificity**
- **Abundant in host**
- **Temporal stability**
- **Geographic continuity**

Microbial Source Tracking Method Classifications

- **Qualitative vs. Quantitative**
- **Phenotypic vs. Genotypic**
- **Library-dependent vs. Library-independent**

Published:

**Simpson, J. M., J. W. Santo Domingo, and D. J. Reasoner. 2002.
Microbial Source Tracking – State of Science. Environ. Sci. & Tech.
36:5279-5288.**

Library Dependent Methods

- **Library = “Fingerprint” database of *E. coli* or fecal enterococci isolates**
- **Requires 1,000s of isolates from water and suspected animal sources**
- **CULTURE-DEPENDENT**


Library Dependent Methods

- **ARA** (antibiotic resistance analysis)
- **CUP** (carbon utilization profiles)
- **PFGE** (pulse field gel electrophoresis)
- **RFLP** (restriction fragment length polymorphism)
- **AFLP** (amplified fragment length polymorphism)
- **RAPD** (random amplified polymorphic DNA)
- **rep-PCR** (repetitive extragenic palindromic)
- **Ribotyping** (RFLP using rDNA probes)

Library Dependent Method Logistics

METHOD	Targets tested	Cultivation	Major Costs	Time Required*
ARA	<ul style="list-style-type: none"> •<i>Escherichia coli</i> •Fecal streptococci •<i>Enterococcus spp.</i> 	<ul style="list-style-type: none"> •Individual •Isolates 	<ul style="list-style-type: none"> •Antibiotics •96-well microplates 	•4-5 days
CUP	<ul style="list-style-type: none"> •<i>Escherichia coli</i> •Fecal streptococci •<i>Enterococcus spp.</i> 	<ul style="list-style-type: none"> •Individual •Isolates 	<ul style="list-style-type: none"> •Microplates with substrates (e.g., Biolog, Phene Plate) 	•2-5 days
rep-PCR	<ul style="list-style-type: none"> •<i>Escherichia coli</i> 	<ul style="list-style-type: none"> •Individual •Isolates 	<ul style="list-style-type: none"> •PCR reagents •PCR disposable •Gel electrophoresis 	•1 day
RAPD	<ul style="list-style-type: none"> •<i>Escherichia coli</i> 	<ul style="list-style-type: none"> •Individual •Isolates 	<ul style="list-style-type: none"> •PCR reagents •PCR disposable •Gel electrophoresis reagents 	•1 day
AFLP	<ul style="list-style-type: none"> •<i>Escherichia coli</i> 	<ul style="list-style-type: none"> •Individual •Isolates 	<ul style="list-style-type: none"> •DNA extraction kit •AFLP kit (\$5 per reaction) 	•5 days
PFGE	<ul style="list-style-type: none"> •<i>Escherichia coli</i> •<i>Enterococcus spp.</i> 	<ul style="list-style-type: none"> •Individual •Isolates 	<ul style="list-style-type: none"> •Plug prep. reagents •Restriction enzymes •Gel electrophoresis reagents 	•2-4 days
Ribotyping	<ul style="list-style-type: none"> •<i>Escherichia coli</i> •Fecal streptococci •<i>Enterococcus spp.</i> 	<ul style="list-style-type: none"> •Individual •Isolates 	<ul style="list-style-type: none"> •DNA purification reagents •Gel electrophoresis reagents •Restriction enzymes •Hybridization/ detection solutions •Labeled gene probe 	•1-3 days

Library Dependent Method: Antibiotic Resistance Analysis



Advantages and Disadvantages of ARA

Advantages

- Easy to type
- Easy to perform
- Easy to interpret
- Inexpensive


Disadvantages

- Transferable trait
- Geographic specific
- Temporal specific
- Culture dependent
- Breaks down in complex watersheds


Library Dependent Method: rep-PCR DNA Fingerprint Patterns

(Dombek et al., 2000)

From this ...



... to this



Advantages and Disadvantages of rep-PCR

Advantages

- Easy to type
- Easy to perform
- Easy to interpret
- Highly reproducible

Disadvantages

- Library dependent
- May be geographic specific
- May be temporal specific
- Culture dependent

Library Independent Methods

- **Phage typing (serotypic or genotypic)**
- **Gene specific PCR**
- **Total Community Analysis**
- **Host-specific PCR**

Library Independent Method Logistics

METHOD	Targets tested	Cultivation	Major Costs	Time Required*
Phage Typing	<ul style="list-style-type: none"> • F+ coliphage 	<ul style="list-style-type: none"> • Individual • Isolates 	<ul style="list-style-type: none"> • Hybridization/ detection solutions • Labeled gene probe • Phage specific antigen 	<ul style="list-style-type: none"> • 1-3 days
Gene Specific PCR	<ul style="list-style-type: none"> • <i>E. coli</i> toxins 	<ul style="list-style-type: none"> • Sample Enrichment 	<ul style="list-style-type: none"> • PCR reagents • PCR disposables 	<ul style="list-style-type: none"> • 2 days
Total Community Analysis	<ul style="list-style-type: none"> • 16S rRNA 	<ul style="list-style-type: none"> • None 	<ul style="list-style-type: none"> • Filtration units • PCR reagents • PCR disposables • DNA sequencing 	<ul style="list-style-type: none"> • 1 month
Host Specific PCR	<ul style="list-style-type: none"> • <i>Bacteroides</i> • <i>Bifidobacteria</i> • <i>Enterococcus</i> • <i>Rhodococcus</i> • <i>F+ coliphage</i> • <i>Enterovirus</i> • <i>Adenovirus</i> 	<ul style="list-style-type: none"> • None 	<ul style="list-style-type: none"> • Filtration units • PCR reagents • PCR disposable 	<ul style="list-style-type: none"> • 6-8 hours

Host-Specific PCR: *Bacteroides* 16S rDNA

Bernard & Field, (2000) AEM 66, 4571-4574

Dick et al, (2005) AEM in press


- **Primer sets that discriminate between human, ruminant, horse, and pig fecal pollution**
- **Target 16S rDNA from fecal *Bacteroides***
- **Successful in fresh and marine waters**

Why *Bacteroides* as a Source Identifier?


- **Only found in feces, rumen, and body cavities**
- **1/3 of fecal flora**
- **Obligate anaerobes**
- **Limited survival in environment**
- **Host-specific variation in animal hosts**

Comparative Sequence Analysis of 16S rRNA

16S rRNA of
Bacteroides moolii



16S rRNA of
Bacteroides horseii



	650	660	670
<i>B. moolii</i>	GCUUGAGUCU	CGUAGAGGGG	GGUAGAAUUC
<i>B. horseii</i>	GCUAGAGUAU	GGGAGAGGAU	GGUAGAAUUC

Advantages of Host-Specific PCR

- Culture independent
- No library required
- Rapid
- Sensitive
- Defined target
- Isolate target in a complex environment
- Automated analysis

Current Limitations of Host-Specific PCR

- PCR inhibition
- Targets only one gene
- Targets only one bacterial group
- Targets are found in low numbers
- Limited number of case studies
- Small target sequence databases
- Current targeted genes have little to do with host/microbe interactions

PART II

US EPA Microbial Source Tracking Guide Document

Office of Research and Development

National Risk Management Research Laboratory

Water Supply & Water Resource Division

Microbial Contaminants Control Branch

Why is a Guide Document Needed?

Recent proliferation of new methods

- Genotypic
- Phenotypic
- Culture-based
- Culture-independent
- Different levels of discrimination

Most useful method depends on circumstances

Content of Microbial Source Tracking Guide Document

I. Introduction

- What is Microbial Source Tracking?
- Definitions of terms

II. Decision Criteria

- When methods should be used
- Importance of sanitary surveys
- Decision tree

Decision Tree

Questions:

- Is the problem adequately defined?
- Has an adequate sanitary survey been conducted?
- How many sources were identified?
- Is the study area of manageable size?
- What is the desired level discrimination?

Content of Microbial Source Tracking Guide Document

III. MST Approaches

- Summary of all current methods
 - *Explanations of how they work*
 - *Summary tables with advantages and disadvantages*
 - *References*

Content of Microbial Source Tracking Guide Document

IV. Data Collection and Analysis

- Design sampling around study objectives
- General principles for sampling
- Library construction and validation
- Spatial and temporal variability
- Similarity measurement methods

Content of Microbial Source Tracking Guide Document

V. Performance Standards

- Universal quality measures
- Method-specific controls
- Method-specific performance criteria

Content of Microbial Source Tracking Guide Document

VI. Assumptions and Limitations


- Characteristics of source identifiers
 - Specificity
 - Distribution in host
 - Geographic range
 - Temporal stability
 - Survival in water

Content of Microbial Source Tracking Guide Document

VII. Applications of Microbial Source Tracking Approaches

- Eight case studies are presented
- A glossary of terms is presented

The Guide Document is Now Available!



Contact Information
Jorge W. Santo Domingo
santodomingo@epa.gov

Office of Research and Development
National Risk Management Research
Laboratory
Water Supply & Water Resource Division
Microbial Contaminants Control Branch

PART III

Current Research

EPA/ORD Source Tracking Group

**Office of Research and Development
National Risk Management Research Laboratory
Water Supply & Water Resource Division
Microbial Contaminants Control Branch**


Current Research

- **Expand Library of 16S rDNA Sequences from Fecal Sources**
- **Validation of *Bacteroides* 16S rDNA Host-Specific PCR Method**
- **Evaluation of Best Management Practices**
- **Discovery of Novel Source Identifiers**

USEPA 16S rDNA Sequence Fecal Microbe Library

- **Fecal Sources (n = 300)**
 - **Domestic animals**
 - **Wildlife**
 - **Humans**


- **16S rDNA sequences**
 - **Bacteroidales (n = 1,000)**
 - **Clostridium (n = 500)**
 - **Enterococci (n = 1,500)**
 - **Bifidobacterium (n = 100)**



Host-Specific 16S rDNA PCR Method Validation: Target Specificity

- **Test host-specific primer sets against fecal library**
 - **Ruminant-specific**
 - **Human-specific**
 - **Pig-specific**
 - **Horse-specific**
- **If cross-specificity observed, then**
 - **Sequence 16S rRNA**
 - **Add to database**


Host-Specific 16S rDNA PCR Method Validation: Spatial Stability



Delaware Project: Best Management Practice Evaluation


(Collaboration with Delaware Department of Natural Resources and Environmental Control)

Method Overview




Progress


- Cows, sediment, & water
- 700 16S rRNA sequences
- Fences installed (Spring 2004)
- Fecal, sediment, & water collection (Current)






Discovery of Novel Source Identifiers



One gene – one group
PHYLOGENETIC Approach




Multiple genes – one bacterial group
GENOMICS Approach



Multiple genes – Multiple bacterial groups
METAGENOMICS Approach



Use of Phylogenetics to Design New Host-Specific 16S rDNA PCR Assays

- Survey evolutionary relationship between *Bacteroides* from different sources
- Visualize relationships with trees
- Sequences from same source can cluster together
- Approach led to a horse-specific PCR assay
(in press, AEM June 2005)



Use of Genomics to Identify Species-Specific DNA Sequences


- Standard for measuring fecal pollution
- Opportunistic pathogens
- Enterococci genomes already sequenced
(*E. faecalis* and *E. faecium*)



Multiple genes – one bacterial group
GENOMICS

Use of Metagenomics to Design Cow-Specific PCR Assays


- Comparison of *Bos taurus* and *Sus scrofa* genome communities
- Access to yet to be cultured microorganism genomes
- Identification of non-16S rRNA host-specific DNA targets



Multiple genes – Multiple bacterial groups
METAGENOMICS

The Future of Source Tracking

- New methods will arise
- Some methods will become obsolete



Indicators

Pathogens

MST

USEPA ORD Support for Microbial Source Tracking

- Guide Document
- Regional Workshops
- Collaboration
- Development of Regional Centers of Excellence

