


Standardized Analytical Methods for Environmental Restoration following Homeland Security Events

Revision 3.1

November 15, 2007



Standardized Analytical Methods for Environmental Restoration following Homeland Security Events

Revision 3.1

November 15, 2007

Prepared by

Computer Sciences Corporation
Alexandria, VA 22304-3540

Prepared under

EPA Contract No. EP-W-06-046

Prepared for

National Homeland Security Research Center
United States Environmental Protection Agency
Office of Research and Development
Cincinnati, OH 45268

Disclaimer

The U.S. Environmental Protection Agency through its Office of Research and Development funded and managed the research described here under Contract EP-W-06-046 with Computer Sciences Corporation (CSC). This document has been subjected to EPA's peer and administrative review and has been approved for publication as an EPA document.

Mention of trade names or commercial products in this document or in the methods referenced in this document does not constitute endorsement or recommendation for use.

Questions concerning this document or its application should be addressed to:

Oba Vincent
National Homeland Security Research Center
Office of Research and Development (163)
U.S. Environmental Protection Agency
26 West Martin Luther King Drive
Cincinnati, OH 45268
(513) 569-7456
vincent.oba@epa.gov

Use of This Document


The information contained in this document represents the latest step in an ongoing National Homeland Security Research Center effort to provide standardized analytical methods for use by laboratories (e.g., EPA contract laboratories) tasked with performing confirmatory analyses of environmental samples following a homeland security incident. Although at this time, some of the methods listed have not been fully validated for a particular analyte (e.g., analytes not explicitly identified in the method) or sample type, the methods are considered to contain the most appropriate currently available techniques. Unless a published method that is listed in this document states specific applicability to the analyte/sample type for which it has been selected, it should be assumed that method testing is needed, and adjustments may be required to accurately account for variations in analyte/sample type characteristics, environmental samples, and target risk levels. Many of the target analytes listed in this document have only recently become an environmental concern. EPA is actively pursuing development and validation of Standard Analytical Protocols (SAPs) based on the methods listed, including optimization of procedures for measuring target compounds. In those cases where method procedures are determined to be insufficient for a particular situation, EPA will provide guidance regarding appropriate actions. This will be an ongoing process as EPA will strive to establish a consistent level of validation for all listed analytes.

Foreword

The U.S. Environmental Protection Agency (EPA) is charged by Congress with protecting the Nation's land, air, and water resources. Under a mandate of national environmental laws, EPA strives to formulate and implement actions leading to a compatible balance between human activities and the ability of natural systems to support and nurture life. To meet this mandate, EPA's research program is providing data and technical support for solving environmental problems today and building a scientific base necessary to manage our ecological resources wisely, understand how pollutants affect our health, and prevent or reduce environmental risks in the future.

The National Homeland Security Research Center (NHSRC) is EPA's center for conducting research to facilitate protection and decontamination of structures and water infrastructure subject to chemical, biological, or radiological (CBR) terror attacks. NHSRC's research is designed to provide appropriate, effective, and validated technologies, methods, and guidance to understand the risks posed by CBR agents and to enhance our ability to detect, contain, and clean up in the event of an incident involving such agents. This document is intended to provide guidance for selecting methods that have a high likelihood of ensuring analytical consistency when laboratories are faced with a large scale environmental restoration crisis. At the same time, the document can be used as a tool to identify analytes that require further methods development and verification to ensure desired performance. NHSRC will also provide direct technical assistance to response personnel in the event of a CBR attack, as well as provide related interagency liaisons.

This publication has been produced as part of the NHSRC's long-term research plan. It is published and made available by EPA's Office of Research and Development to assist the user community and to link researchers with their clients.

A handwritten signature in dark ink, reading "Jonathan G. Herrmann". The signature is fluid and cursive, with the first letters of the first and last names being capitalized and prominent.

Jonathan G. Herrmann, Director
National Homeland Security Research Center

Abbreviations and Acronyms

AEM	<i>Applied and Environmental Microbiology</i>
AGI sampler	All Glass Impinger Sampler
2-Am-DNT	2-Amino-4,6-dinitrotoluene
4-Am-DNT	4-Amino-2,6-dinitrotoluene
Amp-ELISA	Amplified-Enzyme-Linked Immunosorbent Assay
APHA	American Public Health Association
AOAC	AOAC International (formerly the Association of Official Analytical Chemists)
ASTM	ASTM International (formerly the American Society for Testing and Materials)
AWWA	American Water Works Association
BCM	Biosynth Chromogen Medium
BGMK	Buffalo Green Monkey Kidney
BS	Bismuth Sulfite
BLEB	Buffered <i>Listeria</i> Enrichment Broth
BSA	Bovine Serum Albumin
BSL	Biosafety Level
°C	Degrees Celsius
Campy-BAC	<i>Campylobacter-Brucella</i> agar base with sheep blood and antibiotics
CAS RN	Chemical Abstracts Service Registry Number
CBR	Chemical, Biological, or Radiological
CCID	Coordinating Center for Infectious Diseases
CDC	Centers for Disease Control and Prevention
CFR	Code of Federal Regulations
CFSAN	Center for Food Safety and Applied Nutrition
CFU	Colony Forming Unit
CLP	Contract Laboratory Program
CoV	Coronavirus
CPE	Cytopathic Effect
cps	counts per second
CS	Tear gas; Chlorobenzylidene malonitrile
CVAA	Cold Vapor Atomic Absorption or 2-Chlorovinylarsonous acid
CVAFS	Cold Vapor Atomic Fluorescence Spectrometry
DAPI	4',6-diamidino-2-phenylindole
DAS	Diacetoxyscipenol
DAS-HG-HAS	Diacetoxyscipenol Hemiglutarate Human Serum Albumin
DAS-HS-HRP	Diacetoxyscipenol Hemisuccinate Horseradish Peroxidase Conjugate
DBPR	Division of Bioterrorism Preparedness and Response
DHS	U.S. Department of Homeland Security
DIC	Differential Interference Contrast
DIG-ELISA	Digoxigenin Labeled Enzyme-Linked Immunosorbent Assay
DIMP	Diisopropyl methylphosphonate
DNA	Deoxyribonucleic Acid
3,5-DNA	3,5-Dinitroaniline
1,3-DNB	1,3-Dinitrobenzene
2,4-DNPH	2,4-Dinitrophenylhydrazine
2,4-DNT	2,4-Dinitrotoluene
2,6-DNT	2,6-Dinitrotoluene
DoD	U.S. Department of Defense
DOE	U.S. Department of Energy
DOT	U.S. Department of Transportation
DPD	N, N-diethyl- <i>p</i> -phenylenediamine
EA2192	Diisopropylaminoethyl methylthiophosphonate

ECD	Electron Capture Detector
ECL	Electrochemiluminescence
ED	Electron Diffraction or Ethyldichloroarsine
EDEA	N-Ethyldiethanolamine
EDL	Estimated Detection Limit
EDTA	Ethylenediaminetetraacetic acid
EDXA	Energy Dispersive X-ray Analysis
EEB	EHEC Enrichment Broth
EHEC	Enterohemorrhagic <i>Escherichia coli</i>
EIA	Enzyme Immunoassay
ELCD	Electrolytic Conductivity Detector
ELISA	Enzyme-Linked Immunosorbent Assay
EMC	Emission Measurement Center
EMJH	Ellinghausen-McCullough Johnson Harris Formulation
EML	Environmental Measurements Laboratory
EMMI	Environmental Monitoring Methods Index
EMPA	Ethyl methylphosphonic acid
EMSL	Environmental Monitoring and Support Laboratory
EPA	U.S. Environmental Protection Agency
EQL	Estimated Quantitation Limit
ESI	Electrospray Ionization
ETV	Environmental Technology Verification
FA	Fluorescence Assay
FBI	Federal Bureau of Investigation
FDA	Food and Drug Administration
FEMS	Federation of European Microbiological Societies
FGI	Fluorescein derivative of <i>Conus geographus</i> α -conotoxin
FID	Flame Ionization Detector
FITC	Fluorescein Isothiocyanate
FRET	Fluorescence Resonance Energy Transfer
FSIS	Food Safety and Inspection Service
Furadan	Carbofuran
GA	Tabun
GC	Gas Chromatograph or Gas Chromatography
GC-MS	Gas Chromatograph/Mass Spectrometer or Gas Chromatography/Mass Spectrometry
GD	Soman
GE	1-Methylethyl ester ethylphosphonofluoridic acid or Genome Equivalent
Ge(Li)	Germanium (Lithium)
GF	Cyclohexyl sarin
GFAA	Graphite Furnace Atomic Absorption Spectrophotometer or Graphite Furnace Atomic Absorption Spectrophotometry
GITC	Guanidinium Isothiocyanate
GTC	Guanidinium Thiocyanate
HA	Hemagglutinin
HAdV	Human Adenoviruses
HASL	Health and Safety Laboratory, currently known as Environmental Measurements Laboratory (EML)
HAV	Hepatitis A Virus
HCoV	Human Coronavirus
HCT	Human Ileocecal Adenocarcinoma
HEV	Hepatitis E Virus
HECD	Hall Electrolytic Conductivity Detector
HEPA	High Efficiency Particulate Air (Filter)

HD	Sulfur mustard / mustard gas; bis(2-chloroethyl) sulfide
HHA	Hand-held Assay (e.g., immunochromatographic test device)
HHS	Health and Human Services
HMTD	Hexamethylenetriperoxidediamine
HMX	Octahydro-1,3,5,7-tetranitro-1,3,5,7-tetrazocine
HN-1	Nitrogen mustard 1; bis(2-chloroethyl)ethylamine
HN-2	Nitrogen mustard 2; 2,2'-dichloro-N-methyldiethylamine N,N-bis(2-chloroethyl)methylamine
HN-3	Nitrogen mustard 3; tris(2-chloroethyl)amine
HPGe	High Purity Germanium
HPLC	High Performance Liquid Chromatograph or High Performance Liquid Chromatography
HPLC-FL	High Performance Liquid Chromatograph – Fluorescence
HPLC-MS	High Performance Liquid Chromatograph – Mass Spectrometer
HPLC-PDA	High Performance Liquid Chromatography – Photodiode Array
HPLC-vis	High Performance Liquid Chromatography – visible
IC	Ion Chromatograph or Ion Chromatography
ICC	Integrated Cell Culture
ICC/RT-PCR	Integrated Cell Culture/Reverse Transcriptase Polymerase Chain Reaction
ICP	Inductively Coupled Plasma
ICP-AES	Inductively Coupled Plasma – Atomic Emission Spectrometry
ICP-MS	Inductively Coupled Plasma – Mass Spectrometry
ICR	Information Collection Rule or Information Collection Request (depending on context)
ID50	A dose which would be infectious to 50% of the population
IMPA	Isopropyl methylphosphonic acid
IMS	Immunomagnetic Separation
INCHEM	INCHEM is a means of rapid access to internationally peer reviewed information on chemicals commonly used throughout the world, which may also occur as contaminants in the environment and food. It consolidates information from a number of intergovernmental organizations whose goal it is to assist in the sound management of chemicals. http://www.inchem.org/
IO	Inorganic
IRIS	Integrated Risk Information System (U.S. EPA)
ISO	International Organization for Standardization
ISE	Ion Specific Electrode
K-D	Kuderna-Danish
L-1	Lewisite 1; 2-Chlorovinylchloroarsine
L-2	Lewisite 2; bis(2-Chlorovinyl)chloroarsine
L-3	Lewisite 3; tris(2-Chlorovinyl)arsine
LB-M	Lim Benyesh-Melnick
LC	Liquid Chromatograph
LC/APCI-MS	Liquid Chromatography Atmospheric Pressure Chemical Ionization Mass Spectrometry
LC/ESI-MS	Liquid Chromatography Electrospray Ionization Mass Spectrometry
LC-MS-MS	Liquid Chromatography Tandem Mass Spectrometry
LIA	Lysine Iron Agar
LLD	Lower Limit of Detection
LOD	Limit of Detection
LRN	Laboratory Response Network
LSE	Liquid/Solid Extraction
M	Molar
Mab	Monoclonal Antibody
MARLAP	Multi-Agency Radiological Laboratory Analytical Protocols
MDL	Method Detection Limit
MLD	Minimum Lethal Dose

MPA	Methylphosphonic acid
MPN	Most Probable Number
MS	Mass Spectrometer or Mass Spectrometry or Matrix Spike
MSD	Matrix Spike Duplicate
MSRV	Modified Semisolid Rappaport-Vasiliadis
MTBE	Methyl <i>tert</i> -butyl ether
MW	Molecular Weight
NA	Not Applicable
NaI(Tl)	Thallium-Activated Sodium Iodide
NB	Nitrobenzene
NBD chloride	7-Chloro-4-nitrobenzo-2-oxa-1,3-diazol
NCPDCID	National Center for the Prevention, Detection, and Control of Infectious Diseases
NEMI	National Environmental Methods Index
NERL-CI	National Exposure Research Laboratory-Cincinnati
NG	Nitroglycerine
NH ₃ -N	Ammonia Nitrogen
NHSRC	National Homeland Security Research Center
NIOSH	National Institute for Occupational Safety and Health
NIST	National Institute of Standards and Technology
nM	Nanomolar
NNSA	National Nuclear Security Administration
NoV	Norovirus
NPD	Nitrogen-phosphorus Detector
NOS	Not Otherwise Specified
NRC	National Research Council
nS	nano Siemens
2-NT	2-Nitrotoluene
3-NT	3-Nitrotoluene
4-NT	4-Nitrotoluene
NTIS	National Technical Information Service
NTU	Nephelometric Turbidity Units
OAQPS	Office of Air Quality Planning and Standards
OEM	U.S. EPA Office of Emergency Management
ONPG	<i>ortho</i> -Nitrophenyl- β -D-galactopyranoside
ORAU	Oak Ridge Associated Universities
ORIA	U.S. EPA Office of Radiation and Indoor Air
ORISE	Oak Ridge Institute of Science and Education
ORD	U.S. EPA Office of Research and Development
ORF	Open Reading Frame
OSWER	U.S. EPA Office of Solid Waste and Emergency Response
OSHA	Occupational Safety and Health Administration
OW	U.S. EPA Office of Water
OXA	Oxford Medium
PAHs	Polycyclic Aromatic Hydrocarbons
PBS	Phosphate Buffered Saline
PCDDs	Polychlorinated Dibenzo- <i>p</i> -dioxins
PCDFs	Polychlorinated Dibenzofurans
pCi	Picocuries
PCR	Polymerase Chain Reaction
PEL	Permissible Exposure Limit
PETN	Pentaerythritol tetranitrate
PFE	Pressurized Fluid Extraction
PFIB	Perfluoroisobutylene

PID50	50% Pig Infectious Dose
PMPA	Pinacolyl methyl phosphonic acid
PubMed	PubMed is a service of the U.S. National Library of Medicine (www.pubmed.gov), containing citations from scientific journals
PV	Partially Validated
QA	Quality Assurance
QC	Quality Control
qPCR	Qualitative Polymerase Chain Reaction
RDX	Hexahydro-1,3,5-trinitro-1,3,5-triazine
RNA	Ribonucleic Acid
RP-HPLC	Reversed-Phase High Performance Liquid Chromatography
rRNA	Ribosomal Ribonucleic Acid
RTECS	Registry of Toxic Effects of Chemical Substances
RT-PCR	Reverse Transcription-Polymerase Chain Reaction
SAED	Select Area Electron Diffraction
SAM	Standardized Analytical Methods for Environmental Restoration following Homeland Security Events
SAP	Standardized Analytical Protocol
SARS	Severe Acute Respiratory Syndrome
SaV	Sapovirus
SM	<i>Standard Methods for the Examination of Water and Wastewater</i>
SPE	Solid-Phase Extraction
spp.	Species (plural)
STEC	Shiga-toxigenic <i>E. coli</i>
STEL	Short Term Exposure Limit
STX	Saxitoxin
Stx	Shiga Toxin
Stx-1	Shiga Toxin Type 1
Stx-2	Shiga Toxin Type 2
SW	Solid Waste
TATP	Triacetone triperoxide
TBD	To Be Determined
TCBS	Thiosulfate Citrate Bile Salts Sucrose
TC SMAC	Tellurite Cefixime Sorbitol MaConkey Agar
TCLP	Toxicity Characteristic Leaching Procedure
TDG	Thiodiglycol
TEA	Triethanolamine
TEM	Transmission Electron Microscope or Microscopy
Temik	Aldicarb
TETR-PCR	Touchdown Enzyme Time Release-Polymerase Chain Reaction
THF	Tetrahydrofuran
1,3,5-TNB	1,3,5-Trinitrobenzene
2,4,6-TNT	2,4,6-Trinitrotoluene
TOFMS	Time of Flight Mass Spectrometry
TOXNET	Toxicology Data Network
TRF	Time Resolved Fluorescence
TRU	Transuranic
TSB	Tryptic Soy Broth
TTX	Tetrodotoxin
TSAye	Trypticase Soy Agar with yeast extract
TSC	Tryptose-Sulfite Cycloserine
TSI	Triple Sugar Iron
TSP	Thermospray

TTN	Technical Transfer Network
TWA	Total Weighted Average
USDA	U.S. Department of Agriculture
USGS	U.S. Geological Survey
UV	Ultraviolet
VCSB	Voluntary Consensus Standard Body
VE	Phosphonothioic acid, ethyl-, S-(2-(diethylamino)ethyl) O-ethyl ester
VEE	Venezuelan Equine Encephalitis
VM	Phosphonothioic acid, S-(2-(diethylamino)ethyl) O,O-diethyl ester
VOCs	Volatile Organic Compounds
VOA	Volatile Organic Analysis
VX	O-ethyl-S-(2-diisopropylaminoethyl)methylphosphonothiolate
WEF	Water Environment Federation
WHO	World Health Organization
XLD	Xylose Lysine Desoxycholate

Acknowledgments

The contributions of the following persons and organizations to the development of this document are gratefully acknowledged:

United States Environmental Protection Agency

- **Office of Research and Development, National Homeland Security Research Center**
Joan Bursey
Kathy Hall
Romy Lee
Alan Lindquist
Matthew Magnuson
Tonya Nichols
Ben Packard
Eugene Rice
Rob Rothman
Frank Schaefer
Sanjiv Shah
Ramona Sherman
Erin Silvestri
Oba Vincent
Charlena Yoder
- **Office of Research and Development, National Exposure Research Laboratory (NERL), Las Vegas**
Don Betowski (Environmental Sciences Division)
Christian Daughton (Environmental Sciences Division)
Jane Denne (Environmental Sciences Division)
Michael Hiatt (Environmental Sciences Division)
Tammy Jones-Lepp (Environmental Sciences Division)
Stephen Pia (Environmental Sciences Division)
Lee Riddick (Environmental Sciences Division)
Charlita Rosal (Environmental Sciences Division)
Brian Schumacher (Environmental Sciences Division)
Wayne Sovocool (Environmental Sciences Division)
Katrina Varner (Environmental Sciences Division)
John Zimmerman (Environmental Sciences Division)
Kim Rogers (Human Exposure and Atmospheric Sciences Division)
Jeannette Van Emon (Human Exposure and Atmospheric Sciences Division)
- **Office of Research and Development, National Exposure Research Laboratory (NERL), Cincinnati**
James Owens (Microbiological and Chemical Exposure Assessment Research Division)
Gerard Stelma (Microbiological and Chemical Exposure Assessment Research Division)
Ann Grimm (Microbiological and Chemical Exposure Assessment Research Division)
- **Office of Solid Waste and Emergency Response (OSWER)**
Barry Lesnik (Office of Solid Waste)
Michael Johnson (Office of Superfund Remediation and Technology Innovation)
Lawrence Kaelin (Office of Emergency Management, National Decon. Team)
Michael Ottlinger (Office of Emergency Management, National Decon. Team)
Terry Smith (Office of Emergency Management)
Shen Yi-yang (Office of Emergency Management)

- **Office of Air and Radiation, Office of Radiation and Indoor Air**
Lindsey Bender (Radiation Protection Division)
George Dilbeck (National Air and Radiation Environmental Laboratory)
Dennis Farmer (Radiation and Indoor Environments National Laboratory)
John Griggs (National Air and Radiation Environmental Laboratory)
Daniel Mackney (National Air and Radiation Environmental Laboratory)
David Musick (Radiation and Indoor Environments National Laboratory)
- **Office of Water, Office of Ground Water and Drinking Water**
Jessica Pulz (Water Security Division)
- **Office of Prevention, Pesticides, and Toxic Substances, Office of Pesticide Programs**
Yaorong Qian (Biological and Economic Analysis Division)
- **Office of Enforcement and Compliance Assurance, Office of Criminal Enforcement, Forensics and Training**
Larry Strattan (National Enforcement Investigations Center)
- **EPA Regions**
John Curry (Region 3)
Joe Dorsey (Region 3)
Diane Gregg (Region 6)
Adrienne Greenlee (Region 7)
Stephanie Harris (Region 10)
Peggy Knight (Region 10)
Ed O'Neill (Region 6)
Peter Philbrook (Region 1)
Steve Reimer (Region 10)
Dennis Revell (Region 4)
Norman Rodriguez (Region 7)
Sue Warner (Region 3)
Laura Webb (Region 7)
Wayne Whipple (Region 5)
Larry Zintek (Region 5)

United States Centers for Disease Control and Prevention (CDC)

Kevin Ashley (National Institute for Occupational Safety and Health)
Raymond Biagini (National Institute for Occupational Safety and Health)
Erin Black (National Center for Infectious Diseases)
Cheryl Bopp (National Center for Infectious Diseases)
Jay Gee (National Center for Infectious Diseases)
Vince Hill (National Center for Infectious Diseases)
Rudolph Johnson (National Center for Environmental Health)
Richard Kellogg (National Center for Infectious Diseases)
Gene Kennedy (National Institute for Occupational Safety and Health)
Stephen A. Morse (National Center for Infectious Diseases)
John Snawder (National Institute for Occupational Safety and Health)
Richard Wang (National Center for Environmental Health)

United States Department of Agriculture

Jim Trapp
Jim Trout
Paul Zimba

United States Department of Commerce

Peter Moeller (National Oceanic & Atmospheric Administration)

United States Food and Drug Administration

Jennifer Brzezinski

David Craft

Eric Garber

United States Department of Homeland Security

Anna Berne (Environmental Measurements Laboratory)

United States Department of Defense

Johnathan Kiel (U.S. Air Force)

United States Department of the Army

Ronald Swatski (Center for Health Promotion and Preventive Medicine)

Sina Bavari (Medical Research Institute of Infectious Diseases, USAMRIID)

Alan Hewitt (U.S. Army Corps of Engineers)

Mark Poli (Medical Research Institute of Infectious Diseases, USAMRIID)

State Agencies

Rick Bokanyi (Ohio Department of Health)

David Degenhardt (Wisconsin State Lab of Hygiene)

Robert Glowacky (Virginia State Department of General Services)

Ted Haigh (Oregon Department of Environmental Quality)

Rebecca Hoffman (Wisconsin State Lab of Hygiene)

Christopher Retarides (Virginia State Department of General Services)

Tom York (Virginia State Department of General Services)

Sandia National Laboratory

Gary Brown

Hamilton Sundstrand

Richard Trubey

Computer Sciences Corporation

Eric Boring

Danielle Carter

Joan Cuddeback

Melody Jensen

Larry Umbaugh

Standardized Analytical Methods for Environmental Restoration following Homeland Security Events

Revision 3.1

November 15, 2007

Contents

Disclaimer	iii
Use of This Document	iv
Foreword	v
Abbreviations and Acronyms	vii
Acknowledgments	xiii
Section 1.0: Introduction	1
Section 2.0: Background	3
Section 3.0: Scope and Application	7
Section 4.0: Points of Contact	11
Section 5.0: Selected Chemical Methods	13
5.1 General Guidance	14
5.1.1 Standard Operating Procedures for Identifying Chemical Methods	14
5.1.2 General Quality Control (QC) Guidance for Chemical Methods	27
5.1.3 Safety and Waste Management	29
5.2 Method Summaries	29
5.2.1 EPA Method 8: Determination of Sulfuric Acid and Sulfur Dioxide Emissions from Stationary Sources	29
5.2.2 EPA Method 200.8: Determination of Trace Elements in Waters and Wastes by Inductively Coupled Plasma-Mass Spectrometry	30
5.2.3 EPA Method 245.2: Mercury (Automated Cold Vapor Technique)	31
5.2.4 EPA Method 252.2: Osmium (Atomic Absorption, Furnace Technique)	31
5.2.5 EPA Method 300.1, Revision 1.0: Determination of Inorganic Anions in Drinking Water by Ion Chromatography	32
5.2.6 EPA Method 335.4: Determination of Total Cyanide by Semi-Automated Colorimetry	32
5.2.7 EPA Method 350.1: Nitrogen, Ammonia (Colorimetric, Automated Phenate)	33
5.2.8 EPA Method 507: Determination of Nitrogen- and Phosphorus-Containing Pesticides in Water by Gas Chromatography with a Nitrogen-Phosphorus Detector	33
5.2.9 EPA Method 524.2: Measurement of Purgeable Organic Compounds in Water by Capillary Column Gas Chromatography / Mass Spectrometry	34

5.2.10	EPA Method 525.2: Determination of Organic Compounds in Drinking Water by Liquid-Solid Extraction and Capillary Column Gas Chromatography / Mass Spectrometry	34
5.2.11	EPA Method 531.2: Measurement of N-Methylcarbamoyloximes and N-Methylcarbamates in Water by Direct Aqueous Injection HPLC with Postcolumn Derivatization.....	35
5.2.12	EPA Method 549.2: Determination of Diquat and Paraquat in Drinking Water by Liquid-Solid Extraction and High-Performance Liquid Chromatography with Ultraviolet Detection	35
5.2.13	EPA Method 551.1: Determination of Chlorination Disinfection Byproducts, Chlorinated Solvents, and Halogenated Pesticides/Herbicides in Drinking Water by Liquid-Liquid Extraction and Gas Chromatography with Electron-Capture Detection	36
5.2.14	EPA Method 614: The Determination of Organophosphorus Pesticides in Municipal and Industrial Wastewater.....	36
5.2.15	EPA Method 3031 (SW-846): Acid Digestion of Oils for Metals Analysis by Atomic Absorption or ICP Spectrometry	37
5.2.16	EPA Method 3050B (SW-846): Acid Digestion of Sediments, Sludges, and Soils.....	37
5.2.17	EPA Method 3520C (SW-846): Continuous Liquid-Liquid Extraction.....	38
5.2.18	EPA Method 3535A (SW-846): Solid-Phase Extraction	40
5.2.19	EPA Method 3541 (SW-846): Automated Soxhlet Extraction	43
5.2.20	EPA Method 3545A (SW-846): Pressurized Fluid Extraction (PFE).....	44
5.2.21	EPA Method 3571 (SW-846): Extraction of Solid and Aqueous Samples for Chemical Agents	47
5.2.22	EPA Method 3580A (SW-846): Waste Dilution.....	47
5.2.23	EPA Method 3585 (SW-846): Waste Dilution for Volatile Organics.....	49
5.2.24	EPA Method 5030C (SW-846): Purge-and-Trap for Aqueous Samples.....	50
5.2.25	EPA Method 5035A (SW-846): Closed-System Purge-and-Trap and Extraction for Volatile Organics in Soil and Waste Samples.....	51
5.2.26	EPA Method 6010C (SW-846): Inductively Coupled Plasma - Atomic Emission Spectrometry	52
5.2.27	EPA Method 6020A (SW-846): Inductively Coupled Plasma - Mass Spectrometry.....	53
5.2.28	EPA Method 7010 (SW-846): Graphite Furnace Atomic Absorption Spectrophotometry...	54
5.2.29	EPA Method 7470A (SW-846): Mercury in Liquid Wastes (Manual Cold-Vapor Technique)	54
5.2.30	EPA Method 7471B (SW-846): Mercury in Solid or Semisolid Wastes (Manual Cold-Vapor Technique).....	54
5.2.31	EPA Method 7473 (SW-846): Mercury in Solids and Solutions by Thermal Decomposition, Amalgamation, and Atomic Absorption Spectrophotometry	55
5.2.32	EPA Method 7580 (SW-846): White Phosphorus by Solvent Extraction and Gas Chromatography (GC)	55
5.2.33	EPA Method 8015C (SW-846): Nonhalogenated Organics Using GC-FID.....	56
5.2.34	EPA Method 8260C (SW-846): Volatile Organic Compounds by Gas Chromatography-Mass Spectrometry (GC-MS).....	56
5.2.35	EPA Method 8270D (SW-846): Semivolatile Organic Compounds by Gas Chromatography/Mass Spectrometry (GC-MS)	57
5.2.36	EPA Method 8315A (SW-846): Determination of Carbonyl Compounds by High Performance Liquid Chromatography (HPLC).....	60
5.2.37	EPA Method 8316 (SW-846): Acrylamide, Acrylonitrile and Acrolein by High Performance Liquid Chromatography (HPLC)	60
5.2.38	EPA Method 8318A (SW-846): N-Methylcarbamates by High Performance Liquid Chromatography (HPLC).....	61
5.2.39	EPA Method 8321B (SW-846): Solvent-Extractable Nonvolatile Compounds by High Performance Liquid Chromatography-Thermospray-Mass Spectrometry (HPLC-TS-MS) or Ultraviolet (UV) Detection.....	61

5.2.40	EPA Method 8330B (SW-846): Nitroaromatics and Nitramines by High Performance Liquid Chromatography (HPLC)	62
5.2.41	EPA ILM05.3 Cyanide: Analytical Methods for Total Cyanide Analysis.....	63
5.2.42	IO [Inorganic] Compendium Method IO-3.1: Selection, Preparation, and Extraction of Filter Material	64
5.2.43	IO [Inorganic] Compendium Method IO-3.4: Determination of Metals in Ambient Particulate Matter Using Inductively Coupled Plasma (ICP) Spectroscopy	65
5.2.44	IO [Inorganic] Compendium Method IO-3.5: Determination of Metals in Ambient Particulate Matter Using Inductively Coupled Plasma/Mass Spectrometry (ICP-MS)	66
5.2.45	IO [Inorganic] Compendium Method IO-5: Sampling and Analysis for Vapor and Particle Phase Mercury in Ambient Air Utilizing Cold Vapor Atomic Fluorescence Spectrometry (CVAFS)	67
5.2.46	EPA Air Method, Toxic Organics - 10A (TO-10A): Determination of Pesticides and Polychlorinated Biphenyls in Ambient Air Using Low Volume Polyurethane Foam (PUF) Sampling Followed by Gas Chromatographic/Multi-Detector Detection (GC- MD).....	67
5.2.47	EPA Air Method, Toxic Organics - 15 (TO-15): Determination of Volatile Organic Compounds (VOCs) in Air Collected in Specially-Prepared Canisters and Analyzed by Gas Chromatography/Mass Spectrometry (GC-MS)	70
5.2.48	NIOSH Method 1612: Propylene Oxide	71
5.2.49	NIOSH Method 2016: Formaldehyde	71
5.2.50	NIOSH Method 2513: Ethylene Chlorohydrin	71
5.2.51	NIOSH Method 3510: Monomethylhydrazine.....	72
5.2.52	NIOSH Method 5600: Organophosphorus Pesticides.....	72
5.2.53	NIOSH Method 5601: Organonitrogen Pesticides	73
5.2.54	NIOSH Method 6001: Arsine	73
5.2.55	NIOSH Method 6002: Phosphine	74
5.2.56	NIOSH Method 6004: Sulfur Dioxide	74
5.2.57	NIOSH Method 6010: Hydrogen Cyanide.....	74
5.2.58	NIOSH Method 6013: Hydrogen Sulfide.....	75
5.2.59	NIOSH Method 6015: Ammonia	75
5.2.60	NIOSH Method 6402: Phosphorus Trichloride.....	75
5.2.61	NIOSH Method 7903: Acids, Inorganic	76
5.2.62	NIOSH Method 7905: Phosphorus	76
5.2.63	NIOSH Method 7906: Fluorides, Aerosol and Gas	77
5.2.64	NIOSH Method S301-1: Fluoroacetate Anion.....	77
5.2.65	OSHA Method 40: Methylamine	77
5.2.66	OSHA Method 54: Methyl Isocyanate.....	78
5.2.67	OSHA Method 61: Phosgene	78
5.2.68	OSHA Method ID-211: Sodium Azide and Hydrazoic Acid in Workplace Atmospheres ...	79
5.2.69	OSHA Method ID-216SG: Boron Trifluoride (BF ₃)	79
5.2.70	OSHA Method PV2004: Acrylamide	79
5.2.71	OSHA Method PV2103: Chloropicrin.....	80
5.2.72	ASTM Method D5755-03: Standard Test Method for Microvacuum Sampling and Indirect Analysis of Dust by Transmission Electron Microscopy (TEM) for Asbestos Structure Number Surface Loading	80
5.2.73	ASTM Method D6480-99: Standard Test Method for Wipe Sampling of Surfaces, Indirect Preparation, and Analysis for Asbestos Structure Number Concentration by Transmission Electron Microscopy	81
5.2.74	ISO Method 10312:1995: Ambient Air - Determination of Asbestos Fibres - Direct-transfer Transmission Electron Microscopy Method (TEM)	81
5.2.75	Standard Method 4500-NH ₃ B: Nitrogen (Ammonia) Preliminary Distillation Step	82
5.2.76	Standard Method 4500-NH ₃ G: Nitrogen (Ammonia) Automated Phenate Method	82

5.2.77	Standard Method 4500-Cl G: DPD Colorimetric Method	82
5.2.78	Literature Reference for Chlorine (Analyst, 1999. 124: 1853–1857)	83
5.2.79	Literature Reference for Fluoroacetate salts (Analytical Letters, 1994. 27 (14): 2703–2718)	83
5.2.80	Literature Reference for Perfluoroisobutylene (Journal of Chromatography A, 2005. 1098: 156–165)	84
5.2.81	Literature Reference for Sodium Azide (Journal of Forensic Science, 1998. 43(1):200-202)	84
Section 6.0: Selected Radiochemical Methods.....		85
6.1	General Guidance	86
6.1.1	Standard Operating Procedures for Identifying Radiochemical Methods.....	86
6.1.2	General Quality Control (QC) Guidance for Radiochemical Methods	89
6.1.3	Safety and Waste Management	89
6.2	Method Summaries.....	90
6.2.1	EPA Method 111: Determination of Polonium-210 Emissions from Stationary Sources ...	91
6.2.2	EPA Method 900.0: Gross Alpha and Gross Beta Radioactivity in Drinking Water.....	91
6.2.3	EPA Method 901.1: Gamma Emitting Radionuclides in Drinking Water	92
6.2.4	EPA Method 903.0: Alpha-Emitting Radium Isotopes in Drinking Water.....	92
6.2.5	EPA Method 903.1: Radium-226 in Drinking Water - Radon Emanation Technique.....	93
6.2.6	EPA Method 908.0: Uranium in Drinking Water - Radiochemical Method.....	93
6.2.7	EPA Method EMSL-19: Determination of Radium-226 and Radium-228 in Water, Soil, Air and Biological Tissue	94
6.2.8	EPA Method EMSL-33: Isotopic Determination of Plutonium, Uranium, and Thorium in Water, Soil, Air, and Biological Tissue	94
6.2.9	EML HASL-300 Method Am-01-RC: Americium in Soil	95
6.2.10	EML HASL-300 Method Am-02-RC: Americium-241 in Soil-Gamma Spectrometry.....	96
6.2.11	EML HASL-300 Method Am-04-RC: Americium in QAP Water and Air Filters - Eichrom's TRU Resin.....	96
6.2.12	EML HASL-300 Method Ga-01-R: Gamma Radioassay	97
6.2.13	EML HASL-300 Method Po-02-RC: Polonium in Water, Vegetation, Soil, and Air Filters	97
6.2.14	EML HASL-300 Method Pu-12-RC: Plutonium and/or Americium in Soil or Sediments...	98
6.2.15	EML HASL-300 Method Sr-03-RC: Strontium-90 in Environmental Samples	98
6.2.16	FRMAC Method Volume 2, Page 33: Gross Alpha and Beta in Air	99
6.2.17	ORISE Method AP-1: Gross Alpha and Beta for Various Matrices	99
6.2.18	ORISE Method AP-11: Sequential Determination of the Actinides in Environmental Samples Using Total Sample Dissolution and Extraction.....	100
6.2.19	ASTM Method D3084: Standard Practice for Alpha Spectrometry in Water.....	101
6.2.20	ASTM Method D3972: Standard Test Method for Isotopic Uranium in Water by Radiochemistry	101
6.2.21	Standard Method 7110 B: Gross Alpha and Gross Beta Radioactivity (Total, Suspended, and Dissolved).....	102
6.2.22	Standard Method 7120: Gamma-Emitting Radionuclides	103
6.2.23	Standard Method 7500-Ra B: Radium: Precipitation Method	103
6.2.24	Standard Method 7500-Ra C: Radium: Emanation Method	104
6.2.25	Standard Method 7500-Sr B: Total Radioactive Strontium and Strontium-90: Precipitation Method	104
6.2.26	Standard Method 7500-U B: Uranium: Radiochemical Method.....	105
6.2.27	Standard Method 7500-U C: Uranium: Isotopic Method.....	105

Section 7.0: Selected Pathogen Methods	107
7.1 General Guidance	108
7.1.1 Standard Operating Procedures for Identifying Pathogen Methods	108
7.1.2 General Quality Control (QC) Guidance for Pathogen Methods	111
7.1.3 Safety and Waste Management	112
7.2 Method Summaries	113
7.2.1 Laboratory Response Network (LRN)	113
7.2.2 USEPA Manual of Methods for Virology, EPA/600/4-84/013, April 2001	115
7.2.3 USEPA Environmental Regulations and Technology, Control of Pathogens and Vector Attraction in Sewage Sludge EPA/625/R-92/013, July 2003: <i>Baylisascaris procyonis</i>	115
7.2.4 EPA Method 1622: <i>Cryptosporidium</i> in Water by Filtration/IMS/FA	116
7.2.5 EPA Method 1623: <i>Cryptosporidium</i> and <i>Giardia</i> in Water by Filtration/IMS/FA	116
7.2.6 EPA Method 1682: <i>Salmonella</i> spp.	117
7.2.7 U.S. FDA Bacteriological Analytical Manual, Chapter 10, 2003: <i>Listeria monocytogenes</i>	118
7.2.8 Standard Methods 9213 B: <i>Staphylococcus aureus</i>	118
7.2.9 Standard Methods 9260 B: General Qualitative Isolation and Identification Procedures for <i>Salmonella</i>	119
7.2.10 Standard Methods 9260 E: <i>Shigella</i>	119
7.2.11 Standard Methods 9260 F: Pathogenic <i>Escherichia coli</i>	120
7.2.12 Standard Methods 9260 G: <i>Campylobacter jejuni</i>	120
7.2.13 Standard Methods 9260 H: <i>Vibrio cholerae</i>	121
7.2.14 Standard Methods 9260 I: <i>Leptospira</i>	121
7.2.15 Literature Reference for <i>Campylobacter jejuni</i> (Molecular and Cellular Probes. 2006. 20: 269–279)	122
7.2.16 Literature Reference for <i>Chlamydophila psittaci</i> (Journal of Clinical Microbiology. 2000. 38: 1085–1093)	122
7.2.17 Literature Reference for Adenoviruses (Applied and Environmental Microbiology. 2005. 71(6): 3131–3136)	123
7.2.18 Literature Reference for Astroviruses (Canadian Journal of Microbiology. 2004. 50: 269–278)	124
7.2.19 Literature Reference for Noroviruses (Journal of Clinical Microbiology. 2004. 42(10): 4679–4685)	124
7.2.20 Literature Reference for Sapoviruses (Journal of Medical Virology. 2006. 78(10): 1347–1353)	125
7.2.21 Literature Reference for Coronaviruses (SARS) (Journal of Virological Methods. 2004. 122: 29–36)	126
7.2.22 Literature Reference for Hepatitis E Virus (Journal of Virological Methods. 2006. 131(1): 65–71)	126
7.2.23 Literature Reference for Influenza H5N1 (Emerging Infectious Diseases. 2005. 11(8): 1303–1305)	127
7.2.24 Literature Reference for Enteric Viruses (Applied and Environmental Microbiology. 2003. 69(6): 3158–3164)	128
7.2.25 Literature Reference for <i>Cryptosporidium</i> spp. (Applied and Environmental Microbiology. 1999. 65(9): 3936–3941)	128
7.2.26 Literature Reference for <i>Entamoeba histolytica</i> (Journal of Parasitology. 1972. 58(2): 306–310)	129
7.2.27 Literature Reference for <i>Entamoeba histolytica</i> (Journal of Clinical Microbiology. 2005. 43(11): 5491–5497)	129
7.2.28 Literature Reference for <i>Giardia</i> spp. (Transactions of the Royal Society of Tropical Medicine and Hygiene. 1983. 77(4): 487–488)	130

7.2.29	Literature Reference for <i>Toxoplasma gondii</i> (Emerging Infectious Diseases. 2006. 12(2): 326–329)	130
7.2.30	Literature Reference for <i>Toxoplasma gondii</i> (Applied and Environmental Microbiology. 2004. 70(7): 4035–4039).....	131
Section 8.0: Selected Biotoxin Methods		133
8.1	General Guidance	134
8.1.1	Standard Operating Procedures for Identifying Biotoxin Methods.....	134
8.1.2	General Quality Control (QC) Guidance for Biotoxin Methods	137
8.1.3	Safety and Waste Management	138
8.2	Method Summaries.....	139
8.2.1	Laboratory Response Network (LRN)	139
8.2.2	U.S. FDA, Bacteriological Analytical Manual Online, Chapter 17, 2001: Botulinum Neurotoxins	140
8.2.3	U.S. FDA, Bacteriological Analytical Manual Online, Appendix 1, 2001: Rapid Methods for Detecting Foodborne Pathogens	141
8.2.4	AOAC Official Method 991.31: Aflatoxins in Corn, Raw Peanuts, and Peanut Butter	141
8.2.5	AOAC Official Method 993.06: Staphylococcal Enterotoxins in Selected Foods	142
8.2.6	AOAC Official Method 994.08: Aflatoxin in Corn, Almonds, Brazil Nuts, Peanuts, and Pistachio Nuts.....	142
8.2.7	Literature Reference for Abrin (119 th AOAC Annual Meeting & Exposition, 2005, p. 613)	143
8.2.8	Literature Reference for Abrin and Shiga and Shiga-like Toxins (Pharmacological Toxicology. 2001. 88(5): 255-260).....	143
8.2.9	Literature Reference for Abrin and Ricin (Analytical Biochemistry. 2006. 357(2): 200–207)	144
8.2.10	Literature Reference for α -Conotoxin (Biochemical Journal. 1997. 328: 245–250)	144
8.2.11	Literature Reference for α -Conotoxin (Journal of Medicinal Chemistry. 2004. 47(5): 1234–1241)	145
8.2.12	Literature Reference for α -Amanitin, Ricin, T-2 Mycotoxin (Journal of Food Protection. 2005. 68(6): 1294–1301).....	145
8.2.13	Literature Reference for α -Amanitin (Journal of Chromatography. 1991. 563(2): 299–311)	146
8.2.14	Literature Reference for Anatoxin-a (Biomedical Chromatography. 1996. 10: 46–47)	146
8.2.15	Literature Reference for Brevetoxins (Environmental Health Perspectives. 2002. 110(2): 179–185)	147
8.2.16	Literature Reference for Brevetoxins (Toxicon. 2004. 43(4): 455–465)	147
8.2.17	Literature Reference for Cylindrospermopsin (FEMS Microbiology Letters. 2002. 216(2): 159–164)	148
8.2.18	Literature Reference for Diacetoxyscirpenol (DAS) (International Journal of Food Microbiology. 1988. 6(1): 9–17).....	148
8.2.19	Literature Reference for Diacetoxyscirpenol (DAS) and T-2 Mycotoxin (Rapid Communications in Mass Spectrometry. 2006. 20(9): 1422–1428).....	149
8.2.20	Literature Reference for Microcystins (Journal of AOAC International. 2001. 84(4): 1035–1044)	149
8.2.21	Literature Reference for Microcystins (Analyst. 1994. 119(7): 1525–1530).....	150
8.2.22	Literature Reference for Picrotoxin (Journal of Pharmaceutical & Biomedical Analysis. 1989. 7(3): 369–375).....	150
8.2.23	Literature Reference for Saxitoxin (Journal of AOAC International. 1995. 78(2): 528–532)	151

8.2.24 Literature Reference for Shiga and Shiga-like Toxin (Journal of Clinical Microbiology. 2007. 45(10): 3377–3380).....	151
8.2.25 Literature Reference for Tetrodotoxin (Journal of Clinical Laboratory Analysis. 1992. 6(2): 65–72)	152
8.2.26 Literature Reference for Tetrodotoxin (Analytical Biochemistry. 2001. 290(1): 10–17) ...	152
8.2.27 Lateral Flow Immunoassay Kits	153
Section 9.0: Conclusions	155

Appendices

Appendix A: Selected Chemical Methods	A-1
Appendix B: Selected Radiochemical Methods	B-1
Appendix C: Selected Pathogen Methods	C-1
Appendix D: Selected Biotoxin Methods	D-1

Figures

Figure 1-1. Environmental Evaluation Analytical Process Roadmap for Homeland Security Events	2
Figure 2-1. SAM Method Selection Process.....	5

Tables

Table 5-1. Chemical Methods and Corresponding Text Section Numbers.....	14
Table 5-2. Sources of Chemical Methods	27
Table 6-1. Radiochemical Methods and Corresponding Text Section Numbers	86
Table 6-2. Sources of Radiochemical Methods	88
Table 7-1. Pathogen Methods and Corresponding Text Section Numbers	109
Table 7-2. Sources of Pathogen Methods	110
Table 8-1. Biotoxin Methods and Corresponding Text Section Numbers	135
Table 8-2. Sources of Biotoxin Methods	136

Section 1.0: Introduction


In the aftermath of the terrorist attacks of September 11, 2001, and the anthrax attacks in the fall of 2001, federal and state personnel carried out their mission to provide response, recovery, and remediation under trying circumstances, including an unprecedented demand on their capabilities to analyze environmental samples. In reviewing these events, EPA identified several areas where the country could better prepare itself in the event of future terrorist incidents. One of the most important areas identified was the need to improve the nation's laboratory capacity and capability to analyze environmental samples following a homeland security event.

In response, EPA formed the Homeland Security Laboratory Capacity Workgroup to identify and implement opportunities for near-term improvements and to develop recommendations for addressing longer-term laboratory issues. The EPA Homeland Security Laboratory Capacity Workgroup consists of representatives from the Office of Research and Development, Office of Air and Radiation, Office of Water, Office of Solid Waste and Emergency Response, Office of Environmental Information, Office of Pollution Prevention and Toxics, and several EPA regional offices.

A critical area identified by the workgroup was the need for a list of analytical methods to be used by all laboratories when analyzing homeland security event samples and, in particular, when analysis of a large number of samples is required over a short period of time. Having standardized methods would reduce confusion, permit sharing of sample load between laboratories, improve data comparability, simplify the task of outsourcing analytical support to the commercial laboratory sector, and improve the follow-up activities of validating results, evaluating data, and making decisions. To this end, workgroup members formed an Analytical Methods Subteam to address homeland security methods issues.

The Analytical Methods Subteam recognized that widely different analytical methods are required for various phases of environmental sample analyses in support of homeland security preparation and response: (1) ongoing surveillance and monitoring, (2) response and rapid screening for determining whether an event has occurred, (3) preliminary site characterizations to determine the extent and type of contamination, and (4) confirmatory laboratory analyses to plan, implement, and evaluate the effectiveness of site remediation. **Figure 1-1** represents these analytical phases. EPA's *Standardized Analytical Methods for Environmental Restoration following Homeland Security Events* (SAM) provides information for analytical methods to be applied during the "Site Remediation" phase.

Figure 1-1. Environmental Evaluation Analytical Process Roadmap for Homeland Security Events



Section 2.0: Background

SAM identifies a single method per analyte/sample type to ensure a consistent analytical approach across multiple laboratories when analyzing environmental samples following an event. In support of this document, EPA periodically assembles methods experts from within EPA and other federal agencies to review methods and, if necessary, revise the methods listed. Method selection is based on consideration of specific criteria that emphasize method performance and include existing laboratory capabilities, laboratory capacity, method applicability to multiple environmental sample types, and method applicability to multiple SAM analytes. For some analytes, the preferred method is a clear choice; for others, competing criteria make the choice more difficult. Final method selections are based on technical recommendations from the SAM work groups. For analytes where limited laboratory testing/experience exists, such as chemical warfare agents, methods were selected based on their applicability to similar chemicals (e.g., nerve agents and some pesticides). In these cases, laboratory studies to test the ability of the selected method to measure the target analyte(s) are planned. **Figure 2-1** summarizes steps and provides the criteria used during the SAM method selection process. It is important to note that the method selection criteria included in this figure are listed in non-hierarchical order and, in some cases, only a subset of the criteria was considered.

In 2004, EPA's National Homeland Security Research Center (NHSRC) brought together experts from across EPA and its sister agencies to develop a compendium of analytical methods to be used when analyzing environmental samples to address site characterization, remediation and clearance following future homeland security events. Participants included representatives from EPA program offices, EPA regions, EPA national laboratories, Centers for Disease Control and Prevention (CDC), Food and Drug Administration (FDA), Department of Homeland Security (DHS), Federal Bureau of Investigation (FBI), Department of Defense (DoD), Department of Agriculture (USDA), and U.S. Geological Survey (USGS). Methodologies were considered for chemical and biological agents of concern in the types of environmental samples that would be anticipated. The primary objective of this effort was to identify appropriate SAM Analytical Methods Subteam consensus methods that represent a balance between providing existing, documented, determinative techniques and providing consistent and valid analytical results.

A survey of available confirmatory analytical methods for approximately 120 biological and chemical analytes was conducted using existing resources including the following:

- National Environmental Methods Index (NEMI) and NEMI-Chemical, Biological, and Radiological (NEMI-CBR)
- Environmental Monitoring Method Index (EMMI)
- EPA Test Methods Index
- EPA Office of Solid Waste SW-846 Methods
- EPA Microbiology Methods
- National Institute for Occupational Safety and Health (NIOSH) method index
- Occupational Safety and Health Administration (OSHA) method index
- AOAC International
- ASTM International
- International Organization for Standardization (ISO) methods
- Standard Methods for the Examination of Water and Wastewater
- PubMed Literature Database


In September 2004, EPA published *Standardized Analytical Methods for Use During Homeland Security Events, Revision 1.0* (EPA/600/R-04/126), SAM Revision 1.0, which provided a list of analytical and sample preparation methods that were selected for measurement of 82 chemical analytes in aqueous/liquid, solid, oily solid, and air samples, and 27 biological analytes in water, dust, and aerosol

samples. During 2005, SAM was expanded to include radioisotopes, several persistent chemical warfare agent degradation products, a drinking water sample type, methods for determination of the viability of biological organisms, and a separate section for biotoxin analytes. Where necessary, the methods included in SAM Revision 1.0 were updated to reflect more recent or appropriate methodologies. Similar efforts to those used for method selection during development of SAM Revision 1.0 were undertaken to select and include methods for measurement of radioisotopes and chemical warfare agent degradation products in all sample types, for measurement of chemical, biological, and radiochemical analytes in drinking water, and to determine the viability of biological organisms. These additional analytes and the corresponding methods selected were included in SAM Revision 2.0.

During 2006, SAM was revised further to incorporate analytes included on updated federal agency lists, provide additional or more current method listings for target analytes, incorporate explosives into the chemical analytes listing, combine identification and viability methods information for pathogens, and address comments from EPA Science Advisory Board's Homeland Security Advisory Committee¹ to clarify the intended use of the document. These changes were included in SAM Revision 3.0. Following publication of SAM Revision 3.0, NHSRC plans to convene SAM Technical Workgroups at least once per year, to evaluate and, if necessary, update the analytes and methods that are listed in SAM. This version, SAM Revision 3.1, reflects the first of these updates and includes the addition of several chemical analytes, one radionuclide, and one biotoxin, along with corresponding selected methods. NHSRC plans to provide the updated information in a Web-based version of the document that will allow users and other stakeholders to search for specific needs and to submit questions and comments regarding the information.

¹EPA Science Advisory Board's Homeland Security Advisory Committee meeting summary:
http://www.epa.gov/sab/06minutes/hsac_ws_and_sam_01_30-31_06_minutes.pdf

Figure 2-1. SAM Method Selection Process



Section 3.0: Scope and Application

The premise and purpose of this document is to standardize the analytical methods that will be used in cases when multiple laboratories are called on to analyze environmental samples following a homeland security event. The document also is intended as a tool that will be available to assist state and local laboratories in planning for and analyzing environmental samples following a homeland security event. The methods presented in this document should be used to:

- Confirm the identification of contaminants;
- Evaluate the extent of contamination; and
- Evaluate the effectiveness of decontamination.

The list of methods provided is limited to those methods that would be used to determine, to the extent possible within analytical limitations, the presence of chemical, radiochemical, pathogen, and biotoxin analytes of concern and their concentrations in environmental media. The methods include detailed laboratory procedures for confirming the identification of analytes and determining their concentrations in environmental samples. The methods, therefore, are not designed to be used for rapid or immediate response or for conducting an initial evaluation (triage or screening) of suspected material to determine if it poses an immediate danger or should be analyzed in specially designed, highly secure facilities. This document also is not intended to provide information regarding sample collection activities or equipment. Methods for addressing these needs are and will be the subject of other efforts.

Methods are provided in this document as corresponding to specific analyte/sample type combinations that are listed in Appendices A (chemical), B (radiochemical), C (pathogen), and D (biotoxin). Summaries of each method are provided in numerical order by the developing agency, throughout Sections 5.2 (chemical methods), 6.2 (radiochemical methods), 7.2 (pathogen methods), and 8.2 (biotoxin methods).

It is important to note that, in some cases, the methods included in this document have not been fully validated for the analyte/sample type combination(s) for which they have been selected. The information contained in this document represents the latest step in an ongoing National Homeland Security Research Center effort to provide standardized analytical methods for use by laboratories (e.g., EPA contract laboratories) tasked with performing confirmatory analyses on environmental samples following a homeland security incident. Although at this time, some of the methods listed have not been fully validated for a particular analyte (e.g., analytes not explicitly identified in the method) or sample type, the methods are considered to contain the most appropriate currently available techniques. Unless a published method that is listed in this document states specific applicability to the analyte/sample type combination for which it has been selected, it should be assumed that method testing is needed, and adjustments may be required to accurately account for variations in analyte characteristics, environmental samples, and target risk levels. Many of the target analytes listed in this document have only recently become an environmental concern. EPA is actively pursuing development and validation of Standard Analytical Protocols (SAPs) based on the methods listed, including optimization of procedures for measuring target compounds. In those cases where method procedures are determined to be insufficient for a particular situation, EPA will provide guidance regarding appropriate actions. This will be an ongoing process as EPA will strive to establish a consistent level of validation for all listed analytes.

EPA recognizes that specification of a single method may limit laboratory capacity and techniques that may be needed to evaluate difficult samples. In those cases where method procedures are determined to be insufficient for a particular situation, EPA will provide guidance regarding appropriate actions (see list of contacts in Section 4). EPA is developing and validating Standardized Analytical Protocols (SAPs) based on the methods that are listed in this document, where further development and testing are necessary. Once validation is complete, data regarding the resulting method performance and data quality objectives will be available. The SAM document and corresponding SAPs will be reviewed frequently. EPA plans to continue to update the SAM document to address the needs of homeland security, reflect improvements in analytical methodology and new technologies, and incorporate changes in analytes based on needs. EPA also anticipates that addenda may be generated to provide guidance regarding issues that currently are not addressed by this document. Any deviations from the methods referenced in this document should be coordinated with the appropriate point(s) of contact identified in Section 4.

Participants in the chemical, radiochemical, pathogen, and biotoxin work groups, including representatives from the U.S. EPA, CDC, FDA, DHS, FBI, DoD, USDA, and USGS evaluated the suitability of existing methodologies and selected this set of methods for use by EPA laboratories and contract laboratories if called upon in an emergency. EPA recognizes that this advanced selection of such methods may pose potential risks, including the following:

- Selecting technologies that may not be the most cost-effective technologies currently available for addressing the particular situation at hand;
- Selecting methodologies that may not be appropriate for use in responding to a particular emergency because EPA did not anticipate having to analyze for a particular analyte or analyte/sample type combination; and
- Preventing development and adoption of new and better measurement technologies.

To address these potential risks as soon as possible, EPA plans to take several steps. These include the following:

- Developing and specifying measurement quality objectives (i.e., required minimum standards of accuracy (bias and precision) and sensitivity for the analysis of samples that support the data quality needs of the particular stage of the emergency response/recovery process) for all analyte/sample type combinations listed in this document;
- Specifying guidance for ensuring the analytical methods listed provide results that are consistent with and support their intended use;
- Working with other government agencies and the private sector to establish a laboratory network to ensure that laboratories selected to assist EPA and its federal, state, and local partners in responding to homeland security events have the requisite expertise and systems to perform this type of testing; and
- Continuing to work with multiple agencies and stakeholders to update SAM and supporting documents periodically.

EPA recognizes that having data of known and documented quality is critical for public officials to assess accurately the activities that may be needed in remediating a site during and following emergency situations. Data must be of sufficient quality to support decision making. Quality control (QC), however, takes time, and time is often critical in emergency-related activities where there will be tremendous pressure to conduct sampling and analytical operations quickly and efficiently. While reduced levels of QC might be tolerated during the rapid screening stage of emergency response, implementation of analytical methods for risk assessment and site release will require a higher and more appropriate level of QC. Many of the methods listed in this document include QC requirements for collecting and analyzing samples. These QC requirements may or may not be appropriate for addressing emergency response

situations, and may be adjusted as necessary to maximize data and decision quality. Specific QC recommendations for analysis of samples for chemical, radiochemical, pathogen, and biotoxin analytes are provided in each corresponding section of this document (i.e., Sections 5.1.2, 6.1.2, 7.1.2, and 8.1.2, respectively).

Section 4.0: Points of Contact

Questions concerning this document, or the methods identified in this document, should be addressed to the appropriate point(s) of contact identified below. These contacts should be consulted regarding any method deviations or modifications, sample problems or interferences, quality control requirements, or the use of potential alternative methods. As previously indicated, any deviations from the recommended method(s) should be reported immediately to ensure data comparability is maintained when responding to homeland security events.

General	
Oba Vincent - Primary National Homeland Security Research Center U.S. EPA Office of Research and Development (NG16) 26 West Martin Luther King Jr. Drive Cincinnati, OH 45268 (513) 569-7456 vincent.oba@epa.gov	Rob Rothman - Alternate National Homeland Security Research Center U.S. EPA Office of Research and Development (NG16) 26 West Martin Luther King Jr. Drive Cincinnati, OH 45268 (513) 569-7187 rothman.rob@epa.gov
Chemical Methods	
Barry Lesnik - Primary Office of Solid Waste and Emergency Response (5307P) U.S. EPA Headquarters Ariel Rios Building 1200 Pennsylvania Avenue, NW Washington, DC 20460 (703) 308-0476 lesnik.barry@epa.gov	Steve Reimer - Alternate USEPA Region 10 - Manchester Laboratory 7411 Beach Drive East Port Orchard, WA 98366 (360) 871-8718 reimer.steve@epa.gov
Radiochemical Methods	
John Griggs - Primary U.S. EPA Office of Radiation and Indoor Air Environmental Laboratory 540 South Morris Avenue Montgomery, AL 36115-2601 (334) 270-3450 griggs.john@epa.gov	Kathy Hall - Alternate National Homeland Security Research Center U.S. EPA Office of Research and Development (NG16) 26 West Martin Luther King Jr. Drive Cincinnati, OH 45268 (513) 379-5260 hall.kathy@epa.gov
Pathogen Methods	
Alan Lindquist - Primary National Homeland Security Research Center U.S. EPA Office of Research and Development (NG16) 26 West Martin Luther King Jr. Drive Cincinnati, OH 45268 (513) 569-7192 lindquist.alan@epa.gov	Sanjiv Shah - Alternate National Homeland Security Research Center U.S. EPA Office of Research and Development (NG16) 26 West Martin Luther King Jr. Drive Cincinnati, OH 45268 (513) 569-7776 shah.sanjiv@epa.gov
Biotoxins Methods	
Matthew Magnuson National Homeland Security Research Center U.S. EPA Office of Research and Development (NG16) 26 West Martin Luther King Jr. Drive Cincinnati, OH 45268 (513) 569-7321 magnuson.matthew@epa.gov	Sanjiv Shah - Alternate National Homeland Security Research Center U.S. EPA Office of Research and Development (NG16) 26 West Martin Luther King Jr. Drive Cincinnati, OH 45268 (513) 569-7776 shah.sanjiv@epa.gov

Section 5.0: Selected Chemical Methods

A list of methods to be used in analyzing environmental samples for chemical contaminants during remediation activities following a homeland security event is provided in Appendix A. Methods are listed for each analyte and for each sample type that potentially may need to be measured and analyzed when responding to an environmental emergency. Procedures from peer-reviewed journal articles are listed for those analyte-sample type combinations where methods are not available. Once standard procedures are available, the literature references will be replaced.

Please note: This section provides guidance for selecting chemical methods that have a high likelihood of assuring analytical consistency when laboratories are faced with a large scale environmental restoration crisis. Not all methods have been verified for the analyte/sample type combination listed in Appendix A. Please refer to the specified method to identify analyte/sample type combinations that have been verified. Any questions regarding information discussed in this section should be addressed to the appropriate contact(s) listed in Section 4.

Appendix A is sorted alphabetically by analyte and includes the following information:

- **Analyte(s).** The component, contaminant, or constituent of interest.
- **Chemical Abstract Survey Registration Number (CAS RN).** A unique identifier for chemical substances that provides an unambiguous way to identify a chemical or molecular structure when there are many possible systematic, generic, or trivial names.
- **Determinative technique.** An analytical instrument or technique used to determine the quantity and identification of compounds or components in a sample.
- **Method type.** Two method types (sample preparation and determinative) are used to complete sample analysis. In some cases, a single method contains information for both sample preparation and determinative procedures. In most instances, however, two separate methods may need to be used in conjunction.
- **Solid samples.** The recommended method/procedure to identify and measure the analyte of interest in solid phase samples.
- **Non-aqueous liquid/organic solid samples.** The recommended method/procedure to identify and measure the analyte of interest in non-aqueous liquid/organic phase samples. An organic solid sample is a solid that completely dissolves in an organic solvent and leaves no solid residue.
- **Aqueous liquid samples.** The recommended method/procedure to identify and measure the analyte of interest in aqueous liquid phase samples.
- **Drinking water samples.** The recommended method/procedure to identify and measure the analyte of interest in drinking water samples.
- **Air samples.** The recommended method/procedure to identify and measure the analyte of interest in air samples.

5.1 General Guidance

This section provides a general overview of how to identify the appropriate chemical method(s) for a given analyte-sample type combination as well as recommendations for quality control procedures.

For additional information on the properties of the chemicals listed in Appendix A, TOXNET (<http://toxnet.nlm.nih.gov/index.html>), a cluster of databases on toxicology, hazardous chemicals, and related areas maintained by the National Library of Medicine, is an excellent resource. Additional resources include:

- Syracuse Research Corporation's Physprop and Chemfate, part of the Environmental Fate Database supported by EPA. <http://www.syrres.com/esc/databases.htm>
- INCHEM at <http://www.inchem.org/> contains both chemical and toxicity information.
- The Registry of Toxic Effects of Chemical Substances (RTECS) database can be accessed via the National Institute for Occupational Safety and Health (NIOSH) Web site at <http://www.cdc.gov/niosh/rtecs/vz72d288.html#JWIDAW> for toxicity information.
- EPA's Integrated Risk Information System (IRIS): <http://www.epa.gov/iris/> contains toxicity information.
- The Forensic Science and Communications Journal published by the Laboratory Division of the Federal Bureau of Investigation. <http://www.fbi.gov/hq/lab/fsc/current/backissu.htm>
- European Chemicals Bureau Toxicology and Chemical Substances: <http://ecb.jrc.it> and <http://ecb.jrc.it/testing-methods/> containing information regarding European Directive 67/548/EEC and Annex V

Additional research on chemical contaminants is ongoing within EPA, and databases to manage this information are currently under development.

5.1.1 Standard Operating Procedures for Identifying Chemical Methods

To determine the appropriate method that is to be used on an environmental sample, locate the analyte of concern in Appendix A: Chemical Methods under the "Analyte(s)" column. After locating the analyte of concern, continue across the table to identify the appropriate determinative technique (e.g., HPLC, GC-MS), then identify the appropriate sample preparation and determinative method(s) for the sample type of interest (solid, non-aqueous liquid/organic solid, aqueous liquid, drinking water, or air). In some cases, two methods (sample preparation and determinative) are needed to complete sample analysis.

Sections 5.2.1 through 5.2.81 below provide summaries of the sample preparation and determinative methods listed in Appendix A. Once a method has been identified in Appendix A, **Table 5-1** can be used to locate the method summary.

Table 5-1. Chemical Methods and Corresponding Text Section Numbers

Analyte	CASRN	Method	Section
Acrylamide	79-06-1	8316 (EPA SW-846)	5.2.37
Acrylonitrile	107-13-1	PV2004 (OSHA)	5.2.70

Analyte	CASRN	Method	Section
Aldicarb (Temik)	116-06-3	531.2 (EPA OW)	5.2.11
		8318A (EPA SW-846)	5.2.38
		5601 (NIOSH)	5.2.53
Allyl alcohol	107-18-6	3585 (EPA SW-846)	5.2.23
		5030C (EPA SW-846)	5.2.24
		5035A (EPA SW-846)	5.2.25
		8260C (EPA SW-846)	5.2.34
		TO-15 (EPA ORD)	5.2.47
2-Amino-4,6-dinitrotoluene (2-Am-DNT) 4-Amino-2,6-dinitrotoluene (4-Am-DNT) 4-Aminopyridine	35572-78-2 19406-51-0 504-24-5	3535A (EPA SW-846)	5.2.18
		8330B (EPA SW-846)	5.2.40
Ammonia	7664-41-7	350.1 (EPA OW)	5.2.7
		6015 (NIOSH)	5.2.59
		4500-NH ₃ B (SM)	5.2.75
		4500-NH ₃ G (SM)	5.2.76
Ammonium metavanadate Arsenic, Total Arsenic trioxide	7803-55-6	200.8 (EPA OW)	5.2.2
		3031 (EPA SW-846)	5.2.15
		3050B (EPA SW-846)	5.2.16
	7440-38-2	6010C (EPA SW-846)	5.2.26
		6020A (EPA SW-846)	5.2.27
	1327-53-3	IO-3.1 (EPA ORD)	5.2.42
		IO-3.4 (EPA ORD)	5.2.43
		IO-3.5 (EPA ORD)	5.2.44
Arsine	7784-42-1	200.8 (EPA OW)	5.2.2
		3050B (EPA SW-846)	5.2.16
		7010 (EPA SW-846)	5.2.28
		6001 (NIOSH)	5.2.54
Asbestos	1332-21-4	D5755-03 (ASTM)	5.2.72
		D6480-99 (ASTM)	5.2.73
		10312:1995 (ISO)	5.2.74
Boron trifluoride	7637-07-2	ID – 216SG (OSHA)	5.2.68
Brodifacoum Bromadiolone	56073-10-0	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
	28772-56-7	3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8321 (EPA SW-846)	5.2.39
Carbofuran (Furadan)	1563-66-2	531.2 (EPA OW)	5.2.11
		8318A (EPA SW-846)	5.2.38
		5601 (NIOSH)	5.2.53
Carbon disulfide	75-15-0	524.2 (EPA OW)	5.2.9
		3585 (EPA SW-846)	5.2.23
		5030C (EPA SW-846)	5.2.24
		5035A (EPA SW-846)	5.2.25
		8260C (EPA SW-846)	5.2.34
		TO-15 (EPA ORD)	5.2.47

Analyte	CASRN	Method	Section
Chlorfenvinphos Chlorpyrifos	470-90-6 2921-88-2	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Chlorine	7782-50-5	4500-Cl G (SM)	5.2.77
		Analyst, Vol. 124, 1999, pp. 1853-1857	5.2.78
2-Chloroethanol	107-07-3	3585 (EPA SW-846)	5.2.23
		5030C (EPA SW-846)	5.2.24
		5035A (EPA SW-846)	5.2.25
		8260C (EPA SW-846)	5.2.34
		2513 (NIOSH)	5.2.50
3-Chloro-1,2-propanediol	96-24-2	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Chloropicrin	76-06-2	551 (EPA OW)	5.2.13
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		PV2103 (OSHA)	5.2.71
Chlorosarin Chlorosoman	1445-76-7 7040-57-5	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
2-Chlorovinylarsonous acid (CVAA)	85090-33-1	200.8 (EPA OW)	5.2.2
		3031 (EPA SW-846)	5.2.15
		3050B (EPA SW-846)	5.2.16
		6010C (EPA SW-846)	5.2.26
		6020A (EPA SW-846)	5.2.27
		IO-3.1 (EPA ORD)	5.2.42
		IO-3.4 (EPA ORD)	5.2.43
		IO-3.5 (EPA ORD)	5.2.44
Crimidine	535-89-7	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35

Analyte	CASRN	Method	Section
Cyanide, Total	57-12-5	335.4 (EPA OW)	5.2.6
		ILM05.3 CN (EPA CLP)	5.2.41
		6010 (NIOSH)	5.2.57
Cyanogen chloride	506-77-4	3585 (EPA SW-846)	5.2.23
		5030C (EPA SW-846)	5.2.24
		5035A (EPA SW-846)	5.2.25
		8260C (EPA SW-846)	5.2.34
		TO-15 (EPA ORD)	5.2.47
Cyclohexyl sarin (GF)	329-99-7	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
1,2-Dichloroethane	107-06-2	524.2 (EPA OW)	5.2.9
		3585 (EPA SW-846)	5.2.23
		5030C (EPA SW-846)	5.2.24
		5035A (EPA SW-846)	5.2.25
		8260C (EPA SW-846)	5.2.34
		TO-15 (EPA ORD)	5.2.47
Dichlorvos	62-73-7	525.2 (EPA OW)	5.2.10
		3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Dicrotophos	141-66-2	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Diesel Range Organics	NA	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8015C (EPA SW-846)	5.2.33
Diisopropyl methylphosphonate (DIMP) Dimethylphosphite	1445-75-6	3535A (EPA SW-846)	5.2.18
	868-85-9	3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46

Analyte	CASRN	Method	Section
Dimethylphosphoramidic acid	33876-51-6	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8321B (EPA SW-846)	5.2.39
		TO-10A (EPA ORD)	5.2.46
3,5-Dinitroaniline (3,5-DNA) 1,3-Dinitrobenzene (1,3-DNB) 2,4-Dinitrotoluene (2,4-DNT) 2,6-Dinitrotoluene (2,6-DNT)	618-87-1 99-65-0 121-14-2 606-20-2	3535A (EPA SW-846)	5.2.18
		8330B (EPA SW-846)	5.2.40
Diphacinone	82-66-6	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8321B (EPA SW-846)	5.2.39
1,4-Dithiane	505-29-3	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
EA2192 [Diisopropylaminoethyl methyl-thiophosphonate] Ethyl methylphosphonic acid (EMPA)	73207-98-4 1832-53-7	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8321B (EPA SW-846)	5.2.39
		TO-10A (EPA ORD)	5.2.46
Ethylchloroarsine (ED)	598-14-1	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-15 (EPA ORD)	5.2.47
N-Ethyldiethanolamine (EDEA)	139-87-7	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8321B (EPA SW-846)	5.2.39
		TO-10A (EPA ORD)	5.2.46
Ethylene oxide	75-21-8	3585 (EPA SW-846)	5.2.23
		5030C (EPA SW-846)	5.2.24
		5035A (EPA SW-846)	5.2.25
		8260C (EPA SW-846)	5.2.34
		TO-15 (EPA ORD)	5.2.47

Analyte	CASRN	Method	Section
Fenamiphos	22224-92-6	525.2 (EPA OW)	5.2.10
		3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Fentanyl	437-38-7	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8321B (EPA SW-846)	5.2.39
		TO-10A (EPA ORD)	5.2.46
Fluoride	16984-48-8	300.1, Rev 1.0 (EPA OW)	5.2.5
Fluoroacetic acid and fluoroacetate salts	NA	300.1, Rev 1.0 (EPA OW)	5.2.5
		S301-1 (NIOSH)	5.2.64
		Analytical Letters, 1994, 27 (14), 2703-2718	5.2.79
Formaldehyde	50-00-0	8315A (EPA SW-846)	5.2.36
		2016 (NIOSH)	5.2.49
Gasoline Range Organics	NA	3585 (EPA SW-846)	5.2.23
		5030C (EPA SW-846)	5.2.24
		5035A (EPA SW-846)	5.2.25
		8015C (EPA SW-846)	5.2.33
Hexahydro-1,3,5-trinitro-1,3,5-triazine (RDX) Hexamethylenetriperoxidediamine (HMTD)	121-82-4 283-66-9	3535A (EPA SW-846)	5.2.18
		8330B (EPA SW-846)	5.2.40
Hydrogen bromide Hydrogen chloride	10035-10-6 7647-01-0	7903 (NIOSH)	5.2.61
Hydrogen cyanide	74-90-8	335.4 (EPA OW)	5.2.6
		ILM05.3 CN (EPA CLP)	5.2.41
		6010 (NIOSH)	5.2.57
Hydrogen fluoride	7664-39-3	7903 (NIOSH)	5.2.61
Hydrogen sulfide	7783-06-4	6013 (NIOSH)	5.2.58
Isopropyl methylphosphonic acid (IMPA)	1832-54-8	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8321B (EPA SW-846)	5.2.39
		TO-10A (EPA ORD)	5.2.46
Kerosene	64742-81-0	3585 (EPA SW-846)	5.2.23
		5030C (EPA SW-846)	5.2.24
		5035A (EPA SW-846)	5.2.25
		8015C (EPA SW-846)	5.2.33

Analyte	CASRN	Method	Section
Lewisite 1 (L-1) [2-chlorovinylchloroarsine] Lewisite 2 (L-2) [bis(2-chlorovinyl)chloroarsine] Lewisite 3 (L-3) [tris(2-chlorovinyl)arsine] Lewisite oxide	541-25-3	200.8 (EPA OW)	5.2.2
		3031 (EPA SW-846)	5.2.15
	40334-69-8	3050B (EPA SW-846)	5.2.16
		6010C (EPA SW-846)	5.2.26
	40334-70-1	6020A (EPA SW-846)	5.2.27
		IO-3.1 (EPA ORD)	5.2.42
	1306-02-1	IO-3.4 (EPA ORD)	5.2.43
		IO-3.5 (EPA ORD)	5.2.44
Mercury, Total	7439-97-6	245.2 (EPA OW)	5.2.3
		7473 (EPA SW-846)	5.2.31
		IO-5 (EPA ORD)	5.2.45
Methamidophos	10265-92-6	3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Methyomyl	16752-77-5	531.2 (EPA OW)	5.2.11
		8318A (EPA SW-846)	5.2.38
		5601 (NIOSH)	5.2.53
Methyl acrylonitrile	126-98-7	8316 (EPA SW-846)	5.2.37
		PV2004 (OSHA)	5.2.70
Methoxyethylmercuric acetate	151-38-2	245.2 (EPA OW)	5.2.3
		7473 (EPA SW-846)	5.2.31
		IO-5 (EPA ORD)	5.2.45
Methyl fluoroacetate	453-18-9	300.1, Rev 1.0 (EPA OW)	5.2.5
		S301-1 (NIOSH)	5.2.64
		Analytical Letters, 1994, 27 (14), 2703-2718	5.2.79
Methyl hydrazine	60-34-4	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		3510 (NIOSH)	5.2.51
Methyl isocyanate	624-83-9	OSHA 54	5.2.66
Methyl parathion	298-00-0	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Methyl-2,4,6-trinitrophenylnitramine (Tetryl)	479-45-8	3535A (EPA SW-846)	5.2.18
		8330B (EPA SW-846)	5.2.40
Methylamine	74-89-5	OSHA 40	5.2.65

Analyte	CASRN	Method	Section
N-Methyldiethanolamine (MDEA)	105-59-9	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8321B (EPA SW-846)	5.2.39
		TO-10A (EPA ORD)	5.2.46
1-Methylethyl ester ethylphosphonofluoridic acid (GE)	1189-87-3	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Methylphosphonic acid (MPA)	993-13-5	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8321B (EPA SW-846)	5.2.39
		TO-10A (EPA ORD)	5.2.46
Mevinphos	7786-34-7	525.2 (EPA OW)	5.2.10
		3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Mustard, nitrogen (HN-1) [bis(2-chloroethyl)-ethylamine]	538-07-8	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
	51-75-2	3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
	555-77-1	8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Mustard, sulfur (HD) / Mustard gas (H)	505-60-2	3571 (EPA SW-846)	5.2.21
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Nicotine sulfate	54-11-5	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
Nitrobenzene (NB) Nitroglycerine (NG) 2-Nitrotoluene (2-NT) 3-Nitrotoluene (3-NT) 4-Nitrotoluene (4-NT)	98-95-3	3535A (EPA SW-846)	5.2.18
	55-63-0		
	88-72-2	8330B (EPA SW-846)	5.2.40
	99-08-1		
Octahydro-1,3,5,7-tetranitro-1,3,5,7-tetrazocine (HMX)	2691-41-0	3535A (EPA SW-846)	5.2.18
		8330B (EPA SW-846)	5.2.40

Analyte	CASRN	Method	Section
Organophosphate pesticides, NOS	NA	507 (EPA OW)	5.2.8
		614 (EPA OW)	5.2.14
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		5600 (NIOSH)	5.2.52
Osmium tetroxide	20816-12-0	252.2 (EPA OW)	5.2.4
		3050B (EPA SW-846)	5.2.16
		6010C (EPA SW-846)	5.2.26
		IO-3.1 (EPA ORD)	5.2.42
		IO-3.4 (EPA ORD)	5.2.43
Oxamyl	23135-22-0	531.2 (EPA OW)	5.2.11
		8318A (EPA SW-846)	5.2.38
		5601 (NIOSH)	5.2.53
Paraquat	4685-14-7	549.2 (EPA OW)	5.2.12
Parathion	56-38-2	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Pentaerythritol tetranitrate (PETN)	78-11-5	3535A (EPA SW-846)	5.2.18
		8330B (EPA SW-846)	5.2.40
Perfluoroisobutylene (PFIB)	382-21-8	OSHA 61	5.2.67
Phencyclidine	77-10-1	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Phenol	108-95-2	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Phorate	298-02-2	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Phosgene	75-44-5	OSHA 61	5.2.67

Analyte	CASRN	Method	Section
Phosphamidon	13171-21-6	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Phosphine	7803-51-2	6002 (NIOSH)	5.2.55
Phosphorus trichloride	7719-12-2	6402 (NIOSH)	5.2.60
Pinacolyl methyl phosphonic acid (PMPA)	616-52-4	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8321B (EPA SW-846)	5.2.39
		TO-10A (EPA ORD)	5.2.46
Propylene oxide	75-56-9	3585 (EPA SW-846)	5.2.23
		5030C (EPA SW-846)	5.2.24
		5035A (EPA SW-846)	5.2.25
		8260C (EPA SW-846)	5.2.34
		1612 (NIOSH)	5.2.48
R-33 (VR) [methylphosphonothioic acid, S-[2-(diethylamino)ethyl] O-2-methylpropyl ester]	159939-87-4	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Sarin	107-44-8	3571 (EPA SW-846)	5.2.21
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Semivolatile Organic Compounds, NOS	NA	525.2 (EPA OW)	5.2.10
		3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Sodium arsenite	7784-46-5	200.8 (EPA OW)	5.2.2
		3031 (EPA SW-846)	5.2.15
		3050B (EPA SW-846)	5.2.16
		6010C (EPA SW-846)	5.2.26
		6020A (EPA SW-846)	5.2.27
		IO-3.1 (EPA ORD)	5.2.42
		IO-3.4 (EPA ORD)	5.2.43
		IO-3.5 (EPA ORD)	5.2.44

Analyte	CASRN	Method	Section
Sodium azide	26628-22-8	300.1, Rev 1.0 (EPA OW)	5.2.5
		3580A (EPA SW-846)	5.2.22
		ID-211 (OSHA)	5.2.68
		J. of Forensic Science, 1998. 43(1): 200–202	5.2.81
Soman (GD)	96-64-0	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Strychnine	57-24-9	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
Sulfur dioxide	7446-09-5	6004 (NIOSH)	5.2.56
Sulfur trioxide	7446-11-9	Method 8 (EPA)	5.2.1
Tabun (GA)	77-81-6	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Tear gas (CS) [chlorobenzylidene malonitrile]	2698-41-1	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Tetraethyl pyrophosphate	107-49-3	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Tetramethylenedisulfotetramine	80-12-6	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
Thallium sulfate	10031-59-1	200.8 (EPA OW)	5.2.2
		3031 (EPA SW-846)	5.2.15
		3050B (EPA SW-846)	5.2.16
		6010C (EPA SW-846)	5.2.26
		6020A (EPA SW-846)	5.2.27
		IO-3.1 (EPA ORD)	5.2.42
		IO-3.4 (EPA ORD)	5.2.43
		IO-3.5 (EPA ORD)	5.2.44

Analyte	CASRN	Method	Section
Thiodiglycol (TDG)	111-48-8	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
1,4-Thioxane	15980-15-1	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
Titanium tetrachloride	7550-45-0	3050B (EPA SW-846)	5.2.16
		6010C (EPA SW-846)	5.2.26
		6020A (EPA SW-846)	5.2.27
Triacetone triperoxide (TATP)	17088-37-8	3535A (EPA SW-846)	5.2.18
		8330B (EPA SW-846)	5.2.40
Triethanolamine (TEA)	102-71-6	3520C (EPA SW-846)	5.2.17
		3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8321B (EPA SW-846)	5.2.39
		TO-10A (EPA ORD)	5.2.46
Trimethyl phosphite	121-45-9	3535A (EPA SW-846)	5.2.18
		3545A (EPA SW-846)	5.2.20
		3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
1,3,5-Trinitrobenzene (1,3,5-TNB) 2,4,6-Trinitrotoluene (2,4,6-TNT)	99-35-4	3535A (EPA SW-846)	5.2.18
	118-96-7	8330B (EPA SW-846)	5.2.40
Vanadium pentoxide	1314-62-1	200.8 (EPA OW)	5.2.2
		3031 (EPA SW-846)	5.2.15
		3050B (EPA SW-846)	5.2.16
		6010C (EPA SW-846)	5.2.26
		6020A (EPA SW-846)	5.2.27
		IO-3.1 (EPA ORD)	5.2.42
		IO-3.4 (EPA ORD)	5.2.43
		IO-3.5 (EPA ORD)	5.2.44
VE [phosphonothioic acid, ethyl-, S-(2-(diethylamino)ethyl) O-ethyl ester]	21738-25-0	3520C (EPA SW-846)	5.2.17
VG [phosphonothioic acid, S-(2-(diethylamino)ethyl) O,O-diethyl ester]	78-53-5	3535A (EPA SW-846)	5.2.18
		3541 (EPA SW-846)	5.2.19
		3545A (EPA SW-846)	5.2.20
VM [phosphonothioic acid, methyl-,S-(2-(diethylamino)ethyl) O-ethyl ester]	21770-86-5	3580A (EPA SW-846)	5.2.22
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46
VX [O-ethyl-S-(2-diisopropylaminoethyl)methyl-phosphonothiolate]	50782-69-9	3571 (EPA SW-846)	5.2.21
		8270D (EPA SW-846)	5.2.35
		TO-10A (EPA ORD)	5.2.46

Analyte	CASRN	Method	Section
White phosphorus	12185-10-3	7580 (EPA SW-846)	5.2.32
		7905 (NIOSH)	5.2.62
The following analytes should be prepared and/or analyzed by the following methods only if problems (e.g., insufficient recovery, interferences) occur when using the sample preparation/determinative techniques identified for these analytes in Appendix A.			
Allyl alcohol	107-18-6	TO-10A (EPA ORD)	5.2.46
3-Chloro-1,2-propanediol	96-24-2	TO-15 (EPA ORD)	5.2.47
Chlorosarin Chlorosoman	1445-76-7 7040-57-5	TO-15 (EPA ORD)	5.2.47
Crimidine	535-89-7	8321B (EPA SW-846)	5.2.39
Diisopropyl methylphosphonate (DIMP)	1445-75-6	TO-15 (EPA ORD)	5.2.47
Dimethylphosphoramidic acid	33876-51-6	8270D (EPA SW-846)	5.2.35
1,4-Dithiane	505-29-3	3585 (EPA SW-846)	5.2.23
		5030C (EPA SW-846)	5.2.24
		5035A (EPA SW-846)	5.2.25
		8260C (EPA SW-846)	5.2.34
EA2192 [Diisopropylaminoethyl methyl-thiophosphonate]	73207-98-4	8270D (EPA SW-846)	5.2.35
Ethyl methylphosphonic acid (EMPA)	1832-53-7		
Hydrogen fluoride	7664-39-3	7906 (NIOSH)	5.2.63
Isopropyl methylphosphonic acid (IMPA)	1832-54-8	8270D (EPA SW-846)	5.2.35
Mercury, Total	7439-97-6	7470A (EPA SW-846)	5.2.29
		7471B (EPA SW-846)	5.2.30
Methamidophos	10265-92-6	5600 (NIOSH)	5.2.52
Methoxymercuric acetate	151-38-2	7470A (EPA SW-846)	5.2.29
		7471B (EPA SW-846)	5.2.30
1-Methylethyl ester ethylphosphonofluoridic acid (GE)	1189-87-3	TO-15 (EPA ORD)	5.2.47
Methylphosphonic acid (MPA)	993-13-5	8270D (EPA SW-846)	5.2.35
Perfluoroisobutylene (PFIB)	382-21-8	J. Chrom. A, Vol 1098, 2005, pp. 156-165	5.2.80
Pinacolyl methyl phosphonic acid (PMPA)	616-52-4	8270D (EPA SW-846)	5.2.35
Sarin Soman (GD)	107-44-8 96-64-0	TO-15 (EPA ORD)	5.2.47
1,4-Thioxane	15980-15-1	3585 (EPA SW-846)	5.2.23
		5030C (EPA SW-846)	5.2.24
		5035A (EPA SW-846)	5.2.25
		8260C (EPA SW-846)	5.2.34

Method summaries are listed in order of method selection hierarchy (see Figure 2-1), starting with EPA methods followed by methods from other federal agencies, voluntary consensus standard bodies (VCSB), and literature references. Methods are listed in numerical order under each publisher. Where available, a direct link to the full text of the method is provided in the method summary. For additional information on preparation procedures and methods available through consensus standards organizations, please use the contact information provided in **Table 5-2**.

Table 5-2. Sources of Chemical Methods

Name	Publisher	Reference
National Environmental Methods Index (NEMI)	EPA, USGS	http://www.nemi.gov
U.S. EPA Office of Water (OW) Methods	EPA Office of Water	http://www.epa.gov/safewater/methods/sourcalt.html
U.S. EPA SW-846 Methods	EPA Office of Solid Waste and Emergency Response	http://www.epa.gov/epaoswer/hazwaste/test/main.htm
U.S. EPA Office of Research and Development Methods	EPA Office of Research and Development	http://www.epa.gov/nerlcwww/ordmeth.htm
U.S. EPA Air Toxics Methods	EPA Office of Air and Radiation	http://www.epa.gov/ttn/amtic/airtox.html
Occupational Safety and Health Administration Methods	OSHA	http://www.osha.gov
National Institutes for Occupational Safety and Health Methods	NIOSH	http://www.cdc.gov/niosh/nmam/
Standard Methods for the Examination of Water and Wastewater, 21 st Edition, 2005*	American Public Health Association (APHA), American Water Works Association (AWWA), and Water Environment Federation (WEF)	http://www.standardmethods.org
Annual Book of ASTM Standards*	ASTM International	http://www.astm.org
European GESTIS database	HVBG	http://www.hvb.de/e/bia/gestis/analytical_methods/index.html
International Organization for Standardization Methods*	ISO	http://www.iso.org
Official Methods of Analysis of AOAC International*	AOAC International	http://www.aoac.org

* Subscription and/or purchase required.

5.1.2 General Quality Control (QC) Guidance for Chemical Methods

Having analytical data of appropriate quality requires that laboratories: (1) conduct the necessary QC to ensure that measurement systems are in control and operating correctly, (2) properly document results of the analyses, and (3) properly document measurement system evaluation of the analysis-specific QC, including corrective actions. In emergency situations, however, speed and efficiency are also important. Laboratories must be prepared with calibrated instruments, the proper standards, standard analytical procedures, and qualified and trained technicians. Laboratories also must be capable of providing rapid turnaround of sample analyses and data reporting.

The level or amount of QC needed during sample analysis and reporting depends on the intended purpose of the data that are generated (i.e., the decision(s) to be made). The specific needs for data generation should be identified. QC requirements and data quality objectives should be derived based on those needs, and should be applied consistently across laboratories when multiple laboratories are used. For almost all of the chemical warfare agents, most laboratories will not have access to analytical standards for calibration and quality control. Use of these agents is strictly controlled by the DoD and access is

limited. For information regarding purchase and distribution of ultradilute agents, please contact Terry Smith, EPA's Office of Emergency Management, at (202) 564-2908.

A minimum set of analytical QC procedures should be planned and conducted for all chemical testing. Some method-specific QC requirements are described in many of the individual methods that are cited in this document and will be referenced in any Standardized Analytical Protocols (SAPs) developed to address specific analytes and sample types of concern. Individual methods, sampling and analysis protocols, or contractual statements of work also should be consulted to determine any additional QC that may be needed. Analytical QC requirements generally consist of analysis of laboratory control samples to document whether the analytical system is in control, matrix spikes (MS) to identify and quantify measurement system accuracy for the media of concern and at the levels of concern, various blanks as a measure of freedom from contamination, and matrix spike duplicates (MSD) or sample replicates to assess data precision.

In general, for measurement of chemical analytes, appropriate QC includes an initial demonstration of measurement system capability as well as ongoing analysis of standards and other samples to ensure the continued reliability of the analytical results. Examples of appropriate quality control include:

- Demonstration that measurement system is operating properly:
 - ▶ Initial calibration; and
 - ▶ Method blanks.
- Demonstration of analytical method suitability for intended use:
 - ▶ Detection and quantitation limits;
 - ▶ Precision and recovery (verify measurement system has adequate accuracy); and
 - ▶ Analyte/matrix/level of concern-specific QC samples (verify that measurement system has adequate sensitivity at levels of concern).
- Demonstration of continued analytical method reliability:
 - ▶ Matrix spike/matrix spike duplicates (recovery and precision);
 - ▶ QC samples (system accuracy and sensitivity at levels of concern);
 - ▶ Surrogate spikes (where appropriate);
 - ▶ Continuing calibration verification; and
 - ▶ Method blanks.

QC tests should be run as frequently as necessary to ensure the reliability of analytical results. As with the identification of needed QC samples, frequency should be established based on an evaluation of data quality objectives. The type and frequency of QC can be focused over time.

Ensuring data quality also requires that laboratory results are properly evaluated and the results of the data quality evaluation are transmitted to decision makers. This evaluation is as important as the data in ensuring informed and effective decisions. While some degree of data evaluation is necessary in order to be able to confirm data quality, 100% verification and/or validation is neither necessary nor conducive to efficient decision making in emergency situations. The level of such reviews should be determined based on the specific situation being assessed and on the corresponding data quality objectives. In every case, the levels of QC and data review necessary to support decision making should be determined as much in advance of data collection as possible.

Please note: The appropriate point of contact identified in Section 4 should be consulted regarding appropriate quality assurance and quality control (QA/QC) procedures prior to sample analysis. These contacts will consult with the EPA Office of Solid Waste and Emergency Response (OSWER) coordinator responsible for laboratory activities during the specific event to ensure QA/QC procedures are performed consistently across laboratories. OSWER is planning to develop QA/QC guidance for laboratory support. EPA program offices will be responsible for ensuring that the QA/QC practices are implemented.

5.1.3 Safety and Waste Management

It is imperative that safety precautions are used during collection, processing, and analysis of environmental samples. Laboratories should have a documented health and safety plan for handling samples that may contain the target chemical, biological, or radiological contaminants, and laboratory staff should be trained in and implement the safety procedures included in the plan. In addition, many of the methods summarized or cited in Section 5.2 contain some specific requirements, guidance, or information regarding safety precautions that should be followed when handling or processing environmental samples and reagents. These methods also provide information regarding waste management. Other resources that can be consulted for additional information include the following:

- Occupational Health and Safety Administration's (OSHA) standard for Occupational Exposure to Hazardous Chemicals in Laboratories (29 CFR 1910.1450)
- OSHA regulations for hazardous waste operations and emergency response (29 CFR 1910)
- Environmental Protection Agency's standards regulating hazardous waste (40 CFR parts 260 - 270)
- U.S. Department of Transportation (DOT) regulations for transporting hazardous materials (49 CFR Part 172)
- U.S. Department of Health and Human Services, Centers for Disease Control and Prevention's requirements for possession, use, and transfer of select agents and toxins (42 CFR Part 1003)

5.2 Method Summaries

Summaries for the analytical methods listed in Appendix A are provided in Sections 5.2.1 through 5.2.81. These sections contain summary information only, extracted from the selected methods. Each method summary contains a table identifying the contaminants in Appendix A to which the method applies, a brief description of the analytical method, and a link to or source for obtaining a full version of the method. The full version of the method should be consulted prior to sample analysis.

5.2.1 EPA Method 8: Determination of Sulfuric Acid and Sulfur Dioxide Emissions from Stationary Sources

This method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Sulfur Trioxide	7446-11-9

A gas sample is extracted isokinetically. Sulfuric acid (H_2SO_4) and sulfur dioxide (SO_2) are separated, and both fractions are measured separately by the barium-thorin titration method. Sulfur trioxide (SO_3) is measured by the analysis of H_2SO_4 . Possible interfering agents include fluorides, free ammonia, and dimethyl aniline. Collaborative tests have shown that the minimum detectable limits of the method are 0.06 mg/m^3 for H_2SO_4 mist/ SO_3 .

Source: EPA Method 8: Determination Of Sulfuric Acid and Sulfur Dioxide Emissions From Stationary Sources. EPA Emission Measurement Center (EMC) of the Office of Air Quality Planning and Standards (OAQPS). <http://www.epa.gov/ttnemc01/methods/method8.html>

5.2.2 EPA Method 200.8: Determination of Trace Elements in Waters and Wastes by Inductively Coupled Plasma-Mass Spectrometry

This method should be used for **preparation** and **analysis** of aqueous liquid and drinking water samples for the contaminants identified below and listed in Appendix A.

Analyte(s)	CAS RN
Ammonium metavanadate	7803-55-6
Arsenic, Total	740-38-2
Arsenic trioxide	1327-53-3
Arsine	7784-42-1
2-Chlorovinylarsonous acid (CVAA)	85090-33-1
Lewisite 1 (L-1) [2-chlorovinylchloroarsine] ¹	541-25-3
Lewisite 2 (L-2) [bis(2-chlorovinyl)chloroarsine]	40334-69-8
Lewisite 3 (L-3) [tris(2-chlorovinyl)arsine]	40334-70-1
Lewisite Oxide	1306-02-1
Sodium arsenite	7784-46-5
Thallium sulfate	10031-59-1
Vanadium pentoxide	1314-62-1

¹ Laboratory testing is currently underway for speciation of lewisite 1 using GC-MS techniques.

This method will determine metal-containing compounds only as the total metal (e.g., total arsenic). An aliquot of a well-mixed, homogeneous sample is accurately weighed or measured for sample processing. For total recoverable analysis of a sample containing undissolved material, analytes are first solubilized by gentle refluxing with nitric and hydrochloric acids. After cooling, the sample is made up to volume, mixed, and centrifuged or allowed to settle overnight prior to analysis. For determination of dissolved analytes in a filtered aqueous sample aliquot, or for the “direct analysis” total recoverable determination of analytes in drinking water where sample turbidity is < 1 NTU, the sample is made ready for analysis by the addition of nitric acid, and then diluted to a predetermined volume and mixed before analysis. The prepared sample is analyzed using Inductively Coupled Plasma-Mass Spectrometry (ICP-MS). Specific analytes targeted by Method 200.8 are listed in Section 1.1 of the method. The recommended calibration range is 10 to 200 µg/L. Method detection limits for arsenic in aqueous samples have been found to be 1.4 µg/L in scanning mode, and 0.4 µg/L in selected ion monitoring mode.

Source: EPA Method 200.8: Determination of Trace Elements in Waters and Wastes by Inductively Coupled Plasma-Mass Spectrometry, Revision 5.4, 1994.

http://infotrek.er.usgs.gov/pls/apex/f?p=119:38:1701708976536019:::P38_METHOD_ID:4665

5.2.3 EPA Method 245.2: Mercury (Automated Cold Vapor Technique)

This method should be used for **preparation** and **analysis** of drinking water samples for the contaminants identified below and listed in Appendix A.

Analyte(s)	CAS RN
Mercury, Total	7439-97-6
Methoxyethylmercuric acetate	151-38-2

This method will determine methoxyethylmercuric acetate as total mercury. If dissolved mercury is targeted, the sample is filtered prior to acidification. To detect total mercury (inorganic and organic mercury), the sample is treated with potassium permanganate and potassium persulfate to oxidize organic mercury compounds prior to analysis. Inorganic mercury is reduced to the elemental state (using tin sulfate or tin chloride) and aerated from solution. The mercury vapor passes through a cell positioned in the light path of a cold vapor atomic absorption (CVAA) spectrophotometer. The concentration of mercury is measured using the CVAA spectrophotometer. Applicable concentration range is 0.2 to 20.0 µg/L.

Source: EPA Method 245.2: Mercury (Automated Cold Vapor Technique), 1974.

http://infotrek.er.usgs.gov/pls/apex/f?p=119:38:1701708976536019:::P38_METHOD_ID:4822

5.2.4 EPA Method 252.2: Osmium (Atomic Absorption, Furnace Technique)

This method should be used for the **preparation** and **analysis** of aqueous liquid and drinking water samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Osmium tetroxide	20816-12-0

This method will determine osmium tetroxide as osmium. Method 252.2 is not a stand alone method in that sections of the method reference “Methods of Chemical Analysis of Water and Waste”, EPA/600/4-79/020, March 1983 (MCAWW). Samples are prepared according to the “direct aspiration method” (See Section 9.1 of the Atomic Absorption Methods section of MCAWW) except that the addition of sulfuric acid is omitted in the final adjustment. If only dissolved osmium is determined, the sample is filtered before acidification with nitric acid. For total osmium, the sample is digested with nitric and hydrochloric acids and made up to volume. Samples are analyzed according to the “furnace procedure” (see Section 9.3 of the Atomic Absorption Methods section of MCAWW), using a graphite furnace atomic absorption spectrometer. A representative aliquot of sample is placed in the graphite tube in the furnace, evaporated to dryness, chaffed, and atomized. Radiation from an excited element is passed through the vapor containing ground state atoms of the element. The intensity of the transmitted radiation decreases in proportion to the amount of the ground state element in the vapor. A monochromator isolates the characteristic radiation from the hollow cathode lamp and a photosensitive device measures the attenuated transmitted radiation. The optimal applicable concentration range is 50 to 500 µg/L. Detection limits for osmium using this method have been found to be 20 µg/L.

Source: EPA Method 252.2: Osmium (AA, Furnace Technique), 1978.

http://infotrek.er.usgs.gov/pls/apex/f?p=119:38:1701708976536019:::P38_METHOD_ID:5299

5.2.5 EPA Method 300.1, Revision 1.0: Determination of Inorganic Anions in Drinking Water by Ion Chromatography

This method should be used for the **preparation** and **analysis** of aqueous liquid and drinking water samples for the contaminants identified below and listed in Appendix A. It also should be used for **analysis** of solid, non-aqueous liquid/organic solid, and air samples for fluoroacetic acid, fluoroacetate salts, and methyl fluoroacetate when appropriate sample preparation techniques have been applied (refer to Appendix A).

Analyte(s)	CAS RN
Fluoride	16984-48-8
Fluoroacetic acid and fluoroacetate salts	NA
Methyl fluoroacetate	453-18-9

This method was developed for analysis of aqueous samples, and can be adapted for analysis of prepared non-aqueous liquid/organic solid, solid, and air samples when appropriate sample preparation techniques have been applied (see Appendix A). A small volume of an aqueous liquid sample (10 µL or 50 µL) is introduced into an ion chromatograph. The volume selected depends on the concentration of fluoroacetate ion in the sample. The anions of interest are separated and measured, using a system comprising a guard column, analytical column, suppressor device, and conductivity detector. The separator columns and guard columns, as well as eluent conditions, are identical. To achieve comparable detection limits, an ion chromatographic system must use suppressed conductivity detection, be properly maintained, and be capable of yielding a baseline with no more than 5 nS noise/drift per minute of monitored response over the background conductivity. The method detection limit varies depending upon the nature of the sample and the specific instrumentation employed. The estimated calibration range is approximately two orders of magnitude. The detection limit for fluoride in reagent water is 0.009 mg/L.

Source: EPA Method 300.1: Determination of Inorganic Anions in Drinking Water by Ion Chromatography, Revision 1.0, 1997.

http://infotrek.er.usgs.gov/pls/apex/f?p=119:38:1701708976536019:::P38_METHOD_ID:4674

5.2.6 EPA Method 335.4: Determination of Total Cyanide by Semi-Automated Colorimetry

This method should be used for **preparation** and **analysis** of drinking water samples for the contaminants identified below and listed in Appendix A.

Analyte(s)	CAS RN
Cyanide, Total	57-12-5
Hydrogen cyanide	74-90-8

Cyanide is released from cyanide complexes as hydrocyanic acid by manual reflux-distillation, and absorbed in a scrubber containing sodium hydroxide solution. The cyanide ion in the absorbing solution is converted to cyanogen chloride by reaction with chloramine-T, which subsequently reacts with pyridine and barbituric acid to give a red-colored complex. Some interferences, such as aldehydes, nitrate-nitrite, oxidizing agents, thiocyanate, thiosulfate, and sulfide, are eliminated or reduced by distillation. The applicable range is 5 to 500 µg/L.

Source: EPA Method 335.4: Determination of Total Cyanide by Semi-Automated Colorimetry, Revision 1.0, 1993.

http://infotrek.er.usgs.gov/pls/apex/f?p=119:38:1701708976536019:::P38_METHOD_ID:5759

5.2.7 EPA Method 350.1: Nitrogen, Ammonia (Colorimetric, Automated Phenate)

This method should be used for **preparation** and **analysis** of drinking water samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Ammonia	7664-41-7

This method identifies and determines the concentration of ammonia in drinking water samples by spectrophotometry. Samples are buffered at a pH of 9.5 with borate buffer to decrease hydrolysis of cyanates and organic nitrogen compounds, and are distilled into a solution of boric acid. Alkaline phenol and hypochlorite react with ammonia to form indophenol blue that is proportional to the ammonia concentration. The blue color formed is intensified with sodium nitroprusside and measured spectrophotometrically. The working range for ammonia is 0.01 to 2.0 mg/L.

Source: EPA Method 350.1: Nitrogen, Ammonia (Colorimetric, Automated Phenate), Revision 2.0, 1993. http://infotrek.er.usgs.gov/pls/apex/f?p=119:38:1701708976536019:::P38_METHOD_ID:5405

5.2.8 EPA Method 507: Determination of Nitrogen- and Phosphorus-Containing Pesticides in Water by Gas Chromatography with a Nitrogen-Phosphorus Detector

This method should be used for **preparation** and **analysis** of drinking water samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Organophosphate Pesticides, NOS ¹	NA

¹ NOS = Not otherwise specified

This is a gas chromatographic (GC) method applicable to the determination of certain nitrogen and phosphorus containing pesticides in ground water and finished drinking water. A 1-L sample is extracted with methylene chloride by shaking in a separatory funnel or mechanical tumbling in a bottle. The methylene chloride extract is isolated, dried, and concentrated to a volume of 5 mL during a solvent exchange to methyl tert-butyl ether (MTBE). The concentrations of pesticides in the extract are measured using a capillary column GC system equipped with a nitrogen-phosphorus detector (NPD). Specific analytes targeted by Method 507 are listed in Section 1.1 of the method. Resulting estimated detection limits and method detection limits (MDLs) differ depending on the specific pesticide.

Please note: The presence of any organophosphate pesticide must be confirmed by either a secondary GC column or by a mass spectrometer.

Source: EPA Method 507: Determination of Nitrogen- and Phosphorus-Containing Pesticides in Water by Gas Chromatography with a Nitrogen-Phosphorus Detector, Revision 2.1, 1995. http://infotrek.er.usgs.gov/pls/apex/f?p=119:38:1701708976536019:::P38_METHOD_ID:4801

5.2.9 EPA Method 524.2: Measurement of Purgeable Organic Compounds in Water by Capillary Column Gas Chromatography / Mass Spectrometry

This method should be used for **preparation** and **analysis** of drinking water samples for the contaminants identified below and listed in Appendix A.

Analyte(s)	CAS RN
Carbon disulfide	75-15-0
1,2-Dichloroethane	107-06-2

Volatile organic compounds and surrogates with low water solubility are extracted (purged) from the sample matrix by bubbling an inert gas through the aqueous sample. Purged sample components are trapped in a tube containing suitable sorbent materials. When purging is complete, the sorbent tube is heated and backflushed with helium to desorb the trapped sample components into a capillary gas chromatography (GC) column interfaced to a mass spectrometer (MS). The column is temperature programmed to facilitate the separation of the method analytes, which are then detected with the MS. Specific analytes targeted by Method 524.2 are listed in Section 1.1 of the method. Detection levels for carbon disulfide and 1,2-dichloroethane in reagent water have been found to be 0.093 µg/L and 0.02 µg/L, respectively.

Source: EPA Method 524.2: Measurement of Purgeable Organic Compounds in Water by Capillary Column Gas Chromatography/Mass Spectrometry, Revision 4.0, 1992.

http://infotrek.er.usgs.gov/pls/apex/f?p=119:38:1701708976536019:::P38_METHOD_ID:4803

5.2.10 EPA Method 525.2: Determination of Organic Compounds in Drinking Water by Liquid-Solid Extraction and Capillary Column Gas Chromatography / Mass Spectrometry

This method should be used for **preparation** and **analysis** of drinking water samples for the contaminants identified below and listed in Appendix A.

Analyte(s)	CAS RN
Dichlorvos	62-73-7
Fenamiphos	22224-92-6
Mevinphos	7786-34-7
Semivolatile Organic Compounds, NOS ¹	NA

¹ NOS = Not otherwise specified

Organic compounds, internal standards, and surrogates are extracted from a water sample by passing 1 L of sample through a cartridge or disk containing a solid matrix with chemically bonded C₁₈ organic phase (liquid-solid extraction, LSE or solid-phase extraction, SPE). The organic compounds are eluted from the LSE (SPE) cartridge or disk with small quantities of ethyl acetate followed by methylene chloride. The resulting extract is concentrated further by evaporation of some of the solvent. Sample components are separated, identified, and measured by injecting an aliquot of the concentrated extract into a high resolution fused silica capillary column of a gas chromatography/mass spectrometry (GC-MS) system. Specific analytes targeted by Method 525.2 are listed in Section 1.1 of the method. The applicable concentration range for most analytes is 0.1 to 10 µg/L.

Source: EPA Method 525.2: Determination of Organic Compounds in Drinking Water by Liquid-Solid Extraction and Capillary Column Gas Chromatography/Mass Spectrometry, Revision 2.0, 1995.

http://infotrek.er.usgs.gov/pls/apex/f?p=119:38:1701708976536019:::P38_METHOD_ID:4804

5.2.11 EPA Method 531.2: Measurement of N-Methylcarbamoyloximes and N-Methylcarbamates in Water by Direct Aqueous Injection HPLC with Postcolumn Derivatization

This method should be used for **preparation** and **analysis** of drinking water samples for the contaminants identified below and listed in Appendix A.

Analyte(s)	CAS RN
Aldicarb (Temik)	116-06-3
Carbofuran (Furadan)	1563-66-2
Methomyl	16752-77-5
Oxamyl	23135-22-0

An aliquot of sample is measured in a volumetric flask. Samples are preserved, spiked with appropriate surrogates and then filtered. Analytes are chromatographically separated by injecting a sample aliquot (up to 1000 µL) into a high performance liquid chromatographic (HPLC) system equipped with a reverse phase (C₁₈) column. After elution from the column, the analytes are hydrolyzed in a post column reaction to form methylamine, which is in turn reacted to form a fluorescent isoindole that is detected by a fluorescence detector. Analytes also are quantitated using the external standard technique.

Source: EPA Method 531.2: Measurement of N-Methylcarbamoyloximes and N-Methylcarbamates in Water by Direct Aqueous Injection HPLC with Postcolumn Derivatization, Revision 1.0, 2001.

http://infotrek.er.usgs.gov/pls/apex/f?p=119:38:1701708976536019:::P38_METHOD_ID:7392

5.2.12 EPA Method 549.2: Determination of Diquat and Paraquat in Drinking Water by Liquid-Solid Extraction and High-Performance Liquid Chromatography with Ultraviolet Detection

This method should be used for **preparation** and **analysis** of aqueous liquid and drinking water samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Paraquat	4685-14-7

A 250-mL sample is extracted using a C₈ liquid/solid extraction (LSE) cartridge or a C₈ disk that has been specially prepared for the reversed-phase, ion-pair mode. The LSE disk or cartridge is eluted with acidic aqueous solvent to yield the eluate/extract. An ion-pair reagent is added to the eluate/extract. The concentrations of paraquat in the eluate/extract are measured using a high performance liquid chromatography (HPLC) system equipped with a UV absorbance detector. A photodiode array detector is used to provide simultaneous detection and confirmation of the method analytes. The analytical range depends on the sample matrix and the instrumentation used.

Source: EPA Method 549.2: Determination of Diquat and Paraquat in Drinking Water by Liquid-Solid Extraction and High Performance Liquid Chromatography with Ultraviolet Detection, Revision 1.0, 1997.

http://www.epa.gov/nerlcwww/m_549_2.pdf

5.2.13 EPA Method 551.1: Determination of Chlorination Disinfection Byproducts, Chlorinated Solvents, and Halogenated Pesticides/Herbicides in Drinking Water by Liquid-Liquid Extraction and Gas Chromatography with Electron-Capture Detection

This method should be used for **preparation** and **analysis** of drinking water samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Chloropicrin	79-06-2

This is a gas chromatographic/electron capture detection (GC-ECD) method applicable to the determination of halogenated analytes in finished drinking water, drinking water during intermediate stages of treatment, and raw source water. A 50-mL sample aliquot is extracted with 3 mL of methyl tert-butyl ether (MTBE) or 5 mL of pentane. Two μ L of the extract is then injected into a GC equipped with a fused silica capillary column and linearized electron capture detector for separation and analysis. This liquid/liquid extraction technique efficiently extracts a wide boiling range of non-polar and polar organic components of the sample. Thus, confirmation is quite important, particularly at lower analyte concentrations. A confirmatory column is suggested for this purpose. The Estimated Detection Limit (EDL) using MTBE and ammonium chloride-preserved reagent water on a DB-1 column has been found to be 0.014 μ g/L.

Source: EPA Method 551.1: Determination of Chlorination Disinfection Byproducts, Chlorinated Solvents, and Halogenated Pesticides/Herbicides in Drinking Water by Liquid-Liquid Extraction and Gas Chromatography with Electron-Capture Detection.

http://infotrek.er.usgs.gov/pls/apex/f?p=119:38:1701708976536019:::P38_METHOD_ID:4809

5.2.14 EPA Method 614: The Determination of Organophosphorus Pesticides in Municipal and Industrial Wastewater

This method should be used for **preparation** and **analysis** of aqueous liquid samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Organophosphate Pesticides, NOS ¹	NA

¹ NOS = Not otherwise specified

This is a gas chromatographic (GC) method applicable to the determination of organophosphate pesticides in industrial and municipal discharges using a GC with a phosphorus-specific flame photometric detector (FPD) or thermionic bead detector in the nitrogen mode. A measured volume of sample, approximately 1 L, is extracted with 15% methylene chloride in hexane using a separatory funnel. The extract is dried and concentrated to a volume of 10 mL or less. Gas chromatographic conditions are described for the separation and measurement of the compounds in the extract by flame photometric or thermionic bead gas chromatography. Specific analytes targeted by Method 614 are listed in Section 1.1 of the method.

Please note: The presence of any organophosphate pesticide must be confirmed by either a secondary GC column or with a mass spectrometer.

Source: “Methods for the Determination of Nonconventional Pesticides in Municipal and Industrial Wastewater – Vol. I,” United States Environmental Protection Agency, available from National Technical Information Service (NTIS), 5285 Port Royal Road, Springfield, VA 22161. Phone: 800-553-6847.

5.2.15 EPA Method 3031 (SW-846): Acid Digestion of Oils for Metals Analysis by Atomic Absorption or ICP Spectrometry

This method should be used for **preparation** of non-aqueous liquid/organic solid samples for the contaminants identified below and listed in Appendix A. Refer to Appendix A for the appropriate determinative method for these analytes.

Analyte(s)	CAS RN
Ammonium metavanadate	7803-55-6
Arsenic, Total	740-38-2
Arsenic trioxide	1327-53-3
2-Chlorovinylarsonous acid (CVAA)	85090-33-1
Lewisite 1 (L-1) [2-chlorovinylchloroarsine]	541-25-3
Lewisite 2 (L-2) [bis(2-chlorovinyl)chloroarsine]	40334-69-8
Lewisite 3 (L-3) [tris(2-chlorovinyl)arsine]	40334-70-1
Lewisite Oxide	1306-02-1
Sodium arsenite	7784-46-5
Thallium sulfate	10031-59-1
Vanadium pentoxide	1314-62-1

This method is used to prepare samples for the determination of arsenic trioxide, lewisite, lewisite degradation products, and sodium arsenite as total arsenic; thallium sulfate as total thallium; and ammonium metavanadate and vanadium pentoxide as total vanadium. A 0.5-g sample of oil, oil sludge, tar, wax, paint, or paint sludge is mixed with potassium permanganate and sulfuric acid. The mixture is then treated with nitric and hydrochloric acids, filtered, and diluted to volume. Excess manganese may be removed with ammonium hydroxide. Digestates are analyzed by Method 6020A or 6010C (SW-846).

Source: EPA Method 3031 (SW-846): Acid Digestion of Oils for Metals Analysis by Atomic Absorption or ICP Spectrometry, Revision 0, 1996. <http://www.epa.gov/epaoswer/hazwaste/test/pdfs/3031.pdf>

5.2.16 EPA Method 3050B (SW-846): Acid Digestion of Sediments, Sludges, and Soils

This method should be used for **preparation** of solid samples for the contaminants identified below and listed in Appendix A. Refer to Appendix A for the appropriate determinative method for these analytes.

Analyte(s)	CAS RN
Ammonium metavanadate	7803-55-6
Arsenic, Total	740-38-2
Arsenic trioxide	1327-53-3
Arsine	7784-42-1

Analyte(s)	CAS RN
2-Chlorovinylarsonous acid (CVAA)	85090-33-1
Lewisite 1 (L-1) [2-chlorovinylchloroarsine]	541-25-3
Lewisite 2 (L-2) [bis(2-chlorovinyl)chloroarsine]	40334-69-8
Lewisite 3 (L-3) [tris(2-chlorovinyl)arsine]	40334-70-1
Lewisite Oxide	1306-02-1
Osmium tetroxide	20816-12-0
Sodium arsenite	7784-46-5
Thallium sulfate	10031-59-1
Titanium tetrachloride	7550-45-0
Vanadium pentoxide	1314-62-1

This method is used to prepare samples for the determination of arsenic trioxide, arsine, lewisite, lewisite degradation products, and sodium arsenite as total arsenic; thallium sulfate as total thallium; titanium tetrachloride as titanium; osmium tetroxide as osmium; and ammonium metavanadate and vanadium pentoxide as total vanadium. A 1-g to 2-g sample is digested with nitric acid and hydrogen peroxide. Sample volumes are reduced, then brought up to a final volume of 100 mL. Samples are analyzed for total arsenic, total thallium, total titanium, or total vanadium by Method 6010C or 6020A (SW-846); use Method 6010C (SW-846) for total osmium.

Source: EPA Method 3050B (SW-846): Acid Digestion of Sediments, Sludges, and Soils, Revision 2, 1996. <http://www.epa.gov/epaoswer/hazwaste/test/pdfs/3050b.pdf>

5.2.17 EPA Method 3520C (SW-846): Continuous Liquid-Liquid Extraction

This method should be used for **preparation** of aqueous liquid and drinking water samples for the contaminants identified below and listed in Appendix A. Note: For fenamiphos and semivolatile organic compounds only, EPA Method 525.2 (rather than Method 3520C) should be used for preparation of drinking water samples. Refer to Appendix A for the appropriate determinative method for these analytes.

Analyte(s)	CAS RN
Brodifacoum	56073-10-0
Bromadiolone	28772-56-7
Chlorfenvinphos	470-90-6
3-Chloro-1,2-propanediol	96-24-2
Chlorosarin	1445-76-7
Chlorosoman	7040-57-5
Chlorpyrifos	2921-88-2
Cyclohexyl sarin (GF)	329-99-7
Diesel Range Organics	NA

Analyte(s)	CAS RN
Diphacinone	82-66-6
N-Ethyldiethanolamine (EDEA)	139-87-7
Fenamiphos	22224-92-6
Fentanyl	437-38-7
Methyl hydrazine	60-34-4
N-Methyldiethanolamine (MDEA)	105-59-9
1-Methylethyl ester ethylphosphonofluoridic acid (GE)	1189-87-3
Mustard, nitrogen (HN-1) [bis(2-chloroethyl)ethylamine]	538-07-8
Mustard, nitrogen (HN-2) [2,2'-dichloro-N-methyldiethylamine N,N-bis(2-chloroethyl)methylamine]	51-75-2
Mustard, nitrogen (HN-3) [tris(2-chloroethyl)amine]	555-77-1
Parathion	56-38-2
Phenol	108-95-2
Phosphamidon	13171-21-6
R-33 (VR) [methylphosphonothioic acid, S-[2-(diethylamino)ethyl] O-2-methylpropyl ester]	159939-87-4
Semivolatile Organic Compounds, NOS ¹	NA
Tear gas (CS) [chlorobenzylidene malonitrile]	2698-41-1
Tetramethylenedisulfotetramine	80-12-6
Triethanolamine (TEA)	102-71-6
VE [phosphonothioic acid, ethyl-, S-(2-(diethylamino)ethyl) O-ethyl ester]	21738-25-0
VG [phosphonothioic acid, S-(2-(diethylamino)ethyl) O,O-diethyl ester]	78-53-5
VM [phosphonothioic acid, methyl-, S-(2-(diethylamino)ethyl) O-ethyl ester]	21770-86-5

¹ NOS = Not otherwise specified

This method is applicable to the isolation and concentration of water-insoluble and slightly soluble organics in preparation for a variety of chromatographic procedures. A measured volume of sample, usually 1 L, is placed into a continuous liquid-liquid extractor, adjusted, if necessary, to a specific pH and extracted with organic solvent for 18 to 24 hours. The extract is filtered through sodium sulfate to remove residual moisture, concentrated, and exchanged as necessary into a solvent compatible with the cleanup or determinative procedure used for analysis.

Source: EPA Method 3520C (SW-846): Continuous Liquid-Liquid Extraction, Revision 3, 1996.
<http://www.epa.gov/epaoswer/hazwaste/test/pdfs/3520c.pdf>

5.2.18 EPA Method 3535A (SW-846): Solid-Phase Extraction

This method should be used for **preparation** of aqueous liquid and/or drinking water samples for the contaminants identified below and listed in Appendix A. Note: For dichlorvos, fenamiphos, mevinphos, and semivolatile organic compounds only, EPA Method 525.2 (rather than Method 3535A) should be used for preparation of drinking water samples. For polychlorinated biphenyls only, EPA Method 508 (rather than Method 3535A) should be used for preparation of drinking water samples. Refer to Appendix A for the appropriate determinative method for these analytes.

Analyte(s)	CAS RN
2-Amino-4,6-dinitrotoluene (2-Am-DNT)	35572-78-2
4-Amino-2,6-dinitrotoluene (4-Am-DNT)	19406-51-0
4-Aminopyridine	504-24-5
Brodifacoum	56073-10-0
Bromadiolone	28772-56-7
Chlorfenvinphos	470-90-6
3-Chloro-1,2-propanediol	96-24-2
Chlorosarin	1445-76-7
Chlorosoman	7040-57-5
Chlorpyrifos	2921-88-2
Crimidine	535-89-7
Cyclohexyl sarin (GF)	329-99-7
Dichlorvos	62-73-7
Dicrotophos	141-66-2
Diesel Range Organics	NA
Diisopropyl methylphosphonate (DIMP)	1445-75-6
Dimethylphosphite	868-85-9
Dimethylphosphoramidic acid	33876-51-6
3,5-Dinitroaniline (3,5-DNA)	618-87-1
1,3-Dinitrobenzene (1,3-DNB)	99-65-0
2,4-Dinitrotoluene (2,4-DNT)	121-14-2
2,6-Dinitrotoluene (2,6-DNT)	606-20-2
Diphacinone	82-66-6
1,4-Dithiane	505-29-3
EA2192 [Diisopropylaminoethyl methylthiophosphonate]	73207-98-4
Ethyl methylphosphonic acid (EMPA)	1832-53-7
Ethyldichloroarsine (ED)	598-14-1

Analyte(s)	CAS RN
N-Ethyldiethanolamine (EDEA)	139-87-7
Fenamiphos	22224-92-6
Fentanyl	437-38-7
Hexahydro-1,3,5-trinitro-1,3,5-triazine (RDX)	121-82-4
Hexamethylenetriperoxidediamine (HMTD)	283-66-9
Isopropyl methylphosphonic acid (IMPA)	1832-54-8
Methamidophos	10265-92-6
Methyl hydrazine	60-34-4
Methyl parathion	298-00-0
Methyl-2,4,6-trinitrophenylnitramine (Tetryl)	479-45-8
N-Methyldiethanolamine (MDEA)	105-59-9
1-Methylethyl ester ethylphosphonofluoridic acid (GE)	1189-87-3
Methylphosphonic acid (MPA)	993-13-5
Mevinphos	7786-34-7
Mustard, nitrogen (HN-1) [bis(2-chloroethyl)ethylamine]	538-07-8
Mustard, nitrogen (HN-2) [2,2'-dichloro-N-methyldiethylamine N,N-bis(2-chloroethyl)methylamine]	51-75-2
Mustard, nitrogen (HN-3) [tris(2-chloroethyl)amine]	555-77-1
Nicotine sulfate	54-11-5
Nitrobenzene (NB)	98-95-3
Nitroglycerin (NG)	55-63-0
2-Nitrotoluene (2-NT)	88-72-2
3-Nitrotoluene (3-NT)	99-08-1
4-Nitrotoluene (4-NT)	99-99-0
Octahydro-1,3,5,7-tetranitro-1,3,5,7-tetrazocine (HMX)	2691-41-0
Parathion	56-38-2
Pentaerythritol tetranitrate (PETN)	78-11-5
Phencyclidine	77-10-1
Phenol	108-95-2
Phorate	298-02-2
Phosphamidon	13171-21-6
Pinacolyl methyl phosphonic acid (PMPA)	616-52-4

Analyte(s)	CAS RN
R-33 (VR) [methylphosphonothioic acid, S-[2-(diethylamino)ethyl] O-2-methylpropyl ester]	159939-87-4
Semivolatile Organic Compounds, NOS ¹	NA
Soman (GD)	96-64-0
Strychnine	57-24-9
Tabun (GA)	77-81-6
Tear gas (CS) [chlorobenzylidene malonitrile]	2698-41-1
Tetraethyl pyrophosphate	107-49-3
Tetramethylenedisulfotetramine ²	80-12-6
Thiodiglycol (TDG)	111-48-8
1,4-Thioxane	15980-15-1
Triacetone triperoxide (TATP)	17088-37-8
Triethanolamine (TEA)	102-71-6
Trimethyl phosphite	121-45-9
1,3,5-Trinitrobenzene (1,3,5-TNB)	99-35-4
2,4,6-Trinitrotoluene (2,4,6-TNT)	118-96-7
VE [phosphonothioic acid, ethyl-, S-(2-(diethylamino)ethyl) O-ethyl ester]	21738-25-0
VG [phosphonothioic acid, S-(2-(diethylamino)ethyl) O,O-diethyl ester]	78-53-5
VM [phosphonothioic acid, methyl-, S-(2-(diethylamino)ethyl) O-ethyl ester]	21770-86-5

¹ NOS = Not otherwise specified

² This analyte may require SPE extraction using acetone or methyl ethylketone

This method describes a procedure for isolating target organic analytes from aqueous and liquid samples using solid-phase extraction (SPE) media. Sample preparation procedures vary by analyte group. Following any necessary pH adjustment, a measured volume of sample is extracted by passing it through the solid-phase extraction medium (disks or cartridges), which is held in an extraction device designed for vacuum filtration of the sample. Target analytes are eluted from the solid-phase media using an appropriate solvent which is collected in a receiving vessel. The resulting solvent extract is dried using sodium sulfate and concentrated, as needed.

Source: EPA Method 3535A (SW-846): Solid-Phase Extraction (SPE), Revision 1, 1998.
<http://www.epa.gov/epaoswer/hazwaste/test/pdfs/3535a.pdf>

5.2.19 EPA Method 3541 (SW-846): Automated Soxhlet Extraction

This method should be used for **preparation** of solid samples for the contaminants identified below and listed in Appendix A. Refer to Appendix A for the appropriate determinative method for these analytes.

Analyte(s)	CAS RN
Brodifacoum	56073-10-0
Bromadiolone	28772-56-7
Chlorfenvinphos	470-90-6
3-Chloro-1,2-propanediol	96-24-2
Chlorosarin	1445-76-7
Chlorosoman	7040-57-5
Chlorpyrifos	2921-88-2
Cyclohexyl sarin (GF)	329-99-7
Diesel Range Organics	NA
Diphacinone	82-66-6
N-Ethyldiethanolamine (EDEA)	139-87-7
Fenamiphos	22224-92-6
Fentanyl	437-38-7
Methamidophos	10265-92-6
Methyl hydrazine	60-34-4
N-Methyldiethanolamine (MDEA)	105-59-9
1-Methylethyl ester ethylphosphonofluoridic acid (GE)	1189-87-3
Mustard, nitrogen (HN-1) [bis(2-chloroethyl)ethylamine]	538-07-8
Mustard, nitrogen (HN-2) [2,2'-dichloro-N-methyldiethylamine N,N-bis(2-chloroethyl)methylamine]	51-75-2
Mustard, nitrogen (HN-3) [tris(2-chloroethyl)amine]	555-77-1
Organophosphate Pesticides, NOS ¹	NA
Parathion	56-38-2
Phenol	108-95-2
Phosphamidon	13171-21-6
R-33 (VR) [methylphosphonothioic acid, S-[2-(diethylamino)ethyl] O-2-methylpropyl ester]	159939-87-4
Semivolatile Organic Compounds, NOS ¹	NA
Tear gas (CS) [chlorobenzylidene malonitrile]	2698-41-1
Tetramethylenedisulfotetramine	80-12-6

Analyte(s)	CAS RN
Triethanolamine (TEA)	102-71-6
VE [phosphonothioic acid, ethyl-, S-(2-(diethylamino)ethyl) O-ethyl ester]	21738-25-0
VG [phosphonothioic acid, S-(2-(diethylamino)ethyl) O,O-diethyl ester]	78-53-5
VM [phosphonothioic acid, methyl-, S-(2-(diethylamino)ethyl) O-ethyl ester]	21770-86-5

¹ NOS = Not otherwise specified

Approximately 10 g of solid sample is mixed with an equal amount of anhydrous sodium sulfate and placed in an extraction thimble or between two plugs of glass wool. After adding the appropriate surrogate amount, the sample is extracted using an appropriate solvent in an automated Soxhlet extractor. The extract is dried with sodium sulfate to remove residual moisture, concentrated and exchanged, as necessary, into a solvent compatible with the cleanup or determinative procedure for analysis.

Source: EPA Method 3541 (SW-846): Automated Soxhlet Extraction, Revision 0, 1994.
<http://www.epa.gov/epaoswer/hazwaste/test/pdfs/3541.pdf>

5.2.20 EPA Method 3545A (SW-846): Pressurized Fluid Extraction (PFE)

This method should be used for **preparation** of solid samples for the contaminants identified below and listed in Appendix A. Refer to Appendix A for the appropriate determinative method for these analytes.

Analyte(s)	CAS RN
Brodifacoum	56073-10-0
Bromadiolone	28772-56-7
Chlorfenvinphos	470-90-6
3-Chloro-1,2-propanediol	96-24-2
Chloropicrin	76-06-2
Chlorosarin	1445-76-7
Chlorosoman	7040-57-5
Chlorpyrifos	2921-88-2
Crimidine	535-89-7
Cyclohexyl sarin (GF)	329-99-7
Dichlorvos	62-73-7
Dicrotophos	141-66-2
Diesel Range Organics	NA
Diisopropyl methylphosphonate (DIMP)	1445-75-6
Dimethylphosphite	868-85-9

Analyte(s)	CAS RN
Dimethylphosphoramidic acid	33876-51-6
Diphacinone	82-66-6
1,4-Dithiane	505-29-3
EA2192 [Diisopropylaminoethyl methylthiophosphonate]	73207-98-4
Ethyl methylphosphonic acid (EMPA)	1832-53-7
Ethylchloroarsine (ED)	598-14-1
N-Ethyl-diethanolamine (EDEA)	139-87-7
Fenamiphos	22224-92-6
Fentanyl	437-38-7
Isopropyl methylphosphonic acid (IMPA)	1832-54-8
Methamidophos	10265-92-6
Methyl hydrazine	60-34-4
Methyl parathion	298-00-0
N-Methyl-diethanolamine (MDEA)	105-59-9
1-Methylethyl ester ethylphosphonofluoridic acid (GE)	1189-87-3
Methylphosphonic acid (MPA)	993-13-5
Mevinphos	7786-34-7
Mustard, nitrogen (HN-1) [bis(2-chloroethyl)ethylamine]	538-07-8
Mustard, nitrogen (HN-2) [2,2'-dichloro-N-methyldiethylamine N,N-bis(2-chloroethyl)methylamine]	51-75-2
Mustard, nitrogen (HN-3) [tris(2-chloroethyl)amine]	555-77-1
Nicotine sulfate	54-11-5
Organophosphate Pesticides, NOS ¹	NA
Parathion	56-38-2
Phencyclidine	77-10-1
Phenol	108-95-2
Phorate	298-02-2
Phosphamidon	13171-21-6
Pinacolyl methyl phosphonic acid (PMPA)	616-52-4
R-33 (VR) [methylphosphonothioic acid, S-[2-(diethylamino)ethyl] O-2-methylpropyl ester]	159939-87-4

Analyte(s)	CAS RN
Semivolatile Organic Compounds, NOS ¹	NA
Soman (GD)	96-64-0
Strychnine	57-24-9
Tabun (GA)	77-81-6
Tear gas (CS) [chlorobenzylidene malonitrile]	2698-41-1
Tetraethyl pyrophosphate	107-49-3
Tetramethylenedisulfotetramine	80-12-6
Thiodiglycol (TDG)	111-48-8
1,4-Thioxane	15980-15-1
Triethanolamine (TEA)	102-71-6
Trimethyl phosphite	121-45-9
VE [phosphonothioic acid, ethyl-, S-(2-(diethylamino)ethyl) O-ethyl ester]	21738-25-0
VG [phosphonothioic acid, S-(2-(diethylamino)ethyl) O,O-diethyl ester]	78-53-5
VM [phosphonothioic acid, methyl-, S-(2-(diethylamino)ethyl) O-ethyl ester]	21770-86-5

¹ NOS = Not otherwise specified

Approximately 10 to 30 g of soil sample is prepared for extraction either by air drying the sample, or by mixing the sample with anhydrous sodium sulfate or pelletized diatomaceous earth. The sample is then ground and loaded into the extraction cell. The extraction cell containing the sample is heated to the extraction temperature, pressurized with the appropriate solvent system, and extracted for 5 minutes (or as recommended by the instrument manufacturer). The extract may be concentrated, if necessary, and exchanged into a solvent compatible with the cleanup or determinative step being employed. This method has been validated for solid matrices containing 250 to 12,500 µg/kg of semivolatile organic compounds, 250 to 2500 µg/kg of organophosphorus pesticides, 5 to 250 µg/kg of organochlorine pesticides, 50 to 5000 µg/kg of chlorinated herbicides, 1 to 1400 µg/kg of PCBs, and 1 to 2500 ng/kg of polychlorinated dibenzo-*p*-dioxins/polychlorinated dibenzofurans (PCDDs/PCDFs).

Please note: Sodium sulfate can cause clogging, and air drying or diatomaceous earth may be preferred.

Source: EPA Method 3545A (SW-846): Pressurized Fluid Extraction (PFE), Revision 1, 1998.

<http://www.epa.gov/epaoswer/hazwaste/test/pdfs/3545a.pdf>

5.2.21 EPA Method 3571 (SW-846): Extraction of Solid and Aqueous Samples for Chemical Agents

This method should be used for **preparation** of solid, organic solid, nonaqueous liquid, aqueous liquid and drinking water samples for the contaminant identified below and listed in Appendix A. Refer to Appendix A for the appropriate determinative method for these analytes.

Analyte(s)	CAS RN
Mustard, sulfur / Mustard gas (HD)_	505-60-2
Sarin (GB)	107-44-8
VX	50782-69-9

This method provides procedures for sample collection and extraction of the referenced compounds from solids and aqueous samples. A separate extract is required for each agent to be measured. Glacial acetic acid is added as a preservative to samples being assayed for GB and glacial acetic acid/sodium chloride is a preservative for samples assayed for HD. No preservative is added for VX. Samples are extracted with 10% isopropanol in dichloromethane by vortex mixing and filtered, if necessary. An optional water wash is included for VX that back-extracts the compound from heavy organics that could interfere with the assay. An optional column cleanup procedure is described to separate GB from heavy organics, if needed. Solvents are used to elute the extract first through the Carboxprep90 column, then the silica column.

Source: EPA Method 3571: Extraction of Solid and Aqueous Samples for Chemical Agents.

<http://www.epa.gov/epaoswer/hazwaste/test/pdfs/3571.pdf>

5.2.22 EPA Method 3580A (SW-846): Waste Dilution

This method should be used for **preparation** of non-aqueous liquid/organic solid samples for the contaminants identified below and listed in Appendix A. Refer to Appendix A for the appropriate determinative method for these analytes.

Analyte(s)	CAS RN
Brodifacoum	56073-10-0
Bromadiolone	28772-56-7
Chlorfenvinphos	470-90-6
3-Chloro-1,2-propanediol	96-24-2
Chloropicrin	76-06-2
Chlorosarin	1445-76-7
Chlorosoman	7040-57-5
Chlorpyrifos	2921-88-2
Crimidine	535-89-7
Cyclohexyl sarin (GF)	329-99-7
Dichlorvos	62-73-7
Dicrotophos	141-66-2

Analyte(s)	CAS RN
Diesel Range Organics	NA
Diisopropyl methylphosphonate (DIMP)	1445-75-6
Dimethylphosphite	868-85-9
Dimethylphosphoramidic acid	33876-51-6
Diphacinone	82-66-6
1,4-Dithiane	505-29-3
EA2192 [Diisopropylaminoethyl methylthiophosphonate]	73207-98-4
Ethyl methylphosphonic acid (EMPA)	1832-53-7
Ethylchloroarsine (ED)	598-14-1
N-Ethyldiethanolamine (EDEA)	139-87-7
Fenamiphos	22224-92-6
Fentanyl	437-38-7
Isopropyl methylphosphonic acid (IMPA)	1832-54-8
Methamidophos	10265-92-6
Methyl hydrazine	60-34-4
Methyl parathion	298-00-0
N-Methyldiethanolamine (MDEA)	105-59-9
1-Methylethyl ester ethylphosphonofluoridic acid (GE)	1189-87-3
Methylphosphonic acid (MPA)	993-13-5
Mevinphos	7786-34-7
Mustard, nitrogen (HN-1) [bis(2-chloroethyl)ethylamine]	538-07-8
Mustard, nitrogen (HN-2) [2,2'-dichloro-N-methyldiethylamine N,N-bis(2-chloroethyl)methylamine]	51-75-2
Mustard, nitrogen (HN-3) [tris(2-chloroethyl)amine]	555-77-1
Nicotine sulfate	54-11-5
Organophosphate Pesticides, NOS ¹	NA
Parathion	56-38-2
Phencyclidine	77-10-1
Phenol	108-95-2
Phorate	298-02-2
Phosphamidon	13171-21-6

Analyte(s)	CAS RN
Pinacolyl methyl phosphonic acid (PMPA)	616-52-4
R-33 (VR) [methylphosphonothioic acid, S-[2-(diethylamino)ethyl] O-2-methylpropyl ester]	159939-87-4
Semivolatile Organic Compounds, NOS ¹	NA
Soman (GD)	96-64-0
Strychnine	57-24-9
Tabun (GA)	77-81-6
Tear gas (CS) [chlorobenzylidene malonitrile]	2698-41-1
Tetraethyl pyrophosphate	107-49-3
Tetramethylenedisulfotetramine	80-12-6
Thiodiglycol (TDG)	111-48-8
1,4-Thioxane	15980-15-1
Triethanolamine (TEA)	102-71-6
Trimethyl phosphite	121-45-9
VE [phosphonothioic acid, ethyl-, S-(2-(diethylamino)ethyl) O-ethyl ester]	21738-25-0
VG [phosphonothioic acid, S-(2-(diethylamino)ethyl) O,O-diethyl ester]	78-53-5
VM [phosphonothioic acid, methyl-, S-(2-(diethylamino)ethyl) O-ethyl ester]	21770-86-5

¹ NOS = Not otherwise specified

This method describes solvent dilution of a non-aqueous waste sample prior to cleanup and/or analysis. One gram of sample is weighed into a capped tube and the sample is diluted to 10.0 mL with an appropriate solvent. The method is designed for wastes that may contain organic chemicals at a concentration greater than 20,000 mg/kg and that are soluble in the dilution solvent.

Source: EPA Method 3580A (SW-846): Waste Dilution, Revision 1, 1992.

<http://www.epa.gov/epaoswer/hazwaste/test/pdfs/3580a.pdf>

5.2.23 EPA Method 3585 (SW-846): Waste Dilution for Volatile Organics

This method should be used for **preparation** of non-aqueous liquid/organic solid samples for the contaminants identified below and listed in Appendix A. Refer to Appendix A for the appropriate determinative method for these analytes.

Analyte(s)	CAS RN
Allyl alcohol	107-18-6
Carbon disulfide	75-15-0
2-Chloroethanol	107-07-3

Analyte(s)	CAS RN
Cyanogen chloride	506-77-4
1,2-Dichloroethane	107-06-2
Ethylene oxide	75-21-8
Gasoline Range Organics	NA
Kerosene	64742-81-0
Propylene oxide	75-56-9
The following analytes should be prepared by this method (and determined by the corresponding SW-846 Method 8260C) only if problems (e.g., insufficient recovery, interferences) occur when using the sample preparation/determinative techniques identified for these analytes in Appendix A.	
1,4-Dithiane	505-29-3
1,4-Thioxane	15980-15-1

This method describes solvent dilution of a non-aqueous waste sample prior to direct injection analysis. It is designed for use in conjunction with GC or GC-MS analysis of wastes that may contain organic chemicals at a concentration greater than 1 mg/kg and that are soluble in the dilution solvent. Highly contaminated or highly complex samples may be diluted prior to analysis for volatiles using direct injection. One gram of sample is weighed into a capped tube or volumetric flask. The sample is diluted to 2.0 to 10.0 mL with *n*-hexadecane or other appropriate solvent. Diluted samples are injected into the GC or GC-MS for analysis.

Source: EPA Method 3585 (SW-846): Waste Dilution for Volatile Organics, Revision 0, 1996.
<http://www.epa.gov/epaoswer/hazwaste/test/pdfs/3585.pdf>

5.2.24 EPA Method 5030C (SW-846): Purge-and-Trap for Aqueous Samples

This method should be used for **preparation** of aqueous liquid and/or drinking water samples for the contaminants identified below and listed in Appendix A. **Note:** For carbon disulfide, 1,2-dichloroethane, and volatile organic compounds only, EPA Method 524.2 (rather than Method 5030C) should be used for preparation of drinking water samples. Refer to Appendix A for the appropriate determinative method for these analytes.

Analyte(s)	CAS RN
Allyl alcohol	107-18-6
Carbon disulfide	75-15-0
2-Chloroethanol	107-07-3
Cyanogen chloride	506-77-4
1,2-Dichloroethane	107-06-2
Ethylene oxide	75-21-8
Gasoline Range Organics	NA
Kerosene	64742-81-0
Propylene oxide	75-56-9

Analyte(s)	CAS RN
The following analytes should be prepared by this method (and determined by the corresponding SW-846 Method 8260C) only if problems (e.g., insufficient recovery, interferences) occur when using the sample preparation/determinative techniques identified for these analytes in Appendix A.	
1,4-Dithiane	505-29-3
1,4-Thioxane	15980-15-1

This method describes a purge-and-trap procedure for the analysis of volatile organic compounds (VOCs) in aqueous liquid samples and water miscible liquid samples. An inert gas is bubbled through a portion of the aqueous liquid sample at ambient temperature, and the volatile components are transferred from the aqueous liquid phase to the vapor phase. The vapor is swept through a sorbent column where the volatile components are adsorbed. After purging is completed, the sorbent column is heated and backflushed with inert gas to desorb the components onto a gas chromatographic column.

Source: EPA Method 5030C (SW-846): Purge-and-Trap for Aqueous Samples, Revision 3, 2003.
<http://www.epa.gov/epaoswer/hazwaste/test/pdfs/5030c.pdf>

5.2.25 EPA Method 5035A (SW-846): Closed-System Purge-and-Trap and Extraction for Volatile Organics in Soil and Waste Samples

This method should be used for **preparation** of solid samples for the contaminants identified below and listed in Appendix A. Refer to Appendix A for the appropriate determinative method for these analytes.

Analyte(s)	CAS RN
Allyl alcohol	107-18-6
Carbon disulfide	75-15-0
2-Chloroethanol	107-07-3
Cyanogen chloride	506-77-4
1,2-Dichloroethane	107-06-2
Ethylene oxide	75-21-8
Gasoline Range Organics	NA
Kerosene	64742-81-0
Propylene oxide	75-56-9
The following analytes should be prepared by this method (and determined by the corresponding SW-846 Method 8260C) only if problems (e.g., insufficient recovery, interferences) occur when using the sample preparation/determinative techniques identified for these analytes in Appendix A.	
1,4-Dithiane	505-29-3
1,4-Thioxane	15980-15-1

This method describes a closed-system purge-and-trap process for analysis of volatile organic compounds (VOCs) in solid samples containing low levels (0.5 to 200 µg/kg) of VOCs. The method also provides specific procedures for preparation of samples containing high levels (>200 µg/kg) of VOCs. For low-level VOCs, a 5-g sample is collected into a vial that is placed into an autosampler device. Reagent water, surrogates, and internal standards are added automatically, and the vial is heated to 40°C. The

volatiles are purged into an appropriate trap using an inert gas combined with sample agitation. When purging is complete, the trap is heated and backflushed with helium to desorb the trapped sample components into a gas chromatograph for analysis. For high-level VOCs, samples are either collected into a vial that contains a water-miscible organic solvent or a portion of sample is removed from the vial and dispersed in a water-miscible solvent. An aliquot of the solvent is added to reagent water, along with surrogates and internal standards, then purged and analyzed using an appropriate determinative method (e.g., Method 8015C or 8260C (SW-846)).

Source: EPA Method 5035A (SW-846): Closed-System Purge-and-Trap and Extraction for Volatile Organics in Soil and Waste Samples, Draft Revision 1, 2002.

http://www.epa.gov/epaoswer/hazwaste/test/pdfs/5035a_r1.pdf

5.2.26 EPA Method 6010C (SW-846): Inductively Coupled Plasma - Atomic Emission Spectrometry

This method should be used for **analysis** of solid and non-aqueous liquid/organic solid samples for the contaminants identified below and listed in Appendix A. Appropriate sample preparation techniques should be used prior to analysis (refer to Appendix A). Note: Osmium tetroxide and titanium tetrachloride are not of concern in non-aqueous liquid/organic solid samples.

Analyte(s)	CAS RN
Ammonium metavanadate	7803-55-6
Arsenic, Total	740-38-2
Arsenic trioxide	1327-53-3
2-Chlorovinylarsonous acid (CVAA)	85090-33-1
Lewisite 1 (L-1) [2-chlorovinylidichloroarsine] ¹	541-25-3
Lewisite 2 (L-2) [bis(2-chlorovinyl)chloroarsine]	40334-69-8
Lewisite 3 (L-3) [tris(2-chlorovinyl)arsine]	40334-70-1
Lewisite Oxide	1306-02-1
Osmium tetroxide	20816-12-0
Sodium arsenite	7784-46-5
Thallium sulfate	10031-59-1
Titanium tetrachloride	7550-45-0
Vanadium pentoxide	1314-62-1

¹ Laboratory testing is currently underway for speciation of lewisite 1 using GC-MS techniques.

This method determines arsenic trioxide, lewisite, lewisite degradation products, and sodium arsenite as total arsenic; osmium tetroxide as osmium; thallium sulfate as thallium; titanium tetrachloride as titanium; and ammonium metavanadate and vanadium pentoxide as total vanadium. Soil samples (prepared using SW-846 Methods 3050B), and non-aqueous liquid/organic solid samples (prepared using SW-846 Methods 3031) are analyzed by Inductively Coupled Plasma - Atomic Emission Spectrometry (ICP-AES). Detection limits vary with each analyte. Estimated instrument detection limits for arsenic and titanium are 30 µg/L and 5.0 µg/L, respectively. The upper end of the analytical range may be extended by sample dilution.

Source: EPA Method 6010C (SW-846): Inductively Coupled Plasma-Atomic Emission Spectrometry, Revision 3, 2000. <http://www.epa.gov/epaoswer/hazwaste/test/pdfs/6010c.pdf>

5.2.27 EPA Method 6020A (SW-846): Inductively Coupled Plasma - Mass Spectrometry

This method should be used for **analysis** of solid and non-aqueous liquid/organic solid samples for the contaminants identified below and listed in Appendix A. Appropriate sample preparation techniques should be used prior to analysis (refer to Appendix A). Note: Titanium tetrachloride is not of concern in non-aqueous liquid/organic solid samples.

Analyte(s)	CAS RN
Ammonium metavanadate	7803-55-6
Arsenic, Total	740-38-2
Arsenic trioxide	1327-53-3
2-Chlorovinylarsonous acid (CVAA)	85090-33-1
Lewisite 1 (L-1) [2-chlorovinylchloroarsine] ¹	541-25-3
Lewisite 2 (L-2) [bis(2-chlorovinyl)chloroarsine]	40334-69-8
Lewisite 3 (L-3) [tris(2-chlorovinyl)arsine]	40334-70-1
Lewisite oxide	1306-02-1
Sodium arsenite	7784-46-5
Thallium sulfate	10031-59-1
Titanium tetrachloride	7550-45-0
Vanadium pentoxide	1314-62-1

¹ Laboratory testing is currently underway for speciation of lewisite 1 using GC-MS techniques.

This method will determine arsenic trioxide, lewisite, lewisite degradation products, and sodium arsenite as total arsenic. The method also will determine thallium sulfate as total thallium, titanium tetrachloride as titanium, and ammonium metavanadate and vanadium pentoxide as total vanadium. Aqueous samples (prepared using SW-846 Method 5050), soil samples (prepared using SW-846 Methods 3050B or 5050), non-aqueous liquid/organic solid samples (prepared using SW-846 Methods 3050B or 3031), and air filter/particle samples (prepared using IO Method 3.5) are analyzed by Inductively Coupled Plasma - Mass Spectrometry. Instrument detection limits, sensitivities, and linear ranges vary with sample type, instrumentation, and operation conditions. In relatively simple sample types, detection limits will generally be below 0.1 µg/L. Less sensitive elements, such as arsenic, may have detection limits of 1.0 µg/L or higher. The upper end of the analytical range may be extended by sample dilution.

Source: EPA Method 6020A (SW-846): Inductively Coupled Plasma-Mass Spectrometry, Revision 1, 1998. <http://www.epa.gov/epaoswer/hazwaste/test/pdfs/6020a.pdf>

5.2.28 EPA Method 7010 (SW-846): Graphite Furnace Atomic Absorption Spectrophotometry

This method should be used for **analysis** of solid samples for the contaminant identified below and listed in Appendix A. Appropriate sample preparation techniques should be used prior to analysis (refer to Appendix A).

Analyte(s)	CAS RN
Arsine	7784-42-1

This method determines arsine as arsenic in environmental samples. Soil samples (prepared using SW-846 Method 3050B) are analyzed by Graphite Furnace Atomic Absorption Spectrophotometry (GFAA). A representative aliquot of the sample is placed in the graphite tube in the furnace, evaporated to dryness, charred, and atomized. Detection limits vary with each sample type and instrument used. The analytical range may be extended by sample dilution.

Source: EPA Method 7010 (SW-846): Graphite Furnace Atomic Absorption Spectrophotometry, Revision 0, 1998. <http://www.epa.gov/epaoswer/hazwaste/test/pdfs/7010.pdf>

5.2.29 EPA Method 7470A (SW-846): Mercury in Liquid Wastes (Manual Cold-Vapor Technique)

This method is not currently listed for any of the analyte/sample type combinations included in Appendix A. If problems occur when using EPA SW-846 Method 7473, then this method should be used for **preparation** and **analysis** of aqueous liquid samples for the contaminant identified below and listed in Appendix A. (See Footnote 12 of Appendix A.)

Analyte(s)	CAS RN
Mercury, Total	7439-97-6
Methoxyethylmercuric acetate	151-38-2

A 100-mL aqueous sample is digested with acids, permanganate solution, persulfate solution, and heat. The sample is cooled and reduced with hydroxylamine-sodium chloride solution. Just prior to analysis, the sample is treated with Sn(II), reducing the mercury to Hg(0). The reduced sample is sparged and the mercury vapor is analyzed by cold vapor atomic absorption. The detection limit for the method is less than 0.2 µg/L. Chloride and copper are potential interferences.

Source: EPA Method 7470A (SW-846): Mercury in Liquid Waste (Manual Cold-Vapor Technique), Revision 1, 1994. <http://www.epa.gov/epaoswer/hazwaste/test/pdfs/7470a.pdf>

5.2.30 EPA Method 7471B (SW-846): Mercury in Solid or Semisolid Wastes (Manual Cold-Vapor Technique)

This method is not currently listed for any of the analyte/sample type combinations included in Appendix A. If problems occur when using EPA SW-846 Method 7473, then this method should be used for **preparation** and **analysis** of solid samples for the contaminant identified below and listed in Appendix A. (See Footnote 12 of Appendix A.)

Analyte(s)	CAS RN
Mercury, Total	7439-97-6
Methoxyethylmercuric acetate	151-38-2

A 0.5-g to 0.6-g sample is digested with aqua regia, permanganate solution, and heat. The sample is cooled and reduced with hydroxylamine-sodium chloride solution. Just prior to analysis, the sample is treated with Sn(II), reducing the mercury to Hg(0). The reduced sample is sparged and the mercury vapor is analyzed by cold vapor atomic absorption. Chloride and copper are potential interferences.

Source: EPA Method 7471B (SW-846): Mercury in Solid or Semisolid Waste (Manual Cold-Vapor Technique), Revision 2, 1998. <http://www.epa.gov/epaoswer/hazwaste/test/pdfs/7471b.pdf>

5.2.31 EPA Method 7473 (SW-846): Mercury in Solids and Solutions by Thermal Decomposition, Amalgamation, and Atomic Absorption Spectrophotometry

This method should be used for **preparation** and **analysis** of solid and aqueous liquid samples for the contaminants identified below and listed in Appendix A.

Analyte(s)	CAS RN
Mercury, Total ¹	7439-97-6
Methoxyethylmercuric acetate ¹	151-38-2

¹ If equipment is not available, use CVAA Methods 7471B (EPA SW-846) for solid samples and 7470A (EPA SW-846) for aqueous liquid samples.

Controlled heating in an oxygenated decomposition furnace is used to liberate mercury from solid and aqueous samples. The sample is dried and then thermally and chemically decomposed within the furnace. The decomposition products are carried by flowing oxygen to the catalytic section of the furnace, where oxidation is completed and halogens and nitrogen/sulfur oxides are trapped. The remaining decomposition products are then carried to an amalgamator that selectively traps mercury. After the system is flushed with oxygen to remove any remaining gases or decomposition products, the amalgamator is rapidly heated, releasing mercury vapor. Flowing oxygen carries the mercury vapor through absorbance cells positioned in the light path of a single wavelength atomic absorption spectrophotometer. Absorbance (peak height or peak area) is measured at 253.7 nm as a function of mercury concentration. The typical working range for this method is 0.05 – 600 ng. The instrument detection limit (IDL) is 0.01-ng total mercury.

Source: EPA Method 7473 (SW-846): Mercury in Solids and Solutions by Thermal Decomposition, Amalgamation, and Atomic Absorption Spectrophotometry, Revision 0, 1998. <http://www.epa.gov/epaoswer/hazwaste/test/pdfs/7473.pdf>

5.2.32 EPA Method 7580 (SW-846): White Phosphorus by Solvent Extraction and Gas Chromatography (GC)

This method should be used for **preparation** and **analysis** of solid, non-aqueous liquid/organic solid, aqueous liquid, and drinking water samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
White Phosphorus	12185-10-3

Method 7580 may be used to determine the concentration of white phosphorus in soil, sediment, and water samples using solvent extraction and gas chromatography (GC). Water samples are extracted by one of two procedures, depending on the sensitivity required. For the more sensitive procedure, a 500-mL water sample is extracted with 50 mL of diethyl ether. The extract is concentrated by back extraction with reagent water, yielding a final extract volume of approximately 1.0 mL. A 1.0 µL aliquot of this extract is injected into a GC equipped with a nitrogen-phosphorus detector (NPD). This procedure

provides sensitivity on the order of 0.01 µg/L. Wet soil or sediment samples are analyzed by extracting a 40 g wet-weight aliquot of the sample with a mixture of 10.0 mL degassed reagent water and 10.0 mL isooctane. The extraction is performed in a glass jar on a platform shaker for 18 hours. A 1.0 µL aliquot of the extract is analyzed by GC-NPD, providing sensitivity on the order of 1 µg/kg.

Source: EPA Method 7580 (SW-846): White Phosphorus (P₄) by Solvent Extraction and Gas Chromatography, Revision 0, 1996. <http://www.epa.gov/sw-846/pdfs/7580.pdf>

5.2.33 EPA Method 8015C (SW-846): Nonhalogenated Organics Using GC-FID

This method should be used for **analysis** of solid, non-aqueous liquid/organic solid, aqueous liquid, and drinking water samples for the contaminants identified below and listed in Appendix A. Appropriate sample preparation techniques should be used prior to analysis (refer to Appendix A).

Analyte(s)	CAS RN
Diesel Range Organics	NA
Gasoline Range Organics	NA
Kerosene	64742-81-0

This method provides gas chromatographic (GC) conditions for the detection of certain nonhalogenated volatile and semivolatile organic compounds. Depending on the analytes of interest, samples may be introduced into the GC by a variety of techniques including purge-and-trap, direct injection of aqueous liquid samples, and solvent extraction. An appropriate column and temperature program are used in the gas chromatograph to separate the organic compounds. Detection is achieved by a flame ionization detector (FID). The method allows the use of packed or capillary columns for the analysis and confirmation of the non-halogenated individual analytes. The estimated method detection limits vary with each analyte and range between 2 and 48 µg/L for aqueous liquid samples. The method detection limits in other matrices have not been evaluated for this method. The analytical range depends on the target analyte(s) and the instrument used.

Source: EPA Method 8015C (SW-846): Nonhalogenated Organics Using GC/FID, Revision 3, 2000. <http://www.epa.gov/epaoswer/hazwaste/test/pdfs/8015c.pdf>

5.2.34 EPA Method 8260C (SW-846): Volatile Organic Compounds by Gas Chromatography-Mass Spectrometry (GC-MS)

This method should be used for **analysis** of solid, non-aqueous liquid/organic solid, aqueous liquid, and drinking water samples for the contaminants identified below and listed in Appendix A. Appropriate sample preparation techniques should be used prior to analysis (refer to Appendix A). Note: For carbon disulfide and 1,2-dichloroethane only, EPA Method 524.2 (rather than 8260C) should be used for analysis of drinking water samples.

Analyte(s)	CAS RN
Allyl alcohol	107-18-6
Carbon disulfide	75-15-0
2-Chloroethanol	107-07-3
Cyanogen chloride	506-77-4
1,2-Dichloroethane	107-06-2

Analyte(s)	CAS RN
Ethylene oxide	75-21-8
Propylene oxide	75-56-9
The following analytes should be determined by this method (and corresponding sample preparation methods) only if problems (e.g., insufficient recovery, interferences) occur when using the sample preparation/determinative techniques identified for these analytes in Appendix A.	
1,4-Dithiane	505-29-3
1,4-Thioxane	15980-15-1

Volatile compounds are introduced into a gas chromatograph (GC) by purge-and-trap or other procedures (see Section 1.2 in Method 8260C). The analytes can be introduced directly to a wide-bore capillary column or cryofocused on a capillary pre-column before being flash evaporated to a narrow-bore capillary for analysis. Alternatively, the effluent from the trap is sent to an injection port operating in the split mode for injection to a narrow-bore capillary column. The column is temperature-programmed to separate the analytes, which are then detected with a mass spectrometer (MS) interfaced to the GC. Analytes eluted from the capillary column are introduced into the mass spectrometer via a jet separator or a direct connection. The estimated quantitation limit (EQL) for an individual analyte is dependent on the instrument as well as the choice of sample preparation/introduction method. Using standard quadrupole instrumentation and the purge-and-trap technique, estimated quantitation limits are 5 µg/kg (wet weight) for soil/sediment samples and 5 µg/L for ground water. Somewhat lower limits may be achieved using an ion trap mass spectrometer or other instrumentation of improved design. No matter which instrument is used, EQLs will be proportionately higher for sample extracts and samples that require dilution or when a reduced sample size is used to avoid saturation of the detector.

Source: EPA Method 8260C (SW-846): Volatile Organic Compounds by Gas Chromatography-Mass Spectrometry (GC-MS), Revision 3, 2006. <http://www.epa.gov/epaoswer/hazwaste/test/pdfs/8260c.pdf>

5.2.35 EPA Method 8270D (SW-846): Semivolatile Organic Compounds by Gas Chromatography/Mass Spectrometry (GC-MS)

This method should be used for **analysis** of solid, non-aqueous liquid/organic solid, aqueous liquid, and drinking water samples for the contaminants identified below and listed in Appendix A. Appropriate sample preparation techniques should be used prior to analysis (refer to Appendix A). **Note:** For dichlorvos, fenamiphos, mevinphos, and semivolatile organic compounds only, EPA Method 525.2 (rather than Method 8270D) should be used for analysis of drinking water samples. For organophosphate pesticides only, EPA Methods 614 and 507 should be used for analysis of aqueous liquid and drinking water samples, respectively. For chloropicrin only, EPA Method 551.1 should be used for analysis of aqueous liquid and drinking water samples.

Analyte(s)	CAS RN
Chlorfenvinphos	470-90-6
3-Chloro-1,2-propanediol ¹	96-24-2
Chloropicrin ²	76-06-2
Chlorosarin	1445-76-7
Chlorosoman	7040-57-5
Chlorpyrifos	2921-88-2
Crimidine ³	535-89-7

Analyte(s)	CAS RN
Cyclohexyl sarin (GF)	329-99-7
Dichlorvos	62-73-7
Dicrotophos	141-66-2
Diisopropyl methylphosphonate (DIMP)	1445-75-6
Dimethylphosphite	868-85-9
1,4-Dithiane ⁴	505-29-3
Ethylchloroarsine (ED)	598-14-1
Fenamiphos	22224-92-6
Methamidophos	10265-92-6
Methyl hydrazine	60-34-4
Methyl parathion	298-00-0
1-Methylethyl ester ethylphosphonofluoridic acid (GE)	1189-87-3
Mevinphos	7786-34-7
Mustard, nitrogen (HN-1) [bis(2-chloroethyl)ethylamine]	538-07-8
Mustard, nitrogen (HN-2) [2,2'-dichloro-N-methyldiethylamine N,N-bis(2-chloroethyl)methylamine]	51-75-2
Mustard, nitrogen (HN-3) [tris(2-chloroethyl)amine]	555-77-1
Mustard, sulfur / Mustard gas (HD) ⁵	505-60-2
Nicotine sulfate	54-11-5
Organophosphate Pesticides, NOS ⁶	NA
Parathion	56-38-2
Phencyclidine	77-10-1
Phenol	108-95-2
Phorate	298-02-2
Phosphamidon	13171-21-6
R-33 (VR) [methylphosphonothioic acid, S-[2-(diethylamino)ethyl] O-2-methylpropyl ester]	159939-87-4
Sarin (GB) ⁵	107-44-8
Semivolatile Organic Compounds, NOS ⁶	NA
Soman (GD)	96-64-0
Strychnine	57-24-9
Tabun (GA)	77-81-6

Analyte(s)	CAS RN
Tear gas (CS) [chlorobenzylidene malonitrile]	2698-41-1
Tetraethyl pyrophosphate	107-49-3
Tetramethylenedisulfotetramine ^{2,7,8}	80-12-6
Thiodiglycol (TDG)	111-48-8
1,4-Thioxane ⁴	15980-15-1
Trimethyl phosphite	121-45-9
VE [phosphonothioic acid, ethyl-, S-(2-(diethylamino)ethyl) O-ethyl ester]	21738-25-0
VG [phosphonothioic acid, S-(2-(diethylamino)ethyl) O,O-diethyl ester]	78-53-5
VM [phosphonothioic acid, methyl-, S-(2-(diethylamino)ethyl) O-ethyl ester]	21770-86-5
VX [O-ethyl-S-(2-diisopropylaminoethyl)methyl-phosphonothiolate] ⁵	50782-69-9
The following analyte should be determined by this method only if LC-MS (electrospray) procedures are not available to the laboratory. These analytes should be analyzed with GC-MS procedures using derivatization based on SW-846 Method 8270D. Sample preparation methods should remain the same. Both electrospray LC-MS and GC/MS derivatization procedures are currently under development.	
Dimethylphosphoramidic acid	33876-51-6
EA2192 [Diisopropylaminoethyl methylthiophosphonate]	73207-98-4
Ethylmethyl phosphonate (EMPA)	1832-53-7
Isopropyl methylphosphonic acid (IMPA)	1832-54-8
Methylphosphonic acid (MPA)	993-13-5
Pinacolyl methyl phosphonic acid (PMPA)	616-52-4

¹ For this analyte, SW-846 Method 8270D must be modified to include a derivatization step.

² This analyte requires determination using an injection port temperature of less than 200°C.

³ If problems occur when using this method, it is recommended that SW-846 Method 8321B be used.

⁴ If problems occur when using this method, it is recommended that SW-846 Method 8260C and appropriate corresponding sample preparation procedures (i.e., 5035A for solid samples, 3585 for non-aqueous liquid/organic solid samples, and 5030C for aqueous liquid and drinking water samples) be used.

⁵ For this analyte, refer to EPA SW-846 Method 8271 for GC-MS conditions

⁶ NOS = Not otherwise specified

⁷ This analyte may require SIM analyses in order to be determined.

⁸ When analyzing for tetramine, the injection temperature must not exceed 250°C (the decomposition temperature of tetramine).

Samples are prepared for analysis by gas chromatography/mass spectrometry using the appropriate sample preparation and, if necessary, sample cleanup procedures. Semivolatile compounds are introduced into the GC-MS by injecting the sample extract into a gas chromatograph (GC) with a narrow-bore fused-silica capillary column. The GC column is temperature-programmed to separate the analytes, which are then detected with a mass spectrometer (MS) connected to the GC. Analytes eluted from the

capillary column are introduced into the MS. For the determination of 3-chloro-1,2-propanediol, dimethylphosphoramidic acid, EA2192, EMPA, IMPA, MPA, and pinacolyl methyl phosphonic acid, a derivatization step is required prior to injection into the GC-MS. The phosphonic acids require derivatization with a trimethylsilyl agent and 3-chloro-1,2-propanediol requires derivatization with a heptafluorobutyl agent. The estimated detection limits vary with each analyte and range between 10 and 1000 µg/L for aqueous liquid samples and 660 and 3300 µg/kg for soil samples. The analytical range depends on the target analyte(s) and the instrument used.

Source: EPA Method 8270D (SW-846): Semivolatile Organic Compounds by Gas Chromatography/Mass Spectrometry (GC/MS), Revision 4, 1998. <http://www.epa.gov/epaoswer/hazwaste/test/pdfs/8270d.pdf>

5.2.36 EPA Method 8315A (SW-846): Determination of Carbonyl Compounds by High Performance Liquid Chromatography (HPLC)

This method should be used for **preparation** and **analysis** of solid, aqueous liquid, and drinking water samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Formaldehyde	50-00-0

A measured volume of aqueous liquid sample (approximately 100 mL), or an appropriate amount of solids extract (approximately 25 g), is buffered to pH 3 and derivatized with 2,4-dinitrophenylhydrazine (DNPH). Using the appropriate extraction technique, the derivatives are extracted using methylene chloride and the extracts are exchanged with acetonitrile prior to HPLC analysis. HPLC conditions are described permitting the separation and measurement of various carbonyl compounds in the extract by absorbance detection at 360 nm. If formaldehyde is the only analyte of interest, the aqueous liquid sample and/or solid sample extract should be buffered to pH 5.0 to minimize the formation of artifact formaldehyde. The method detection limit for formaldehyde varies depending on sample conditions and instrumentation, but is approximately 6.2 µg/L for aqueous liquid samples.

Source: EPA Method 8315A (SW-846): Determination of Carbonyl Compounds by High Performance Liquid Chromatography (HPLC), Revision 1, 1996.
<http://www.epa.gov/epaoswer/hazwaste/test/pdfs/8315a.pdf>

5.2.37 EPA Method 8316 (SW-846): Acrylamide, Acrylonitrile and Acrolein by High Performance Liquid Chromatography (HPLC)

This method should be used for **preparation** and/or **analysis** of solid, non-aqueous liquid/organic solid, aqueous liquid, and drinking water samples for the contaminants identified below and listed in Appendix A. *Note:* Solid and non-aqueous liquid/organic solid samples are extracted with water prior to analysis.

Analyte(s)	CAS RN
Acrylamide	79-06-1
Acrylonitrile	107-13-1
Methyl acrylonitrile	126-98-7

Samples are analyzed by High Performance Liquid Chromatography (HPLC). A 200-µL aliquot is injected onto a C₁₈ reverse-phase column, and compounds in the effluent are detected with an ultraviolet (UV) detector. Solid and non-aqueous liquid/organic solid samples should be water extracted prior to injection. Aqueous liquid and drinking water samples can be directly injected. Acrylamide has an MDL of 10 µg/L; acrylonitrile has an MDL of 20 µg/L.

Source: EPA Method 8316 (SW-846): Acrylamide, Acrylonitrile and Acrolein by High Performance Liquid Chromatography (HPLC). <http://www.epa.gov/epaoswer/hazwaste/test/pdfs/8316.pdf>

5.2.38 EPA Method 8318A (SW-846): N-Methylcarbamates by High Performance Liquid Chromatography (HPLC)

This method should be used for **preparation** and **analysis** of solid, non-aqueous liquid/organic solid, and aqueous liquid samples for the contaminants identified below and listed in Appendix A.

Analyte(s)	CAS RN
Aldicarb (Temik)	116-06-3
Carbofuran (Furadan)	1563-66-2
Methomyl	16752-77-5
Oxamyl	23135-22-0

N-methylcarbamates are extracted from aqueous liquid samples with methylene chloride, and from soils, oily solid waste, and oils with acetonitrile. The extract solvent is exchanged to methanol/ethylene glycol, and the extract is cleaned using a C₁₈ cartridge, filtered, and eluted on a C₁₈ analytical column. After separation, the target analytes are hydrolyzed and derivatized post-column, then quantified fluorometrically. The sensitivity of the method usually depends on the level of interferences present, rather than on instrument conditions. Waste samples with a high level of extractable fluorescent compounds are expected to yield significantly higher detection limits. The estimated method detection limits vary with each analyte and range between 1.7 to 9.4 µg/L for aqueous liquid samples and 10 to 50 µg/kg for soil samples.

Source: EPA Method 8318A (SW-846): N-Methylcarbamates by High Performance Liquid Chromatography (HPLC), Revision 1, 2000. <http://www.epa.gov/epaoswer/hazwaste/test/pdfs/8318a.pdf>

5.2.39 EPA Method 8321B (SW-846): Solvent-Extractable Nonvolatile Compounds by High Performance Liquid Chromatography-Thermospray-Mass Spectrometry (HPLC-TS-MS) or Ultraviolet (UV) Detection

This method should be used for **analysis** of solid, non-aqueous liquid/organic solid, aqueous liquid, and drinking water samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Brodifacoum	56073-10-0
Bromadiolone	28772-56-7
Dimethylphosphoramidic acid ¹	33876-51-6
Diphacinone	82-66-6
EA2192 [Diisopropylaminoethyl methylthiophosphonate] ¹	73207-98-4
Ethylmethyl phosphonate (EMPA) ¹	1832-53-7
N-Ethyldiethanolamine (EDEA)	139-87-7
Fentanyl	437-38-7
Isopropyl methylphosphonic acid (IMPA) ¹	1832-54-8

Analyte(s)	CAS RN
N-Methyldiethanolamine (MDEA)	105-59-9
Methylphosphonic acid (MPA) ¹	993-13-5
Pinacolyl methyl phosphonic acid (PMPA) ¹	616-52-4
Triethanolamine (TEA)	102-71-6
The following analyte should be determined by this method only if problems (e.g., insufficient recovery, interferences) occur when using SW-846 Method 8270D. Sample preparation methods should remain the same as those listed in Appendix A.	
Crimidine ²	535-89-7

¹ LC-MS (electrospray) procedures are preferred for these analytes; however, if this technique is not available to the laboratory, GC-MS procedures using derivatization based on SW-846 Method 8270D may be used. Sample preparation methods should remain the same. Both electrospray LC-MS and GC/MS derivatization procedures are currently under development.

² This analyte needs to be determined using a wavelength of 230 nm.

This method provides reversed-phase high performance liquid chromatographic (RP-HPLC), thermospray (TSP) mass spectrometric (MS), and ultraviolet (UV) conditions for detection of the target analytes. Sample extracts can be analyzed by direct injection into the thermospray or onto a liquid chromatographic-thermospray interface. A gradient elution program is used to separate the compounds. Primary analysis may be performed by UV detection; however, positive results should be confirmed by TSP-MS. Quantitative analysis may be performed by either TSP-MS or UV detection, using either an external or internal standard approach. TSP-MS detection may be performed in either a negative ionization (discharge electrode) mode or a positive ionization mode, with a single quadrupole mass spectrometer. The use of MS-MS techniques is an option. The analytical range and detection limits vary depending on the target analyte and instrument used.

Source: EPA Method 8321B (SW-846): Solvent-Extractable Nonvolatile Compounds by High Performance Liquid Chromatography-Thermospray-Mass Spectrometry (HPLC-TSP-MS) or Ultraviolet (UV) Detection, Revision 2, 1998. <http://www.epa.gov/epaoswer/hazwaste/test/pdfs/8321b.pdf>

5.2.40 EPA Method 8330B (SW-846): Nitroaromatics and Nitramines by High Performance Liquid Chromatography (HPLC)

This method should be used for **preparation** and/or **analysis** of solid, non-aqueous liquid/organic solid, aqueous liquid, and drinking water samples for the contaminants identified below and listed in Appendix A. Note: Aqueous liquid and drinking water samples are prepared using Methods 3535A or 8330B prior to analysis.

Analyte(s)	CAS RN
2-Amino-4,6-dinitrotoluene (2-Am-DNT)	35572-78-2
4-Amino-2,6-dinitrotoluene (4-Am-DNT)	19406-51-0
4-Aminopyridine	504-24-5
3,5-Dinitroaniline (3,5-DNA)	618-87-1
1,3-Dinitrobenzene (1,3-DNB)	99-65-0
2,4-Dinitrotoluene (2,4-DNT)	121-14-2

Analyte(s)	CAS RN
2,6-Dinitrotoluene (2,6-DNT)	606-20-2
Hexahydro-1,3,5-trinitro-1,3,5-triazine (RDX)	121-82-4
Hexamethylenetriperoxidediamine (HMTD)	283-66-9
Methyl-2,4,6-trinitrophenylnitramine (Tetryl)	479-45-8
Nitrobenzene (NB)	98-95-3
Nitroglycerin (NG)	55-63-0
2-Nitrotoluene (2-NT)	88-72-2
3-Nitrotoluene (3-NT)	99-08-1
4-Nitrotoluene (4-NT)	99-99-0
Octahydro-1,3,5,7-tetranitro-1,3,5,7-tetrazocine (HMX)	2691-41-0
Pentaerythritol tetranitrate (PETN)	78-11-5
Triacetone triperoxide (TATP)	17088-37-8
1,3,5-Trinitrobenzene (1,3,5-TNB)	99-35-4
2,4,6-Trinitrotoluene (2,4,6-TNT)	118-96-7

This method is intended for the trace analysis of explosives and propellant residues by high performance liquid chromatography (HPLC) using a dual wavelength ultraviolet (UV) detector in a water, soil, or sediment matrix. All of the compounds listed in this method are either used in the manufacture of explosives or propellants, or they are the degradation products of compounds used for that purpose. Samples are prepared for analysis by HPLC-UV using the appropriate sample preparation technique (solid phase extraction by 3535A or solvent extraction by 8330B) and, if necessary, sample cleanup procedures. Method 8330 provides a salting-out extraction procedure for low concentration (ppt or ng/L) of explosives residues in surface or ground water. Direct injection of diluted and filtered water samples can be used for water samples of higher concentration. Soil and sediment samples are extracted using acetonitrile in an ultrasonic bath, filtered and chromatographed. The detection limits, ranges, and interferences depend on the target compound.

Source: EPA Method 8330B (SW-846): Nitroaromatics, Nitramines, and Nitrate Esters by High Performance Liquid Chromatography (HPLC), Revision 2, 2006.

<http://www.epa.gov/epaoswer/hazwaste/test/pdfs/8330b.pdf>

5.2.41 EPA ILM05.3 Cyanide: Analytical Methods for Total Cyanide Analysis

This method should be used for **preparation** and **analysis** of solid and aqueous liquid samples for the contaminants identified below and listed in Appendix A. Note: Hydrogen cyanide is not of concern in solid samples.

Analyte(s)	CAS RN
Cyanide, Total	57-12-5
Hydrogen cyanide	74-90-8

The method allows for either large volume (500-mL aqueous liquid samples or 1-g to 5-g solid samples mixed with 500 mL of reagent water) or medium volume (50-mL aqueous liquid samples or 1-g solid samples mixed with 50 mL of reagent water) sample preparation. Aqueous liquid samples are tested for sulfides and oxidizing agents prior to preparation. Sulfides are removed with cadmium carbonate or lead carbonate. Samples are treated with sulfuric acid and magnesium chloride and distilled into a sodium hydroxide solution. The solution is treated with color agents and the cyanide determined as an ion complex by visible spectrophotometry. The method quantitation limits are 10 µg/L or 2.5 mg/kg. Surfactants may interfere with the distillation procedure.

Source: EPA ILM05.3: Exhibit D – Part D: Analytical Methods for Total Cyanide Analysis.

<http://www.epa.gov/superfund/programs/clp/download/ilm/ilm53d.pdf>

5.2.42 IO [Inorganic] Compendium Method IO-3.1: Selection, Preparation, and Extraction of Filter Material

This method should be used for **preparation** of air samples for the contaminants identified below and listed in Appendix A. Refer to Appendix A for the appropriate determinative method for these analytes.

Analyte(s)	CAS RN
Ammonium metavanadate	7803-55-6
Arsenic, Total	740-38-2
Arsenic trioxide	1327-53-3
2-Chlorovinylarsonous acid (CVAA)	85090-33-1
Lewisite 1 (L-1) [2-chlorovinylchloroarsine] ¹	541-25-3
Lewisite 2 (L-2) [bis(2-chlorovinyl)chloroarsine]	40334-69-8
Lewisite 3 (L-3) [tris(2-chlorovinyl)arsine]	40334-70-1
Lewisite Oxide	1306-02-1
Osmium tetroxide	20816-12-0
Sodium arsenite	7784-46-5
Thallium sulfate	10031-59-1
Vanadium pentoxide	1314-62-1

¹ Laboratory testing is currently underway for speciation of lewisite 1 using GC-MS techniques.

This method supports determination of arsenic trioxide, lewisite compounds, lewisite oxide, CVAA, and sodium arsenite as total arsenic. Thallium sulfate is determined as total thallium and ammonium metavanadate and vanadium pentoxide are determined as total vanadium. A subsample (one-ninth of the overall filter) is obtained by cutting a strip from the filter used to collect the sample. The filter strip is extracted using a hydrochloric/nitric acid mix and microwave or hotplate heating. The extract is filtered, worked up to 20 mL, and analyzed using either Method IO-3.4 or Method IO-3.5.

Source: IO Compendium Method IO-3.1: Compendium of Methods for the Determination of Inorganic Compounds in Ambient Air: Selection, Preparation and Extraction of Filter Material, 1999.

<http://www.epa.gov/ttn/amtic/files/ambient/inorganic/mthd-3-1.pdf>

5.2.43 IO [Inorganic] Compendium Method IO-3.4: Determination of Metals in Ambient Particulate Matter Using Inductively Coupled Plasma (ICP) Spectroscopy

This method should be used for **analysis** of air samples for the contaminants identified below and listed in Appendix A.

Analyte(s)	CAS RN
Ammonium metavanadate	7803-55-6
Arsenic, Total	740-38-2
Arsenic trioxide	1327-53-3
2-Chlorovinylarsonous acid (CVAA)	85090-33-1
Lewisite 1 (L-1) [2-chlorovinylchloroarsine] ¹	541-25-3
Lewisite 2 (L-2) [bis(2-chlorovinyl)chloroarsine]	40334-69-8
Lewisite 3 (L-3) [tris(2-chlorovinyl)arsine]	40334-70-1
Lewisite Oxide	1306-02-1
Osmium tetroxide	20816-12-0
Sodium arsenite	7784-46-5
Thallium sulfate	10031-59-1
Vanadium pentoxide	1314-62-1

¹ Laboratory testing is currently underway for speciation of lewisite 1 using GC-MS techniques.

This method determines arsenic trioxide, lewisite, lewisite degradation products, and sodium arsenite as total arsenic. Osmium tetroxide is determined as total osmium, thallium sulfate is determined as total thallium, and ammonium metavanadate and vanadium pentoxide are determined as total vanadium. Ambient air is sampled by high-volume filters using Method IO-2.1 (a sampling method) and the filters are extracted by Method IO-3.1. The extracts are analyzed by Inductively Coupled Plasma - Atomic Emission Spectroscopy (ICP-AES) or Inductively Coupled Plasma - Mass Spectrometry (ICP-MS) (see Method IO-3.5 in Section 5.2.41). Detection limits, ranges, and interference corrections are dependent on the analyte and the instrument used.

Source: IO Compendium Method IO-3.4: Compendium of Methods for the Determination of Inorganic Compounds in Ambient Air: Determination of Metals in Ambient Particulate Matter Using Inductively Coupled Plasma (ICP) Spectroscopy, 1999.

<http://www.epa.gov/ttn/amtic/files/ambient/inorganic/mthd-3-4.pdf>

IO Compendium Method IO-2.1: Compendium of Methods for the Determination of Inorganic Compounds in Ambient Air: Sampling of Ambient Air for Total Suspended Particulate Matter (SPM) and PM₁₀ Using High Volume (HV) Sampler, 1999.

<http://www.epa.gov/ttn/amtic/files/ambient/inorganic/mthd-2-1.pdf>

5.2.44 IO [Inorganic] Compendium Method IO-3.5: Determination of Metals in Ambient Particulate Matter Using Inductively Coupled Plasma/Mass Spectrometry (ICP-MS)

This method should be used for **analysis** of air samples for the contaminants identified below and listed in Appendix A.

Analyte(s)	CAS RN
Ammonium metavanadate	7803-55-6
Arsenic, Total	740-38-2
Arsenic trioxide	1327-53-3
2-Chlorovinylarsonous acid (CVAA)	85090-33-1
Lewisite 1 (L-1) [2-chlorovinylchloroarsine] ¹	541-25-3
Lewisite 2 (L-2) [bis(2-chlorovinyl)chloroarsine]	40334-69-8
Lewisite 3 (L-3) [tris(2-chlorovinyl)arsine]	40334-70-1
Lewisite Oxide	1306-02-1
Sodium arsenite	7784-46-5
Thallium sulfate	10031-59-1
Vanadium pentoxide	1314-62-1

¹ Laboratory testing is currently underway for speciation of lewisite 1 using GC-MS techniques.

This method determines arsenic trioxide, lewisite, lewisite degradation products, and sodium arsenite as total arsenic. Thallium sulfate is determined as total thallium and ammonium metavanadate and vanadium pentoxide are determined as total vanadium. Ambient air is sampled by high-volume filters using Method IO-2.1 (a sampling method). The filters are extracted by Method IO-3.1 and the extracts analyzed by Inductively Coupled Plasma/Mass Spectrometry (ICP/MS) or Inductively Coupled Plasma - Atomic Emission Spectroscopy (ICP-AES) (see Method IO-3.4 in Section 5.2.40). Detection limits, ranges, and interference corrections are dependent on the analyte and the instrument used.

Source: IO Compendium Method IO-3.5: Compendium of Methods for the Determination of Inorganic Compounds in Ambient Air: Determination of Metals in Ambient Particulate Matter Using Inductively Coupled Plasma/Mass Spectrometry (ICP/MS), 1999.

<http://www.epa.gov/ttn/amtic/files/ambient/inorganic/mthd-3-5.pdf>

IO Compendium Method IO-2.1: Compendium of Methods for the Determination of Inorganic Compounds in Ambient Air: Sampling of Ambient Air for Total Suspended Particulate Matter (SPM) and PM₁₀ Using High Volume (HV) Sampler, 1999.

<http://www.epa.gov/ttn/amtic/files/ambient/inorganic/mthd-2-1.pdf>

5.2.45 IO [Inorganic] Compendium Method IO-5: Sampling and Analysis for Vapor and Particle Phase Mercury in Ambient Air Utilizing Cold Vapor Atomic Fluorescence Spectrometry (CVAFS)

This method should be used for **preparation** and **analysis** of air samples for the contaminants identified below and listed in Appendix A.

Analyte(s)	CAS RN
Mercury, Total	7439-97-6
Methoxyethylmercuric acetate	151-38-2

Vapor phase mercury is collected using gold-coated glass bead traps at a flow rate of 0.3 L/min. The traps are directly desorbed onto a second (analytical) trap. The mercury desorbed from the analytical trap is determined by Atomic Fluorescence Spectrometry. Particulate mercury is sampled on glass-fiber filters at a flow rate of 30 L/min. The filters are extracted with nitric acid and microwave heating. The extract is oxidized with BrCl, then reduced with stannous chloride and purged from solution onto a gold-coated glass bead trap. This trap is desorbed onto a second trap, the second trap is desorbed, and the mercury is determined by Atomic Fluorescence Spectrometry. The detection limits are 30 pg/m³ for particulate mercury and 45 pg/m³ for vapor phase mercury. Detection limits, analytical range, and interferences are dependent on the instrument used. There are no known positive interferences at 253.7 nm wavelength. Water vapor will cause a negative interference.

Source: IO Compendium Method IO-5: Compendium of Methods for the Determination of Inorganic Compounds in Ambient Air: Sampling and Analysis for Vapor and Particle Phase Mercury in Ambient Air Utilizing Cold Vapor Atomic Fluorescence Spectrometry (CVAFS).

<http://www.epa.gov/ttn/amtic/files/ambient/inorganic/mthd-5.pdf>

5.2.46 EPA Air Method, Toxic Organics - 10A (TO-10A): Determination of Pesticides and Polychlorinated Biphenyls in Ambient Air Using Low Volume Polyurethane Foam (PUF) Sampling Followed by Gas Chromatographic/Multi-Detector Detection (GC-MD)

This method should be used for **preparation** and **analysis** of air samples for the contaminants identified below and listed in Appendix A.

Analyte(s)	CAS RN
Chlorfenvinphos	470-90-6
Chlorpyrifos	2921-88-2
3-Chloro-1,2-propanediol ^{1,2}	96-24-2
Chlorosarin ¹	1445-76-7
Chlorosoman ¹	7040-57-5
Cyclohexyl sarin (GF)	329-99-7
Dichlorvos	62-73-7
Dicrotophos	141-66-2
Diisopropyl methylphosphonate (DIMP) ¹	1445-75-6
Dimethylphosphite	868-85-9
Dimethylphosphoramidic acid ²	33876-51-6

Analyte(s)	CAS RN
EA2192 [Diisopropylaminoethyl methylthiophosphonate] ²	73207-98-4
Ethyl methylphosphonic acid (EMPA) ²	1832-53-7
N-Ethyldiethanolamine (EDEA)	139-87-7
Fenamiphos	22224-92-6
Fentanyl	437-38-7
Isopropyl methylphosphonic acid (IMPA) ²	1832-54-8
Methamidophos ³	10265-92-6
Methyl parathion	298-00-0
N-Methyldiethanolamine (MDEA)	105-59-9
1-Methylethyl ester ethylphosphonofluoridic acid (GE)	1189-87-3
Methylphosphonic acid (MPA) ²	993-13-5
Mevinphos	7786-34-7
Mustard, nitrogen (HN-1) [bis(2-chloroethyl)ethylamine]	538-07-8
Mustard, nitrogen (HN-2) [2,2'-dichloro-N-methyldiethylamine N,N-bis(2-chloroethyl)methylamine]	51-75-2
Mustard, nitrogen (HN-3) [tris(2-chloroethyl)amine]	555-77-1
Mustard, sulfur / Mustard gas (HD)	505-60-2
Parathion	56-38-2
Phencyclidine	77-10-1
Phenol	108-95-2
Phorate	298-02-2
Phosphamidon	13171-21-6
Pinacolyl methyl phosphonic acid (PMPA) ²	616-52-4
R-33 (VR) [methylphosphonothioic acid, S-[2-(diethylamino)ethyl] O-2-methylpropyl ester]	159939-87-4
Sarin (GB) ¹	107-44-8
Semivolatile Organic Compounds, NOS ⁴	NA
Soman (GD) ¹	96-64-0
Tabun (GA)	77-81-6
Tear gas (CS) [chlorobenzylidene malonitrile]	2698-41-1
Tetraethyl pyrophosphate	107-49-3

Analyte(s)	CAS RN
Tetramethylenedisulfotetramine	80-12-6
Thiodiglycol (TDG)	111-48-8
Triethanolamine (TEA)	102-71-6
Trimethyl phosphite	121-45-9
VE [phosphonothioic acid, ethyl-, S-(2-(diethylamino)ethyl) O-ethyl ester]	21738-25-0
VG [phosphonothioic acid, S-(2-(diethylamino)ethyl) O,O-diethyl ester]	78-53-5
VM [phosphonothioic acid, methyl-, S-(2-(diethylamino)ethyl) O-ethyl ester]	21770-86-5
VX [O-ethyl-S-(2-diisopropylaminoethyl)methyl-phosphonothiolate]	50782-69-9
The following analyte should be determined by this method only if problems (e.g., insufficient recovery, interferences) occur when using Method TO-15.	
Allyl alcohol	107-18-6

¹ If problems occur when using this method, it is recommended that the canister Method TO-15 be used.

² For this analyte, Method TO-10A must be modified to include a derivatization step.

³ If problems occur when using this method, it is recommended that NIOSH Method 5600 be used.

⁴ NOS = Not otherwise specified

A low-volume (1 to 5 L/minute) sample collection rate is used to collect vapors on a sorbent cartridge containing polyurethane foam (PUF) or PUF in combination with another solid sorbent. Airborne particles also are collected, but the sampling efficiency is not known. Pesticides and other chemicals are extracted from the sorbent cartridge with 5% diethyl ether in hexane and determined by gas chromatography coupled with an electron capture detector (ECD). A nitrogen-phosphorous detector (NPD), flame photometric detector (FPD), Hall electrolytic conductivity detector (HECD), or mass spectrometer (MS) also may be used. Dimethylphosphoramidic acid, EA2192, EMPA, IMPA, MPA, and PMPA require derivatization with a trimethylsilyl agent prior to injection into the GC. This method is applicable to multicomponent atmospheres, 0.001 to 50 µg/m³ concentrations, and 4 to 24-hour sampling periods. The limit of detection will depend on the specific compounds measured, the concentration level, and the degree of specificity required. If analytical difficulties are noted, the canister procedures described in Method TO-15 may be appropriate.

Source: EPA Air Method, Toxic Organics-10A (TO-10A): Compendium of Methods for the Determination of Inorganic Compounds in Ambient Air: Determination of Pesticides and Polychlorinated Biphenyls in Ambient Air Using Low Volume Polyurethane Foam (PUF) Sampling Followed by Gas Chromatographic/Multi-Detector Detection (GC/MD), 1999.

<http://www.epa.gov/ttnamtl1/files/ambient/airtox/to-10ar.pdf>

5.2.47 EPA Air Method, Toxic Organics - 15 (TO-15): Determination of Volatile Organic Compounds (VOCs) in Air Collected in Specially-Prepared Canisters and Analyzed by Gas Chromatography/Mass Spectrometry (GC-MS)

This method should be used for **preparation** and **analysis** of air samples for the contaminants identified below and listed in Appendix A.

Analyte(s)	CAS RN
Allyl alcohol ¹	107-18-6
Carbon disulfide	75-15-0
Cyanogen chloride	506-77-4
1,2-Dichloroethane	107-06-2
Ethylchloroarsine (ED)	598-14-1
Ethylene oxide	75-21-8
The following analytes should be determined by this method only if problems (e.g., insufficient recovery, interferences) occur when using Method TO-10A.	
3-Chloro-1,2-propanediol ²	96-24-2
Chlorosarin	1445-76-7
Chlorosoman	7040-57-5
Diisopropyl methylphosphonate (DIMP)	1445-75-6
1-Methylethyl ester ethylphosphonofluoridic acid (GE)	1189-87-3
Sarin (GB)	107-44-8
Soman (GD)	96-64-0

¹ If problems occur when using this method, it is recommended that Method TO-10A be used.

² For this analyte, Method TO-15 must be modified to include a derivatization step.

The atmosphere is sampled by introduction of air into a specially prepared stainless steel canister (specially electropolished or silica-coated). A sample of air is drawn through a sampling train comprising components that regulate the rate and duration of sampling into the pre-evacuated and passivated canister. Grab samples also may be collected. After the air sample is collected, the canister valve is closed, an identification tag is attached to the canister, and the canister is transported to the laboratory for analysis. To analyze the sample, a known volume of sample is directed from the canister through a solid multisorbent concentrator. Recovery of less volatile compounds may require heating the canister. After the concentration and drying steps are completed, VOCs are thermally desorbed, entrained in a carrier gas stream, and then focused in a small volume by trapping on a cryo-focusing (ultra-low temperature) trap or small volume multisorbent trap. The sample is then released by thermal desorption and analyzed by Gas Chromatography/Mass Spectrometry (GC-MS). This method applies to ambient concentrations of VOCs above 0.5 ppbv and typically requires Volatile Organic Compounds (VOC) enrichment by concentrating up to 1 L of a sample volume; however, when using current technologies, quantifications of approximately 100 pptv have been achieved with 0.5-L sample volumes.

Source: EPA Air Method, Toxic Organics-15 (TO-15): Compendium of Methods for the Determination of Inorganic Compounds in Ambient Air: Determination of Volatile Organic Compounds (VOCs) in Air Collected in Specially-Prepared Canisters and Analyzed by Gas Chromatography/Mass Spectrometry (GC/MS), 1999.

<http://www.epa.gov/ttn/amtic/files/ambient/airtox/to-15r.pdf>

5.2.48 NIOSH Method 1612: Propylene Oxide

This method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Propylene oxide	75-56-9

A sample tube containing coconut shell charcoal is used for sample collection with a flow rate of 0.01 to 0.2 L/min. One milliliter of carbon disulfide (CS₂) is added to the vial and allowed to sit for 30 minutes prior to analysis with occasional agitation. No interferences have been found. The working range is between 8 to 295 ppm for air samples of 5 L.

Source: NIOSH Method 1612: Propylene Oxide, Issue 2, 1994.

<http://www.cdc.gov/niosh/nmam/pdfs/1612.pdf>

5.2.49 NIOSH Method 2016: Formaldehyde

This method should be used for **preparation** and **analysis** of samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Formaldehyde	50-00-0

This method can be used for the determination of formaldehyde using high performance liquid chromatography (HPLC) with an ultraviolet (UV) detector. Air is sampled onto a cartridge containing silica gel coated with 2,4-dinitrophenyl hydrazine, at a rate of 0.03 to 1.5 L/min. The cartridge is extracted with 10 mL of acetonitrile and analyzed by HPLC-UV at a wavelength of 360 nm. The working range is 0.015 to 2.5 mg/m³ (0.012 to 2.0 ppm) for a 15-L sample. The detection limit for formaldehyde is 0.07 µg/sample. Ozone has been observed to consume the 2,4-dinitrophenylhydrazine (2,4-DNPH) reagent and to degrade the formaldehyde derivative. Ketones and other aldehydes can react with 2,4-DNPH; the derivatives produced, however, are separated chromatographically from the formaldehyde derivative.

Source: NIOSH Method 2016: Formaldehyde, Issue 2, 2003.

<http://www.cdc.gov/niosh/nmam/pdfs/2016.pdf>

5.2.50 NIOSH Method 2513: Ethylene Chlorohydrin

This method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
2-Chloroethanol (ethylene chlorohydrin)	107-07-3

Samples are drawn into a tube containing petroleum charcoal at a rate of 0.01 to 0.2 L/min and transferred into vials containing eluent (carbon disulfide, 2-propanol, and *n*-pentadiene as an internal standard). Vials must sit for 30 minutes prior to analysis by gas chromatography – flame ionization detection (GC-FID). No interferences have been identified. Humidity may decrease the breakthrough volume during sample collection. The working range of the method is 0.5 to 15 ppm for a 20-L air sample.

Source: NIOSH Method 2513: Ethylene Chlorohydrin, Issue 2, 1994.

<http://www.cdc.gov/niosh/nmam/pdfs/2513.pdf>

5.2.51 NIOSH Method 3510: Monomethylhydrazine

This method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Methyl hydrazine (monomethylhydrazine)	60-34-4

Samples are collected into a bubbler containing HCl using a flow rate of 0.5 to 1.5 L/min. Samples are then mixed with phosphomolybdic acid solution and transferred to a large test tube for spectrophotometric analysis. Positive interferences that have been noted include stannous ion, ferrous ion, zinc, sulfur dioxide, and hydrogen sulfide. Negative interferences may occur by oxidation of the monomethylhydrazine by halogens, oxygen (especially in the presence of copper (I) ions) and hydrogen dioxide. The working range of the method is 0.027 to 2.7 ppm for a 20-L air sample.

Source: NIOSH Method 3510: Monomethylhydrazine, Issue 1, 1994.

<http://www.cdc.gov/niosh/nmam/pdfs/3510.pdf>

5.2.52 NIOSH Method 5600: Organophosphorus Pesticides

This method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Organophosphate Pesticides, NOS ¹	NA
The following analyte should be determined by this method only if problems (e.g., insufficient recovery, interferences) occur when using Method TO-10A.	
Methamidophos	10265-92-6

¹ NOS = Not otherwise specified

This method is used for the detection of organophosphorus pesticides using gas chromatography (GC) with a flame photometric detector (FPD). The method also may be applicable to the determination of other organophosphorus compounds after evaluation for desorption efficiency, sample capacity, sample stability, and precision and accuracy. The working range for each analyte is provided in Table 5 of the method. These ranges cover from 0.1 to 2 times the OSHA Permissible Exposure Limits (PELs) (see Table 5 of the method). The method also is applicable to Short Term Exposure Limit (STEL) measurements using 12-L samples. The detection limit depends on the compound being analyzed. Several organophosphates may co-elute with either target analytes or internal standards causing integration errors. These include other pesticides, and the following: tributyl phosphate, tris-(2-butoxy ethyl) phosphate, tricresyl phosphate, and triphenyl phosphate.

Source: NIOSH Method 5600: Organophosphorus Pesticides, Issue 1, 1994.

<http://www.cdc.gov/niosh/nmam/pdfs/5600.pdf>

5.2.53 NIOSH Method 5601: Organonitrogen Pesticides

This method should be used for **preparation** and **analysis** of air samples for the contaminants identified below and listed in Appendix A.

Analyte(s)	CAS RN
Aldicarb	116-06-3
Carbofuran	1563-66-2
Methomyl	16752-77-5
Oxamyl	23135-22-0

This method can be used for the determination of organonitrogen pesticides using high performance liquid chromatography (HPLC) with an ultraviolet (UV) detector. The method also may be applicable to the determination of other organonitrogen compounds and to a broad range of pesticides having UV chromophores, e.g., acetanilides, acid herbicides, organophosphates, phenols, pyrethroids, sulfonyl ureas, sulfonamides, triazines, and uracil pesticides. The detection limit for aldicarb is 1.2 µg per sample and 0.6 µg per sample for carbofuran, methomyl, and oxamyl. The working ranges for aldicarb, carbofuran, and oxamyl range from 0.5 to 10 times the OSHA Permissible Exposure Limit (PEL) (see Table 2 of the method). Because of the broad response of the UV detector at shorter wavelengths, there are many potential interferences. Those tested include solvents (chloroform and toluene), antioxidants (BHT), plasticizers (dialkyl phthalates), nitrogen compounds (nicotine and caffeine), impurities in HPLC reagents (e.g., in triethylamine), other pesticides (2,4-D, atrazine, parathion, etc.), and pesticide hydrolysis products (1-naphthol). Confirmation techniques are recommended when analyte identity is uncertain.

Source: NIOSH Method 5601: Organonitrogen Pesticides, Issue 1, 1998.

<http://www.cdc.gov/niosh/nmam/pdfs/5601.pdf>

5.2.54 NIOSH Method 6001: Arsine

This method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Arsine	7784-42-1

Arsine is determined as arsenic. A 0.1 to 10-L volume of air is drawn through a sorbent tube containing activated charcoal. The sorbent is extracted with a nitric acid solution, and arsenic is determined by graphite furnace atomic absorption. The working range of the method is 0.001 to 0.2 mg/m³ for a 10-L sample. The method is subject to interferences from other arsenic compounds.

Source: NIOSH Method 6001: Arsine, Issue 2, 1994. <http://www.cdc.gov/niosh/nmam/pdfs/6001.pdf>

5.2.55 NIOSH Method 6002: Phosphine

This method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Phosphine	7803-51-2

In this method, phosphine is determined as phosphate. A volume of 1 to 16 L of air are drawn through a sorbent tube containing silica gel coated with $\text{Hg}(\text{CN})_2$. The sorbent is extracted with a potassium permanganate/ sulfuric acid solution and washed with reagent water. Following treatment with the color agent and extraction into organic solvent, phosphate is determined by visible spectrometry. The working range of the method is 0.02 to 0.9 mg/m^3 for a 16-L sample. The method is subject to interferences from phosphorus trichloride, phosphorus pentachloride, and organic phosphorus compounds.

Source: NIOSH Method 6002: Phosphine, Issue 2, 1994. <http://www.cdc.gov/niosh/nmam/pdfs/6002.pdf>

5.2.56 NIOSH Method 6004: Sulfur Dioxide

This method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Sulfur dioxide	7446-09-5

In this method, sulfur dioxide is determined as sulfite plus sulfate. A volume of 40 to 200 L of air is drawn through a sodium carbonate-treated filter that is preceded by a 0.8 μm filter to remove particulates and sulfuric acid. The treated filter is extracted with a carbonate/bicarbonate solution and the extract analyzed by ion chromatography for sulfite and sulfate. The sulfur dioxide is present as sulfite on the filter; however, because sulfite oxidizes to sulfate, both ions must be determined and the results summed. The working range of the method is 0.5 to 20 mg/m^3 for a 100-L sample. The method is subject to interference from sulfur trioxide in dry conditions.

Source: NIOSH Method 6004: Sulfur Dioxide, Issue 2, 1994.
<http://www.cdc.gov/niosh/nmam/pdfs/6004.pdf>

5.2.57 NIOSH Method 6010: Hydrogen Cyanide

This method should be used for **preparation** and **analysis** of air samples for the contaminants identified below and listed in Appendix A.

Analyte(s)	CAS RN
Cyanide, Total	57-12-5
Hydrogen cyanide	74-90-8

Hydrogen cyanide is determined as a cyanide ion complex by this method. A volume of 2 to 90 L of air is drawn through a soda lime sorbent tube. A glass-fiber filter is used to remove particulate cyanides prior to the sorbent tube. Cyanide is extracted from the sorbent with reagent water treated with NaOH. The extract is pH adjusted with HCl, oxidized with N-chlorosuccinimide/succinimide, and treated with the coupling-color agent (barbituric acid/pyridine). The cyanide ion is determined by visible spectrophotometry. The working range of the method is 3 to 260 mg/m^3 for a 3-L sample. The method is subject to interference from high concentrations of hydrogen sulfide.

Source: NIOSH Method 6010: Hydrogen Cyanide, Issue 2, 1994.

<http://www.cdc.gov/niosh/nmam/pdfs/6010.pdf>

5.2.58 NIOSH Method 6013: Hydrogen Sulfide

This method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Hydrogen sulfide	7783-06-4

Hydrogen sulfide is determined as sulfate by this method. A volume of 1.2 to 40 L of air is drawn through charcoal sorbent. A prefilter is used to remove particulates. The sorbent portions are extracted with an ammonium hydroxide/hydrogen peroxide solution and the extract is analyzed for sulfate by ion chromatography. The working range of the method is 0.9 to 20 mg/m³ for a 20-L sample. The method is subject to interference from sulfur dioxide.

Source: NIOSH Method 6013: Hydrogen Sulfide, Issue 1, 1994.

<http://www.cdc.gov/niosh/nmam/pdfs/6013.pdf>

5.2.59 NIOSH Method 6015: Ammonia

This method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Ammonia	7664-41-7

Ammonia is determined as indophenol blue by this method. A volume of 0.1 to 96 L of air is drawn through a sulfuric acid-treated silica gel sorbent. A prefilter is used to remove particulates. The sorbent is extracted with reagent water, the pH adjusted, and reagents are added to generate the indophenol blue compound in the presence of ammonium. The extract is analyzed by visible spectrophotometry. The working range of the method is 0.15 to 300 mg/m³ for a 10-L sample. Twice the recommended sample volume should be collected in order to achieve an action level of 70 µg/m³. No interferences have been identified.

Source: NIOSH Method 6015: Ammonia, Issue 2, 1994. <http://www.cdc.gov/niosh/nmam/pdfs/6015.pdf>

5.2.60 NIOSH Method 6402: Phosphorus Trichloride

This method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Phosphorus trichloride	7719-12-2

In this method, phosphorus trichloride is determined as phosphate. A volume of 11 to 100 L of air is drawn through a bubbler containing reagent water. The resulting H₃PO₃ solution is oxidized to H₃PO₄ and color agents are added. The solution is analyzed by visible spectrophotometry. The working range of the method is 1.2 to 80 mg/m³ for a 25-L sample. Phosphorus (V) compounds do not interfere. The sample solutions are stable to oxidation by air during sampling.

Source: NIOSH Method 6402: Phosphorus Trichloride, Issue 2, 1994.

<http://www.cdc.gov/niosh/nmam/pdfs/6402.pdf>

5.2.61 NIOSH Method 7903: Acids, Inorganic

This method should be used for **preparation** and **analysis** of air samples for the contaminants identified below and listed in Appendix A.

Analyte(s)	CAS RN
Hydrogen bromide	10035-10-6
Hydrogen chloride	7647-01-0
Hydrogen fluoride ¹	7664-39-3

¹ If problems occur when using this method, it is recommended that NIOSH 7906 be used.

Acids are analyzed as bromide, chloride, and fluoride, respectively, by this method. A volume of 3 to 100 L of air is drawn through a silica gel sorbent. The sorbent portions are extracted with a buffered carbonate/bicarbonate solution and the extract is analyzed by ion chromatography. The working range of this method is 0.01 to 5 mg/m³ for a 50-L sample. Particulate salts of the acids are an interference (trapped on the glass wool filter plug in the sorbent tube). Chlorine and bromine are also interferences. Acetate, formate, and propionate interferences may be reduced by use of a weaker eluent.

Source: NIOSH Method 7903: Acids, Inorganic, Issue 2, 1994.

<http://www.cdc.gov/niosh/nmam/pdfs/7903.pdf>

5.2.62 NIOSH Method 7905: Phosphorus

This method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
White phosphorus	12185-10-3

This method identifies and determines the concentration of white phosphorus in air by using a gas chromatography/flame photometric detector (GC-FPD). Samples prepared by this method can be analyzed alternatively by gas chromatography/nitrogen phosphorus detection (GC-NPD). Five to 100 L of air are drawn through a GC solid sorbent tube, and the sorbent is extracted (desorbed) with xylene. Phosphorus is determined by GC-FPD or GC-NPD. The working range for samples analyzed by GC-FPD is 0.056 to 0.24 mg/m³ for a 12-L sample. The limit of detection (LOD) for samples analyzed by GC-FPD is 0.005 µg per sample. The method is applicable to vapor-phase phosphorus only; if particulate phosphorus is expected, a filter could be used in the sampling train.

Source: NIOSH Method 7905: Phosphorus, Issue 2, 1994.

<http://www.cdc.gov/niosh/nmam/pdfs/7905.pdf>

5.2.63 NIOSH Method 7906: Fluorides, Aerosol and Gas

This method is not currently listed for any of the analyte/sample type combinations included in Appendix A. If problems occur when using NIOSH Method 7903 for the analysis of hydrogen fluoride, then this method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A. (See Footnote 9 of Appendix A.)

Analyte(s)	CAS RN
Hydrogen fluoride ¹	7664-39-3

¹ This analyte should be determined by this method (7906) **only** if problems occur when using NIOSH Method 7903.

Hydrogen fluoride is determined as fluoride ion by this method. A volume of 1 to 800 L of air is drawn through a 0.8- μ m cellulose ester membrane (to trap particulate fluorides) and a cellulose pad treated with sodium carbonate (to trap gaseous fluoride). The pad is extracted with reagent water and the extract is analyzed for fluoride by ion chromatography. The working range of the method is 0.04 to 8 mg/m³ for 250-L samples. If other aerosols are present, gaseous fluoride may be slightly underestimated due to adsorption onto or reaction with particles, with concurrent overestimation of particulate/gaseous fluoride ratio.

Source: NIOSH Method 7906: Fluorides, Aerosol and Gas by IC, Issue 1, 1994.

<http://www.cdc.gov/niosh/nmam/pdfs/7906.pdf>

5.2.64 NIOSH Method S301-1: Fluoroacetate Anion

This method should be used for **preparation** of air samples for the contaminant identified below and listed in Appendix A. Note: EPA Method 300.1 Rev 1.0 should be used as the determinative method (refer to Appendix A).

Analyte(s)	CAS RN
Fluoroacetic acid and fluoroacetate salts	NA
Methyl fluoroacetate	453-18-9

This method was developed specifically for sodium fluoroacetate, but also may be applicable to other fluoroacetate salts. The method determines fluoroacetate salts as fluoroacetate anion. A known volume of air (e.g., 480 L was used in validation of this method) is drawn through a cellulose ester membrane filter to collect sodium fluoroacetate. Sodium fluoroacetate is extracted from the filter with 5 mL of deionized water, and the resulting sample is analyzed by ion chromatography using electrolytic conductivity detection. The analytical range of this method is estimated to be 0.01 to 0.16 mg/m³. The detection limit is estimated to be 20 ng of sodium fluoroacetate per injection, corresponding to a 100- μ L aliquot of a 0.2- μ g/mL standard.

Source: NIOSH Method S301-1: Sodium Fluoroacetate, 1977. <http://www.cdc.gov/niosh/pdfs/s301.pdf>

5.2.65 OSHA Method 40: Methylamine

This method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Methylamine	74-89-5

This method is used for detection of methylamine using high performance liquid chromatography (HPLC) with a fluorescence (FL) or visible (vis) detector. Samples are collected by drawing 10-L volumes of air

at a rate of 0.2 L/min through standard size sampling tubes containing XAD-7 resin coated with 10% NBD chloride by weight. Samples are desorbed with 5% (w/v) 7-chloro-4-nitrobenzo-2-oxa-1,3-diazole (NBD chloride) in tetrahydrofuran (with a small amount of sodium bicarbonate present), heated in a hot water bath, and analyzed by HPLC-FL or HPLC-vis. The detection limit of the overall procedure is 0.35 µg per sample (28 ppb or 35 µg/m³). Quantitation limits of 28 ppb (35 µg/m³) have been reliably achieved. This is the smallest amount of methylamine that can be quantified within the requirements of a recovery of at least 75% and a precision (standard deviation of 1.96) of ±25% or better.

Source: OSHA Method 40: Methylamine, 1982.

<http://www.osha.gov/dts/sltc/methods/organic/org040/org040.html>

5.2.66 OSHA Method 54: Methyl Isocyanate

This method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Methyl isocyanate	624-83-9

This method determines the concentration of methyl isocyanate in air by using high performance liquid chromatography (HPLC) with a fluorescence or ultraviolet (UV) detector. Samples are collected by drawing a known volume of air through XAD-7 tubes coated with 0.3 mg of 1-(2-pyridyl)piperazine (1-2PP). Samples are desorbed with acetonitrile (ACN) and analyzed by HPLC using a fluorescence or UV detector.

Source: OSHA Method 54: Methyl Isocyanate (MIC), 1985.

<http://www.osha.gov/dts/sltc/methods/organic/org054/org054.html>

5.2.67 OSHA Method 61: Phosgene

This method should be used for **preparation** and **analysis** of air samples for the contaminants identified below and listed in Appendix A.

Analyte(s)	CAS RN
Perfluoroisobutylene (PFIB) ¹	382-21-8
Phosgene	75-44-5

¹ If problems occur when using this method, it is recommended that a method based on the following journal article be used: J. Chrom. A, 1098: (2005) 156–165.

This method determines the concentration of phosgene in air by using gas chromatography with a nitrogen selective detector. Air samples are collected by drawing known volumes of air through sampling tubes containing XAD-2 adsorbent that has been coated with 2-(hydroxymethyl)piperidine. The samples are desorbed with toluene and then analyzed by gas chromatography using a nitrogen selective detector.

Source: OSHA Method 61: Phosgene, 1986.

<http://www.osha.gov/dts/sltc/methods/organic/org061/org061.html>

5.2.68 OSHA Method ID-211: Sodium Azide and Hydrazoic Acid in Workplace Atmospheres

This method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Sodium azide	26628-22-8

This method describes sample collection and analysis of airborne azides [as sodium azide (NaN_3) and hydrazoic acid (HN_3)]. Particulate sodium azide (NaN_3) is collected on a PVC filter or in the glass wool plug of the sampling tube. Gaseous HN_3 is collected and converted to NaN_3 by the impregnated silica gel (ISG) sorbent within the sampling tube. The collected azide on either media is desorbed in a weak buffer solution, and the resultant anion (N_3^-) is analyzed by Ion Chromatography (IC) using a variable wavelength UV detector at 210 nm. A gravimetric conversion is used to calculate the amount of NaN_3 or HN_3 collected. The detection limit was found to be 0.001 ppm as HN_3 or 0.003 mg/m^3 as NaN_3 for a 5-L air sample. The quantitation limit was found to be 0.004 ppm as HN_3 or 0.011 mg/m^3 as NaN_3 for a 5-L air sample.

Source: OSHA Method ID-211: Sodium Azide and Hydrazoic Acid in Workplace Atmospheres.
<http://www.osha.gov/dts/sltc/methods/inorganic/id211/id211.html>

5.2.69 OSHA Method ID-216SG: Boron Trifluoride (BF_3)

This method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Boron trifluoride	7637-07-2

Boron trifluoride is determined as fluoroborate by this method. A volume of 30 to 480 L of air is drawn through a bubbler containing 0.1 M ammonium fluoride. The solution is diluted and analyzed with a fluoroborate ion specific electrode (ISE). The detection limit is 10 μg in a 30-L sample.

Source: OSHA Method ID-216SG: Boron Trifluoride (BF_3), 1989.
<http://www.osha.gov/dts/sltc/methods/partial/id216sg/id216sg.html>

5.2.70 OSHA Method PV2004: Acrylamide

This method should be used for **preparation** and **analysis** of air samples for the contaminants identified below and listed in Appendix A.

Analyte(s)	CAS RN
Acrylamide	79-06-1
Acrylonitrile	107-13-1
Methylacrylonitrile	126-98-7

This method determines the concentration of acrylamide in air by using high performance liquid chromatography (HPLC) with an ultraviolet (UV) detector. Samples are collected by drawing known

volumes of air through OSHA versatile sampler (OVS-7) tubes, each containing a glass fiber filter and two sections of XAD-7 adsorbent. Samples are desorbed with a solution of 5% methanol/95% water and analyzed by HPLC using a UV detector. The detection limit was found to be 0.7 µg/mL (0.006 mg/m³ for a 1-mL desorption volume or 0.029 mg/m³ for a 5-mL desorption volume based on a 120-L air volume). Applicable working ranges for a 1-mL and 5-mL desorption volume are 0.017 – 1.5 mg/m³ and 0.083 – 7.5 mg/m³, respectively.

Source: OSHA Method PV2004: Acrylamide.

<http://www.osha.gov/dts/sltc/methods/partial/pv2004/2004.html>

5.2.71 OSHA Method PV2103: Chloropicrin

This method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Chloropicrin	79-06-2

This method determines the concentration of chloropicrin in air by gas chromatography/electron capture detection (GC-ECD). Samples are collected by drawing a known volume of air through two XAD-4 tubes in series. Samples are desorbed with ethyl acetate and analyzed by GC-ECD. Chloropicrin is light sensitive, and samples should be protected from light. The detection limit of the analytical procedure is 0.01 ng, with a 1-µL injection volume. This is the smallest amount that could be detected under normal operating conditions. The working range is 0.25 to 0.99 ppmv.

Source: OSHA Method PV2103: Chloropicrin. <http://www.osha.gov/dts/sltc/methods/partial/t-pv2103-01-9111-ch/t-pv2103-01-9111-ch.html>

5.2.72 ASTM Method D5755-03: Standard Test Method for Microvacuum Sampling and Indirect Analysis of Dust by Transmission Electron Microscopy (TEM) for Asbestos Structure Number Surface Loading

This method should be used for **preparation** and **analysis** of solid samples (e.g., soft surfaces-microvac) for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Asbestos	1332-21-4

This method describes procedures to identify asbestos in dust and provide an estimate of the surface loading of asbestos reported as the number of asbestos structures per unit area of sampled surface. The sample is collected by vacuuming a known surface area with a standard 25- or 37-mm air sampling cassette using a plastic tube that is attached to the inlet orifice, which acts as a nozzle. The sample is transferred from inside the cassette to an aqueous suspension of known volume. Aliquots of the suspension are then filtered through a membrane, and a section of the membrane is prepared and transferred to a transmission electron microscopy (TEM) grid using a direct transfer method. The asbestiform structures are identified, sized, and counted by TEM, using select area electron diffraction (SAED) and energy dispersive X-ray analysis (EDXA) at a magnification of 15,000 to 20,000X.

Source: ASTM Method D5755-03: Standard Test Method for Microvacuum Sampling and Indirect Analysis of Dust by Transmission Electron Microscopy for Asbestos Structure Number Surface Loading, 2003. http://www.astm.org/cgi-bin/SoftCart.exe/STORE/filtrexx40.cgi?U+mystore+tavs3076+-L+D5755:03+/usr6/htdocs/astm.org/DATABASE.CART/REDLINE_PAGES/D5755.htm

5.2.73 ASTM Method D6480-99: Standard Test Method for Wipe Sampling of Surfaces, Indirect Preparation, and Analysis for Asbestos Structure Number Concentration by Transmission Electron Microscopy

This method should be used for **preparation** and **analysis** of solid samples (e.g., hard surfaces-wipes) for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Asbestos	1332-21-4

This method describes a procedure to identify asbestos in samples wiped from surfaces and to provide an estimate of the concentration of asbestos reported as the number of asbestos structures per unit area of sampled surface. A sample is collected by wiping a surface of known area with a wipe material. The sample is transferred from the wipe material to an aqueous suspension of known volume. Aliquots of the suspension are then filtered through a membrane filter, and a section of the membrane filter is prepared and transferred to a TEM grid, using the direct transfer method. The asbestiform structures are identified, sized, and counted by TEM, using electron diffraction (ED) and energy dispersive X-ray analysis (EDXA) at a magnification from 15,000 to 20,000X.

Source: ASTM Method D6480-99: Standard Test Method for Wipe Sampling of Surfaces, Indirect Preparation, and Analysis for Asbestos Structure Number Concentration by Transmission Electron Microscopy. <http://www.astm.org/cgi-bin/SoftCart.exe/STORE/filtrexx40.cgi?U+mystore+tavs3076+-L+D6480:99+/usr6/htdocs/astm.org/DATABASE.CART/HISTORICAL/D6480-99.htm>

5.2.74 ISO Method 10312:1995: Ambient Air - Determination of Asbestos Fibres - Direct-transfer Transmission Electron Microscopy Method (TEM)

This method should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Asbestos	1332-21-4

This method determines the type(s) of asbestos fibers present, but cannot discriminate between individual fibers of the asbestos and non-asbestos analogues of the same amphibole mineral. The method is defined for polycarbonate capillan/pore filters or cellulose ester (either mixed esters of cellulose or cellulose nitrate) filters through which a known volume of air has been drawn. The method is suitable for determination of asbestos in both exterior and building atmospheres. The range of concentrations that can be determined is 50 structures/mm² to 7,000 structures/mm² on the filter. In a 4000-L air sample with approximately 10 pg/m³ (typical of clean or rural atmospheres), an analytical sensitivity of 0.5 structure/L can be obtained. This is equivalent to a detection limit of 1.8 structure/L when an area of 0.195 mm of the TEM specimen is examined.

Source: ISO Method 10312: 1995: Ambient Air—Determination of Asbestos Fibres—Direct Transfer Transmission Electron Microscopy Method, 2005. <http://www.iso.org/iso/en/CatalogueDetailPage.CatalogueDetail?CSNUMBER=18358&ICS1=13&ICS2=40&ICS3=20>

5.2.75 Standard Method 4500-NH₃ B: Nitrogen (Ammonia) Preliminary Distillation Step

This method should be used for **preparation** of aqueous/liquid samples for the contaminant identified below and listed in Appendix A. Note: Standard Method 4500-NH₃ G should be used as the determinative method (refer to Appendix A).

Analyte(s)	CAS RN
Ammonia	7664-41-7

A 0.5- to 1-L sample is dechlorinated, buffered, adjusted to pH 9.5, and distilled into a sulfuric acid solution. The distillate is brought up to volume, neutralized with sodium hydroxide, and analyzed by Method 4500-NH₃ G.

Source: American Public Health Association, American Water Works Association, and Water Environment Federation. 2005. *Standard Methods for the Examination of Water and Wastewater*. 21st Edition. <http://www.standardmethods.org/>

5.2.76 Standard Method 4500-NH₃ G: Nitrogen (Ammonia) Automated Phenate Method

This method should be used for **analysis** of aqueous liquid samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Ammonia	7664-41-7

Ammonia is determined as indophenol blue by this method. A portion of the neutralized sample distillate (from procedure 4500-NH₃ B) is run through a manifold. The ammonium in the distillate reacts with solutions of disodium ethylenediaminetetraacetic acid (EDTA), sodium phenate, sodium hypochlorite, and sodium nitroprusside. The resulting indophenol blue is detected by colorimetry in a flow cell. The range of the method is 0.02 to 2.0 mg/L.

Source: American Public Health Association, American Water Works Association, and Water Environment Federation. 2005. *Standard Methods for the Examination of Water and Wastewater*. 21st Edition. <http://www.standardmethods.org/>

5.2.77 Standard Method 4500-Cl G: DPD Colorimetric Method

This method should be used for **preparation** and **analysis** of aqueous liquid and drinking water samples for the contaminant identified below and listed in Appendix A. It also should be used for analysis of air samples when appropriate sample preparation techniques have been applied (refer to Appendix A).

Analyte(s)	CAS RN
Chlorine	7782-50-5

A portion of aqueous liquid sample is buffered and reacted with N,N-diethyl-*p*-phenylenediamine (DPD) color agent. The resulting free chlorine is determined by colorimetry. If total chlorine (including chloramines and nitrogen trichloride) is to be determined, KI crystals are added. Results for chromate and manganese are blank corrected using thioacetamide solution. The method can detect 10 µg/L chlorine. Organic contaminants and strong oxidizers may cause interference.

Source: American Public Health Association, American Water Works Association, and Water Environment Federation. 2005. *Standard Methods for the Examination of Water and Wastewater*. 21st Edition. <http://www.standardmethods.org/>

5.2.78 Literature Reference for Chlorine (Analyst, 1999. 124: 1853–1857)

This procedure should be used for **preparation** of air samples for the contaminant identified below and listed in Appendix A. Note: Standard Method 4500-Cl G should be used as the determinative method (refer to Appendix A).

Analyte(s)	CAS RN
Chlorine	7782-50-5

A procedure is described for determination of total combined gas-phase active chlorine (i.e., Cl₂, HOCl and chloramines) and is based on a sulfonamide-functionalized silica gel sorbent. For determination of the collected chlorine, a modified version of the N,N-diethyl-*p*-phenylenediamine (DPD) colorimetric procedure is used, which yielded a detection limit of 0.1 mg of chlorine. At flow rates ranging from 31 to 294 mL/min, the collection efficiency was >90% based on breakthrough analysis. Recovery of chlorine spikes from 0.05-g aliquots of the sorbent was not quantitative (~60%) but was reproducible; the recovery is accounted for in samples by adding weighed amounts of sorbent to the standards.

Source: Johnson, B.J., Emerson, D.W., Song, L., Floyd, J., and Tadepalli, B. “Determination of active chlorine in air by bonded phase sorbent collection and spectrophotometric analysis,” *Analyst*, 124: 1853–1857 (1999).

5.2.79 Literature Reference for Fluoroacetate salts (Analytical Letters, 1994. 27 (14): 2703–2718)

The initial portion of this procedure (ultrasonic extraction) should be used for **preparation** of solid and non-aqueous liquid/organic solid samples for the contaminant identified below and listed in Appendix A. Note: EPA Method 300.1, Revision 1.0 should be used as the determinative method (refer to Appendix A).

Analyte(s)	CAS RN
Fluoroacetic acid and fluoroacetate salts	NA
Methyl fluoroacetate	453-18-9

Sodium fluoroacetate is determined at sub-microgram per gram concentrations in small (~1 g) soil samples. Samples are ultrasonically extracted with water, filtered, and analyzed by Method 300.1.

Source: Tomkins, B.A., “Screening-Procedure for Sodium Fluoroacetate (Compound 1080) at Sub-Microgram/Gram Concentrations in Soils,” *Analytical Letters*. 27(14), 2703–2718 (1994).

5.2.80 Literature Reference for Perfluoroisobutylene (Journal of Chromatography A, 2005. 1098: 156–165)

This procedure is not currently listed for any of the analyte/sample type combinations included in Appendix A. If problems occur when using OSHA Method 61 for the analysis of perfluoroisobutylene, then this procedure should be used for **preparation** and **analysis** of air samples for the contaminant identified below and listed in Appendix A. (See Footnote 15 of Appendix A.)

Analyte(s)	CAS RN
Perfluoroisobutylene (PFIB) ¹	382-21-8

¹ This analyte should be determined using this article **only** if problems occur when using OSHA Method 61.

This procedure is for the identification and measurement of phosgene and perfluoroisobutylene in air using gas chromatography/mass spectrometry (GC-MS), gas chromatography/nitrogen phosphorus detector (GC-NPD), or gas chromatography/flame photometric detector (GC-FPD). This procedure is to be used in the event that OSHA Method 61 is problematic for the determination of PFIB. Air samples are collected by drawing known volumes of air through sampling tubes containing Tenax TA sorbent coated with 2-aminothiophenol (ATP) or 3,4-dimercaptotoluene (DMT). The derivatized samples are thermally desorbed and analyzed with one of the GC techniques. The 2-aminothiophenol derivative can be analyzed by either of the three techniques, whereas the 3,4-dimercaptotoluene is not suitable for analysis by GC-NPD or GC-FPD. Limits of detection for PFIB-ATP, and PFIB-DMT, using 10-L air samples (typical sampling volume) by GC-MS analyses were determined to be 2 and 19 ng/m³ respectively.

Source: Muir, B., Cooper, D.B., Carric, W.A., Timperley, C.M., Slater, B.J., and S. Quick, “Analysis of Chemical Warfare Agents III. Use of bis-Nucleophiles in the Trace Level Determination of Phosgene and Perfluoroisobutylene,” *Journal of Chromatography A*, 1098, 156–165 (2005).

5.2.81 Literature Reference for Sodium Azide (Journal of Forensic Science, 1998. 43(1):200-202)

This procedure should be used for **preparation** of solid, non-aqueous liquid/organic solid, aqueous liquid, and drinking water samples for the contaminant identified below and listed in Appendix A.

Analyte(s)	CAS RN
Sodium azide	26628-22-8

Sample aliquots are acidified, converting azide to volatile hydrazoic acid (HN₃), which is then trapped in sodium hydroxide. The resulting aliquots are analyzed by ion chromatography, using a sodium tetraborate eluent and suppressed conductivity detection. This method can routinely quantify to at least 100 µg/L, and the detection limit is estimated to be 30 µg/L. Water extraction, filtration, and acidification steps should be used for the preparation of solid samples. Filtration and acidification steps should be used for preparation of aqueous liquid and drinking water samples. The acidification step should be used with EPA SW-846 Method 3580A for preparation of non-aqueous liquid/organic solid samples.

Please note: The procedure described above has been developed for the analysis of sodium azide in blood samples.

Source: Kruszyna R, Smith RP, Kruszyna H. “Determining Sodium Azide Concentration in the Blood by Ion Chromatography,” *Journal of Forensic Science*. 43(1):200-202 (1998).

Section 6.0: Selected Radiochemical Methods

A list of analytical methods to be used in analyzing environmental samples for radiochemical contaminants during homeland security events is provided in Appendix B. Methods are listed for each isotope and for each sample type that potentially may need to be measured and analyzed when responding to an environmental emergency.

Please note: This section provides guidance for selecting radiochemical methods that have a high likelihood of assuring analytical consistency when laboratories are faced with a large scale environmental restoration crisis. Not all methods have been verified for the analyte/sample type combination listed in Appendix B. Please refer to the specified method to identify analyte/sample type combinations that have been verified. Any questions regarding information discussed in this section should be addressed to the appropriate contact(s) listed in Section 4.

Appendix B is sorted alphabetically by analyte and includes the following information:

- **Analyte(s).** The radionuclide(s) or contaminant(s) of interest.
- **Chemical Abstract Survey Registration Number (CAS RN).** A unique identifier for chemical substances that provides an unambiguous way to identify a chemical or molecular structure when there are many possible systematic, generic, or trivial names. In this section (Section 6.0) and Appendix B, the CAS RNs correspond to the specific radionuclide identified.
- **Determinative technique.** An analytical instrument or technique used for qualitative and confirmatory determination of compounds or components in a sample.
- **Drinking water sample methods.** The recommended methods/procedures for sample preparation and analysis to measure the analyte of interest in drinking water samples. Methods have been identified for qualitative and confirmatory determination.
- **Aqueous and liquid phase sample methods.** The recommended methods/procedures for sample preparation and analysis to measure the analyte of interest in aqueous and/or non-aqueous liquid phase samples. Methods have been identified for qualitative and confirmatory determination.
- **Soil and sediment phase sample methods.** The recommended methods/procedures for sample preparation and analysis to measure the analyte of interest in soil and sediment samples. Methods have been identified for qualitative and confirmatory determination.
- **Surface wipe sample methods.** The recommended methods/procedures for sample preparation and analysis to measure the analyte of interest in surface wipe samples. Methods have been identified for qualitative and confirmatory determination.
- **Air filter sample methods.** The recommended methods/procedures for sample preparation and analysis to measure the analyte of interest in air filter samples. Methods have been identified for qualitative and confirmatory determination.
- **Qualitative determination method identifier.** A unique identifier or number assigned to an analytical method by the method publisher. The identified method is intended to determine the presence of a radiological element or isotope. These methods are less precise than confirmatory methods, and are used when greater sample throughput and more rapid reporting of results is required.
- **Confirmatory method identifier.** A unique identifier or number assigned to an analytical method by the method publisher. The identified method is for measurement of the activity from a particular radioisotope per unit of mass, volume, or area sampled.

6.1 General Guidance

The guidance summarized in this section provides a general overview of how to identify the appropriate radiochemical method(s) for a given analyte-sample type combination as well as recommendations for quality control procedures.

For additional information on the properties of the radionuclides listed in Appendix B, TOXNET (<http://toxnet.nlm.nih.gov/index.html>), a cluster of databases on toxicology, hazardous chemicals, and related areas maintained by the National Library of Medicine, is an excellent resource. EPA's Radiation Information (<http://www.epa.gov/radiation/radionuclides/index.html>) and the Multi-Agency Radiological Laboratory Analytical Protocols Manual (MARLAP) (<http://www.epa.gov/radiation/marlap/manual.html>) Web sites provide some additional information pertaining to radionuclides of interest and radiochemical methods.

6.1.1 Standard Operating Procedures for Identifying Radiochemical Methods

To determine the appropriate method that is to be used on an environmental sample, locate the analyte of concern in Appendix B: Radiochemical Methods under the "Analyte Class" or "Analyte(s)" column. After locating the analyte of concern, continue across the table to identify the appropriate determinative technique (e.g., alpha or beta spectrometry), then identify the appropriate qualitative and/or confirmatory method for the sample type of interest (drinking water, aqueous and liquid phase, soil and sediment, surface wipes, and air filters) for the particular analyte.

Sections 6.2.1 through 6.2.27 below provide summaries of the qualitative and confirmatory methods listed in Appendix B. Once a method has been identified in Appendix B, **Table 6-1** can be used to locate the method summary.

Table 6-1. Radiochemical Methods and Corresponding Text Section Numbers

Analyte / Analyte Class	CASRN	Method	Section
Gross Alpha	NA	900.0 (EPA)	6.2.2
		FRMAC, Vol 2, pg. 33	6.2.16
Gross Beta	NA	AP1 (ORISE)	6.2.17
		7110 B (SM)	6.2.21
Gamma	NA	901.1 (EPA)	6.2.3
		Ga-01-R (HASL-300)	6.2.12
Americium-241	14596-10-2	Am-01-RC (HASL-300)	6.2.9
		Am-02-RC (HASL-300)	6.2.10
		Am-04-RC (HASL-300)	6.2.11
		Pu-12-RC (HASL-300)	6.2.14
		AP11 (ORISE)	6.2.18
		D3084 (ASTM)	6.2.19
Californium-252	13981-17-4	Am-01-RC (HASL-300)	6.2.9
		Am-04-RC (HASL-300)	6.2.11
		Pu-12-RC (HASL-300)	6.2.14
		AP11 (ORISE)	6.2.18
		D3084 (ASTM)	6.2.19

Analyte / Analyte Class	CASRN	Method	Section
Cesium-137	10045-97-3	901.1 (EPA)	6.2.3
		Ga-01-R (HASL-300)	6.2.12
Cobalt-60	10198-40-0	7120 (SM)	6.2.22
Curium-244	13981-15-2	Am-01-RC (HASL-300)	6.2.9
		Am-04-RC (HASL-300)	6.2.11
		Pu-12-RC (HASL-300)	6.2.14
		AP11 (ORISE)	6.2.18
		D3084 (ASTM)	6.2.19
Europium-154	15585-10-1	901.1 (EPA)	6.2.3
Iridium-192	14694-69-0	Ga-01-R (HASL-300)	6.2.12
		7120 (SM)	6.2.22
Plutonium-238	13981-16-3	EMSL-33 (EPA)	6.2.8
Plutonium-239	15117-48-3	AP11 (ORISE)	6.2.18
		D3084 (ASTM)	6.2.19
Polonium-210	13981-52-7	Method 111 (EPA)	6.2.1
		Po-02-RC (HASL-300)	6.2.13
Radium-226	13982-63-3	903.0 (EPA)	6.2.4
		903.1 (EPA)	6.2.5
		EMSL-19 (EPA)	6.2.7
		D3084 (ASTM)	6.2.19
		7500-Ra B (SM)	6.2.23
		7500-Ra C (SM)	6.2.24
Ruthenium-103	13968-53-1	901.1 (EPA)	6.2.3
Ruthenium-106	13967-48-1	Ga-01-R (HASL-300)	6.2.12
Selenium-75	14265-71-5	7120 (SM)	6.2.22
Strontium-90	10098-97-2	Sr-03-RC (HASL-300)	6.2.15
		7500-Sr B (SM)	6.2.25
Uranium-238	7440-61-1	908.0 (EPA)	6.2.6
		EMSL-33 (EPA)	6.2.8
		AP11 (ORISE)	6.2.18
		D3084 (ASTM)	6.2.19
		D3972 (ASTM)	6.2.20
		7500-U B (SM)	6.2.26
		7500-U C (SM)	6.2.27

The method summaries are listed in order of method selection hierarchy (see Figure 2-1), starting with EPA methods followed by methods from other federal agencies and voluntary consensus standard bodies (VCSB). Methods are listed in numerical order under each publisher. Where available, a direct link to the full text of the selected analytical method is provided in the method summary. For additional information regarding sample preparation and analysis procedures and on methods available through consensus standards organizations, please use the contact information provided in **Table 6-2**.

Table 6-2. Sources of Radiochemical Methods

Name	Publisher	Reference
National Environmental Methods Index (NEMI)	U.S. Environmental Protection Agency (USEPA), United States Geological Survey (USGS)	http://www.nemi.gov
CFR Promulgated Test Methods (TM)	U.S. Environmental Protection Agency (USEPA), Technical Transfer Network (TTN) Emission Measurement Center (EMC)	http://www.epa.gov/ttn/emc/promgate.html
Prescribed Procedures for Measurement of Radioactivity in Drinking Water (EPA-600 4-80-032, August 1980)	U.S. Environmental Protection Agency (USEPA), Office of Research and Development (ORD), Environmental Monitoring and Support Laboratory (EMSL)	Available from National Technical Information Service (NTIS)*. NTIS, U.S. Department of Commerce, 5285 Port Royal Road, Springfield, VA 22161, (703) 605-6000.
Radiochemical Analytical Procedures for Analysis of Environmental Samples, March 1978. EMSL-LV-0539-17	United States Environmental Protection Agency (USEPA) Environmental Monitoring and Support Laboratory (EMSL)	Available from National Technical Information Service (NTIS)*. NTIS, U.S. Department of Commerce, 5285 Port Royal Road, Springfield, VA 22161, (703) 605-6000.
EML Procedures Manual, HASL-300, 28 th Edition, February, 1997	U.S. Department of Energy (DOE), Environmental Measurements Laboratory (EML) / Now, U.S. Department of Homeland Security (DHS)	http://www.eml.st.dhs.gov/publications/procman.cfm Also available from National Technical Information Service (NTIS)*. NTIS, U.S. Department of Commerce, 5285 Port Royal Road, Springfield, VA 22161, (703) 605-6000.
Federal Radiological Monitoring and Assessment Center (FRMAC) Laboratory Manual	United States Department of Energy (DOE) National Nuclear Security Administration (NNSA)	http://www.nv.doe.gov/nationalsecurity/homelandsecurity/frmac/manuals.aspx
Oak Ridge Institute for Science and Education (ORISE) Laboratory Procedures Manual	Oak Ridge Institute of Science and Education (ORISE) Independent Environmental Assessment and Verification	http://orise.ornl.gov/ieav/survey-projects/lab-manual.htm
Annual Book of ASTM Standards, Vol. 11.02*	American Society for Testing and Materials (ASTM) International	http://www.astm.org
Standard Methods for the Examination of Water and Wastewater, 21 st Edition, 2005*	American Public Health Association (APHA), American Water Works Association (AWWA), and Water Environment Federation (WEF)	http://www.standardmethods.org

* Subscription and/or purchase required.

6.1.2 General Quality Control (QC) Guidance for Radiochemical Methods

Having data of known and documented quality is critical for public officials to assess accurately the activities that may be needed in responding to emergency situations. Having such data requires that laboratories: (1) conduct the necessary QC to ensure that measurement systems are in control and operating correctly, (2) properly document results of the analyses, and (3) properly document measurement system evaluation of the analysis-specific QC. Ensuring data quality also requires that laboratory results are properly evaluated and the results of the data quality evaluation are transmitted to decision makers.

The level or amount of QC needed often depends on the intended purpose of the data that are generated. Various levels of QC may be required if the data are generated during contaminant presence/absence qualitative determinations versus confirmatory analyses. The specific needs for data generation should be identified. Quality control requirements and data quality objectives should be derived based on those needs, and should be applied consistently across laboratories when multiple laboratories are used. For example, during rapid sample screening analyses, minimal QC samples (e.g., blanks, duplicates) and documentation might be required to ensure data quality. Implementation of the analytical methods for evaluation of environmental samples during site assessment through site clearance, such as those identified in this document, might require increased QC.

Some method-specific QC requirements are described in many of the individual methods that are cited in this manual. QC requirements will be referenced in standardized analytical protocols developed to address specific analytes and sample types of concern. Individual methods, sampling and analysis protocols, or contractual statements of work also should be consulted to determine any additional QC that may be needed.

QC samples are required to assess the precision, bias, and reliability of sample results. All QC results are tracked on control charts and reviewed for acceptability and trends in analysis or instrument operation. Quality control parameters are measured as required per method at the prescribed frequency. QC of laboratory analyses using radiochemical methods includes ongoing analysis of QC samples and tracking QC parameters including, but not limited to the following:

- Method blanks;
- Calibration checks;
- Sample and sample duplicates;
- Laboratory control sample recoveries for samples that are not chemically prepared; or
- Matrix spike and matrix spike duplicate recoveries for samples that are chemically prepared; and
- Tracer and/or carrier yield.

Please note: The appropriate point of contact identified in Section 4 should be consulted regarding appropriate QA/QC procedures prior to sample analysis. These contacts will consult with the EPA coordinator responsible for laboratory activities during the specific event to ensure QA/QC procedures are performed consistently across laboratories. EPA program offices will be responsible for ensuring that the QA/QC practices are implemented.

6.1.3 Safety and Waste Management

It is imperative that safety precautions be used during collection, processing, and analysis of environmental samples. Laboratories should have a documented health and safety plan for handling samples that may contain target chemical, biological, or radiological contaminants, and laboratory staff should be trained in and implement the safety procedures included in the plan. In addition, many of the methods summarized or cited in Section 6.2 contain specific requirements, guidance, or information regarding safety precautions that should be followed when handling or processing environmental samples

and reagents. These methods also provide information regarding waste management. Other resources that can be consulted for additional information include the following:

- Occupational Health and Safety Administration’s Standard for Occupational Exposure to Hazardous Chemicals in Laboratories (29 CFR 1910.1450)
- Environmental Protection Agency’s Standards Regulating Hazardous Waste (40 CFR parts 260–270)
- Standards for Protection Against Radiation (10 CFR part 20)
- U.S. Department of Energy (DOE). Order O 435.1: Radioactive Waste Management. July 1, 1999. Available at: www.directives.doe.gov/pdfs/doe/doetext/neword/435/o4351.html
- U.S. Department of Energy (DOE). M 435.1-1. *Radioactive Waste Management Manual*. Office of Environmental Management. July 9, 1999. Available at: <http://www.directives.doe.gov/pdfs/doe/doetext/neword/435/m4351-1.html>
- U.S. Department of Energy (DOE). *Compendium of EPA-Approved Analytical Methods for Measuring Radionuclides in Drinking Water*. Prepared by the Office of Environmental Policy and Assistance Air, Water and Radiation Division (EH-412). June 1998
- U.S. Environmental Protection Agency (EPA). 1996. *Profile and Management Options for EPA Laboratory Generated Mixed Waste*. Office of Radiation and Indoor Air, Washington, DC
- EPA 402-R-96-015. August. Profile and Management Options for EPA Laboratory Generated Mixed Waste. Available at: <http://epa.gov/radiation/docs/mixed-waste/402-r-96-015.pdf>
- U.S. Environmental Protection Agency (EPA). 2001. Changes to 40 CFR 266 (Storage, Treatment, Transportation, and Disposal of Mixed Waste), *Federal Register* 66:27217-27266, May 16
- U.S. Environmental Protection Agency (EPA). 2002. *RCRA Orientation Manual*. Office of Solid Waste, Washington, DC. EPA530-R-02-016. 259 pp. Available at: <http://www.epa.gov/epaoswer/general/orientat/>
- Waste Management in a Radioanalytical Laboratory, Chapter 17 Multi-Agency Radiological Laboratory Analytical Protocols (MARLAP) Manual, July 2004
- National Research Council. 1995. *Prudent Practices in the Laboratory; Handling and Disposal of Chemicals*, National Academy Press, Washington, DC
- National Council on Radiation Protection and Measurements (NCRP). 2002. *Risk-Based Classification of Radioactive and Hazardous Chemical Wastes*, 7910 Woodmont Avenue, Suite 400, Bethesda, MD 20814-3095
- U.S. Nuclear Regulatory Commission/U.S. Environmental Protection Agency (NRC/EPA). 1995. Low-Level Mixed Waste Storage Guidance, *Federal Register* 60:40204-40211, August 7

6.2 Method Summaries

Summaries for the analytical methods listed in Appendix B are provided in Sections 6.2.1 through 6.2.27. These summaries contain information that has been extracted from the selected methods. Each method summary contains a table identifying the contaminants in Appendix B to which the method applies, a brief description of the analytical method, and a link to the full version of the method or a source for obtaining a full version of the method. The full version of the method should be consulted prior to sample analysis.

6.2.1 EPA Method 111: Determination of Polonium-210 Emissions from Stationary Sources

This method should be used for **qualitative** and **confirmatory analysis** of surface wipes and air filters for the contaminant identified below and listed in Appendix B.

Analyte(s)	CAS RN
Polonium-210	13981-52-7

This method covers the determination of polonium-210 in particulate matter samples collected from stationary sources such as exhaust stacks. Polonium-210 in the sample is put in solution, deposited on a metal disc, and the radioactive disintegration rate measured. Polonium in acid solution spontaneously deposits on surface metals that are more electropositive than polonium.

Source: 40 CFR 61 Appendix B; National Emission Standards for Air Pollutants; Appendix B to Part 61 - Test Methods, U.S. Environmental Protection Agency (EPA). Emission Measurement Center (EMC). February 2000. Prepared by the Office of Air Quality Planning and Standards (OAQPS), Research Triangle Park, North Carolina, 27711. Also at: <http://www.epa.gov/ttn/emc/promgate.html>

6.2.2 EPA Method 900.0: Gross Alpha and Gross Beta Radioactivity in Drinking Water

This method should be used for **gross alpha** and **gross beta determination** in drinking water samples. The method provides an indication of the presence of alpha and beta emitters, including the following SAM analytes:

• Americium-241	(CAS RN 14596-10-2)	Alpha emitter
• Californium-252	(CAS RN 13981-17-4)	Alpha emitter
• Cesium-137	(CAS RN 10045-97-3)	Beta emitter
• Cobalt-60	(CAS RN 10198-40-0)	Beta emitter
• Curium-244	(CAS RN 13981-15-2)	Alpha emitter
• Europium-154	(CAS RN 15585-10-1)	Beta emitter
• Iridium-192	(CAS RN 14694-69-0)	Beta emitter
• Plutonium-238	(CAS RN 13981-16-3)	Alpha emitter
• Plutonium-239	(CAS RN 15117-48-3)	Alpha emitter
• Polonium-210	(CAS RN 13981-52-7)	Alpha emitter
• Radium-226	(CAS RN 13982-63-3)	Alpha emitter
• Ruthenium-103	(CAS RN 13968-53-1)	Beta emitter
• Ruthenium-106	(CAS RN 13967-48-1)	Beta emitter
• Strontium-90	(CAS RN 10098-97-2)	Beta emitter
• Uranium-238	(CAS RN 7440-16-1)	Alpha emitter

An aliquot of a preserved drinking water sample is evaporated to a small volume (3 to 5 mL) and transferred quantitatively to a tared 2-inch planchet. The aliquot volume is determined based on a maximum total solids content of 100 mg. The sample aliquot is evaporated to dryness in the planchet to a constant weight, cooled, and counted using a gas proportional or scintillation counting system. The counting system is calibrated with thorium-230 for gross alpha, or with cesium-137 for gross beta analysis. A traceable standards-based absorption curve must be developed for each calibration nuclide (Th-230 or Cs-137) based on a range of total solids content in the 2-inch planchet from 0 to 100 mg. The results are corrected for the absorption factor based total solids on the planchet.

Source: *Prescribed Procedures for Measurement of Radioactivity in Drinking Water*, National Exposure Risk Laboratory-Cincinnati (NERL-CI), EPA/600/4/80/032, August 1980, available from National Technical Information Service (NTIS), 5285 Port Royal Road, Springfield, VA 22161. Phone: 800-553-6847. Also at http://infotrek.er.usgs.gov/pls/apex/f?p=119:38:1908830694668232:::P38_METHOD_ID:4730

6.2.3 EPA Method 901.1: Gamma Emitting Radionuclides in Drinking Water

This method should be used for **qualitative** and **confirmatory analysis** of drinking water samples for the contaminants identified below and listed in Appendix B.

Analyte(s)	CAS RN
Cesium-137 ¹	10045-97-3
Cobalt-60	10198-40-0
Europium-154	15585-10-1
Iridium-192	14694-69-0
Ruthenium-103	13968-53-1
Ruthenium-106 ¹	13967-48-1
Selenium-75	14265-71-5

¹ The method will measure decay products of these isotopes

This method is applicable for analysis of water samples that contain radionuclides that emit gamma photons with energies ranging from approximately 60 to 2000 keV. The method uses gamma spectroscopy for measurement of gamma photons emitted from radionuclides without separating them from the sample matrix. A homogeneous aliquot of water is placed into a standard geometry (normally a Marinelli beaker) for gamma counting, typically using a High Purity Germanium (HPGe) detector. Detectors such as Germanium (Lithium) (Ge(Li)) or thallium-activated sodium iodide (NaI(Tl)) also can be used. Sample aliquots are counted long enough to meet the required sensitivity of measurement. To reduce adsorbance of radionuclides on the walls of the counting container, the sample is acidified at collection time. Due to its lower resolution, significant interference can occur using the NaI(Tl) detector when counting a sample containing radionuclides that emit gamma photons of similar energies. When using this method, shielding is needed to reduce background interference. Detection limits are dependent on sample volume, geometry (physical shape), and counting time.

Source: Prescribed Procedures for Measurement of Radioactivity in Drinking Water, National Exposure Risk Laboratory-Cincinnati (NERL-CI), EPA/600/4/80/032, August 1980, available from National Technical Information Service (NTIS), 5285 Port Royal Road, Springfield, VA 22161. Phone: 800-553-6847.

6.2.4 EPA Method 903.0: Alpha-Emitting Radium Isotopes in Drinking Water

This method should be used for **qualitative determination** in drinking water samples for the contaminant identified below and listed in Appendix B.

Analyte(s)	CAS RN
Radium-226 ¹	13982-63-3

¹ The method will measure decay products of these isotopes

This method covers measurement of the total soluble alpha emitting radioisotopes of radium, namely radium-223, radium-224 and radium-226 in drinking water. The method does not give an accurate measurement of radium-226 content in the sample when other alpha emitters are present. If radium-223 and radium-224 are present, the results can be used to provide a gross determination of radium-226. When the total radium alpha activity of a drinking water sample is greater than 5 pCi/L, use of Method 903.1 (Radium-226 in Drinking Water) is preferred. Radium in the water sample is collected by coprecipitation

with barium and lead sulfate, and purified by re-precipitation from ethylenediaminetetraacetic acid (EDTA) solution. Citric acid is added to ensure that complete interchange occurs before the first precipitation step. The final barium sulfate precipitate is alpha counted to determine the total disintegration rate of the radium isotopes. By making a correction for the ingrowth of radon and its alpha emitting progeny for the elapsed time after separation, one can determine radium activity in the sample. Presence of significant natural barium in the sample can result in a falsely high yield. Based on a 1000-mL sample and 100-minute counting time, the minimum detectable level for this method is 0.5 pCi/L.

Source: Prescribed Procedures for Measurement of Radioactivity in Drinking Water, National Exposure Risk Laboratory-Cincinnati (NERL-CI), EPA/600/4/80/032, August 1980, available from National Technical Information Service (NTIS), 5285 Port Royal Road, Springfield, VA 22161. Phone: 800-553-6847. Also at http://infotrek.er.usgs.gov/pls/apex/f?p=119:38:1908830694668232:::P38_METHOD_ID:9183

6.2.5 EPA Method 903.1: Radium-226 in Drinking Water - Radon Emanation Technique

This method should be used for **confirmatory analysis** of drinking water samples for the contaminant identified below and listed in Appendix B.

Analyte(s)	CAS RN
Radium-226 ¹	13982-63-3

¹ The method will measure decay products of these isotopes

This method is specific for radium-226, and is based on the emanation and scintillation counting of radon-222, a daughter product of radium-226. Radium-226 is concentrated and separated from the water sample by coprecipitation on barium sulfate. The precipitate is dissolved in EDTA reagent, placed in a sealed bubbler and stored for ingrowth of radon-222. After ingrowth, the gas is purged into a scintillation cell. When the short-lived radon-222 daughters are in equilibrium with the parent (after ~4h), the scintillation cell is counted for activity. The absolute measurement of radium-226 is effected by calibrating the scintillation cell system with a standard solution of the nuclide. There are no radioactive interferences in this method. Based on a 1000-mL sample and 100-minute counting time, the minimum detectable level for this method is 0.5 pCi/L.

Source: Prescribed Procedures for Measurement of Radioactivity in Drinking Water, National Exposure Risk Laboratory-Cincinnati (NERL-CI), EPA/600/4/80/032, August 1980, available from National Technical Information Service (NTIS), 5285 Port Royal Road, Springfield, VA 22161. Phone: 800-553-6847. Also at http://infotrek.er.usgs.gov/pls/apex/f?p=119:38:1701708976536019:::P38_METHOD_ID:4732

6.2.6 EPA Method 908.0: Uranium in Drinking Water - Radiochemical Method

This method should be used for **qualitative determination** in drinking water samples for the contaminant identified below and listed in Appendix B.

Analyte(s)	CAS RN
Uranium-238 ¹	7440-61-1

¹ If it is suspected that the sample exists in refractory form (i.e., non-digestable or dissolvable material after normal digestion methods) or if there is a matrix interference problem, use ORISE Method AP11

This method measures total uranium alpha activity of a sample, without doing an isotopic uranium analysis. The sample is acidified with hydrochloric acid and boiled to eliminate carbonate and bicarbonate ions. Uranium is coprecipitated with ferric hydroxide and separated from the sample. The uranium is then separated from other radionuclides that were carried down with the ferric hydroxide by dissolving the hydroxide precipitate in hydrochloric acid, putting the solution through an anion exchange

column, washing the column with hydrochloric acid, and finally eluting the uranium with hydrochloric acid. The uranium eluate is evaporated and the uranium chemical form is converted to nitrate. The residue is transferred to a stainless steel planchet, dried, flamed, and counted for alpha particle activity. Since uranium is a naturally occurring radionuclide, reagents must be checked for uranium contamination by analyzing a complete reagent blank by the same procedure as used for the samples. Based on a 1000-mL sample and 100-minute counting time in a single laboratory study, the minimum detectable level for this method is 1.0 pCi/L.

Source: Prescribed Procedures for Measurement of Radioactivity in Drinking Water, National Exposure Risk Laboratory-Cincinnati (NERL-CI), EPA/600/4/80/032, August 1980, available from National Technical Information Service (NTIS), 5285 Port Royal Road, Springfield, VA 22161. Phone: 800-553-6847. Also at http://infotrek.er.usgs.gov/pls/apex/f?p=119:38:1908830694668232:::P38_METHOD_ID:9184

6.2.7 EPA Method EMSL-19: Determination of Radium-226 and Radium-228 in Water, Soil, Air and Biological Tissue

This method should be used for **confirmatory analysis** of soil/sediment, surface wipe, and air filter samples for the contaminant identified below and listed in Appendix B.

Analyte(s)	CAS RN
Radium-226 ¹	13982-63-3

¹ The method will measure decay products of these isotopes

Following acid digestion and filtration of soil, sediment, surface wipe, or air filter samples, radium is precipitated with barium sulfate. Barium-radium-sulfate is dissolved in a pentasodium diethylenetriamine-pentaacetate solution and transferred to an emanation tube. The radon is allowed to come to equilibrium for approximately 30 days. Radium-226 decays by alpha emission to radon-222. Radon-222 is separated and collected from the liquid by a de-emanation technique. The radon is counted by alpha scintillation 4.5 hours after de-emanation, at which time the short-lived progeny have reached 97% of equilibrium. An applicable measurement range has not been determined; however, samples that contain 0.1 pCi of Radium-226 have been analyzed.

Source: Radiochemical Analytical Procedures for Analysis of Environmental Samples, United States Environmental Protection Agency, Environmental Monitoring and Support Laboratory (EMSL), March 1979, available from National Technical Information Service (NTIS), 5285 Port Royal Road, Springfield, VA 22161. Phone: 800-553-6847.

6.2.8 EPA Method EMSL-33: Isotopic Determination of Plutonium, Uranium, and Thorium in Water, Soil, Air, and Biological Tissue

This method should be used for **confirmatory analysis** of drinking water, aqueous/liquid, soil/sediment, surface wipe, and/or air filter samples for the contaminants identified below and listed in Appendix B.

Analyte(s)	CAS RN
Plutonium-238 ¹	13981-16-3
Plutonium-239 ¹	15117-48-3
Uranium-238 ¹	7440-61-1

¹ If it is suspected that the sample exists in refractory form (i.e., non-digestable or dissolvable material after normal digestion methods) or if there is a matrix interference problem, use ORISE Method AP11

This method is appropriate for the analysis of isotopic plutonium, uranium, and thorium, together or individually, by alpha spectrometry. Plutonium-236, uranium-232, and thorium-234 tracer standards are added for the determination of chemical yields. Samples are decomposed by nitric-hydrofluoric acid digestion or ignition to assure that all of the plutonium is dissolved and chemically separated from the sample by coprecipitation with sodium and ammonium hydroxide, anion exchange, and electrodeposition. The residues are dissolved in dilute nitric acid and successive sodium and ammonium hydroxide precipitations are performed in the presence of boric acid to remove fluoride and soluble salts. The hydroxide precipitate is dissolved, the solution is pH-adjusted with hydrochloric acid, and plutonium and uranium are adsorbed on an anion exchange column, separating them from thorium. Plutonium is eluted with hydrobromic acid. The actinides are electrodeposited on stainless steel discs from an ammonium sulfate solution and subsequently counted by alpha spectrometry. This method is designed to detect environmental levels of activity as low as 0.02 pCi per sample. To avoid possible cross-contamination, sample aliquot activities should be limited to 25 pCi or less.

Source: Radiochemical Analytical Procedures for Analysis of Environmental Samples, United States Environmental Protection Agency, Environmental Monitoring and Support Laboratory (EMSL), March 1979, available from National Technical Information Service (NTIS), 5285 Port Royal Road, Springfield, VA 22161. Phone: 800-553-6847.

6.2.9 EML HASL-300 Method Am-01-RC: Americium in Soil

This method should be used for **confirmatory analysis** of soil/sediment samples for the contaminants identified below and listed in Appendix B.

Analyte(s)	CAS RN
Americium-241 ¹	14596-10-2
Californium-252 ¹	13981-17-4
Curium-244 ¹	13981-15-2

¹If it is suspected that the sample exists in refractory form (i.e., non-digestible or dissolvable material after normal digestion methods) or if there is a matrix interference problem, use ORISE Method AP11

This method uses alpha spectrometry for determination of americium-241 in soil, and also can be applied for determination of californium. Americium is leached from soil with nitric acid and hydrochloric acid. Americium-243 is added as a tracer to determine chemical yield. The soil is processed through the plutonium separation steps using ion exchange resin according to Method Pu-11-RC. Americium is collected with a calcium oxalate precipitation and finally isolated and purified by ion exchange. Californium-252 and curium-244 are eluted with americium as americium is stripped off the column. After source preparation by microprecipitation, americium-241, californium-252, and curium-244 are determined by alpha spectrometry analysis. The counting period chosen depends on the sensitivity required of the measurement and the degree of uncertainty in the result that is acceptable. The lower limit of detection (LLD) for americium-241 is 0.5 mBq when counted for 1000 minutes. In cases where less than 100 g of sample is available, use of Pu-12-RC is recommended.

Source: “Am-01-RC, Pu-11-RC, and Pu-12-RC.” EML Procedures Manual, HASL-300, 28th Edition, Environmental Measurements Laboratory (EML), Department of Energy (EML is currently part of the U.S. Department of Homeland Security), February 1997. Web: <http://www.eml.st.dhs.gov/publications/procman.cfm>

6.2.10 EML HASL-300 Method Am-02-RC: Americium-241 in Soil-Gamma Spectrometry

This method should be used for **qualitative determination in** soil/sediment samples for the contaminant identified below and listed in Appendix B.

Analyte(s)	CAS RN
Americium-241 ¹	14596-10-2

¹If it is suspected that the sample exists in refractory form (i.e., non-digestable or dissolvable material after normal digestion methods) or if there is a matrix interference problem, use ORISE Method AP11

This method uses gamma spectrometry for determination of americium-241 in soil. Americium-241 decays with the emission of a gamma ray at 59.5 keV with a decay frequency (abundance or yield) of 35.9%. The sample is placed into an appropriately sized standard geometry (normally a Marinelli beaker) after drying and grinding the sample for homogenization. Gamma-ray attenuation corrections are required if the calibration source and the sample are in a different matrix or are of different densities. The lower limit of detection (LLD) for 600 to 800 g of soil in a Marinelli beaker is 0.74 mBq for a 1000-minute count.

Source: *EML Procedures Manual*, HASL-300, 28th Edition, Environmental Measurements Laboratory (EML), Department of Energy (EML is currently part of the U.S. Department of Homeland Security), February 1997. Web: <http://www.eml.st.dhs.gov/publications/procman.cfm>

6.2.11 EML HASL-300 Method Am-04-RC: Americium in QAP Water and Air Filters - Eichrom's TRU Resin

This method should be used for **confirmatory analysis** of drinking water and aqueous/liquid samples for the contaminants identified below and listed in Appendix B.

Analyte(s)	CAS RN
Americium-241 ¹	14596-10-2
Californium-252 ¹	13981-17-4
Curium-244 ¹	13981-15-2

¹If it is suspected that the sample exists in refractory form (i.e., non-digestable or dissolvable material after normal digestion methods) or if there is a matrix interference problem, use ORISE Method AP11

This method is specific to measurement of americium isotopes in samples that do not contain lanthanides, but also can be used for measurement of californium and curium. The method uses microprecipitation and determination by alpha spectrometry. Americium-243 is added to the sample to determine chemical yield. The sample is processed through separation steps using ion exchange resins. The eluate from the ion exchange column containing americium (and all other ions, except plutonium) is evaporated, redissolved, and loaded onto a Transuranic (TRU) Resin extraction column. Americium (and curium and californium, if present) is separated and purified on the column and finally stripped with dilute nitric acid stripping solution. Microprecipitation is used to prepare for alpha spectrometry. The method involves sample preparation steps from EML HASL-300 Method Pu-10-RC for water samples. The lower limit of detection (LLD) for total americium is 0.3 mBq when counted for 1000 minutes.

Source: *EML Procedures Manual*, HASL-300, 28th Edition, Environmental Measurements Laboratory (EML), Department of Energy (EML is currently part of the U.S. Department of Homeland Security), February 1997. Web: <http://www.eml.st.dhs.gov/publications/procman.cfm>

6.2.12 EML HASL-300 Method Ga-01-R: Gamma Radioassay

This method should be used for **qualitative** and/or **confirmatory analysis** of soil/sediment, surface wipes, and/or air filter samples for the contaminants identified below and listed in Appendix B.

Analyte(s)	CAS RN
Cesium-137 ¹	10045-97-3
Cobalt-60	10198-40-0
Europium-154	15585-10-1
Iridium-192	14694-69-0
Ruthenium-103	13968-53-1
Ruthenium-106 ¹	13967-48-1
Selenium-75	14265-71-5

¹ The method will measure decay products of these isotopes

This method uses gamma spectroscopy for the measurement of gamma photons emitted from radionuclides without separating them from the sample matrix. Samples are placed into a standard geometry for gamma counting, typically using a High Purity Germanium (HPGe) detector. Detectors such as Germanium (Lithium) (Ge(Li)) or thallium-activated sodium iodide (NaI(Tl)) also can be used. The sample is placed into a standard geometry (physical shape) for gamma counting. Soil samples and sludge are placed into an appropriately sized Marinelli beaker after drying and grinding the sample for homogenization. Air filters and surface wipes can be counted directly or pressed into a planchet and counted. Samples are counted long enough to meet the required sensitivity of measurement. For typical counting systems and sample types, activity levels of approximately 40 Bq are measured, and sensitivities as low as 0.002 Bq can be achieved for many nuclides. Because of electronic limitations, count rates higher than 2000 counts per second (cps) should be avoided. High activity samples may be diluted, reduced in size, or moved away from the detector (a limited distance) to reduce the count rate and allow for analysis. The method is applicable for analysis of samples that contain radionuclides emitting gamma photons with energies ranging from approximately >40 keV for Germanium (Lithium) (Ge(Li)) and 100 keV for thallium-activated sodium iodide (NaI(Tl)) detectors.

Source: *EML Procedures Manual*, HASL-300, 28th Edition, Environmental Measurements Laboratory (EML), Department of Energy (EML is currently part of the U.S. Department of Homeland Security), February, 1997. Web: <http://www.eml.st.dhs.gov/publications/procman.cfm>

6.2.13 EML HASL-300 Method Po-02-RC: Polonium in Water, Vegetation, Soil, and Air Filters

This method should be used for **qualitative** and **confirmatory analysis** of drinking water, aqueous/liquid phases and soil/sediment samples for the contaminant identified below and listed in Appendix B.

Analyte(s)	CAS RN
Polonium-210	1-3981-52-7

This method uses alpha spectrometry for determination of polonium in water, vegetation, soil, and air filter samples. Polonium equilibrated with Po-208 or Po-209 tracer is isolated from most other elements by coprecipitation with lead sulfide. The sulfide precipitate is dissolved in weak hydrochloric acid solution. Polonium is quantitatively deposited on a nickel disc, and the plated disc is counted on an alpha spectrometer to measure chemical yield and activity of the sample. The solution from the deposition may

be retained and analyzed for Pb-210. When counted for 1000 minutes, the lower level of detection (LLD) for polonium is 1.0 mBq for water and 1.3 mBq for vegetation, soil and filters.

Source: *EML Procedures Manual*, HASL-300, 28th Edition, Environmental Measurements Laboratory (EML), Department of Energy (EML is currently part of the U.S. Department of Homeland Security), February, 1997. Web: <http://www.eml.st.dhs.gov/publications/procman.cfm>

6.2.14 EML HASL-300 Method Pu-12-RC: Plutonium and/or Americium in Soil or Sediments

This method is not currently listed for any of the analyte/sample type combinations included in Appendix B. In cases where only small sample volumes (≤ 100 g) will be analyzed, this method could be used for **confirmatory analysis** of soil/sediment samples for the contaminants identified below and listed in Appendix B.

Analyte(s)	CAS RN
Americium-241 ¹	14596-10-2
Californium-252 ¹	13981-17-4
Curium-244 ¹	13981-15-2

¹If it is suspected that the sample exists in refractory form (i.e., non-digestable or dissolvable material after normal digestion methods) or if there is a matrix interference problem, use ORISE Method AP11

A sample of soil of up to 100 g in size is equilibrated with Am-243 tracer. Contaminant isotopes are leached with nitric and hydrochloric acid. Plutonium is removed by ion exchange. The eluent from the plutonium separation is saved for determination of americium, curium, and californium. Americium, curium, and californium are collected with a calcium oxalate coprecipitation, isolated and purified by extraction chromatography. Microprecipitation is used to prepare the sample for analysis by alpha spectrometry of americium, curium, and californium. The lower limit of detection for Americium is 0.5 mBq when counted for 1000 minutes.

Source: *EML Procedures Manual*, HASL-300, 28th Edition, Environmental Measurements Laboratory (EML), Department of Energy (EML is currently part of the U.S. Department of Homeland Security), February, 1997. Web: <http://www.eml.st.dhs.gov/publications/procman.cfm>

6.2.15 EML HASL-300 Method Sr-03-RC: Strontium-90 in Environmental Samples

This method should be used for **qualitative** and **confirmatory analysis** of soil/sediment, surface wipes, and air filter samples for the contaminant identified below and listed in Appendix B.

Analyte(s)	CAS RN
Strontium-90 ¹	10098-97-2

¹ The method will measure decay products of these isotopes

Strontium is separated from calcium, other fission products, and natural radioactive elements. Fuming nitric acid separations remove the calcium and most other interfering ions. Radium, lead and barium are removed with barium chromate. Traces of other fission products are scavenged with iron hydroxide. After strontium-90 and yttrium-90 equilibrium has been attained, yttrium-90 is precipitated as the hydroxide and converted to oxalate for counting on a low-background gas proportional beta counter. Chemical yield is determined with strontium-85 tracer by counting in a gamma well detector.

Source: *EML Procedures Manual*, HASL-300, 28th Edition, Environmental Measurements Laboratory (EML), Department of Energy (EML is currently part of the U.S. Department of Homeland Security), February, 1997. Web: <http://www.eml.st.dhs.gov/publications/procman.cfm>

6.2.16 FRMAC Method Volume 2, Page 33: Gross Alpha and Beta in Air

This method should be used for **gross alpha** and **gross beta determination** in air filters, and can also be used as a direct counting of surface wipes. The method provides an indication of the presence of alpha and beta emitters, including the following SAM analytes:

• Americium-241	(CAS RN 14596-10-2)	Alpha emitter
• Californium-252	(CAS RN 13981-17-4)	Alpha emitter
• Cesium-137	(CAS RN 10045-97-3)	Beta emitter
• Cobalt-60	(CAS RN 10198-40-0)	Beta emitter
• Curium-244	(CAS RN 13981-15-2)	Alpha emitter
• Europium-154	(CAS RN 15585-10-1)	Beta emitter
• Iridium-192	(CAS RN 14694-69-0)	Beta emitter
• Plutonium-238	(CAS RN 13981-16-3)	Alpha emitter
• Plutonium-239	(CAS RN 15117-48-3)	Alpha emitter
• Polonium-210	(CAS RN 13981-52-7)	Alpha emitter
• Radium-226	(CAS RN 13982-63-3)	Alpha emitter
• Ruthenium-103	(CAS RN 13968-53-1)	Beta emitter
• Ruthenium-106	(CAS RN 13967-48-1)	Beta emitter
• Strontium-90	(CAS RN 10098-97-2)	Beta emitter
• Uranium-238	(CAS RN 7440-16-1)	Alpha emitter

This method allows for measurement of gross alpha and gross beta radiation in air samples. The method also can be applied for the analysis of surface swipes. A thin-window gas-flow proportional counter is used for counting gross alpha and beta radioactivity. The method supplies an approximation of the alpha and beta activity present in the air or the removable surface activity dependent on the sample type. For this application, the procedure requires the use of thorium-230 for alpha counting efficiency and cesium-137 for beta counting efficiency in the calibration of the detector. An air filter or swipe sample is placed onto a planchet, then counted for alpha and beta radioactivity. Activity is reported in activity units per volume of air sampled, as units of activity per surface area sampled, or as total units of activity in cases where sample collection information is not available.

Source: *FRMAC Monitoring and Analysis Manual – Sample Preparation and Analysis - Volume 2*, Federal Radiological Monitoring and Assessment Center (FRMAC), DOE/NV/1178-181 Vol. 2, UC-707, August 1998. Web: <http://www.nv.doe.gov/nationalsecurity/homelandsecurity/frmac/manuals.aspx>

6.2.17 ORISE Method AP-1: Gross Alpha and Beta for Various Matrices

This method should be used for **gross alpha** and **gross beta determination** in soil/sediment samples for the contaminants identified below and listed in Appendix B. This method provides an indication of the presence of alpha and beta emitters, including the following SAM analytes:

• Americium-241	(CAS RN 14596-10-2)	Alpha emitter
• Californium-252	(CAS RN 13981-17-4)	Alpha emitter
• Cesium-137	(CAS RN 10045-97-3)	Beta emitter
• Cobalt-60	(CAS RN 10198-40-0)	Beta emitter
• Curium-244	(CAS RN 13981-15-2)	Alpha emitter
• Europium-154	(CAS RN 15585-10-1)	Beta emitter
• Iridium-192	(CAS RN 14694-69-0)	Beta emitter
• Plutonium-238	(CAS RN 13981-16-3)	Alpha emitter

• Plutonium-239	(CAS RN 15117-48-3)	Alpha emitter
• Polonium-210	(CAS RN 13981-52-7)	Alpha emitter
• Radium-226	(CAS RN 13982-63-3)	Alpha emitter
• Ruthenium-103	(CAS RN 13968-53-1)	Beta emitter
• Ruthenium-106	(CAS RN 13967-48-1)	Beta emitter
• Strontium-90	(CAS RN 10098-97-2)	Beta emitter
• Uranium-238	(CAS RN 7440-16-1)	Alpha emitter

This method covers the measurement of gross alpha and gross beta in various sample types. This procedure provides screening measurements to indicate whether specific chemical analyses are required for water, soil, vegetation, and other solids. Liquid samples are acidified, concentrated, dried in a planchet, and counted in a low-background proportional counter. Solid samples are dried and processed to provide homogeneity, and a known quantity is transferred to a planchet and counted in a low-background proportional counter.

Source: *Laboratory Procedures Manual for the Environmental Survey and Site Assessment Program*; May 2001; Oak Ridge Institute for Science and Education (ORISE), Oak Ridge Associated Universities (ORAU). Web: <http://orise.orau.gov/ieav/survey-projects/pubs/lab-manual/19ap1.pdf>

6.2.18 ORISE Method AP-11: Sequential Determination of the Actinides in Environmental Samples Using Total Sample Dissolution and Extraction

This method is not currently listed for any of the analyte/sample type combinations included in Appendix B. If it is suspected that a sample exists in a refractory form (i.e., non-digestable or dissolvable material after normal digestion methods) or if there is a matrix interference problem, this method should be routinely used for **confirmatory analysis** and, in the event of refractory radioactive material, should be used for both **qualitative determination** and **confirmatory analysis** of drinking water, aqueous/liquid, soil/sediment, surface wipes, and air filter samples for the contaminants identified below and listed in Appendix B.

Analyte(s)	CAS RN
Americium-241	14596-10-2
Californium-252	13981-17-4
Curium-244	13981-15-2
Plutonium-238	13981-16-3
Plutonium-239	15117-48-3
Uranium-238	7440-61-1

This method covers the measurement of the actinides americium-241, californium-252, curium-244, plutonium-238 and -239, and uranium-238 in water, soil, and other solids. The method is used if there is an indication that the material is in a highly refractory or “non-dissolvable” form. Solid and unfiltered aqueous samples are completely dissolved by a combination of potassium hydrogen fluoride and pyrosulfate fusions. Filtered aqueous samples are evaporated to dryness followed by a pyrosulfate fusion. The fusion cake is dissolved and, for analyses requiring uranium only, two barium sulfate precipitations are performed and the uranium is separated using EDTA. For all other analyses, one barium sulfate precipitation is performed and all alpha emitters are coprecipitated on barium sulfate. The barium sulfate is dissolved and the actinides are separated by extraction chromatography. An optional section is presented for the separation of americium from the lanthanides. All actinides are coprecipitated on cerium fluoride and counted with an alpha spectrometer system.

Source: *Laboratory Procedures Manual for the Environmental Survey and Site Assessment Program*; May 2001; Oak Ridge Institute for Science and Education (ORISE), Oak Ridge Associated Universities (ORAU). Web: <http://orise.orau.gov/ieav/survey-projects/pubs/lab-manual/ap11.pdf>

6.2.19 ASTM Method D3084: Standard Practice for Alpha Spectrometry in Water

This method should be used for **qualitative determination analysis** of drinking water, aqueous/liquid, soil and sediment, surface wipes, and air filter samples for the contaminants identified below and listed in Appendix B.

Analyte(s)	CAS RN
Americium-241 ¹	14596-10-2
Californium-252 ¹	13981-17-4
Curium-244 ¹	13981-15-2
Plutonium-238 ¹	13981-16-3
Plutonium-239 ¹	15117-48-3
Radium-226 ²	13982-63-3
Uranium-238 ¹	7440-61-1

¹ If it is suspected that the sample exists in refractory form (i.e., non-digestible or dissolvable material after normal digestion methods) or if there is a matrix interference problem, use ORISE Method AP11

² The method will measure decay products of these isotopes

This standard practice covers the process that is required to obtain well-resolved alpha spectra from water samples and discusses the associated problems. This practice is typically followed with specific chemical separations and mounting techniques that are included in referenced methods. A chemical procedure is required to isolate and purify the radionuclides (see ASTM Methods D3865, *Standard Test Method for Plutonium in Water* and D3972, *Standard Test Method for Isotopic Uranium in Water by Radiochemistry*), and a radioactive tracer is added to determine yield. A source is prepared by employing electrodeposition, microprecipitation, or evaporation (depositing the solution onto a stainless steel or platinum disc). Electrodeposition and microprecipitation are preferred. The source's radioactivity is then measured in an alpha spectrometer according to manufacturer's operating instructions. The counting period chosen depends on the sensitivity required of the measurement and the degree of uncertainty in the result that is acceptable.

Source: *Annual Book of ASTM Standards*, Vol. 11.02, American Society for Testing and Materials (ASTM), 1996, ASTM International, 100 Barr Harbor Drive West, Conshohocken, PA 19428. Phone: 610-832-9500. Web: <http://www.astm.org>. Use method number when ordering.

6.2.20 ASTM Method D3972: Standard Test Method for Isotopic Uranium in Water by Radiochemistry

This method should be used for **confirmatory analysis** of drinking water samples for the contaminant identified below and listed in Appendix B.

Analyte(s)	CAS RN
Uranium-238 ¹	7440-61-1

¹ If it is suspected that the sample exists in refractory form (i.e., non-digestible or dissolvable material after normal digestion methods) or if there is a matrix interference problem, use ORISE Method AP11

This method covers the determination of uranium isotopes in water by means of chemical separations and alpha spectrometry analysis. Uranium is chemically separated from a water sample by coprecipitation with ferrous hydroxide followed by anion exchange, and electrodeposition. When suspended matter is present, an acid dissolution step is added to ensure that all of the uranium dissolves. The sample is acidified, and uranium-232 is added as an isotopic tracer to determine chemical yield. Uranium is coprecipitated from the sample with ferrous hydroxide. This precipitate is dissolved in concentrated hydrochloric acid, or is subjected to acid dissolution with concentrated nitric and hydrofluoric acids, if the hydrochloric acid fails to dissolve the precipitate. Uranium is separated from other radionuclides by adsorption on anion exchange resins, followed by elution with hydrochloric acid. The uranium is finally electrodeposited onto a stainless steel disc and counted using alpha spectrometry.

Source: Annual Book of ASTM Standards, Vol. 11.02, American Society for Testing and Materials (ASTM), 2002, ASTM International, 100 Barr Harbor Drive West, Conshohocken, PA 19428. Phone: 610-832-9500. Web: <http://www.astm.org>. Use method number when ordering.

6.2.21 Standard Method 7110 B: Gross Alpha and Gross Beta Radioactivity (Total, Suspended, and Dissolved)

This method should be used for **gross alpha** and **gross beta determination** in aqueous/liquid samples for the contaminants identified below and listed in Appendix B. The method provides an indication of the presence of alpha and beta emitters, including the following SAM analytes:

• Americium-241	(CAS RN 14596-10-2)	Alpha emitter
• Californium-252	(CAS RN 13981-17-4)	Alpha emitter
• Cesium-137	(CAS RN 10045-97-3)	Beta emitter
• Cobalt-60	(CAS RN 10198-40-0)	Beta emitter
• Curium-244	(CAS RN 13981-15-2)	Alpha emitter
• Europium-154	(CAS RN 15585-10-1)	Beta emitter
• Iridium-192	(CAS RN 14694-69-0)	Beta emitter
• Plutonium-238	(CAS RN 13981-16-3)	Alpha emitter
• Plutonium-239	(CAS RN 15117-48-3)	Alpha emitter
• Polonium-210	(CAS RN 13981-52-7)	Alpha emitter
• Radium-226	(CAS RN 13982-63-3)	Alpha emitter
• Ruthenium-103	(CAS RN 13968-53-1)	Beta emitter
• Ruthenium-106	(CAS RN 13967-48-1)	Beta emitter
• Strontium-90	(CAS RN 10098-97-2)	Beta emitter
• Uranium-238	(CAS RN 7440-16-1)	Alpha emitter

This method allows for measurement of gross alpha and gross beta radiation in water samples. This method recommends using a thin-window gas-flow proportional counter for counting gross alpha and beta radioactivity. An internal proportional or Geiger counter may also be used. An aliquot of sample is evaporated to a small volume and transferred to a tared counting pan. The sample residue is dried to constant weight, cooled, and reweighed to determine dry residue weight, then counted for alpha and beta radioactivity.

Source: Standard Methods for Examination of Water and Wastewater, 21st Edition, American Public Health Association (APHA), American Water Works Association (AWWA), and Water Environment Federation (WEF), 2008. Web: <http://www.standardmethods.org/>

6.2.22 Standard Method 7120: Gamma-Emitting Radionuclides

This method should be used for **qualitative** and **confirmatory analysis** of aqueous/liquid samples for the contaminants identified below and listed in Appendix B and with gamma photon energies ranging from approximately 60 to 2000 keV.

Analyte(s)	CAS RN
Cesium-137 ¹	10045-97-3
Cobalt-60	10198-40-0
Europium-154	15585-10-1
Iridium-192	14694-69-0
Ruthenium-103	13968-53-1
Ruthenium-106 ¹	13967-48-1
Selenium-75	14265-71-5

¹ The method will measure decay products of these isotopes

The method uses gamma spectroscopy using either germanium (Ge) detectors or thallium-activated sodium iodide (NaI(Tl)) crystals for the measurement of gamma photons emitted from radionuclides present in water. The method can be used for qualitative and confirmatory determinations with Ge detectors or semi-qualitative and semi-quantitative determinations (using NaI(Tl) detectors). Exact confirmation using NaI is possible for single nuclides or when the gamma emissions are limited to a few well-separated energies. A homogeneous water sample is placed into a standard geometry (normally a Marinelli beaker) for gamma counting. Sample portions are counted long enough to meet the required sensitivity of measurement. A standard containing a mixture of gamma energies from approximately 100 to 2000 keV is used for energy calibration.

Source: *Standard Methods for Examination of Water and Wastewater*, 21st Edition, American Public Health Association (APHA), American Water Works Association (AWWA), and Water Environment Federation (WEF), 2005. Web: <http://www.standardmethods.org/>

6.2.23 Standard Method 7500-Ra B: Radium: Precipitation Method

This method should be used for **qualitative determination** in aqueous/liquid samples for the contaminant identified below and listed in Appendix B.

Analyte(s)	CAS RN
Radium-226 ¹	13982-63-3

¹ The method will measure decay products of these isotopes

This method is for determination of all alpha-emitting radium isotopes by alpha decay analysis. Lead and barium carriers are added to the sample containing alkaline citrate, then sulfuric acid is added to precipitate radium, barium, and lead as sulfates. The precipitate is purified by washing with nitric acid, dissolving in alkaline ethylenediaminetetraacetic acid (EDTA), and re-precipitating as radium-barium sulfate after pH adjustment to 4.5. This slightly acidic EDTA keeps other naturally occurring alpha-emitters and the lead carrier in solution. Radium-223, -224, and -226 are identified by the rate of ingrowth of their daughter products in barium sulfate precipitate. The results are corrected by the rate of ingrowth to determine radium activity. This method involves alpha counting by a gas-flow internal proportional counter, scintillation counter, or thin end-window gas-flow proportional counter.

Source: *Standard Methods for Examination of Water and Wastewater*, 21st Edition, American Public Health Association (APHA), American Water Works Association (AWWA), and Water Environment Federation (WEF), 2005. Web: <http://www.standardmethods.org/>

6.2.24 Standard Method 7500-Ra C: Radium: Emanation Method

This method should be used for **confirmatory analysis** of aqueous/liquid samples for the contaminant identified below and listed in Appendix B.

Analyte(s)	CAS RN
Radium-226 ¹	13982-63-3

¹ The method will measure decay products of these isotopes

This method is for determination of radium-226 by alpha counting. Radium in water is concentrated and separated from sample solids by coprecipitation with a relatively large amount of barium as the sulfate. The precipitate is treated to remove silicates, if present, and to decompose insoluble radium compounds, fumed with phosphoric acid to remove sulfite, and dissolved in hydrochloric acid. The completely dissolved radium is placed in a bubbler, which is then closed and stored for a period of several days to 4 weeks for ingrowth of radon. The bubbler is connected to an evacuation system and the radon gas is removed from the liquid by aeration and helium, dried with a desiccant, and collected in a counting cell. Four hours after radon collection, the cell is counted. The activity of the radon is equal to the radium concentration. The minimum detectable concentration depends on counter characteristics, background-counting rate of scintillation cell, cell efficiency, length of counting period, and contamination of apparatus and environment by radium-226. Without reagent purification, the overall reagent blank (excluding background) should be between 0.03 and 0.05 pCi radium-226, which may be considered the minimum detectable amount under routine conditions.

Source: *Standard Methods for Examination of Water and Wastewater*, 21st Edition, American Public Health Association (APHA), American Water Works Association (AWWA), and Water Environment Federation (WEF), 2005. Web: <http://www.standardmethods.org/>

6.2.25 Standard Method 7500-Sr B: Total Radioactive Strontium and Strontium-90: Precipitation Method

This method should be used for **qualitative** and **confirmatory analysis** of drinking water and aqueous/liquid samples for the contaminant identified below and listed in Appendix B.

Analyte(s)	CAS RN
Strontium-90 ¹	10098-97-2

¹ The method will measure decay products of these isotopes

A known amount of inactive strontium ions, in the form of strontium nitrate, is added as a “carrier.” The carrier, alkaline earths, and rare earths are precipitated as the carbonate to concentrate the radiostrontium. The carrier, along with the radionuclides of strontium, is separated from other radioactive elements and inactive sample solids by precipitation as strontium nitrate using fuming nitric acid solution. The carrier and radionuclides of strontium are precipitated as strontium carbonate, which is dried, weighed to determine recovery of carrier, and measured for radioactivity. The activity of the final precipitate is due to radioactive strontium only, because all other radioactive elements have been removed. Because it is impossible to separate the isotopes of strontium-89 and strontium-90 by any chemical procedure, the amount of strontium-90 is determined by separating and measuring the activity of yttrium-90, its daughter product. This method involves beta counting by a gas-flow internal proportional counter or thin end-window low-background proportional counter. A correction is applied to compensate for loss of carriers and activity during the various purification steps.

Source: *Standard Methods for Examination of Water and Wastewater*, 21st Edition, American Public Health Association (APHA), American Water Works Association (AWWA), and Water Environment Federation (WEF), 2005. Web: <http://www.standardmethods.org/>

6.2.26 Standard Method 7500-U B: Uranium: Radiochemical Method

This method should be used for **qualitative determination** in aqueous/liquid samples for the contaminant identified below and listed in Appendix B.

Analyte(s)	CAS RN
Uranium-238 ¹	7440-61-1

¹ If it is suspected that the sample exists in refractory form (i.e., non-digestible or dissolvable material after normal digestion methods) or if there is a matrix interference problem, use ORISE Method AP11

The sample is acidified with hydrochloric or nitric acid and boiled to eliminate carbonate and bicarbonate ions. Uranium is coprecipitated with ferric hydroxide and subsequently separated. The ferric hydroxide is dissolved, passed through an anion-exchange column, and washed with acid, and the uranium is eluted with dilute hydrochloric acid. The acid eluate is evaporated to near dryness, the residual salt is converted to nitrate, and the alpha activity is counted by a gas-flow proportional counter or alpha scintillation counter.

Source: *Standard Methods for Examination of Water and Wastewater*, 21st Edition, American Public Health Association (APHA), American Water Works Association (AWWA), and Water Environment Federation (WEF), 2005. Web: <http://www.standardmethods.org/>

6.2.27 Standard Method 7500-U C: Uranium: Isotopic Method

This method should be used for **confirmatory analysis** of aqueous/liquid samples for the contaminant identified below and listed in Appendix B.

Analyte(s)	CAS RN
Uranium-238 ¹	7440-61-1

¹ If it is suspected that the sample exists in refractory form (i.e., non-digestible or dissolvable material after normal digestion methods) or if there is a matrix interference problem, use ORISE Method AP11

This method is a radiochemical procedure for determination of the isotopic content of uranium alpha activity; it is consistent with determining the differences among naturally occurring, depleted, and enriched uranium. The sample is acidified with hydrochloric or nitric acid, and uranium-232 is added as an isotopic tracer. Uranium is coprecipitated with ferric hydroxide and subsequently separated. The ferric hydroxide is dissolved, passed through an anion-exchange column, and washed with acid, and the uranium is eluted with dilute hydrochloric acid. The acid eluate is evaporated to near dryness, and the residual salt is converted to nitrate and electrodeposited onto a stainless steel disc and counted by alpha spectrometry.

Source: *Standard Methods for Examination of Water and Wastewater*, 21st Edition, American Public Health Association (APHA), American Water Works Association (AWWA), and Water Environment Federation (WEF), 2005. Web: <http://www.standardmethods.org/>

Section 7.0: Selected Pathogen Methods

A list of the most appropriate methods currently available for use in analyzing environmental samples for pathogens is provided in Appendix C. This list represents an initial effort towards the goal of providing standardized analytical procedures. These methods should be used to support remediation activities (site assessment through clearance) following a homeland security event. The purpose of this section is to provide summary information regarding the procedures listed. Methods are listed for each pathogen that may need to be measured and analyzed following an event. Appendix C is sorted alphabetically within pathogen categories (i.e., bacteria, viruses, protozoa, and helminths).

Protocols from peer-reviewed journal articles are listed for pathogens where standardized methods are not currently available. Future steps include the development and verification of standardized methods. The literature references will be replaced as standardized, verified protocols become available.

Pathogens that are categorized as Biosafety Level 4 (BSL-4), such as hemorrhagic fever viruses and smallpox, will be handled only by reference laboratories with BSL-4 capability and are not included in this document. Pathogens that are categorized as Biosafety Level-3 (BSL-3) select agents will be handled by appropriate Laboratory Response Network (LRN) facilities and are discussed in Section 7.2.1, with the exception of influenza H5N1 virus and *Chlamydophila psittaci*, which are analyzed using SAM procedures. All other pathogens are to be handled as Biosafety Level-2 (BSL-2) using SAM procedures. Pathogens that are considered to be solely of agricultural concern (i.e., animal and plant pathogens) also are not currently included. Such pathogens may be considered for possible inclusion in future document revisions. Some of the pathogens addressed by this document are commonly found in the environment, and the methods listed in Appendix C assume that analyses will be used to evaluate contamination levels that are above background conditions. If possible, an investigation of initial background levels, as well as appropriate controls for background levels, should be performed. It is anticipated that additional site- or event-specific information also will be required to support the analytical results.

Selection of methods from Appendix C should be based on specific data and information needs, including whether there is a need to determine the presence of a pathogen, the viability of a pathogen, or both presence and viability. Although culture-based methods have been selected for many of the pathogens, molecular techniques such as polymerase chain reaction (PCR) will likely be used for viruses because of the difficulty and time required to propagate these agents in host cell cultures. It should be noted that, although molecular techniques are appropriate for evaluation of the presence of a pathogen, these techniques have inherent limitations with regard to the determination of the viability or infectivity of the pathogen. In situations where viability determinations are required (e.g., evaluation of the efficacy of disinfection), viability procedures should be used. Viability procedures are listed for each pathogen where available.

Appendix C includes the following information:

- **Pathogen(s).** A specific causative agent (e.g., viruses, bacteria) of disease.
- **Viability.** Ability to grow and/or develop. If the analytical technique can assess viability, it is denoted by a “yes” in Appendix C. If a technique is not able to assess viability, it is denoted by a “no.”
- **Analytical technique.** An analytical instrument or procedure used to determine the identity, quantity, and/or viability of a pathogen.
- **Analytical method.** The unique identifier or number assigned to an analytical method by the method publisher.
- **Solid (soil, powder).** The recommended method/procedure for the pathogen of interest in solid samples such as soil and powders.

- **Particulate (swabs, wipes, filters).** The recommended method/procedure to measure the pathogen of interest in particulate sample collection devices such as swabs, wipes and high efficiency particulate air (HEPA) filters.
- **Liquid/water (filter, grab).** The recommended method/procedure for the pathogen of interest in liquid or aqueous samples that have been collected using a filter or as a grab sample.
- **Drinking water (filter, grab).** The recommended method/procedure for the pathogen of interest in potable water samples that have been collected using a filter or as a grab sample.
- **Aerosol (growth media, filter, liquid).** The recommended method/procedure for the pathogen of interest in air sample collectors such as growth media, filters, or liquid.

Some of the methods in Appendix C include multiple analytical techniques by inference. The analytical technique listed in Appendix C is intended to be a description of the predominant technique that is required to provide the data quality parameter (viability or detection and identification). This description does not preclude the use of other techniques that are within or referenced by the method. For example, a viability test listed as “culture” may include immunochemical or PCR based assays for the identification of the culture colony isolates.

Several of the methods listed in Appendix C also include options such as the potential for use of multiple cell culture media for primary isolation, allowance for selection of a defined subset of a larger number of biochemical tests for biochemical testing, or use of alternative devices for sample concentration. The method may provide guidance as to which options should be used under particular circumstances, or this may be left to the discretion of the laboratory.

7.1 General Guidance

This section provides a general overview of how to identify the appropriate pathogen method(s) for a given pathogen as well as recommendations for quality control procedures.

For additional information on the properties of the pathogens listed in Appendix C, TOXNET (<http://toxnet.nlm.nih.gov/index.html>), a cluster of databases on toxicology, hazardous chemicals, and related areas maintained by the National Library of Medicine, is an excellent resource. Also informative are CDC’s Emergency Preparedness and Response Web site (<http://www.bt.cdc.gov/>) and the FDA Center for Food Safety and Applied Nutrition (CFSAN) “Bad Bug Book” (<http://www.cfsan.fda.gov/~mow/intro.html>). Further research on pathogens is ongoing within EPA, and databases to manage this information are currently under development.

7.1.1 Standard Operating Procedures for Identifying Pathogen Methods

To determine the appropriate method that is to be used on an environmental sample, locate the pathogen in Appendix C: Pathogen Methods under the “Pathogen(s)” column. After locating the pathogen, continue across the table and select an analytical technique that either assesses viability or does not (i.e., yes or no, respectively). After an analytical technique has been chosen (e.g., culture, polymerase chain reaction (PCR), immunoassay), select the analytical method applicable to the sample type of interest (solid, particulate, liquid/drinking water or aerosol).

Sections 7.2.1 through 7.2.30 below provide summaries of the analytical methods listed in Appendix C. Once a method has been identified in Appendix C, **Table 7-1** can be used to locate the method summary.

Table 7-1. Pathogen Methods and Corresponding Text Section Numbers

Note: Where procedures are listed in this table for both viability assessment and detection, viability procedures are listed first.

Pathogen	Method	Section
Bacteria		
<i>Campylobacter jejuni</i>	SM 9260 G	7.2.12
	Molecular and Cellular Probes. 2006. 20: 269–279	7.2.15
<i>Chlamydomphila psittaci</i>	Journal of Clinical Microbiology. 2000. 38: 1085–1093	7.2.16
<i>Escherichia coli</i> O157:H7	SM 9260 F	7.2.11
<i>Leptospira</i> spp.	SM 9260 I	7.2.14
<i>Listeria monocytogenes</i>	FDA/Bacteriological Analytical Manual Chapter 10, 2003	7.2.7
<i>Salmonella</i> spp. (Method not applicable for <i>Salmonella typhi</i>)	Method 1682	7.2.6
<i>Salmonella typhi</i>	SM 9260 B	7.2.9
<i>Shigella</i> spp.	SM 9260 E	7.2.10
<i>Staphylococcus aureus</i>	SM 9213 B	7.2.8
<i>Vibrio cholerae</i> O1 and O139	SM 9260 H	7.2.13
Viruses		
Adenoviruses: A-F	Applied and Environmental Microbiology. 2005. 71(6): 3131–3136	7.2.17
Astroviruses	Canadian Journal of Microbiology. 2004, 50: 268–278	7.2.18
Caliciviruses: Noroviruses	Journal of Clinical Microbiology. 2004, 42(10): 4679–4658	7.2.19
Caliciviruses: Sapoviruses	Journal of Medical Virology, 2006, 78 (10): 1347–1353	7.2.20
Coronaviruses: SARS	Journal of Virological Methods. 2004. 122: 29–36	7.2.21
Hepatitis E Virus	Journal of Virological Methods. 2006. 131(1): 65–71	7.2.22
Influenza H5N1 Virus	Emerging Infectious Diseases. 2005. 11(8): 1303–1305	7.2.23
Picornaviruses: Enteroviruses	USEPA Manual of Methods for Virology EPA/600/4–84/013, 2001	7.2.2
	Applied and Environmental Microbiology. 2003. 69(6): 3158–3164	7.2.24
Picornaviruses: Hepatitis A Virus	Applied and Environmental Microbiology. 2003. 69(6): 3158–3164	7.2.24
Reoviruses: Rotaviruses	Applied and Environmental Microbiology. 2003. 69(6): 3158–3164	7.2.24
Protozoa		
<i>Cryptosporidium</i> spp.	Method 1622	7.2.4
	Method 1623	7.2.5
	Applied and Environmental Microbiology. 1999. 65(9): 3936–3941	7.2.25

<i>Entamoeba histolytica</i>	Journal of Parasitology. 1972. 58(2): 306–310	7.2.26
	Journal of Clinical Microbiology. 2005. 43(11): 5491–5497	7.2.27
<i>Giardia</i> spp.	Method 1623	7.2.5
	Transactions of the Royal Society of Tropical Medicine and Hygiene. 1983. 77(4): 487–488	7.2.28
<i>Toxoplasma gondii</i>	Emerging Infectious Diseases. 2006. 12(2): 326–329	7.2.29
	Applied and Environmental Microbiology. 2004. 70(7): 4035–4039	7.2.30
Helminths		
<i>Baylisascaris procyonis</i>	EPA/625/R92/013	7.2.3

Method summaries are listed in order of method selection hierarchy (see Figure 2-1), starting with EPA methods followed by methods from other federal agencies, voluntary consensus standard bodies (VCSB), and journal articles. Methods are listed in numerical order under each publisher. For additional information regarding sample preparation and analysis procedures available through consensus standards organizations, please use the contact information provided in **Table 7-2**.

Table 7-2. Sources of Pathogen Methods

Name	Publisher	Reference
National Environmental Methods Index (NEMI)	U.S. EPA, U.S. Geographical Society (USGS)	http://www.nemi.gov
U.S. EPA Microbiology Methods	U.S. EPA	http://www.epa.gov/microbes/
ICR Microbial Laboratory Manual	U.S. EPA Office of Research and Development	http://www.epa.gov/nerlcwww/icrmicro.pdf
USEPA Manual of Methods for Virology	U.S. EPA	http://www.epa.gov/nerlcwww/about.htm
Environmental Regulations and Technology: Control of Pathogens and Vector Attraction in Sewage and Sludge	U.S. EPA, National Risk Management Research Laboratory (NRMRL)	http://www.epa.gov/nrmrl/pubs/625r92013/625r92013.htm
USDA/FSIS Microbiology Laboratory Guidebook	U.S. Department of Agriculture (USDA) Food Safety and Inspection Service	http://www.fsis.usda.gov/Science/Microbiological_Lab_Guidebook/index.asp
Bacteriological Analytical Manual	U.S. Food and Drug Administration (FDA), Center for Food Safety and Applied Nutrition, (CFSAN)	http://www.cfsan.fda.gov/~ebam/bam-toc.html
Occupational Safety and Health Administration Methods	Occupational Health and Safety Administration (OSHA)	http://www.osha.gov
National Institutes for Occupational Safety and Health Methods	National Institute for Occupational Safety and Health (NIOSH)	http://www.cdc.gov/niosh/nmam/

Name	Publisher	Reference
Standard Methods for the Examination of Water and Wastewater, 21 st Edition, 2005*	American Public Health Association (APHA), American Water Works Association (AWWA), and Water Environment Federation (WEF)	http://www.standardmethods.org
Annual Book of ASTM Standards*	ASTM International	http://www.astm.org
Applied and Environmental Microbiology*	American Society for Microbiology	http://aem.asm.org/
Journal of Clinical Microbiology*	American Society for Microbiology	http://jcm.asm.org/
Clinical Microbiology Procedures Handbook, 2 nd Edition, 2004*	American Society for Microbiology	http://estore.asm.org/viewItemDetails.asp?ItemID=323
Molecular and Cellular Probes*	Elsevier	http://www.elsevier.com
Canadian Journal of Microbiology*	NRC Research Press	http://pubs.nrc-cnrc.gc.ca/
Journal of Medical Virology*	Wiley InterScience	http://www3.interscience.wiley.com/cgi-bin/home
Journal of Virological Methods*	Elsevier	http://www.elsevier.com
Emerging Infectious Diseases	Centers for Disease Control	http://www.cdc.gov/ncidod/EID/
Journal of Parasitology*	American Society of Parasitologists	http://www.bioone.org
Transactions of the Royal Society of Tropical Medicine and Hygiene*	The Royal Society of Tropical Medicine and Hygiene	http://www.rstmh.org/
Diagnostic Procedures in Veterinary Bacteriology and Mycology	Academic Press	http://www.pubmedcentral.nih.gov/articlerender.fcgi?artid=1481267

* Subscription and/or purchase required.

7.1.2 General Quality Control (QC) Guidance for Pathogen Methods

Generation of analytical data of known and documented quality is a critical factor in the accurate assessment of and appropriate response to emergency situations. The generation of data of sufficient quality requires that analytical laboratories: (1) have appropriately trained personnel, (2) acquire and maintain required supplies, equipment, and reagents, (3) conduct the appropriate QC procedures to ensure that all measurement systems are in control and operating properly, (4) properly document all analytical results, and (5) properly document analytical QC procedures and corrective actions.

The level or amount of QC needed depends on the intended purpose of the data generated. Various levels of QC may be required if the data are generated for presence/absence determinations versus quantitative results. Specific data needs should be identified, and QC requirements should be based on those needs, and should be applied consistently across laboratories when multiple laboratories are used. The individual methods listed, sampling and analytical protocols, or contractual statements of work should be consulted to determine if additional QC procedures are required.

Method-specific QC requirements are described in many of the methods cited in this manual and will be included in protocols developed to address specific pathogen/sample type combinations of concern. In general, analytical QC requirements for pathogen methods include an initial demonstration of

measurement system capability, as well as the capability of the laboratory and the analyst to perform the method with the required precision and accuracy.

Ongoing analysis of control samples to ensure the continued reliability of the analytical results should also be performed. At a minimum, the following QC analyses should be conducted on an ongoing basis:

- Media and reagent sterility checks;
- Positive and negative controls;
- Method blanks;
- Reference matrix spikes to evaluate initial and ongoing method/analyst performance, if available;
- Matrix spikes to evaluate method performance in the sample type of interest;
- Matrix spike duplicates (MSD) and/or sample replicates to assess method precision; and
- Instrument calibration checks and temperature controls.

QC procedures should be performed as frequently as necessary to ensure the reliability of analytical results.

Please note: The appropriate point of contact identified in Section 4 should be consulted regarding appropriate QA/QC procedures prior to sample analysis. These contacts will consult with the EPA OSWER coordinator responsible for laboratory activities during the specific event to ensure QA/QC procedures are performed consistently across laboratories. OSWER is planning to develop QA/QC guidance for laboratory support. EPA program offices will be responsible for ensuring that the QA/QC practices are implemented.

7.1.3 Safety and Waste Management

It is imperative that safety precautions be used during collection, processing, and analysis of environmental samples. Laboratories should have a documented health and safety plan for handling samples that may contain target chemical, biological, or radiological contaminants, and laboratory staff should be trained in and implement the safety procedures included in the plan. Pathogens in samples taken from areas contaminated as the result of a homeland security event may be more hazardous than naturally occurring pathogens of the same genus and species. The pathogens may have been manufactured or treated in such a manner as to enhance dispersion or virulence characteristics. These conditions may warrant special handling for samples arising from intentional contamination incidents. A laboratory must be made aware of these potential circumstances, and should carefully consider implementing additional safety measures before agreeing to accept these samples.

In addition, many of the methods listed in Appendix C and summarized or cited in Section 7.2 contain specific requirements, guidance, or information regarding safety precautions that should be followed when handling or processing environmental samples and reagents. These methods also provide information regarding waste management. Other resources that can be consulted for additional information include the following:

- Environmental Protection Agency's standards regulating hazardous waste (40 CFR parts 260–270), found at http://ecfr.gpoaccess.gov/cgi/t/text/text-idx?sid=cac9da30cd241fa70d461e4a917eb75e&c=ecfr&tpl=/ecfrbrowse/Title40/40tab_02.tpl
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition, found at <http://www.cdc.gov/OD/ohs/biosfty/bmbl4/bmbl4toc.htm>
- “Laboratory Security and Emergency Response Guidance for Laboratories Working with Select Agents,” *Morbidity and Mortality Weekly Report*, Vol. 51, No. RR-19, 1–6, December 6, 2002, found at <http://www.cdc.gov/mmwr/pdf/rr/rr5119.pdf>.
- Microbiology Biosafety for Level A Laboratories, found at <http://www.bt.cdc.gov/documents/PPTResponse/table3bbiosafety.pdf>

- OSHA Standards for Hazardous Waste Operations and Emergency Response (29 CFR 1910.120) found at http://www.osha.gov/pls/oshaweb/owadisp.show_document?p_table=STANDARDS&p_id=9765
- OSHA Standards for Occupational Exposure to Hazardous Chemicals in Laboratories (29 CFR 1910.1450) found at http://www.osha.gov/pls/oshaweb/owadisp.show_document?p_table=STANDARDS&p_id=10106
- OSHA Standards for Respiratory Protection (29 CFR 1910.134) found at http://www.osha.gov/pls/oshaweb/owadisp.show_document?p_id=12716&p_table=STANDARDS
- DOT Hazardous Materials Shipment and Packaging (49 CFR 171–180) http://ecfr.gpoaccess.gov/cgi/t/text/text-idx?sid=585c275ee19254ba07625d8c92fe925f&c=ecfr&tpl=/ecfrbrowse/Title49/49cfrv2_02.tpl
- Select Agent Rules and Regulations (42 CFR 73 and 9 CFR 121) found at http://www.access.gpo.gov/nara/cfr/waisidx_03/42cfr73_03.html and http://ecfr.gpoaccess.gov/cgi/t/text/text-idx?c=ecfr&tpl=/ecfrbrowse/Title09/9cfr121_main_02.tpl

7.2 Method Summaries

Summaries of the analytical methods listed in Appendix C are provided in Sections 7.2.1 through 7.2.30. Each summary contains a table identifying the pathogen(s) and sample type to which the method applies, a brief description of the method, performance data (if available), and a link to, or source for, obtaining a full version of the method.

7.2.1 Laboratory Response Network (LRN)

The agents identified below and listed in Appendix C are included in the U.S. Department of Health and Human Services/U.S. Department of Agriculture (HHS/USDA) select agent list and should be analyzed in accordance with appropriate LRN protocols.

The LRN was created in accordance with Presidential Directive 39, which established terrorism preparedness responsibilities for federal agencies. The LRN is primarily a national network of local, state, federal, military, food, agricultural, veterinary, and environmental laboratories; however, additional LRN laboratories are located in strategic international locations. The Centers for Disease Control and Prevention (CDC) provides technical and scientific support to member laboratories as well as secure access to standardized procedures and reagents for rapid (within 4–6 hours) presumptive detection of biothreat agents and emerging infectious disease agents. These rapid presumptive assays are part of sample type/analyte specific algorithms of assays which lead to a confirmed result. The algorithm for a confirmed result is often a combination of one or more presumptive positive results from a rapid assay in combination with a positive result from one of the “gold standard” methods, such as culture. The standardized procedures, reagents, and agent-specific algorithms are considered to be sensitive and are available only to LRN member laboratories. Thus, these procedures are not available to the general public and are not discussed in this document.

It is important to note that, in some cases, the procedures may not be fully developed or validated for each environmental sample type/pathogen combination listed in Appendix C, nor are all LRN member laboratories necessarily capable of analyzing all of the sample type/pathogen combinations.

Pathogen(s)	Agent Category
<i>Bacillus anthracis</i> [Anthrax]	Bacteria
<i>Brucella</i> spp. [Brucellosis]	Bacteria

Pathogen(s)	Agent Category
<i>Burkholderia mallei</i> [Glanders]	Bacteria
<i>Burkholderia pseudomallei</i> [Meliodiosis]	Bacteria
<i>Coxiella burnetii</i> [Q-fever]	Bacteria
<i>Francisella tularensis</i> [Tularemia]	Bacteria
<i>Yersinia pestis</i> [Plague]	Bacteria
Orthopoxviruses (Variola, Vaccinia, Cowpox, Monkeypox, Camelpox, Ectromelia, and Gerbilpox)	Viruses

Please note: Not all methods have been verified for the pathogen/sample type combination listed in Appendix C. Please refer to the agent-specific method to identify the pathogen/sample type combinations that have been verified. Any questions regarding information discussed in this section should be referred to the appropriate contact(s) listed in Section 4.

For additional information on the LRN, including selection of a laboratory capable of receiving and processing the specified sample type/pathogen, please use the contact information provided below or visit <http://www.bt.cdc.gov/lrn/>.

Centers for Disease Control and Prevention

Laboratory Response Branch

Division of Bioterrorism Preparedness and Response (DBPR)

National Center for the Prevention, Detection, and Control of Infectious Diseases (NCPDCID)

Coordinating Center for Infectious Diseases (CCID)

Centers for Disease Control and Prevention (CDC)

1600 Clifton Road NE, Mailstop C-18

Atlanta, GA 30333

Telephone: (404) 639-2790

E-mail: lrn@cdc.gov

Local public health laboratories, private laboratories, and commercial laboratories with questions about the LRN should contact their state public health laboratory director or the Association of Public Health Laboratories (contact information provided below).

Association of Public Health Laboratories

8515 Georgia Avenue, Suite 700

Silver Spring, MD 20910 Telephone: (240) 485-2745

Fax: (240) 485-2700

Web site: www.aphl.org

E-mail: info@aphl.org

7.2.2 USEPA Manual of Methods for Virology, EPA/600/4-84/013, April 2001

This method should be used for **detection** and **viability** assessment of picornaviruses (enteroviruses) in solid, particulate, liquid, and water samples. Note: This manual also describes procedures for preparation of samples for non-LRN viruses [i.e., adenovirus, astrovirus, norovirus, sapovirus, coronavirus (SARS), hepatitis E virus, influenza H5N1 virus, picornaviruses (enterovirus and hepatitis A virus), and reovirus (rotavirus)].

Pathogen(s)	Agent Category
Picornaviruses: Enteroviruses	Viruses

This manual describes procedures for determining the infectivity of enteroviruses, including a neutralization test used to identify these viruses. The test uses reference-typing sera directed against isolated waterborne viruses, and consists of simultaneously inoculating virus and antiserum into a microtiter plate, incubating the virus-antibody mixture for two hours, adding a suspension of host cells to the mixture, incubating the host cells-virus-antibody mixture for three days, and then examining the cells daily for five more days for the absence or presence of cytopathic effect (CPE). The test uses the Lim Benyesh-Melnick (LB-M) antiserum pools, which consist of 61 equine antisera, including LB-M antiserum pools A–H for the identification of 41 enteroviruses. Chapters 7 and 14 in this manual describe procedures for the collection and preparation of virus samples. Sample preparation procedures described include concentration and processing of waterborne viruses by positively charged 1MDS cartridge filters and flocculation. These general procedures can be used for many viruses and may be adapted for analysis of particulate, liquid, water and aerosol samples. (See Footnote 4 of Appendix C.)

Please note: These procedures have been applied to analysis of water samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank.

Source: USEPA Manual of Methods for Virology EPA/600/4-84/013. April 2001. *Identification of Enteroviruses*. Chapter 12. (<http://www.epa.gov/microbes/chapt12.pdf>). For information regarding this manual, please contact Ann Grimm (see Section 4 for contact information).

7.2.3 USEPA Environmental Regulations and Technology, Control of Pathogens and Vector Attraction in Sewage Sludge EPA/625/R-92/013, July 2003: *Baylisascaris procyonis*

These protocols should be used for **detection** and **viability** assessment of *Baylisascaris procyonis* in solid, particulate, liquid, and water samples.

Pathogen(s)	Agent Category
<i>Baylisascaris procyonis</i>	Helminths

The protocol describes procedures for analysis of solid and wastewater samples and may be adapted to analysis of particulate, liquid, water, and aerosol samples. Samples are processed by blending with buffered water containing a surfactant. The blend is screened to remove large particles, the solids in the screened portion are allowed to settle out, and the supernatant is decanted. The sediment is subjected to density gradient centrifugation using magnesium sulfate. This flotation procedure yields a layer likely to contain *Ascaris* and other parasite ova, if present in the sample. Small particulates are removed by a second screening on a small mesh size screen. The resulting concentrate is incubated until control *Ascaris* eggs are fully embryonated. The concentrate is then microscopically examined for the categories of *Ascaris* ova on a counting chamber.

Please note: This protocol has been applied to analysis of wastewater, sludge, and compost samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank.

Source: USEPA Environmental Regulations and Technology: Control of Pathogens and Vector Attractions in Sewage Sludge, “Appendix I: Test Method for Detecting, Enumerating, and Determining the Viability of Ascaris Ova in Sludge.” EPA/625/R-92/013. July 2003.
(<http://www.epa.gov/nrmrl/pubs/625r92013/625r92013.htm>).

7.2.4 EPA Method 1622: *Cryptosporidium* in Water by Filtration/IMS/FA

This method should be used for **detection** of *Cryptosporidium* spp. in solid, particulate, and drinking water samples, but is not suitable for determining viability.

Pathogen(s)	Agent Category
<i>Cryptosporidium</i> spp. [Cryptosporidiosis]	Protozoa

This method describes procedures for analysis of drinking water samples and may be adapted for analysis of solid and particulate samples. A water sample is filtered and the oocysts and extraneous materials are retained on the filter. Materials on the filter are eluted, the eluate is centrifuged to pellet the oocysts, and the supernatant fluid is aspirated. The oocysts are magnetized by attachment of magnetic beads conjugated to anti-*Cryptosporidium* antibodies. The magnetized oocysts are separated from the extraneous materials using a magnet, and the extraneous materials are discarded. The magnetic bead complex is then detached from the oocysts. The oocysts are stained on well slides with fluorescently labeled monoclonal antibodies and 4',6-diamidino-2-phenylindole (DAPI). The stained sample is examined using fluorescence and differential interference contrast (DIC) microscopy. Qualitative analysis is performed by scanning each slide well for objects that meet the size, shape, and fluorescence characteristics of *Cryptosporidium* oocysts. Quantitative analysis is performed by counting the total number of objects on the slide confirmed as oocysts. This method is not intended to determine viability of the oocysts.

Please note: This method has been applied to analysis of drinking water samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, matrix spike and matrix spike duplicate (MS/MSD), and blank.

*Source: USEPA. 2005. *Cryptosporidium* in Water by Filtration/IMS/FA. United States Environmental Protection Agency, Washington, D.C.* (<http://www.epa.gov/nerlcwww/1622de05.pdf>)

7.2.5 EPA Method 1623: *Cryptosporidium* and *Giardia* in Water by Filtration/IMS/FA

This method should be used for the **detection** of *Cryptosporidium* spp. in solid, particulate, and liquid/water samples and *Giardia* spp. in solid, particulate, and liquid/water samples, but is not suitable for determining viability.

Pathogen(s)	Agent Category
<i>Cryptosporidium</i> spp.	Protozoa
<i>Giardia</i> spp.	Protozoa

This method describes procedures for analysis of wastewater samples and may be adapted for assessment of solid, particulate, and liquid samples. A water sample is filtered and the oocysts and cysts and extraneous materials are retained on the filter. Materials on the filter are eluted, the eluate is centrifuged

to pellet the oocysts and cysts, and the supernatant fluid is aspirated. The oocysts and cysts are magnetized by attachment of magnetic beads conjugated to anti-*Cryptosporidium* and anti-*Giardia* antibodies. The magnetized oocysts and cysts are separated from the extraneous materials using a magnet, and the extraneous materials are discarded. The magnetic bead complex is then detached from the oocysts and cysts. The oocysts and cysts are stained on well slides with fluorescently labeled monoclonal antibodies and 4',6-diamidino-2-phenylindole (DAPI). The stained sample is examined using fluorescence and differential interference contrast (DIC) microscopy. Qualitative analysis is performed by scanning each slide well for objects that meet the size, shape, and fluorescence characteristics of *Cryptosporidium* oocysts and *Giardia* cysts. Quantitative analysis is performed by counting the total number of objects on the slide confirmed as oocysts or cysts. This method is not intended to determine viability of the parasites.

Please note: This method has been applied to analysis of wastewater samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, matrix spike and matrix spike duplicate (MS/MSD), and blank.

Source: USEPA. 2001. *Cryptosporidium* and *Giardia* in Water by Filtration/IMS/FA. United States Environmental Protection Agency, Washington, D.C. (<http://www.epa.gov/microbes/1623de05.pdf>)

7.2.6 EPA Method 1682: *Salmonella* spp.

This method should be used for **detection** and **viability** assessment of *Samonella* spp. (except *S. typhi*) in solid, particulate, and liquid/water samples.

Pathogen(s)	Agent Category
<i>Salmonella</i> spp. (Method not applicable to <i>Salmonella typhi</i>)	Bacteria

This method describes procedures for analysis of solid samples (biosolids) and may be adapted for assessment of water, liquid, particulate, and aerosol samples. Prepared samples are inoculated into tubes of tryptic soy broth (TSB) and incubated for 24 hours. Positive (turbid) tubes are spotted onto plates of modified semisolid Rappaport-Vassiliadis (MSRV) medium and incubated at 42°C for 16 to 18 hours. The MSRV medium uses novobiocin and malachite green to inhibit non-*Salmonella* species, while allowing most *Salmonella* species to grow. MSRV and the elevated incubation temperature (42°C) are inhibitory for *S. typhi*. Presumptive colonies are isolated on xylose-lysine desoxycholate (XLD) agar and confirmed using lysine-iron (LIA), triple sugar iron (TSI) agar, and urea broth, followed by serological typing using polyvalent O antisera.

Please note: This procedure has been applied to analysis of solid (biosolids) samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank.

Source: USEPA. 2006. *Salmonella* in Sewage Sludge (Biosolids) by Modified Semisolid Rappaport-Vassiliadis (MSRV) Medium. United States Environmental Protection Agency, Washington, D.C. (<http://www.epa.gov/waterscience/methods/method/biological/1682.pdf>)

7.2.7 U.S. FDA Bacteriological Analytical Manual, Chapter 10, 2003: *Listeria monocytogenes*

This method should be used for **detection** and **viability** assessment of *Listeria monocytogenes* in solid, particulate, aerosol, liquid and water samples.

Pathogen(s)	Agent Category
<i>Listeria monocytogenes</i>	Bacteria

Procedures are described for analysis of food samples and may be adapted for assessment of solid, particulate, aerosol, liquid, and water samples. Prepared samples are incubated for 4 hours in Buffered *Listeria* Enrichment Broth (BLEB) without selective agents. Cycloheximide is added and incubation continued. At 24 and 48 hours, BLEB cultures are streaked onto esculin-containing selective isolation agar (i.e., Oxford Medium, OXA) and incubated for an additional 24 to 48 hours. A *L. monocytogenes*-*L. ivanovii* differential selective agar, such as Biosynth Chromogen Medium (BCM), is streaked at 48 hours. Presumptive *Listeria* colonies are black with a black halo on esculin-containing media and blue on BCM plates. Isolated colonies are streaked onto Trypticase soy agar with yeast extract (TSAye), incubated for 24 to 48 hours, and examined for morphological and biochemical characteristics. *L. monocytogenes* is a rod-shaped Gram-positive, motile bacterium. It is catalase positive, ramnose-positive, and mannitol- and xylose-negative. Purified isolates may be rapidly identified using commercially available biochemical typing kits. Confirmation is performed with commercially available sera.

Please note: This method has been applied to analysis of food samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank.

Source: U.S. Food & Drug Administration, Center for Food Safety & Applied Nutrition. 2001. *Bacteriological Analytical Manual Online*. Chapter 10. (<http://www.cfsan.fda.gov/~ebam/bam-16.html>)

7.2.8 Standard Methods 9213 B: *Staphylococcus aureus*

This method should be used for **detection** and **viability** assessment of *Staphylococcus aureus* in solid, particulate, aerosol, liquid, and water samples.

Pathogen(s)	Agent Category
<i>Staphylococcus aureus</i>	Bacteria

Procedures are described for analysis of water samples and may be adapted for assessment of solid, liquid, particulate, and aerosol samples. Prepared samples are inoculated into tubes of M-staphylococcus broth and incubated for 24 hours. Positive (turbid) tubes are streaked onto plates of Baird-Parker agar and incubated for 48 hours. Presumptive *S. aureus* colonies are tested for mannitol fermentation by the addition of a drop of bromthymol blue, a pH indicator. Isolated colonies are examined for morphological and biochemical characteristics. *S. aureus* is a Gram-positive coccus. Biochemical characterizations include catalase-positive, coagulase-positive, fermentation of mannitol, and anaerobic fermentation of glucose.

Please note: This method has been applied to analysis of water samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank.

Source: American Public Health Association, American Water Works Association, and Water Environment Federation. 2005. *Standard Methods for the Examination of Water and Wastewater*. 21st Edition. (<http://www.standardmethods.org/>)

7.2.9 Standard Methods 9260 B: General Qualitative Isolation and Identification Procedures for *Salmonella*

This method should be used for **detection** and **viability** assessment of *Salmonella typhi* in solid, particulate, liquid, and water samples.

Pathogen(s)	Agent Category
<i>Salmonella typhi</i> [Typhoid fever]	Bacteria

This method describes procedures for analysis of water samples and may be adapted for assessment of solid, particulate, and liquid samples. Concentrated samples are enriched in either selenite cystine, selenite-F, or tetrathionate broths and incubated at 35°C to 37°C for up to 5 days. An aliquot from each turbid tube is streaked onto bismuth sulfite (BS) plates and incubated at 35°C to 37°C for 24 to 48 hours. Presumptive positive colonies are then subjected to biochemical characterization. Confirmation is through serological testing using polyvalent O and Vi antiserum.

Please note: This method has been applied to analysis of water samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank.

Source: American Public Health Association, American Water Works Association, and Water Environment Federation. 2005. *Standard Methods for the Examination of Water and Wastewater*. 21st Edition. (<http://www.standardmethods.org/>)

7.2.10 Standard Methods 9260 E: *Shigella*

This method should be used for **detection** and **viability** assessment of *Shigella* species in solid, particulate, liquid, and water samples.

Pathogen(s)	Agent Category
<i>Shigella</i> spp. [Shigellosis]	Bacteria

This method describes procedures for analysis of water samples and may be adapted for assessment of solid, particulate, and liquid samples. This method contains two options for sample concentration: membrane filtration (liquid samples) and centrifugation (liquid and solid samples) for analyses. Both options include inoculation of an enrichment medium (Selenite F broth). Isolation of the target pathogen is achieved by plating onto Xylose Lysine Deoxycholate (XLD) and/or MacConkey agar. Biochemical identification consists of inoculating Triple Sugar Iron (TSI) and Lysine Iron Agar (LIA) slants. Confirmation is performed by slide agglutination tests using polyvalent antisera.

Please note: This method has been applied to analysis of water and solid samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank.

Source: American Public Health Association, American Water Works Association, and Water Environment Federation. 2005. *Standard Methods for the Examination of Water and Wastewater*. 21st Edition. (<http://www.standardmethods.org/>)

7.2.11 Standard Methods 9260 F: Pathogenic *Escherichia coli*

This method should be used for **detection** and **viability** assessment of *Escherichia coli* O157:H7 in solid, particulate, liquid, and water samples.

Pathogen(s)	Agent Category
<i>Escherichia coli</i> (<i>E. coli</i>) O157:H7	Bacteria

This method describes procedures for analysis of water samples and may be adapted for assessment of solid, particulate, and liquid samples. The method allows for two options, one being a modification of SM 9221B followed by plating and biochemical identification. The second option, modification of a food method, allows for the analysis of large sample volumes. A 200-mL water sample is centrifuged, resuspended in *E. coli* enrichment broth (EEB) and incubated for 6 hours. Tellurite Cefixime Sorbitol-MacConkey (TC SMAC) plates are inoculated with the enriched EEB culture, and incubated for up to 24 hours. Colorless colonies on TC SMAC are tested for indole production. Presumptive positive colonies are then subjected to biochemical characterization. Confirmation is through serological testing.

Please note: This method has been applied to analysis of water samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank.

Source: American Public Health Association, American Water Works Association, and Water Environment Federation. 2005. *Standard Methods for the Examination of Water and Wastewater*. 21st Edition. (<http://www.standardmethods.org/>)

7.2.12 Standard Methods 9260 G: *Campylobacter jejuni*

This method should be used for **detection** and **viability** assessment of *Campylobacter jejuni* in solid, particulate, liquid, and water samples.

Pathogen(s)	Agent Category
<i>Campylobacter jejuni</i>	Bacteria

This method describes procedures for analysis of water samples and may be adapted for assessment of solid, particulate, and liquid samples. Water samples (1 to several liter volumes) are filtered using a cellulose nitrate membrane filter. Filters are placed face down on Skirrow's medium and incubated for 24 hours at 42°C in a microaerobic environment. Alternatively, samples are enriched in *Campylobacter* broth supplemented with antibiotics and lysed horse blood in a microaerobic environment at 37°C for 4 hours, then at 42°C for 24 to 48 hours prior to streaking on Skirrow's medium. Identification is made by culture examination, microscopy, motility testing, and biochemical testing. Biochemical tests include oxidase, catalase, nitrite and nitrate reduction, H₂S production, and hippurate hydrolysis. Confirmation is performed using commercially available rapid serological test kits. Skirrow's and other selective media containing antibiotics (trimethoprim, vancomycin, polymyxin) may prevent the growth of injured organisms.

Please note: This method has been applied to analysis of water samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank.

Source: American Public Health Association, American Water Works Association, and Water Environment Federation. 2005. *Standard Methods for the Examination of Water and Wastewater*. 21st Edition. (<http://www.standardmethods.org/>)

7.2.13 Standard Methods 9260 H: *Vibrio cholerae*

This method should be used for **detection** and **viability** assessment of *Vibrio cholerae* in solid, particulate, liquid, and water samples.

Pathogen(s)	Agent Category
<i>Vibrio cholerae</i> [Cholera]	Bacteria

This method describes procedures for analysis of water samples and may be adapted for assessment of solid, particulate, and liquid samples. Samples are enriched in alkaline peptone broth and incubated for up to 8 hours. Thiosulfate-citrate-bile salts-sucrose (TCBS) agar plates are inoculated with the incubated broth and incubated for 24 hours. Yellow sucrose-fermenting colonies are presumptive for *V. cholera* and are plated on tryptic soy agar with 0.5% NaCl. Presumptive positive colonies are then subjected to biochemical characterization. Confirmation is performed using slide agglutination assays for serological identification.

Please note: This method has been applied to analysis of water samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank.

Source: American Public Health Association, American Water Works Association, and Water Environment Federation. 2005. *Standard Methods for the Examination of Water and Wastewater*. 21st Edition. (<http://www.standardmethods.org/>)

7.2.14 Standard Methods 9260 I: *Leptospira*

This method should be used for **detection** and **viability** assessment of *Leptospira interrogans* in solid, particulate, liquid, and water samples.

Pathogen(s)	Agent Category
<i>Leptospira interrogans</i> [Leptospirosis]	Bacteria

This method describes procedures for analysis of water samples and may be adapted for assessment of solid, particulate, and liquid samples. Filter samples through a 0.22- μ m filter, retaining the filtrate as inoculum. If the sample is turbid, a succession of filters of decreasing pore size may be used prior to the 0.22- μ m filter. A tube of *Leptospira* Medium Base (Ellinghausen-McCullough Johnson Harris formulation (EMJH), supplemented with bovine serum albumin (BSA) and Tween®), is inoculated with the sample and incubated at 30°C for up to 6 weeks. Alternatively, the sample is added directly to the media, incubated overnight, passed through a 0.22- μ m membrane filter, and incubation continued for up to 6 weeks. Cultures are examined by darkfield microscopy for motile leptospires. Confirmation is performed by microscopic agglutination test using reference antisera.

Please note: This method has been applied to analysis of water samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank.

Source: American Public Health Association, American Water Works Association, and Water Environment Federation. 2005. *Standard Methods for the Examination of Water and Wastewater*. 21st Edition. (<http://www.standardmethods.org/>)

7.2.15 Literature Reference for *Campylobacter jejuni* (Molecular and Cellular Probes. 2006. 20: 269–279)

These procedures should be used for **detection** of *Campylobacter jejuni* in solid, particulate, aerosol, liquid, and water samples, but are not suitable for determining viability.

Pathogen(s)	Agent Category
<i>Campylobacter jejuni</i>	Bacteria

Procedures are described for analysis of clinical samples and may be adapted for assessment of solid, particulate, aerosol, liquid, and water samples. The method uses real-time quantitative polymerase chain reaction (PCR) for identification of *C. jejuni*. A high through-put method using guanidinium thiocyanate (GTC) and glass beads is used for extraction of deoxyribonucleic acid (DNA). *C. jejuni* uses the CJ F primer, CJ R primer, and CJ probe for real-time quantitative PCR. Results are evaluated against standard curves made with a 10-fold dilution series of *C. jejuni* NCTC 11168 DNA in a background of cecum DNA.

Please note: These procedures have been applied to analysis of clinical samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank. PCR quality control checks should be performed according to *EPA Draft Quality Assurance/Quality Control Guidance for Laboratories Performing PCR Analyses on Environmental Samples* document (www.epa.gov/nerlcwww/qa_qc_pcr10_04.pdf) or consult the point of contact identified in Section 4.

Source: Skanseng, B., Kaldhusdal, M., and Rudi, K. 2006. “Comparison of Chicken Gut Colonization by the Pathogens *Campylobacter jejuni* and *Clostridium perfringens* by the Real-time Quantitative PCR.” *Molecular and Cellular Probes*. 20: 269–279.

7.2.16 Literature Reference for *Chlamydomphila psittaci* (Journal of Clinical Microbiology. 2000. 38: 1085–1093)

These procedures should be used for **detection** and **viability** assessment of *Chlamydomphila psittaci* in solid, particulate, aerosol, liquid and water samples.

Pathogen(s)	Agent Category
<i>Chlamydomphila psittaci</i> (formerly known as <i>Chlamydia psittaci</i>)	Bacteria

Procedures are described for analysis of clinical samples and may be adapted for assessment of solid, particulate, aerosol, liquid, and water samples. The method uses touchdown enzyme time release-polymerase chain reaction (TETR-PCR) for detection and identification of *Chlamydomphila psittaci*. Deoxyribonucleic acid (DNA) is extracted from chlamydia cultures (Buffalo Green Monkey Kidney (BGMK) cells), minced clinical tissues, and respiratory samples by mixing with a chelating resin, heating first to 56°C for 15 to 30 minutes, then 100°C for 8 to 10 minutes. Primer sets specific for *C. psittaci* are designed based on the DNA sequences of the 16S ribosomal ribonucleic acid (rRNA) and 16S-23S rRNA genes (CPS 100/101). PCR products are separated by electrophoresis in 12% polyacrylamide gels with Tris-borate–ethylenediaminetetraacetic acid (EDTA) buffer and visualized with ethidium bromide. *Chlamydomphila psittaci* samples are to be assessed at BioSafety Level-3 (BSL-3) laboratories.

Please note: These procedures have been applied to analysis of clinical (vaginal and respiratory) samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank. PCR quality control checks should be performed according to *EPA Draft Quality Assurance/Quality Control Guidance for Laboratories Performing PCR Analyses on Environmental Samples* document (www.epa.gov/nerlcwww/qa_qc_pcr10_04.pdf) or consult the point of contact identified in Section 4.

Source: Madico, G., Quinn, T.C., Boman, J., and Gaydos, C.A. 2000. “Touchdown Enzyme Time Release-PCR for Detection and Identification of *Chlamydia trachomatis*, *C. pneumoniae*, and *C. psittaci* Using the 16S and 16S-23S Spacer rRNA Genes.” *Journal of Clinical Microbiology*. 38(3): 1085–1093. <http://www.pubmedcentral.nih.gov/picrender.fcgi?artid=86346&blobtype=pdf>

7.2.17 Literature Reference for Adenoviruses (Applied and Environmental Microbiology. 2005. 71(6): 3131–3136)

These procedures should be used for **detection** and **viability** assessment of Adenoviruses, enteric and non-enteric (A–F) in solid, particulate, aerosol, liquid, and water samples. Note: Procedures described in the USEPA Manual of Methods for Virology (EPA/600/4-84/013, April 2001) can be adapted for the preparation of samples.

Pathogen(s)	Agent Category
Adenoviruses: Subgroups A–F	Viruses

Procedures are described for analysis of cell culture lysates and may be adapted for assessment of solid, particulate, aerosol, liquid, and water samples. The detection procedure uses a broadly reactive fluorogenic 5' nuclease (TaqMan®) quantitative real-time polymerase chain reaction (PCR) assay for the detection of all six species (A–F) of human adenoviruses (HAdV). Sensitive detection and discrimination of adenovirus F species (AdV40 and AdV41) can be achieved by using a real-time fluorescence resonance energy transfer (FRET)-based PCR assay.

For the viability assessment of adenovirus 40 and 41, given that they can be difficult to grow in culture, cell lines such as G293 (Journal of Medical Virology, 11(3): 215–231) or CaCo-2 (Journal of Medical Virology. 1994. 44(3): 310–315) may be considered when these viruses are suspected to be present. As detection of adenovirus in environmental samples can be difficult, additional methods such as described in *Effect of Adenovirus Resistance on UV Disinfection Experiments: A Report on the State of Adenovirus Science* (J. AWWA. 2006. 98(6):93–106) also may be useful.

Please note: These procedures have been applied to analysis of cell culture lysates. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank. PCR quality control checks should be performed according to *EPA Draft Quality Assurance/Quality Control Guidance for Laboratories Performing PCR Analyses on Environmental Samples* document (www.epa.gov/nerlcwww/qa_qc_pcr10_04.pdf) or consult the point of contact identified in Section 4.

Source: Jothikumar, N., Cromeans, T. L., Hill, V. R., Lu, X., Sobsey, M., and Erdman, D. D. 2005. “Quantitative Real-Time PCR Assays for Detection of Human Adenoviruses and Identification of Serotypes 40 and 41.” *Applied and Environmental Microbiology*. 71 (6): 3131–3136. <http://aem.asm.org/cgi/reprint/71/6/3131>

7.2.18 Literature Reference for Astroviruses (Canadian Journal of Microbiology. 2004. 50: 269–278)

These procedures should be used for **detection** and **viability** assessment of Astroviruses in solid, particulate, aerosol, liquid, and water samples. Note: Procedures described in the USEPA Manual of Methods for Virology (EPA/600/4-84/013, April 2001) can be adapted for the preparation of samples.

Pathogen(s)	Agent Category
Astroviruses	Viruses

Procedures are described for analysis of clinical samples and may be adapted for assessment of solid, particulate, aerosol, liquid, and water samples. The method detects eight astrovirus serotypes. The method is a reverse transcription-polymerase chain reaction (RT-PCR) procedure optimized for use in a real-time PCR assay and can be integrated with sample-cell culture (CaCo-2 cells) to enhance sensitivity (ICC/RT-PCR). Water samples are collected by filtration (1MDS filter), and viruses are eluted using a beef extract solution (1.5%, pH 9.5). Viruses are concentrated using celite adsorption (pH 4.0), filtration, and celite-elution with sodium phosphate (0.15 M, pH 9.0), followed by further concentration and processing to remove inhibitors (ultracentrifugation, solvent extraction, and MW-exclusion filtration). Concentrated samples are analyzed directly or indirectly (following cell culture) by a two-step RT-PCR (RT followed by PCR) assay using astrovirus-specific primer sets. Detection of amplicons is by gel electrophoresis with subsequent confirmation by hybridization (dot-blot) using digoxigenin-labeled internal (nested) probes or by real-time detection using fluorogenic probes.

Please note: These procedures have been applied to analysis of clinical (fecal) samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank. PCR quality control checks should be performed according to *EPA Draft Quality Assurance/Quality Control Guidance for Laboratories Performing PCR Analyses on Environmental Samples* document (www.epa.gov/nerlcwww/qa_qc_pcr10_04.pdf) or consult the point of contact identified in Section 4.

Source: Grimm, A. C., Cashdollar, J. L., Williams, F. P., and Fout, G. S. 2004. “Development of an Astrovirus RT-PCR Detection Assay for Use with Conventional, Real-Time, and Integrated Cell Culture/RT-PCR.” *Canadian Journal of Microbiology*. 50(4): 269–268.

7.2.19 Literature Reference for Noroviruses (Journal of Clinical Microbiology. 2004. 42(10): 4679–4685)

These procedures should be used for the **detection** assessment of noroviruses in solid, particulate, aerosol, liquid, and water samples, but are not suitable for determining viability. Note: Procedures described in the USEPA Manual of Methods for Virology (EPA/600/4-84/013, April 2001) can be adapted for the preparation of samples.

Pathogen(s)	Agent Category
Caliciviruses: Noroviruses	Viruses

Procedures are described for analysis of clinical samples and may be adapted for assessment of solid, particulate, aerosol, liquid, and water samples. This method is an assay for the detection and quantitation of norovirus using LightCycler® real-time reverse transcription-polymerase chain reaction (real-time LC RT-PCR) technology. Viral ribonucleic acid (RNA) is extracted using either a commercial kit or a silica-based method. For Norovirus G-I, primers based on the capsid gene sequence are used, and for Norovirus G-II, primers based on the polymerase gene sequence are used. A SYBR® Green I system is used in the reaction for visualization. External standard curves for the quantification of norovirus are established using RNA transcripts from strains S5 and S19, corresponding to G-I/4 and G-II/12, respectively.

Please note: These procedures have been applied to analysis of clinical (fecal) samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank. PCR quality control checks should be performed according to *EPA Draft Quality Assurance/Quality Control Guidance for Laboratories Performing PCR Analyses on Environmental Samples* document (www.epa.gov/nerlcwww/qa_qc_pcr10_04.pdf) or consult the point of contact identified in Section 4.

Source: Pang, X., Lee, B., Chui, L., Preiksaitis, J.K., and Monroe, S.S. 2004. "Evaluation and Validation of Real-Time Reverse Transcription-PCR Assay Using the LightCycler System for Detection and Quantitation of Norovirus." *Journal of Clinical Microbiology*. 42(10): 4679–4685.
<http://jcm.asm.org/cgi/reprint/42/10/4679>

7.2.20 Literature Reference for Sapoviruses (Journal of Medical Virology. 2006. 78(10): 1347–1353)

These procedures should be used for **detection** of sapoviruses in solid, particulate, aerosol, liquid, and water samples, but are not suitable for determining viability. **Note:** Procedures described in the USEPA Manual of Methods for Virology (EPA/600/4-84/013, April 2001) can be adapted for the preparation of samples.

Pathogen(s)	Agent Category
Caliciviruses: Sapoviruses	Viruses

Procedures are described for analysis of clinical samples and may be adapted for assessment of solid, particulate, aerosol, liquid, and water samples. The method is a TaqMan®-based real-time reverse transcriptase polymerase chain reaction (RT-PCR) assay that has the ability to detect four of the five distinct sapovirus (SaV) genogroups (GI–GV). Sets of primers, based on the multiple alignment of 27 gene sequences for the polymerase-capsid junction in open reading frame 1 (ORF1), are used to detect human SaV GI, GII, GIV, and GV sequences in a single tube. Sensitivity using control plasmids range from 2.5×10^1 to 2.5×10^7 copies per tube. No cross-reactivity is observed against other enteric viruses, including norovirus (NoV), rotavirus, astrovirus, and adenovirus.

Please note: These procedures have been applied to analysis of clinical (fecal) samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank. PCR quality control checks should be performed according to *EPA Draft Quality Assurance/Quality Control Guidance for Laboratories Performing PCR Analyses on Environmental Samples* document (www.epa.gov/nerlcwww/qa_qc_pcr10_04.pdf) or consult the point of contact identified in Section 4.

Source: Oka T, Katayama K, Hansman GS, Kageyama T, Ogawa S, Wu FT, White PA, and Takeda N. 2006. "Detection of Human Sapovirus by Real-time Reverse Transcription-Polymerase Chain Reaction." *Journal of Medical Virology*. 78(10): 1347–1353.

7.2.21 Literature Reference for Coronaviruses (SARS) (Journal of Virological Methods. 2004. 122: 29–36)

These procedures should be used for **detection** of severe acute respiratory syndrome-associated human coronavirus (SARS-HCoV) in solid, particulate, aerosol, liquid, and water samples, but are not suitable for determining viability. Note: Procedures described in the USEPA Manual of Methods for Virology (EPA/600/4-84/013, April 2001) can be adapted for the preparation of samples.

Pathogen(s)	Agent Category
Coronaviruses: SARS-associated human coronavirus	Viruses

Procedures are described for analysis of clinical samples and may be adapted for assessment of solid, particulate, aerosol, liquid, and water samples. This method uses a conventional single-tube reverse transcription-polymerase chain reaction (RT-PCR) procedure based on consensus primer sequences targeting conserved regions of coronavirus genome sequences. End-point amplicon analysis is by electrophoresis and subsequent visualization. The assay can detect the SARS-HCoV as well as several other human respiratory coronaviruses (HCoV-OC43 and HCoV-229E). Species identification is provided by sequencing the amplicon, although rapid screening can be performed by restriction enzyme analysis.

Please note: These procedures have been applied to analysis of clinical (throat, lung, blood, stool, urine) samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank. PCR quality control checks should be performed according to EPA Draft Quality Assurance/Quality Control Guidance for Laboratories Performing PCR Analyses on Environmental Samples document (www.epa.gov/nerlcwww/qa_qc_pcr10_04.pdf) or consult the point of contact identified in Section 4.

Source: Adachi, D., Johnson, G., Draker, R., Ayers, M., Mazzulli, T., Talbot, P. J., and Tellier, R. 2004. “Comprehensive Detection and Identification of Human Coronaviruses, Including the SARS-associated Coronavirus, with a Single RT-PCR Assay.” *Journal of Virological Methods*. 122: 29–36.

7.2.22 Literature Reference for Hepatitis E Virus (Journal of Virological Methods. 2006. 131(1): 65–71)

These procedures should be used for the **detection** of Hepatitis E virus in solid, particulate, aerosol, liquid, and water samples, but are not suitable for determining viability. Note: Procedures described in the USEPA Manual of Methods for Virology (EPA/600/4-84/013, April 2001) can be adapted for the preparation of samples.

Pathogen(s)	Agent Category
Hepatitis E virus (HEV)	Viruses

Procedures are described for analysis of spiked water samples and may be adapted for assessment of solid, particulate, aerosol, and liquid samples. The method uses a TaqMan® real-time reverse transcriptase-polymerase chain reaction (RT-PCR) assay in order to detect and quantitate all four major HEV strains that may be present in clinical and environmental samples. Primers and probes are based on the multiple sequence alignments of 27 gene sequences for the ORF3 region. Thirteen HEV isolates representing genotypes 1–4 are used to standardize the real-time RT-PCR assay. The assay can detect as few as four genome equivalent (GE) copies of HEV plasmid DNA and 1.2 50% pig infectious dose (PID50) of swine HEV. Concentrations of swine HEV from 1.2 to 120 PID50 spiked into a surface water sample can also be detected.

Please note: These procedures have been applied to analysis of spiked water. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative

control, and blank. PCR quality control checks should be performed according to *EPA Draft Quality Assurance/Quality Control Guidance for Laboratories Performing PCR Analyses on Environmental Samples* document (www.epa.gov/nerlcwww/qa_qc_pcr10_04.pdf) or consult the point of contact identified in Section 4.

Source: Jothikumar, N., Cromeans, T.L., Robertson, B.H., Meng, X.J., and Hill, V.R. 2006. “A Broadly Reactive One-step Real-time RT-PCR Assay for Rapid and Sensitive Detection of Hepatitis E Virus.” *Journal of Virological Methods*, Vol. 131(1): 65–71.

7.2.23 Literature Reference for Influenza H5N1 (Emerging Infectious Diseases. 2005. 11(8): 1303–1305)

These procedures should be used for the **detection** of Influenza H5N1 virus in solid, particulate, aerosol, liquid, and water samples, but are not suitable for determining viability. **Note:** Procedures described in the USEPA Manual of Methods for Virology (EPA/600/4-84/013, April 2001) can be adapted for the preparation of samples.

Pathogen(s)	Agent Category
Influenza H5N1 virus	Viruses

Procedures are described for analysis of clinical samples and may be adapted for assessment of solid, particulate, aerosol, liquid, and water samples. This is a two-step, real-time reverse transcriptase-polymerase chain reaction (RT-PCR) multiplex assay. It employs a mixture of two sets of primers and dual-labeled fluorescent probes that specifically target two different regions of the HA gene of H5N1. Viral ribonucleic acid (RNA) is extracted using a commercial viral RNA extraction kit and reverse transcribed with random hexamers. Five microliters of the complementary deoxyribonucleic acid (cDNA) is used for PCR. At the end of each annealing step, the fluorescent signal of each reaction is measured at a wavelength of 530 nm with the fluorimeter. The assay is specific for the H5 subtype. Influenza H5N1 virus samples are to be handled at BioSafety Level-3 (BSL-3) laboratories.

Please note: These procedures have been applied to analysis of clinical (respiratory) samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank. PCR quality control checks should be performed according to *EPA Draft Quality Assurance/Quality Control Guidance for Laboratories Performing PCR Analyses on Environmental Samples* document (www.epa.gov/nerlcwww/qa_qc_pcr10_04.pdf) or consult the point of contact identified in Section 4.

Source: Ng, E.K.O., Cheng, P.K.C., Ng, A.Y.Y., Hoang, T.L., and Lim, W.W.L. 2005. “Influenza A H5N1 Detection.” *Emerging Infectious Diseases*. 11(8): 1303–1305. <http://www.cdc.gov/ncidod/EID/vol11no08/pdfs/04-1317.pdf>

7.2.24 Literature Reference for Enteric Viruses (Applied and Environmental Microbiology. 2003. 69(6): 3158–3164)

These procedures should be used for **detection** of enterovirus and hepatitis A virus (HAV) in solid, particulate, aerosol, liquid, and water samples, and but are not suitable for determining viability. These procedures should also be used for **detection** and **viability** assessment of rotavirus (Group A) in solid, particulate, aerosol, liquid, and water samples, although rotavirus is not likely to be viable in aerosol samples. Note: Procedures described in the USEPA Manual of Methods for Virology (EPA/600/4-84/013, April 2001) can be adapted for the preparation of samples.

Pathogen(s)	Agent Category
Picornaviruses: Enteroviruses	Viruses
Picornaviruses: Hepatitis A virus (HAV)	Viruses
Reoviruses: Rotaviruses (Group A)	Viruses

Procedures are described for analysis of water samples and may be adapted for assessment of solid, particulate, aerosol, and liquid samples. The method is used to detect human enteric viruses (enteroviruses, HAV, rotavirus) in ground water samples. It is a multiplex reverse-transcription polymerase chain reaction (RT-PCR) procedure optimized for the simultaneous detection of enteroviruses, HAV, reoviruses, and rotaviruses. Water samples are collected by filtration and viruses are eluted using a beef extract solution (1.5%, pH 9.5). Viruses are concentrated using celite adsorption (pH 4.0), filtration, and celite-elution with sodium phosphate (0.15 M, pH 9.0), followed by further concentration and processing to remove inhibitors (ultracentrifugation, solvent extraction, and molecular weight (MW)-exclusion filtration). Concentrated samples are analyzed by a two-step multiplex RT-PCR using virus-specific primer sets. Detection of amplicons is by gel electrophoresis with subsequent confirmation by hybridization (dot-blot) using digoxigenin-labeled internal (nested) probes.

Please note: These procedures have been applied to analysis of ground water. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank. PCR quality control checks should be performed according to *EPA Draft Quality Assurance/Quality Control Guidance for Laboratories Performing PCR Analyses on Environmental Samples* document (www.epa.gov/nerlcwww/qa_qc_pcr10_04.pdf) or consult the point of contact identified in Section 4.

Source: Fout, G. S., Martinson, B. C., Moyer, M. W. N., and Dahling, D. R. 2003. “A Multiplex Reverse Transcription-PCR Method for Detection of Human Enteric Viruses in Groundwater.” *Applied and Environmental Microbiology*. 69(6): 3158–3164. (<http://aem.asm.org/cgi/reprint/69/6/3158.pdf>)

7.2.25 Literature Reference for *Cryptosporidium* spp. (Applied and Environmental Microbiology. 1999. 65(9): 3936–3941)

These procedures should be used for **detection** and **viability** assessment of *Cryptosporidium* spp. in solid, particulate, liquid, and water samples.

Pathogen(s)	Agent Category
<i>Cryptosporidium</i> spp. [Cryptosporidiosis]	Protozoa

Procedures are described for analysis of animal samples and may be adapted for assessment of solid, particulate, liquid, and water samples. A cell culture infectivity assay capable of detecting infectious oocysts is used to quantify viable oocysts through sporozoite invasion and clustering of foci. Oocysts diluted in a standard 5- or 10-fold multiple dilution format are inoculated onto human ileocecal adenocarcinoma (HCT-8) cell monolayers. After incubation for 48 hours, anti-sporozoite polyclonal

antibody is used to detect sporozoite invasion, and microscopy is used to confirm replication (life stages present). Levels of infection and clustering are used to determine the most probable number (MPN) of infectious oocysts in the stock suspension. For oocysts less than 30 days of age, the correlation between the initial oocyst inoculum and the MPN calculation is 0.9726. The relationship between the oocyst inoculum and the MPN diverge as the oocysts age. The 50% infective dose (ID50) in the cell culture system is approximately 10 oocysts.

Please note: These procedures have been applied to analysis of animal samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank.

Source: Slifko, T.R., Huffman, D.E., and Rose, J.B. 1999. “A Most-Probable-Number Assay for Enumeration of Infectious *Cryptosporidium parvum* Oocysts.” *Applied and Environmental Microbiology*. 65(9): 3936–3941. <http://aem.asm.org/cgi/reprint/65/9/3936>

7.2.26 Literature Reference for *Entamoeba histolytica* (Journal of Parasitology. 1972. 58(2): 306–310)

These procedures should be used for **detection** and **viability** assessment of *Entamoeba histolytica* in solid, particulate, liquid, and water samples.

Pathogen(s)	Agent Category
<i>Entamoeba histolytica</i>	Protozoa

Procedures are described for analysis of clinical samples and may be adapted for assessment of solid, particulate, liquid, and water samples. *Entamoeba histolytica* cysts are placed in a modified TP-S-1 medium and incubated for 10 hours. Live amoebae excyst through a rupture in the cyst wall, whereas non-viable amoebae remain encysted. Microscopic examination of an aliquot of the incubated excystation culture allows calculation of the percent of empty (live) cysts and full (dead) cysts in a population.

Please note: These procedures have been applied to analysis of clinical (fecal) samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank.

Source: Stringert, R.P. 1972. “New Bioassay System for Evaluating Percent Survival of *Entamoeba histolytica* Cysts.” *The Journal of Parasitology*. 53(2): 306–310.

7.2.27 Literature Reference for *Entamoeba histolytica* (Journal of Clinical Microbiology. 2005. 43(11): 5491–5497)

These procedures should be used for **detection** of *Entamoeba histolytica* in solid, particulate, liquid, and water samples, but are not suitable for determining viability.

Pathogen(s)	Agent Category
<i>Entamoeba histolytica</i>	Protozoa

Procedures are described for analysis of clinical samples and may be adapted for assessment of solid, particulate, liquid, and water samples. The method is a real-time polymerase chain reaction (PCR) assay that targets the 18S rRNA gene sequence of *E. histolytica*. Deoxyribonucleic acid (DNA) is extracted using cell disruption and a commercial DNA extraction kit, with a second commercial kit used to remove potential PCR inhibitors. TaqMan® real-time PCR is performed on the purified product. The limit of

detection is 1 (\pm 4) cells per mL of sample within 4 hours. The method differentiates between *E. histolytica* and *E. dispar*.

Please note: These procedures have been applied to analysis of clinical (fecal and liver abscess) samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank. PCR quality control checks should be performed according to *EPA Draft Quality Assurance/Quality Control Guidance for Laboratories Performing PCR Analyses on Environmental Samples* document (www.epa.gov/nerlcwww/qa_qc_pcr10_04.pdf) or consult the point of contact identified in Section 4.

Source: Qvarnstrom, Y., James, C., Xayavong, M., Holloway, B.P., Visvesvara, G.S., Sriram, R., and da Silva, A.J. 2005. "Comparison of Real-time PCR Protocols for Differential Laboratory Diagnosis of Amebiasis." *Journal of Clinical Microbiology*. 43(11): 5491–5497.
<http://www.pubmedcentral.nih.gov/picrender.fcgi?artid=1287814&blobtype=pdf>

7.2.28 Literature Reference for *Giardia spp.* (Transactions of the Royal Society of Tropical Medicine and Hygiene. 1983. 77(4): 487–488)

These procedures should be used for **detection** and **viability** assessment of *Giardia spp.* in solid, particulate, liquid, and water samples.

Pathogen(s)	Agent Category
<i>Giardia spp.</i>	Protozoa

Procedures are described for analysis of cell culture samples and may be adapted for assessment of solid, particulate, liquid, and water samples. TYI-S-33 medium supplemented with bovine bile and additional cysteine is used to isolate and culture *Giardia lamblia*. *G. lamblia* is incubated for intervals of 72 and 96 hours at 36°C in borosilicate glass tubes. The cells form a dense, adherent monolayer on the surface of the glass.

Please note: These procedures have been applied to analysis of cell culture samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank.

Source: Keister, D. 1983. "Axenic Culture of *Giardia lamblia* in TYI-S-33 Medium Supplemented with Bile." *Transactions of the Royal Society of Tropical Medicine and Hygiene*. 77(4): 487–488.

7.2.29 Literature Reference for *Toxoplasma gondii* (Emerging Infectious Diseases. 2006. 12(2): 326–329)

These procedures should be used for **detection** and **viability** assessment of *Toxoplasma gondii* in solid, particulate, liquid, and water samples.

Pathogen(s)	Agent Category
<i>Toxoplasma gondii</i> [Toxoplasmosis]	Protozoa

Procedures are described for analysis of water samples and may be adapted for assessment of solid, particulate, liquid, and water samples. Water samples are filtered through fluoropore membrane filters and concentrated by centrifugation. The filters can be assayed by any of three methods. The first method involves performing a bioassay in *T. gondii*-seronegative chickens. Serum samples are tested by enzyme-linked immunosorbent assay (ELISA) and/or modified agglutination test until seroconversion, with the organs from seropositive animals examined microscopically for *T. gondii*. Mice are injected with brain

and heart tissue of seropositive chickens, with parasites found in the lungs of mice being confirmatory for *T. gondii*. The second method is a similar bioassay with pigs and cats. For the third assay, deoxyribonucleic acid (DNA) is extracted from the fluoropore membranes for polymerase chain reaction (PCR) identification of isolates.

Please note: These procedures have been applied to analysis of water samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank.

Source: de Moura, L., Bahia-Oliveira, L.M.G., Wada, M.Y., Jones, J.L., Tuboi, S.H., Carmo, E.H., Ramalho, W.M., Camargo, N.J., Trevisan, R., Graca, R.M.T., da Silva, A.J., Moura, I., Dubey, J.P., and Garrett, D.O. 2006. “Waterborne Toxoplasmosis, Brazil, from Field to Gene.” *Emerging Infectious Diseases*. 12(2): 326–329. <http://www.cdc.gov/ncidod/EID/vol12no02/04-1115.htm>

7.2.30 Literature Reference for *Toxoplasma gondii* (Applied and Environmental Microbiology. 2004. 70(7): 4035–4039)

These procedures should be used for **detection** of *Toxoplasma gondii* in solid, particulate, liquid, and water samples, but are not suitable for determining viability.

Pathogen(s)	Agent Category
<i>Toxoplasma gondii</i> [Toxoplasmosis]	Protozoa

Procedures are described for analysis of water samples and may be adapted for assessment of solid, particulate, and liquid samples. The method uses a fluorogenic 5' nuclease (TaqMan®) real-time polymerase chain reaction (PCR) assay for the detection of *T. gondii* oocyst deoxyribonucleic acid (DNA) using gene-specific (B1 gene) primers and probe. Water samples (10 to 100 L) are filtered to concentrate oocysts. Filters are eluted and recovered oocysts are further purified and concentrated by differential flotation and centrifugation. Final sample pellets are split and subjected to PCR detection and mouse bioassay. In experimental seeding assays, a parasite density of 10 oocysts/L is detectable in 100% of the cases, and a density of 1 oocyst/L is observable in 60% of the cases.

Please note: These procedures have been applied to analysis of water samples. Further research is needed to develop and standardize the protocols for other sample types. At a minimum, the following quality control checks should be performed and evaluated before using this protocol: positive control, negative control, and blank. PCR quality control checks should be performed according to *EPA Draft Quality Assurance/Quality Control Guidance for Laboratories Performing PCR Analyses on Environmental Samples* document (www.epa.gov/nerlcwww/qa_qc_pcr10_04.pdf) or consult the point of contact identified in Section 4.

Source: Villena, I., Aubert, D., Gomis, P., Ferte, H., Ingland, J-C., Denise-Bisiaux, H., Dondon, J-M., Pisano, E., Ortis, N., and Pinon, J.M. 2004. “Evaluation of a Strategy for *Toxoplasma gondii* Oocyst Detection in Water.” *Applied and Environmental Microbiology*. 70(7): 4035–4039. <http://www.pubmedcentral.nih.gov/picrender.fcgi?artid=444816&blobtype=pdf>

Section 8.0: Selected Biotoxin Methods

A list of methods or procedures to be used in analyzing environmental samples for biotoxin contaminants is provided in Appendix D. These methods should be used to support remediation activities (site assessment through clearance) following a homeland security event. Procedures have been compiled for each biotoxin that may need to be identified and/or quantified following a contamination incident. Analytical procedures are not currently available for all the analyte-sample type combinations included in this document. Future research needs include identification of additional methods and development and validation of the procedures listed. Appendix D is sorted alphabetically by analyte, within each of two analyte types (i.e., protein and small molecule).

Please note: This section provides guidance for selecting biotoxin methods that have a high likelihood of assuring analytical consistency when laboratories are faced with a large-scale environmental restoration crisis. Not all methods have been verified for the analyte/sample type combination listed in Appendix D. Please refer to the specified method to identify analyte/sample type combinations that have been verified. Any questions regarding information discussed in this section should be addressed to the appropriate contact(s) listed in Section 4.

Appendix D provides the following information:

- **Analyte(s).** The compound or compound(s) of interest.
- **Analysis type.** Tests are either for presumptive identification, confirmatory identification, or biological activity determination.
- **Analytical Technique.** An analytical instrument or technique used to determine the quantity and identification of compounds or components in a sample.
- **Analytical Method.** The recommended method or procedure, and the corresponding publisher.
- **Aerosol (filter/cassette or liquid impinger).** The recommended method/procedure to measure the analyte of interest in air sample collection media such as filter cassettes and liquid impingers.
- **Solid (soil, powder).** The recommended method/procedure to measure the analyte of interest in solid samples such as soil and powders.
- **Particulate (swabs, wipes, filters).** The recommended method/procedure to measure the analyte of interest in particulate sample collection media such as swabs, wipes and high efficiency particulate air (HEPA) filters.
- **Liquid/drinking water.** The recommended method/procedure to measure the analyte of interest in liquid and drinking water samples.

Some of the biotoxins addressed in this document are commonly found in the environment, and the methods listed in Appendix D assume that analysis will be used to evaluate contamination levels that are above background conditions. If possible, an investigation of initial background levels, as well as controls for background levels, should be performed.

Procedures listed in Appendix D for protein biotoxins are intended to address presumptive, confirmatory, and biological activity determinations. Because the confirmatory procedures listed for the small molecule biotoxins involve a determination of intact compound structure (an indication of biological activity capability), only presumptive and confirmatory methods are listed for these biotoxins. In terms of this document, presumptive methods, or methods that support a reasonable basis for accurate results, should be used in situations that require a large number of samples to be processed. Most of the presumptive methods listed in Appendix D use the immunoassay technique and are designed for large scale sample processing. The confirmatory method, or the method that corroborates the presumptive results, should be used on the smaller subset of samples for which presumptive analysis indicates the presence of the biotoxin. A variety of techniques are listed in Appendix D as confirmatory and generally are more time

consuming and expensive. The use of these terms in this document is not intended to redefine LRN usage of these terms. The terms presumptive and confirmatory as used by the LRN are described in Section 8.2.1. If it is necessary to determine the biological activity of a toxin, either an assay (for proteins) or a technique such as HPLC that determines whether the structure of the biotoxin is intact and likely to be biologically active (for small molecules) may be used. Biological activity analysis should be applied on an as-needed basis following analysis with the confirmatory technique.

EPA's NHSRC is working on a sample collection document that is intended as a companion to SAM. This sample collection document will provide information regarding sampling container/media, preservation, holding time, sample size, and shipping and is intended to complement the laboratory analytical methods that are the focus of the SAM document.

8.1 General Guidance

This section provides a general overview of how to identify the appropriate method(s) for a given biotoxin as well as recommendations for quality control procedures.

For additional information on the properties of the biotoxins listed in Appendix D, TOXNET (<http://toxnet.nlm.nih.gov/index.html>), a cluster of databases on toxicology, hazardous chemicals, and related areas maintained by the National Library of Medicine, is an excellent resource.

Additional resources include:

- A U.S. Army Medical Research Institute of Infectious Diseases' document at <http://www.usamriid.army.mil/education/defensetox/toxdefbook.pdf> contains information regarding sample collection, toxin analysis and identification, as well as decontamination and water treatment.
- Public Health – Select Agents and Toxins, 42 CFR Part 73, found at http://www.cdc.gov/od/sap/pdfs/42_cfr_73_final_rule.pdf
- The U.S. Centers for Disease Control has information regarding toxins, including 42 CFR Part 1003 regulations for possession, use, and transfer of select agents and toxins, on the following Web site: <http://www.cdc.gov/od/sap/sap/toxinamt.htm>
- Syracuse Research Corporation's Physprop and Chemfate, part of the Environmental Fate Database supported by EPA. See <http://www.syrres.com/esc/databases.htm>.
- INCHEM at <http://www.inchem.org/> contains both chemical and toxicity information.
- The Registry of Toxic Effects of Chemical Substances (RTECS) database can be accessed via the NIOSH Web site at <http://www.cdc.gov/niosh/rtecs/vz72d288.html#JWIDAW> for toxicity information.
- EPA's Integrated Risk Information System (IRIS): <http://www.epa.gov/iris/> contains toxicity information.
- The Forensic Science and Communications Journal published by the Laboratory Division of the Federal Bureau of Investigation. See <http://www.fbi.gov/hq/lab/fsc/current/backissu.htm>.

Additional research on biotoxin contaminants is ongoing within EPA.

8.1.1 Standard Operating Procedures for Identifying Biotoxin Methods

To determine the appropriate method that is to be used on an environmental sample, locate the biotoxin of concern in Appendix D: Biotoxin Methods under the "Analyte(s)" column. After locating the biotoxin, continue across the table and identify the appropriate analysis type. After an analysis type has been chosen, find the analytical technique (e.g., immunoassay) and analytical method applicable to the sample

type of interest (solid, particulate, liquid/drinking water, or aerosol) corresponding to that particular analyte.

Sections 8.2.1 through 8.2.27 below provide summaries of the analytical methods listed in Appendix D. Once a method has been identified in Appendix D, **Table 8-1** can be used to locate the method summary.

Table 8-1. Biotoxin Methods and Corresponding Text Section Numbers

Biotoxin	Method	Section
Proteins		
Abrin	119 th AOAC Annual Meeting & Exposition, 2005, p.613	8.2.7
	Pharmacology & Toxicology. 2001. 88(5): 255–260	8.2.8
	Analytical Biochemistry. 2006. 357(2): 200–207	8.2.9
Botulinum neurotoxins ¹	U.S. FDA, Bacteriological Analytical Manual Online, January 2001, Chapter 17, <i>Clostridium botulinum</i>	8.2.2
	Lateral Flow Immunoassay Kits	8.2.27
α -Conotoxin	Biochemical Journal. 1997. 328: 245–250	8.2.10
	Journal of Medicinal Chemistry. 2004. 47(5): 1234–1241	8.2.11
Ricin ¹	Analytical Biochemistry. 2006. 357(2): 200–207	8.2.9
	Journal of Food Protection. 2005. 68(6): 1294–1301	8.2.12
	Lateral Flow Immunoassay Kits	8.2.27
Shiga and Shiga-like toxins (Stx, Stx-1, Stx-2)	U.S. FDA, Bacteriological Analytical Manual Online, January 2001, Appendix 1, Rapid Methods for Detecting Foodborne Pathogens	8.2.3
	Pharmacology & Toxicology. 2001. 88(5): 255–260	8.2.8
	Journal of Clinical Microbiology. 2007. 45(10) 3377–3380	8.2.24
Staphylococcal enterotoxins (SEB) ¹	AOAC Official Method 993.06	8.2.5
Staphylococcal enterotoxins (SEA, SEC)	AOAC Official Method 993.06	8.2.5
Small Molecules		
Aflatoxin (Type B1)	AOAC Official Method 991.31	8.2.4
	AOAC Official Method 994.08	8.2.6
α -Amanitin	Journal of Food Protection. 2005. 68(6):1294–1301	8.2.12
	Journal of Chromatography. 1991. 563(2): 299–311	8.2.13
Anatoxin-a	Biomedical Chromatography. 1996. 10: 46–47	8.2.14
Brevetoxins (B form)	Environmental Health Perspectives. 2002. 110(2): 179–185	8.2.15
	Toxicon. 2004. 43(4): 455–465	8.2.16
Cylindrospermopsin	FEMS Microbiology Letters. 2002. 216: 159–164	8.2.17
Diacetoxyscirpenol (DAS)	International Journal of Food Microbiology. 1988. 6(1): 9–17	8.2.18
	Rapid Communications in Mass Spectrometry. 2006. 20(9): 1422–1428	8.2.19
Microcystins (Principal isoforms: LR, YR, RR, LW)	Journal of AOAC International. 2001. 84(4): 1035–1044	8.2.20
	Analyst. 1994. 119(7): 1525–1530	8.2.21
Picrotoxin	Journal of Pharmaceutical and Biomedical Analysis. 1989. 7(3): 369–375	8.2.22
Saxitoxin (STX, NEOSAX, GTX, dcGTX, dcSTX)	Journal of AOAC International. 1995. 78: 528–532	8.2.23

Biotoxin	Method	Section
T-2 Mycotoxin	Journal of Food Protection. 2005. 68(6): 1294–1301	8.2.12
	Rapid Communications in Mass Spectrometry. 2006. 20(9): 1422–1428	8.2.19
Tetrodotoxin	Journal of Clinical Laboratory Analysis. 1992. 6: 65–72	8.2.25
	Analytical Biochemistry. 2001. 290: 10-17	8.2.26

¹For solid, particulate, liquid, and water sample types, LRN methods are used. See Section 8.2.1.

Method summaries are listed in order of method selection hierarchy (see Figure 2-1), starting with EPA methods followed by methods from other federal agencies, voluntary consensus standard bodies (VCSB), and journal articles. Methods are listed in numerical order under each publisher. Where available, a direct link to the full text of the method is provided with the method summary. For additional information on sample preparation procedures and methods available through consensus standards organizations, please use the contact information provided in **Table 8-2**.

Table 8-2. Sources of Biotoxin Methods

Name	Publisher	Reference
U.S. FDA, Bacteriological Analytical Manual Online, Chapter 17	U.S. Food and Drug Administration	http://www.cfsan.fda.gov/~ebam/bam-17.html
U.S. FDA, Bacteriological Analytical Manual Online, Appendix 1	U.S. Food and Drug Administration	http://www.cfsan.fda.gov/~ebam/bam-a1.html
Official Methods of Analysis of AOAC International*	AOAC International	http://www.aoac.org
Pharmacology & Toxicology*	Blackwell Synergy	http://www.blackwell-synergy.com/
Analytical Biochemistry*	Science Direct	http://www.sciencedirect.com/
Biochemical Journal	Biochemical Journal	http://www.biochemj.org/
Journal of Medicinal Chemistry*	American Chemical Society	http://www.acs.org/
Journal of Food Protection*	International Association for Food Protection	http://www.foodprotection.org/
Journal of Chromatography*	Elsevier Science Publishers	http://www.elsevier.com/
Biomedical Chromatography*	John Wiley And Sons Ltd	http://www.wiley.com/
Environmental Health Perspectives*	National Institute of Environmental Health Sciences	http://www.niehs.nih.gov/
Toxicon*	Elsevier Science Publishers	http://www.elsevier.com/
FEMS Microbiology Letters*	Blackwell Publishing	http://www.blackwellpublishing.com/
International Journal of Food Microbiology*	Elsevier Science Publishers	http://www.elsevier.com/
Rapid Communications in Mass Spectrometry *	John Wiley And Sons Ltd.	http://www.wiley.com/
Journal of AOAC International*	AOAC International	http://www.aoac.org

Name	Publisher	Reference
Analyst*	Royal Society of Chemistry	http://www.rsc.org/
Journal of Pharmaceutical and Biomedical Analysis*	Elsevier Science Publishers	http://www.elsevier.com/
Journal of Clinical Microbiology	American Society for Microbiology	http://www.asm.org/
Journal of Clinical Laboratory Analysis*	John Wiley And Sons Ltd.	http://www.wiley.com/
Lateral Flow Immunoassay Kits	Environmental Technology Verification (ETV) Program	http://www.epa.gov/etv/

* Subscription and/or purchase required.

8.1.2 General Quality Control (QC) Guidance for Biotoxin Methods

Having data of known and documented quality is critical for public officials to determine the activities that may be needed during environmental remediation following emergency events. Having such data requires that laboratories: (1) conduct the necessary QC to ensure that measurement systems are in control and operating properly, (2) properly document results of the analyses, and (3) properly document measurement system evaluation of the analysis-specific QC. Ensuring data quality also requires that laboratory results are properly evaluated and the results of the data quality evaluation are transmitted to decision makers.

The level or amount of QC needed often depends on the intended purpose of the data that are generated. Various levels of QC may be required if the data are generated during presence/absence determinations versus confirmatory analyses. The specific needs for data generation should be identified, and quality control requirements and data quality objectives should be derived based on those needs, and should be applied consistently across laboratories when multiple laboratories are used. For example, during rapid sample screening, minimal QC samples (e.g., blanks, duplicates) and documentation might be required to ensure data quality. Sample analyses for environmental evaluation during site assessment through site clearance, such as those identified in this document, might require increased QC (e.g., demonstrations of method sensitivity, precision, and accuracy).

While method-specific QC requirements may be included in many of the procedures that are cited in this document, and will be referenced in any standardized analytical protocols developed to address specific analytes and sample types of concern, the following describes a minimum set of QC samples and procedures that should be conducted for all analyses. Individual methods, sampling and analysis protocols, or contractual statements of work also should be consulted to determine any additional QC that may be needed. QC tests should be run as frequently as necessary to ensure the reliability of analytical results. In general, sufficient QC includes an initial demonstration of measurement system capability as well as ongoing assessments to ensure the continued reliability of the analytical results. Examples of sufficient quality control for the presumptive tests listed in Appendix D include:

- Method blanks;
- Positive test samples / negative test samples;
- Calibration check samples;
- Use of test kits and reagents prior to expiration; and
- Accurate temperature controls.

Examples of sufficient quality control for the confirmatory tests listed in Appendix D include:

- Demonstration that the measurement system is operating properly
 - Initial calibration
 - Method blanks
- Demonstration of measurement system suitability for intended use
 - Precision and recovery (verify measurement system has adequate accuracy)
 - Analyte/sample type/level of concern-specific QC samples (verify that measurement system has adequate sensitivity at levels of concern)
- Demonstration of continued measurement system reliability
 - Matrix spike/matrix spike duplicates (recovery and precision)
 - QC samples (system accuracy and sensitivity at levels of concern)
 - Continuing calibration verification
 - Method blanks

Please note: The appropriate point of contact identified in Section 4 should be consulted regarding appropriate QA/QC procedures prior to sample analysis. These contacts will consult with the EPA OSWER coordinator responsible for laboratory activities during the specific event to ensure QA/QC procedures are performed consistently across laboratories. OSWER is planning to develop QA/QC guidance for laboratory support. EPA program offices will be responsible for ensuring that the QA/QC practices are implemented.

8.1.3 Safety and Waste Management

It is imperative that safety precautions be used during collection, processing, and analysis of environmental samples. Laboratories should have a documented health and safety plan for handling samples that may contain target chemical, biological, or radiological contaminants, and laboratory staff should be trained in and implement the safety procedures included in the plan. In addition, many of the methods summarized or cited in Section 8.2 contain some specific requirements, guidance, or information regarding safety precautions that should be followed when handling or processing environmental samples and reagents. These methods also provide information regarding waste management. Other resources that can be consulted for additional information include the following:

- Occupational Health and Safety Administration's (OSHA) standard for Occupational Exposure to Hazardous Chemicals in Laboratories (29 CFR 1910.1450);
- OSHA regulations for hazardous waste operations and emergency response (29 CFR 1910.120);
- Environmental Protection Agency's standards regulating hazardous waste (40 CFR Parts 260 – 270);
- U.S. Department of Transportation (DOT) regulations for transporting hazardous materials (49 CFR Part 172);
- U.S. Department of Health and Human Services, Centers for Disease Control and Prevention's requirements for possession, use, and transfer of select agents and toxins (42 CFR Part 1003);
- American Society for Microbiology, Biological Safety: Principles and Practices, 4th Ed. (<http://estore.asm.org/>); and
- American Biological Safety Association, Risk Group Classifications for Infectious Agents (<http://www.absa.org/resriskgroup.html>).
- Select Agent Rules and Regulations (42 CFR Part 73 and 9 CFR Part 121) found at http://www.access.gpo.gov/nara/cfr/waisidx_03/42cfr73_03.html and http://ecfr.gpoaccess.gov/cgi/t/text/text-idx?c=ecfr&tpl=/ecfrbrowse/Title09/9cfr121_main_02.tpl

8.2 Method Summaries

Summaries for the analytical methods listed in Appendix D are provided in Sections 8.2.1 through 8.2.27. These sections contain summary information only, extracted from the selected methods. The full version of the method should be consulted prior to sample analysis.

Each summary contains a table identifying the biotoxin(s) and sample type to which the method applies, a brief description of the method, performance data (if available), and a link to or source for obtaining a full version of the method.

8.2.1 Laboratory Response Network (LRN)

The agents and sample types identified below and listed in Appendix D are included in the U.S. Department of Health and Human Services/U.S. Department of Agriculture (HHS/USDA) select agent list and should be analyzed in accordance with the appropriate LRN protocols.

The LRN was created in accordance with Presidential Directive 39, which established terrorism preparedness responsibilities for federal agencies. The LRN is primarily a national network of local, state, federal, military, food, agricultural, veterinary, and environmental laboratories; however, additional LRN laboratories are located in strategic international locations. The Centers for Disease Control and Prevention (CDC) provides technical and scientific support to member laboratories as well as secure access to standardized procedures and reagents for rapid (within 4 to 6 hours) presumptive detection of biothreat agents and emerging infectious disease agents. These rapid presumptive assays are part of agent-specific algorithms of assays which lead to a confirmed result. The algorithm for a confirmed result is often a combination of one or more presumptive positive results from a rapid assay in combination with a positive result from one of the “gold standard” methods, such as culture. The standardized procedures, reagents, and agent-specific algorithms are considered to be sensitive and are only available to LRN member laboratories. Thus, these procedures are not available to the general public and are not discussed in this document.

It is important to note that, in some cases, the procedures may not be fully developed or validated for each environmental sample type/analyte combination listed in Appendix D, nor are all LRN member laboratories necessarily capable of analyzing all of the sample type/analyte combinations.

Analyte(s)	Sample Type	CAS RN
Botulinum neurotoxins (Serotypes A, B, E, F)	Solid, Particulate, Liquid/Water	NA
Ricin toxin	Solid, Particulate, Liquid/Water	9009-86-3
Staphylococcal enterotoxin B (SEB)	Solid, Particulate, Liquid/Water	39424-53-8

Please note: Not all methods have been verified for the biotoxin/sample type combination listed in Appendix D. Please refer to the agent-specific method to identify the biotoxin/sample type combinations that have been validated. Any questions regarding information discussed in this section should be referred to the appropriate contact(s) listed in Section 4.

For additional information on the LRN, including selection of a laboratory capable of receiving and processing the specified sample type/analyte, please use the contact information provided below or visit <http://www.bt.cdc.gov/lrn/>.

Centers for Disease Control and Prevention

Laboratory Response Branch

Division of Bioterrorism Preparedness and Response (DBPR)

National Center for Prevention, Detection, and Control of Infectious Diseases (NCPDCID)

Coordinating Center for Infectious Diseases (CCID)

Centers for Disease Control and Prevention (CDC)

1600 Clifton Road NE, Mailstop C-18

Atlanta, GA 30333

Telephone: (404) 639-2790

E-mail: lrn@cdc.gov

Local public health laboratories, private, and commercial laboratories with questions about the LRN should contact their state public health laboratory director or the Association of Public Health Laboratories (contact information provided below).

Association of Public Health Laboratories

8515 Georgia Avenue, Suite 700

Silver Spring, MD, 20910

Telephone: (240) 485-2745

Fax: (240) 485-2700

Web site: www.aphl.orgE-mail: info@aphl.org**8.2.2 U.S. FDA, Bacteriological Analytical Manual Online, Chapter 17, 2001: Botulinum Neurotoxins**

These methods should be used for the **confirmatory** and **biological activity** analysis of Types A, B, E, and F botulinum neurotoxins in aerosol samples.

Analyte(s)	Agent Category	CAS RN
Botulinum neurotoxins (Serotypes A, B, E, F)	Protein	NA

An amplified-enzyme-linked immunosorbent assay (amp-ELISA) and a digoxigenin-labeled enzyme-linked immunosorbent assay (DIG-ELISA) are described for the detection of Types A, B, E, and F botulinum neurotoxins in food products. The amp-ELISA method uses goat anti-A or E, rabbit anti-B, or horse anti-F serum to capture the toxins in a 96-well plate, and a corresponding biotinylated goat antitoxin to detect the toxin. Visualization is with streptavidin-alkaline phosphatase. The DIG-ELISA method is a modification of the amp-ELISA method, with digoxigenin-labeled antitoxin IgG's substituted for the streptavidin-alkaline phosphatase. Toxin can be detected at approximately 10 minimum lethal doses (MLD)/mL (0.12 to 0.25 ng/mL). Very toxic samples (greater than 10,000 MLD/mL) may give a positive absorbance for more than one toxin type. Further dilution of the sample will remove cross-reactivity.

The mouse bioassay detects biologically active toxin using a three part approach: toxin screening; toxin titer; and finally, toxin neutralization using monovalent antitoxins. Samples are prepared by centrifugation for suspended solids under refrigeration, or solids are extracted with an equal volume of pH 6.2 gel-phosphate buffer and then centrifuged. Toxins of nonproteolytic types may need trypsin activation to be detected. Serial dilutions of untreated and trypsin-treated sample fluids are injected in separate pairs of mice intraperitoneally (i.p.). Mice are also injected with heated untreated, undiluted sample. Death of mice, along with symptoms of botulism, confirms presence of botulinum toxin. After calculation of an MLD, dilute monovalent antitoxins types A, B, E, and F are injected into mice 30 minutes to 1 hour before challenging them with the i.p. injection of each dilution that gave the highest MLD from the toxic preparation.

Please note: This method has been applied to analysis of food products. Further research is needed to develop and standardize the protocols for other sample types.

Source: U.S. Food & Drug Administration, Center for Food Safety & Applied Nutrition. 2001. *Bacteriological Analytical Manual Online*. Chapter 17. (<http://www.cfsan.fda.gov/~ebam/bam-17.html>)

8.2.3 U.S. FDA, Bacteriological Analytical Manual Online, Appendix 1, 2001: Rapid Methods for Detecting Foodborne Pathogens

This method should be used for the **presumptive** analysis of Shiga (Stx) and Shiga-like (Stx-1 and Stx-2) biotoxins in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Shiga and Shiga-like toxins (Stx, Stx-1, Stx-2)	Protein	75757-64-1

Shiga toxin (Stx) is produced by *Shigella dysenteriae* and Shiga-like toxins (Stx-1 and Stx-2) are produced by various Shiga-toxigenic *E. coli* (STEC). An enzyme-linked immunosorbent assay (ELISA) is described for the detection of these toxins. Diluted samples are added to microwells coated with an anti-Shiga toxin capture antibody. After incubation at room temperature, a wash is performed to remove unbound material. A second anti-Shiga toxin antibody is added for detection and incubation continued at room temperature. A wash is performed to remove unbound antibody. Enzyme conjugated anti-IgG visualization antibody is added and the plate incubated then rinsed. Substrate is added, and after incubation to develop the color, stop solution is added. The results are interpreted spectrophotometrically.

Please note: This method has not been applied to analysis of environmental samples. Further research is needed to develop and standardize the protocols for these sample types.

Source: U.S. Food & Drug Administration, Center for Food Safety & Applied Nutrition. 2001. *Bacteriological Analytical Manual Online*. Appendix 1. (<http://www.cfsan.fda.gov/~ebam/bam-a1.html>)

8.2.4 AOAC Official Method 991.31: Aflatoxins in Corn, Raw Peanuts, and Peanut Butter

This method should be used for **presumptive** analysis of aflatoxin Type B1 in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Aflatoxin (Type B1)	Small Molecule	1402-68-2

This method is for the detection of aflatoxins in agricultural products. The sample is extracted with methanol-water (7 + 3), filtered, diluted with water, and applied to an affinity column containing monoclonal antibody specific for aflatoxins B1, B2, G1, and G2. Antibody-bound aflatoxins are removed from the column with methanol. Total aflatoxins are quantified by fluorescence measurement after reaction with bromine solution. Individual aflatoxins are quantitated by LC with fluorescence detection and postcolumn iodine derivatization. Method performance was characterized using various commodities (e.g., corn) at aflatoxin levels over a range of 10 to 30 ng/g. This method was originally designed for the analysis of aflatoxins (B₁, B₂, G₁, and G₂) in samples where cleanup was necessary to remove food components, such as fats and proteins; the cleanup procedure may not be necessary for analysis of water samples.

Please note: This method has been applied to analysis of agricultural products. Further research is needed to develop and standardize the protocols for other sample types.

Source: AOAC International. 1994. *Official Methods of Analysis of AOAC International*. 16th Edition, 4th Revision; Vol II.

8.2.5 AOAC Official Method 993.06: Staphylococcal Enterotoxins in Selected Foods

This method should be used for **presumptive** analysis of staphylococcal enterotoxins Type B in aerosol samples, and Types A and C in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Staphylococcal enterotoxins (SEB, SEA, SEC)	Protein	39424-53-8, 37337-57-8, 39424-54-9

This method is an enzyme immunoassay (EIA) using a mixture of high-affinity capture for identification of toxin(s) in food samples. Samples are prepared by dilution in Tris buffer, centrifugation, and filtration of the supernatant through a syringe, with adjustment to a final pH of 7.0 to 8.0. Samples are incubated in 96-well plates with the mixture of antibodies conjugated to horseradish peroxidase, and visualized with an acetic peroxidase substrate. Assay results are determined visually or using a microtiter plate reader. Test is considered positive for staphylococcal enterotoxins if absorbance is >0.200 and is considered negative if absorbance is ≤0.200. Specific toxin serotypes are not differentiated. This method detects from 1.3 to 3.3 ng/mL staphylococcal enterotoxin in extracts prepared from food containing 4 to 10 ng/mL staphylococcal enterotoxin.

Please note: This method has been applied to analysis of food products. Further research is needed to develop and standardize the protocols for other sample types.

Source: AOAC International. 1994. *Official Methods of Analysis of AOAC International*. 16th Edition, 4th Revision; Vol I.

8.2.6 AOAC Official Method 994.08: Aflatoxin in Corn, Almonds, Brazil Nuts, Peanuts, and Pistachio Nuts

This method should be used for **confirmatory** analysis of aflatoxin Type B1 in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Aflatoxin (Type B1)	Small Molecule	1402-68-2

This method is for the identification of aflatoxins in agricultural products. Samples are extracted using an acetonitrile-water (9 + 1) solution. Sample extracts are then run through a multifunctional cleanup column. The purified extract and standards are derivatized with trifluoroacetic acid, and then analyzed using a high performance liquid chromatography (HPLC) system with a fluorescence detector. Specific aflatoxins can be identified by their retention time and quantified using standard curves. Method performance was characterized using various commodities (e.g., corn) at aflatoxin levels over a range of 5 to 30 ng/g. This method was originally designed for the analysis of aflatoxins (B₁, B₂, G₁, and G₂) in commodities where cleanup was necessary to remove other food components, such as fats and proteins; the cleanup procedure may not be necessary for water analyses. Coupling the procedures, or a modification of the procedures, included in this method with an immunoassay and/or viability test (where available) will provide more information regarding specificity and toxicity.

Please note: This method has been applied to analysis of agricultural products. Further research is needed to develop and standardize the protocols for other sample types.

Source: AOAC International. 1998. *Official Methods of Analysis of AOAC International*. 16th Edition, 4th Revision; Vol II.

8.2.7 Literature Reference for Abrin (119th AOAC Annual Meeting & Exposition, 2005, p. 613)

These procedures should be used for **presumptive** analysis of abrin in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Abrin	Protein	1393-62-0

Procedures are described for using mouse monoclonal antibodies and rabbit derived polyclonal antibodies prepared against a mixture of abrin isozymes for three separate enzyme-linked immunosorbent assay (ELISA) and electrochemiluminescence (ECL)- based assays in food products. The three assays vary by use of antibody combination: (1) polyclonal (capture)/polyclonal (detection) ELISA, (2) polyclonal/monoclonal ELISA, and (3) polyclonal/monoclonal ECL assay. The limits of detection (LOD) with purified Abrin C and various abrin extracts in buffer are between 0.1 and 0.5 ng/mL for all three assays. The LOD for abrin dissolved into food products ranges from 0.1 to 0.5 ng/mL, using the ECL assay. The LOD for abrin dissolved into food products for the ELISA assays range between 1 and 4 ng/mL, depending on the assay configuration. In all cases, the LODs are considerably less than the concentration at which abrin may pose a health concern.

Please note: These procedures have been applied to analysis of food products. Further research is needed to develop and standardize the protocols for other sample types.

Source: Garber, E.A., Aldrich, J.L., Wang, J., Brewer, V.A., O'Brien, T.W., and Sigal, G. 2005. "Detection of Abrin Foods Using ELISA and Electrochemiluminescence (ECL) Technologies." 119th AOAC Annual Meeting & Exposition. p. 613.

8.2.8 Literature Reference for Abrin and Shiga and Shiga-like Toxins (Pharmacological Toxicology. 2001. 88(5): 255-260)

These procedures should be used for **confirmatory** analysis of abrin and Shiga and Shiga-like toxins in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Abrin	Protein	1393-62-0
Shiga and Shiga-like Toxins (Stx, Stx-1, Stx-2)	Protein	75757-64-1

Procedures are described for measuring the biological activity of ribosome-inactivating proteins using a microtiter plate format for detection of abrin in phosphate buffered saline (PBS). Nuclease-treated rabbit reticulocyte lysate containing luciferase messenger ribonucleic acid (mRNA) is used to measure toxin activity via inhibition of protein synthesis. The relative biological activity is determined by comparing luminescence levels in treated samples versus those of untreated controls. The amount of luciferase translated, as measured by luminescence, is inversely proportional to the toxin concentration. Linear dose response curves are generated for abrin, with a 50% inhibition of translation at 0.5 nM. Coupling this procedure, or a modification of this procedure, with an immunoassay will provide more information regarding the specificity and toxicity of the target biotoxin.

Please note: These procedures have been applied to analysis of abrin in PBS. Further research is needed to develop and standardize the protocols for other sample types.

Source: Hale, M.L. 2001. “Microtiter-based Assay for Evaluating the Biological Activity of Ribosome-inactivation Proteins.” *Pharmacological Toxicology*. 88(5): 255–260.

8.2.9 Literature Reference for Abrin and Ricin (Analytical Biochemistry. 2006. 357(2): 200–207)

These procedures should be used for **biological activity** analysis of abrin and ricin in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Abrin	Protein	1393-62-0
Ricin	Protein	9009-86-3

This assay is an N-glycosidase enzyme activity assay for the detection of purified abrin and ricin in a reaction buffer consisting of nuclease free water with Triton X-100. Synthetic biotinylated RNA substrates are cleaved by the combined actions of the toxin and a chemical agent, N,N'-dimethylethyldiamine. Annealing of the product with a ruthenylated oligodeoxynucleotide results in the capture of ruthenium chelate onto magnetic beads. Inclusion of the monoclonal antibody (Mab) 9C3 is used to enhance the N-glycosidase activity. Commercially available electrochemiluminescence (ECL)-based reagents are used to detect the product. Polyclonal/monoclonal antibodies are commercially available as an ELISA test kit. The activity assay exhibits similar limits of detection just below signal with 0.1 ng/ml of ricin; the ECL response was linear as the ricin concentration increased by two orders of magnitude. Abrin II similarly cleaves the biotinylated RNA substrate, with the N-glycosidase activity of the toxin enhanced to a greater degree by Mab 9C3, enabling differentiation from ricin.

Please note: These procedures have been applied to analysis of purified toxins in reaction buffer. Further research is needed to develop and standardize the protocols for other sample types.

Source: Keener, W.K., Rivera, V.R., Young, C.C., and Poli, M.A. 2006. “An Activity-dependent assay for Ricin and Related RNA N-glycosidases Based on Electrochemiluminescence.” *Analytical Biochemistry*. 357(2): 200–207.

8.2.10 Literature Reference for α -Conotoxin (Biochemical Journal. 1997. 328: 245–250)

These procedures should be used for **presumptive** and **biological activity** analysis of α -conotoxin in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
α -Conotoxin	Protein	156467-85-5

A biologically active fluorescein derivative of *Conus geographus* α -conotoxin (FGI) is used in solution-phase-binding assays with purified *Torpedo californica* nicotinic acetylcholine receptor (nAChR) and monoclonal antibodies (Mabs) to detect the toxin in laboratory samples. FGI binding to *T. californica* nAChR or Mabs is determined by spin-column at room temperature. For competitive ligand-displacement assays, FGI is premixed with various dilutions of unlabelled ligands and then incubated with nAChR or mAb. Fluorescence is measured in ratio mode using cuvettes with excitation and emission monochromators set at gamma = 490 nm and gamma = 525 nm respectively. The binding of FGI to nAChR or antibody had apparent dissociation constants of 10 to 100 nM.

Please note: These procedures have been applied to analysis of purified toxin in phosphate buffer. Further research is needed to develop and standardize the protocols for other sample types.

Source: Ashcom, J.D., and Stiles, B.G. 1997. “Characterization of α -Conotoxin Interactions with the Nicotinic Acetylcholine Receptor and Monoclonal Antibodies.” *Biochemical Journal*. 328: 245–250. <http://www.biochemj.org/bj/328/0245/3280245.pdf>

8.2.11 Literature Reference for α -Conotoxin (Journal of Medicinal Chemistry. 2004. 47(5): 1234–1241)

These procedures should be used for **confirmatory** analysis of α -conotoxin in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
α -Conotoxin	Protein	156467-85-5

Procedures are discussed for the detection of peptides within the α -conotoxin molecular mass range using high performance liquid chromatography-mass spectrometry (HPLC-MS). A crude extract of the sample is made using 30% acetonitrile/water acidified with 0.1% trifluoroacetic acid (TFA), with the insoluble portion of the sample removed by centrifugation. A portion of the sample extract is fractionated by size-exclusion chromatography in order to prepare a sample containing small peptides in the range of 1000 to 2500 Da. Chromatography conditions are elution with 30% acetonitrile / 0.048% TFA at a flow rate of 0.5 mL/minute, with detection at 214 nm. Three sulfated alpha-conotoxins (AnIA, AnIB, and AnIC) can be identified by LC-MS. Peptides can be quantified by reversed-phase high performance liquid chromatography (RP-HPLC) using an external reference standard for each peptide.

Please note: These procedures have been applied to analysis of *Conus anemone* venom duct material. Further research is needed to develop and standardize the protocols for other sample types.

Source: Loughnan, M.L., Nicke, A., Jones, A., Adams, D.J., Alewood, P.F., and Lewis, R.J. 2004. “Chemical and Functional Identification and Characterization of Novel Sulfated Alpha-conotoxins from the Cone Snail *Conus anemone*.” *Journal of Medicinal Chemistry*. 47(2): 1234–1241.

8.2.12 Literature Reference for α -Amanitin, Ricin, T-2 Mycotoxin (Journal of Food Protection. 2005. 68(6): 1294–1301)

These procedures should be used for **presumptive** analysis of α -amanitin and T-2 toxin in aerosol, solid, particulate, liquid, and water samples and for **confirmatory** analysis of ricin in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
α -Amanitin	Small Molecule	23109-05-9
Ricin	Protein	9009-86-3
T-2 Mycotoxin	Protein	21259-20-1

Commercially available enzyme-linked immunosorbent assays (ELISAs) are described and assessed for detection of ricin, amanitin, and T-2 toxin at levels below those described as a health concern in food samples. Solid food samples are prepared by washing the sample with sodium phosphate buffer followed by dilution with phosphate-buffered saline. Liquid beverage samples are prepared by dilution in sodium phosphate buffer. Amanitin samples are similarly prepared using water instead of buffer, and T-2 toxin samples are similarly prepared using 35% methanol in water instead of buffer. The prepared samples are used with commercially obtained ELISA kits according to the manufacturer’s directions, except for the

incorporation of an eight-point calibration curve and reading the plates at both 405 and 650 nm after 26 minutes of incubation at 37°C. This assay detects ricin in food products at 0.01 µg/mL with acceptable background levels. Amanitin can be detected in food products at 1 µg/g with the ELISA. Background responses occurred, but at less than the equivalent of 0.5 ppm for amanitin. The ELISA kit will successfully detect T-2 toxin at targeted levels of 0.2 µg/g. The ELISA kit successfully detects T-2 toxin at targeted levels of 0.2 µg/g, the immunoassay for T-2 toxin, however, shows significant background responses and varies up to 0.1 ppm.

Please note: These procedures have been applied to analysis of food products. Further research is needed to develop and standardize the protocols for other sample types.

Source: Garber, E.A., Eppley, R.M., Stack, M.E., McLaughlin, M.A., and Park, D.L. 2005. “Feasibility of Immunodiagnostic Devices for the Detection of Ricin, Amanitin, and T-2 Toxin in Food.” *Journal of Food Protection*. 68(6): 1294–1301.

8.2.13 Literature Reference for α -Amanitin (Journal of Chromatography. 1991. 563(2): 299–311)

These procedures should be used for **confirmatory** analysis of α -amanitin in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
α -Amanitin	Small Molecule	23109-05-9

Procedures are described for the selective determination in human plasma of α -amanitin using high-performance liquid chromatography (HPLC) with amperometric detection. After extraction of plasma with disposable C₁₈ silica cartridges, the extracts are separated by isocratic reversed-phase chromatography using a macroporous polystyrene-divinylbenzene column and a mobile phase of 0.05 M phosphate buffer-acetonitrile (91:9) at pH 9.5. Amperometric detection is performed by applying an oxidation potential as low as +350 mV (vs. Ag/AgCl) to a glassy carbon electrode, in a thin-layer flow-cell. The linear range for alpha-amanitin is 3 to 200 ng/mL, and the relative limit of detection in plasma is 2 ng/mL at a signal-to-noise ratio of 2. The intra-assay precision has been evaluated at levels of 10 and 200 ng/mL.

Please note: These procedures have been applied to analysis of plasma. Further research is needed to develop and standardize the protocols for other sample types.

Source: Tagliaro, F., Schiavon, G., Bontempelli, G., Carli, G., and Marigo, M. 1991. “Improved High-performance Liquid Chromatographic Determination with Amperometric Detection of Alpha-amanitin in Human Plasma Based on its Voltammetric Study.” *Journal of Chromatography*. 563(2): 299–311.

8.2.14 Literature Reference for Anatoxin-a (Biomedical Chromatography. 1996. 10: 46–47)

These procedures should be used for **confirmatory** analysis of anatoxin-a in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Anatoxin-a	Small Molecule	64285-06-9

Procedures are described for high-performance liquid chromatography (HPLC) analysis with fluorimetric detection of anatoxin-a in water samples after derivatization with 7-fluoro-4-nitro-2,1,3-benzoxadiazole (NBD-F). Samples are extracted at pH 7 with solid phase extraction (SPE) using a weak cation exchanger. The toxin is eluted with methanol containing 0.2% trifluoroacetic acid. Samples are

evaporated, reconstituted with acetonitrile, and re-evaporated prior to derivatization. This procedure detects anatoxin-a at concentrations of 0.1 µg/L with a good linear calibration.

Please note: These procedures have been applied to analysis of water. Further research is needed to develop and standardize the protocols for other sample types.

Source: James, K.J., and Sherlock, I.R. 1996. “Determination of the Cyanobacterial Neurotoxin, Anatoxin-a, by Derivatization Using 7-Fluoro-4-Nitro-2,1,3-Benzoxadiazole (NBD-F) and HPLC Analysis with Fluorimetric Detection.” *Biomedical Chromatography*. 10: 46–47.

8.2.15 Literature Reference for Brevetoxins (Environmental Health Perspectives. 2002. 110(2): 179–185)

These procedures should be used for **presumptive** analysis of brevetoxins in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Brevetoxins (B form)	Small Molecule	98225-48-0

Procedures are described for a competitive enzyme-linked immunosorbent assay (ELISA) used to detect brevetoxins in shellfish. The assay uses goat anti-brevetoxin antibodies in combination with a three-step signal amplification process: (1) secondary biotinylated antibodies; (2) streptavidine-horseradish peroxidase conjugate; and (3) chromogenic enzyme substrate. Sample preparation for liquids is dilution in phosphate buffered saline (PBS). Sample preparation for solids (oysters) is homogenization in PBS, or extraction in acetone. The working range for the assay is 0.2 to 2 ng/mL for diluted and undiluted liquid samples, and 0.2 to 2 ng/mL for solid samples, corresponding to 0.8 to 8 µg brevetoxins/100 g shellfish. The method has been compared to the mouse bioassay and is equivalent in sensitivity.

Please note: These procedures have been applied to analysis of shellfish. Further research is needed to develop and standardize the protocols for other sample types.

Source: Naar, J., Bourdelais, A., Tomas, C., Kubanek, J., Whitney, P.L., Flewelling, L., Steidinger, K., Lancaster, J., and Badan, D.G. 2002. “A Competitive ELISA to Detect Brevetoxins from *Karenia brevis* (Formerly *Gymnodinium breve*) in Seawater, Shellfish, and Mammalian Body Fluid.” *Environmental Health Perspectives*. 110(2): 179–185.

<http://www.pubmedcentral.nih.gov/picrender.fcgi?artid=1240733&blobtype=pdf>

8.2.16 Literature Reference for Brevetoxins (Toxicon. 2004. 43(4): 455–465)

These procedures should be used for **confirmatory** analysis of brevetoxins in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Brevetoxins (B form)	Small Molecule	98225-48-0

Shellfish sample homogenates are extracted with acetone, and centrifuged. The supernatants are combined, evaporated, and re-solubilized in 80% methanol. Following a wash with 95% n-hexane, the methanolic layer is evaporated, and the residue re-solubilized in 25% methanol and applied to a C₁₈ SPE column. Analytes are eluted with 100% methanol, evaporated, and re-solubilized in methanol for analysis. Analysis of prepared samples is performed using liquid chromatography tandem mass spectrometers (LC-MS-MS) with a mobile phase of water and acetonitrile with acetic acid. Analytes are detected by mass spectrometer with electrospray ionization (ESI) interface. Brevetoxins are extensively metabolized, with

many sub-forms. This method describes multiple LC-MS-ESI profiles for metabolites of brevetoxins from oysters.

Please note: These procedures have been applied to analysis of shellfish. Further research is needed to develop and standardize the protocols for other sample types.

Source: Wang, Z., Plakas, S.M., El Said, K.R., Jester, E.L., Granade, H.R., and Dickey, R.W. 2004. “LC/MS Analysis of Brevetoxin Metabolites in the Eastern Oyster (*Crassostrea virginica*).” *Toxicon*. 43(4): 455–465.

8.2.17 Literature Reference for Cylindrospermopsin (FEMS Microbiology Letters. 2002. 216(2): 159–164)

These procedures should be used for **confirmatory** analysis of cylindrospermopsin in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Cylindrospermopsin	Small Molecule	143545-90-8

Cylindrospermopsin is detected using high performance liquid chromatography with photodiode array detector (HPLC-PDA) in environmental waters. The suggested solvent range for cylindrospermopsin is below 50% methanol and 30% acetonitrile. Complex samples (culture medium) are purified on a C₁₈ column with a linear gradient of 1 to 12% (v/v) methanol/water over 24 minutes at 40°C, with monitoring at 262 nm. The use of C₁₈ columns for environmental waters is suggested for removal of the large number of organic compounds that may be present. This method detects and recovers cylindrospermopsin from spiked environmental water samples at 1 µg/L.

Please note: These procedures have been applied to analysis of water. Further research is needed to develop and standardize the protocols for other sample types.

Source: Metcalf, J.S., Beattie, K.A., Saker, M.L., and Codd, G.A. 2002. “Effects of Organic Solvents on the High Performance Liquid Chromatographic Analysis of the Cyanobacterial Toxin Cylindrospermopsin and Its Recovery from Environmental Eutrophic Waters by Solid Phase Extraction.” *FEMS Microbiology Letters*. 216(2): 159–164.

8.2.18 Literature Reference for Diacetoxyscirpenol (DAS) (International Journal of Food Microbiology. 1988. 6(1): 9–17)

These procedures should be used for **presumptive** analysis of diacetoxyscirpenol (DAS) in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Diacetoxyscirpenol (DAS)	Small Molecule	2270-40-8

An enzyme-linked immunosorbent assay (ELISA) is used for the detection of diacetoxyscirpenol (DAS) in food samples. Antibodies against DAS are obtained after immunization of rabbits with DAS-hemiglutamate-human serum albumin (DAS-HG-HSA), and a DAS-hemisuccinate-horseradish peroxidase-conjugate (DAS-HS-HRP) is prepared by an ester method for use as enzyme-labeled toxin in the competitive assay. The detection limit for DAS using this assay is approximately 10 pg/mL. The relative cross reactivities of the assay are 597.5, 5.2, 100.0, 2.5, and 1.5% for 3 alpha-acetyl-DAS, DAS, T-2 toxin, neosolaniol, and 15-acetoxyscirpenol, respectively.

Please note: These procedures have been applied to analysis of food. Further research is needed to develop and standardize the protocols for other sample types.

Source: Klaffer, U., Martlbauer, E., and Terplan, G. 1988. “Development of a Sensitive Enzyme-linked Immunosorbent Assay for the Detection of Diacetoxyscirpenol.” *International Journal of Food Microbiology*. 6(1): 9–17.

8.2.19 Literature Reference for Diacetoxyscirpenol (DAS) and T-2 Mycotoxin (Rapid Communications in Mass Spectrometry. 2006. 20(9): 1422–1428)

These procedures should be used for **confirmatory** analysis of diacetoxyscirpenol (DAS) and T-2 mycotoxin in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Diacetoxyscirpenol (DAS)	Small Molecule	2270-40-8
T-2 Mycotoxin	Small Molecule	21259-20-1

A liquid chromatography/atmospheric pressure chemical ionization mass spectrometry (LC/APCI-MS) procedure based on time-of-flight mass spectrometry (TOFMS), with a real-time reference mass correction, is used for simultaneous determination of *Fusarium* mycotoxins (to include diacetoxyscirpenol (DAS) and T-2 mycotoxin) in foodstuffs. Mycotoxin samples are extracted with acetonitrile/water (85:15) and centrifuged, and the supernatant is applied to a column for cleanup. Prepared samples are separated by liquid chromatography with an aqueous mobile phase of ammonium acetate and methanol detection is provided in exact mass chromatograms with a mass window of 0.03 Th. The limits of detection range from 0.1 to 6.1 ng/g in analyzed foodstuffs.

Please note: These procedures have been applied to analysis of food. Further research is needed to develop and standardize the protocols for other sample types.

Source: Tanaka, H., Takino, M., Sugita-Konishi, Y., and Tanaka, T. 2006. “Development of Liquid Chromatography/Time-of-flight Mass Spectrometric Method for the Simultaneous Determination of Trichothecenes, Zearalenone, and Aflatoxins in Foodstuffs.” *Rapid Communications in Mass Spectrometry*. 20(9): 1422–1428.

8.2.20 Literature Reference for Microcystins (Journal of AOAC International. 2001. 84(4): 1035–1044)

These procedures should be used for **presumptive** analysis of microcystins in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Microcystins (Principal Isoforms: LR, YR, RR, LW)	Small Molecule	77238-39-2

Enzyme-linked immunosorbent assay (ELISA) and protein phosphatase inhibition assay are used to detect microcystins in algae products. Solid samples are prepared by homogenization in methanol (75% in water), with centrifugation to remove solids. Immunoassays are performed on the prepared samples using a commercially available ELISA test kit as described by the manufacturer. Samples are quantitated by comparison with a microcystin LR standard curve. Quantitation with the colorimetric protein phosphatase inhibition assay is based on a comparison with a microcystin LR standard curve. ELISA and phosphatase assay results agree over a concentration range of 0.5 to 35 µg/g. Neither assay is specific for a particular isoform.

Please note: These procedures have been applied to analysis of algae products. Further research is needed to develop and standardize the protocols for other sample types.

Source: Lawrence, J.F., Niedzwiadek, B., Menard, C., Lau, B.P., Lewis, D., Kuper-Goodman, T., Carbone, S., and Holmes, C. 2001. “Comparison of Liquid Chromatography/Mass Spectrometry, ELISA, and Phosphatase Assay for the Determination of Microcystins in Blue-green Algae Products.” *Journal of AOAC International*. 84(4): 1035–1044.

8.2.21 Literature Reference for Microcystins (Analyst. 1994. 119(7): 1525–1530)

These procedures should be used for **confirmatory** analysis of microcystins in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Microcystins (Principal Isoforms: LR, YR, RR, LW)	Small Molecule	77238-39-2

Procedures are discussed to test the presence of microcystin-LR, -LY, -LW, -LF, and -RR in treated and untreated water samples. Cyanobacterial cells are separated from the water by filtration through 110-mm GF/C discs. The cellular components collected on the discs are extracted three times with methanol; the collected extraction fluids are combined and dried. The residue is resuspended in methanol and analyzed by photodiode-array high-performance liquid chromatography (PDA-HPLC). The liquid portion of the filtered water sample is subjected to trace enrichment using a C₁₈ solid-phase extraction cartridge, followed by identification and determination by PDA-HPLC. This procedure can detect microcystin concentrations as low as 250 ng/L and is the basis of the World Health Organization (WHO) method for the detection of microcystins.

Please note: These procedures have been applied to analysis of water. Further research is needed to develop and standardize the protocols for other sample types.

Source: Lawton, L.A., Edwards, C., and Codd, G.A. 1994. “Extraction and High-performance Liquid Chromatographic Method for the Determination of Microcystins in Raw and Untreated Waters.” *Analyst*. 119(7): 1525–1530.

8.2.22 Literature Reference for Picrotoxin (Journal of Pharmaceutical & Biomedical Analysis. 1989. 7(3): 369–375)

These procedures should be used for **confirmatory** analysis of picrotoxin in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Picrotoxin	Small Molecule	124-87-8

Procedures are described for quantification of the two components of picrotoxin (picrotin and picrotoxinin) in serum samples. Serum samples are prepared by washing with *n*-hexane, followed by extraction with chloroform. The chloroform is evaporated and the sample is reconstituted in acetonitrile-1 mM ammonium acetate buffer (pH 6.4) 34:66 (v/v) for assay by reversed-phase high-performance liquid chromatography (HPLC). The effluent is monitored at 200 nm, and quantification is based on peak-height ratio of analyte to the internal standard. A linear response is obtained for both analytes (picrotin and picrotoxinin) in the range 0.2 to 20 µg/mL.

Please note: These procedures have been applied to analysis of serum. Further research is needed to develop and standardize the protocols for other sample types.

Source: Soto-Otero, R., Mendez-Alvarez, E., Sierra-Paredes, G., Galan-Valiente, J., Aguilar-Veiga, E., and Sierra-Marcuno, G. 1989. “Simultaneous Determination of the Two Components of Picrotoxin in Serum by Reversed-phase High-performance Liquid Chromatography with Application to a Pharmacokinetic Study in Rats.” *Journal of Pharmaceutical & Biomedical Analysis*. 7(3): 369–375.

8.2.23 Literature Reference for Saxitoxin (Journal of AOAC International. 1995. 78(2): 528–532)

These procedures should be used for **confirmatory** analysis of saxitoxin in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Saxitoxin (STX, NEOSAX, GTX, dcGTX, dcSTX, STX)	Small Molecule	35523-89-8

Procedures are described to detect multiple analogues of saxitoxin in shellfish using ion-interaction chromatography on a silica-based reversed-phase (C₈) column with postcolumn periodate oxidation and fluorescence detection. Toxin groups of different net charges are determined separately by isocratic elution using either sodium 1-heptanesulfonate in ammonium phosphate (GTX-1, GTX-6, dcGTX2, dcGTX3) or sodium 1-heptanesulfonate in ammonium phosphate and acetonitrile (STX, neoSTX, dcSTX). For biological matrices, a cleanup procedure using a C₁₈ solid-phase extraction cartridge is effective in preventing false peaks. High sensitivity with detection limits ranging from 20 to 110 fmol are achieved as a result of reduced band broadening and optimized reaction conditions. This method, when applied to low-toxicity shellfish, gives higher values than the standard mouse bioassay.

Please note: These procedures have been applied to analysis of shellfish. Further research is needed to develop and standardize the protocols for other sample types.

Source: Oshima, Y. 1995. “Postcolumn Derivatization Liquid Chromatographic Method for Paralytic Shellfish Toxins.” *Journal of AOAC International*. 78(2): 528–532.

8.2.24 Literature Reference for Shiga and Shiga-like Toxin (Journal of Clinical Microbiology. 2007. 45(10): 3377–3380)

These procedures should be used for **presumptive** analysis of Shiga and Shiga-like toxins in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Shiga and Shiga-like Toxins (Stx, Stx-1, Stx-2)	Protein	75757-64-1

Procedures are described for a rapid optical immunoassay for the detection of Shiga toxin Types 1 (Stx-1) and 2 (Stx-2) using a commercially available kit. Fecal samples (742 specimens) are assayed for Shiga toxins with and without enrichment of the specimens in broth. Duplicate assays are applied using either the rapid optical immunoassay or a commercially available ELISA kit. Samples producing positive results by immunoassay are confirmed by Vero cell cytotoxicity assay. Sensitivities of 96.8% are achieved for direct stool sample assays.

Please note: These procedures have not been applied to analysis of environmental samples. Further research is needed to develop and standardize the protocols for these sample types.

Source: Teel, L.D., Daly, J.A., Jerris, R.C., Maul, D., Svanas, G., O'Brien, A.D., Park, C.H. "Rapid Detection of Shiga Toxin-Producing *Escherichia coli* by Optical Immunoassay." *Journal of Clinical Microbiology*. 45(10): 3377–3380.

8.2.25 Literature Reference for Tetrodotoxin (Journal of Clinical Laboratory Analysis. 1992. 6(2): 65–72)

These procedures should be used for **presumptive** analysis of tetrodotoxin (TTX) in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Tetrodotoxin	Small Molecule	9014-39-5

Procedures are described for a competitive inhibition enzyme immunoassay (CIEIA) for TTX in biological samples. An anti-TTX monoclonal antibody (Mab), designated T20G10, is directly labeled with alkaline phosphatase for use in the assay. Sensitivities of 6 to 7 ng/mL (IC 50) and 2 to 3 ng/mL (IC20) are achieved.

Please note: These procedures have been applied to analysis of biological samples. Further research is needed to develop and standardize the protocols for other sample types.

Source: Raybould, T.J., Bignami, G.S., Inouye, L.K., Simpson, S.B., Byrnes, J.B., Grothaus, P.G., and Vann, D.C. 1992. "A Monoclonal Antibody-based Immunoassay for Detecting Tetrodotoxin in Biological Samples." *Journal of Clinical Laboratory Analysis*. 6(2): 65–72.

8.2.26 Literature Reference for Tetrodotoxin (Analytical Biochemistry. 2001. 290(1): 10–17)

These procedures should be used for **confirmatory** analysis of tetrodotoxin (TTX) in aerosol, solid, particulate, liquid, and water samples.

Analyte(s)	Agent Category	CAS RN
Tetrodotoxin	Small Molecule	9014-39-5

Procedures are described for liquid chromatography/electrospray ionization mass spectrometry (LC/ESI-MS) analysis of TTXs in tissue samples from puffer fish and newts by a combination of chromatography on a reversed-phase column with long carbon chains (C30) and with the mobile phase containing an ion pair reagent (ammonium heptafluorobutyrate). The relationship between the amount of applied standard TTX and its peak area on the mass chromatogram (m/z 320) shows good linearity over a range of 50 to 1000 pmol. The detection limit of TTX in the selected ion monitoring (SIM) mode is estimated to be 0.7 pmol, with a signal to noise ratio of 2:1.

Please note: These procedures have been applied to analysis of tissue samples. Further research is needed to develop and standardize the protocols for other sample types.

Source: Shoji, Y., Yotsu-Yamashita, M., Miyazawa, T., and Yasumoto, T. 2001. "Electrospray Ionization Mass Spectrometry of Tetrodotoxin and its Analogs: Liquid Chromatography/Mass Spectrometry, Tandem Mass Spectrometry, and Liquid Chromatography/Tandem Mass Spectrometry." *Analytical Biochemistry*. 290(1): 10–17.

8.2.27 Lateral Flow Immunoassay Kits

These procedures should be used for **presumptive** analysis of botulinum neurotoxins Types A and B and ricin in aerosol samples.

Analyte(s)	Agent Category	CASRN
Botulinum neurotoxins (Types A, B)	Protein	NA
Ricin	Protein	9009-86-3

Test strips are self-contained, qualitative assays for screening environmental samples for the presence of botulinum toxin and ricin. After the sample is collected, it is transferred onto the test strip where dye-labeled antibodies detect trace amounts of the contaminant, as indicated by the presence of two bands in the test result window. After 15 minutes, the results are read visually. Botulinum neurotoxin Type A can be detected at 5 mg/L and Type B at 4 mg/L, 33% of the time. Ricin toxin can be detected at 20 mg/L, with no cross-reactivity certain substances (i.e., lectin from soybeans).

An alternative lateral flow immunochromatographic device also can be used. This device uses two antibodies in combination to specifically detect target antigen in solution. When a sufficient amount of target antigen is present, the colloidal gold label accumulates in the sample window on a test strip, forming a visible reddish-brown colored line. The presence of two bands indicates a positive reading. Botulinum neurotoxin Type A can be detected at 0.01 mg/L and Type B at 0.5 mg/L, with no false negatives detected when interferents are present. Ricin toxin can be detected at 0.035 mg/L, with 88% accuracy.

These two lateral flow immunoassay kits have been evaluated by the U.S. EPA Environmental Technology Verification (ETV) Program (http://www.epa.gov/etv/pdfs/vrvs/01_vr_badd.pdf and http://www.epa.gov/etv/pdfs/vrvs/01_vr_biothreat.pdf) for the detection of botulinum neurotoxins Types A and B and ricin. Information regarding the evaluation of test strips can be accessed at these sites.

Please note: These procedures have been applied to analysis of toxin in buffer or water. Further research is needed to develop and standardize the protocols for other sample types.

Source: Environmental Technology Verification. 2006. <http://www.epa.gov/etv/>

Section 9.0: Conclusions

Methods listed in Appendix A (chemical methods), Appendix B (radiochemical methods), Appendix C (pathogen methods), and Appendix D (biotoxin methods) are recommended for use in assessment of the extent of contamination and the effectiveness of decontamination following a homeland security event.

The primary objective of this document is not necessarily to identify the “best” analytical methods, but rather to identify appropriate methods that represent a balance between providing existing, documented determinative techniques and providing consistent and valid analytical results. The method selected for each analyte/sample type combination was deemed the most general, appropriate, and broadly applicable of available methods. This is a living guidance document for use by EPA and EPA-contracted laboratories tasked with analysis of environmental samples following a homeland security event. Recommended methods are subject to change based on procedure testing and advances in technology.

Any questions concerning the information in this document should be directed to the appropriate point(s) of contact listed in Section 4.

Appendix A: Selected Chemical Methods

Appendix A: Selected Chemical Methods

Analyte(s)	CAS RN	Determinative Technique	Method Type	Solid Samples	Non-aqueous Liquid/Organic Solid Samples ¹	Aqueous Liquid Samples	Drinking Water Samples	Air Samples
Acrylamide	79-06-1	HPLC	Sample Prep	Water extraction	Water extraction	8316 (EPA SW-846)	8316 (EPA SW-846)	PV2004 (OSHA)
			Determinative	8316 (EPA SW-846)	8316 (EPA SW-846)			
Acrylonitrile	107-13-1	HPLC	Sample Prep	Water extraction	Water extraction	8316 (EPA SW-846)	8316 (EPA SW-846)	PV2004 (OSHA)
			Determinative	8316 (EPA SW-846)	8316 (EPA SW-846)			
Aldicarb (Temik)	116-06-3	HPLC	Sample Prep	8318A (EPA SW-846)	8318A (EPA SW-846)	8318A (EPA SW-846)	531.2 (EPA OW)	5601 (NIOSH)
			Determinative					
Allyl alcohol	107-18-6	GC-MS	Sample Prep	5035A (EPA SW-846)	3585 (EPA SW-846)	5030C (EPA SW-846)	5030C (EPA SW-846)	TO-15 ² (EPA ORD)
			Determinative	8260C (EPA SW-846)	8260C (EPA SW-846)	8260C (EPA SW-846)	8260C (EPA SW-846)	
2-Amino-4,6-dinitrotoluene (2-Am-DNT)	35572-78-2	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	
4-Amino-2,6-dinitrotoluene (4-Am-DNT)	19406-51-0	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	
4-Aminopyridine	504-24-5	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	
Ammonia	7664-41-7	Spectrophotometry	Sample Prep	Not of concern	Not of concern	4500- NH ₃ B (SM)	350.1 (EPA OW)	6015 (NIOSH)
			Determinative			4500- NH ₃ G (SM)		
Ammonium metavanadate (analyze for total vanadium)	7803-55-6	ICP-MS / ICP-AES	Sample Prep	3050B (EPA SW-846)	3031 (EPA SW-846)	200.8 (EPA OW)	200.8 (EPA OW)	IO-3.1 (EPA ORD)
			Determinative	6020A/6010C (EPA SW-846)	6020A/6010C (EPA SW-846)			IO-3.4/IO-3.5 (EPA ORD)
Arsenic, Total	7440-38-2	ICP-MS / ICP-AES	Sample Prep	3050B (EPA SW-846)	3031 (EPA SW-846)	200.8 (EPA OW)	200.8 (EPA OW)	IO-3.1 (EPA ORD)
			Determinative	6020A/6010C (EPA SW-846)	6020A/6010C (EPA SW-846)			IO-3.4/IO-3.5 (EPA ORD)
Arsenic trioxide (analyze for total arsenic)	1327-53-3	ICP-MS / ICP-AES	Sample Prep	3050B (EPA SW-846)	3031 (EPA SW-846)	200.8 (EPA OW)	200.8 (EPA OW)	IO-3.1 (EPA ORD)
			Determinative	6020A/6010C (EPA SW-846)	6020A/6010C (EPA SW-846)			IO-3.4/IO-3.5 (EPA ORD)

Analyte(s)	CAS RN	Determinative Technique	Method Type	Solid Samples	Non-aqueous Liquid/Organic Solid Samples ¹	Aqueous Liquid Samples	Drinking Water Samples	Air Samples
Arsine	7784-42-1	GFAA	Sample Prep	3050B (EPA SW-846)	Not of concern	200.8 (EPA OW)	200.8 (EPA OW)	6001 (NIOSH)
			Determinative	7010 (EPA SW-846)				
Asbestos	1332-21-4	TEM	Sample Prep	D5755-03 (soft surfaces- microvac) or D6480-99 (hard surfaces-wipes) (ASTM)	Not of concern	Not of concern	Not of concern	10312:1995 (ISO)
			Determinative					
Boron trifluoride	7637-07-2	ISE	Sample Prep	Not of concern	Not of concern	Not of concern	Not of concern	ID - 216SG (OSHA)
			Determinative					
Brodifacoum	56073-10-0	HPLC	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	Not of concern
			Determinative	8321B (EPA SW-846)	8321B (EPA SW-846)	8321B (EPA SW-846)	8321B (EPA SW-846)	
Bromadiolone	28772-56-7	HPLC	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	Not of concern
			Determinative	8321B (EPA SW-846)	8321B (EPA SW-846)	8321B (EPA SW-846)	8321B (EPA SW-846)	
Carbofuran (Furadan)	1563-66-2	HPLC	Sample Prep	8318A (EPA SW-846)	8318A (EPA SW-846)	8318A (EPA SW-846)	531.2 (EPA OW)	5601 (NIOSH)
			Determinative					
Carbon disulfide	75-15-0	GC-MS	Sample Prep	5035A (EPA SW-846)	3585 (EPA SW-846)	5030C (EPA SW-846)	524.2 (EPA OW)	TO-15 (EPA ORD)
			Determinative	8260C (EPA SW-846)	8260C (EPA SW-846)	8260C (EPA SW-846)		
Chlorfenvinphos	470-90-6	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Chlorine	7782-50-5	Spectrophotometry	Sample Prep	Not of concern	Not of concern	4500-Cl G (SM)	4500-Cl G (SM)	Analyst, Vol. 124, 1999, pp. 1853-1857
			Determinative					4500-Cl G (SM)
2-Chloroethanol	107-07-3	GC-MS / GC-FID	Sample Prep	5035A (EPA SW-846)	3585 (EPA SW-846)	5030C (EPA SW-846)	5030C (EPA SW-846)	2513 (NIOSH)
			Determinative	8260C (EPA SW-846)	8260C (EPA SW-846)	8260C (EPA SW-846)	8260C (EPA SW-846)	
3-Chloro-1,2-propanediol	96-24-2	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A ⁴ (EPA ORD)
			Determinative	8270D ³ (EPA SW-846)	8270D ³ (EPA SW-846)	8270D ³ (EPA SW-846)	8270D ³ (EPA SW-846)	

Analyte(s)	CAS RN	Determinative Technique	Method Type	Solid Samples	Non-aqueous Liquid/Organic Solid Samples ¹	Aqueous Liquid Samples	Drinking Water Samples	Air Samples
Chloropicrin	76-06-2	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	551.1 (EPA OW)	551.1 (EPA OW)	PV2103 (OSHA)
			Determinative	8270D ⁵ (EPA SW-846)	8270D ⁵ (EPA SW-846)			
Chlorosarin	1445-76-7	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A ⁴ (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Chlorosoman	7040-57-5	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A ⁴ (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
2-Chlorovinylarsonous acid (CVAA) (degradation product of Lewisite)	85090-33-1	ICP-MS / ICP-AES	Sample Prep	3050B (EPA SW-846)	3031 (EPA SW-846)	200.8 (EPA OW)	200.8 (EPA OW)	IO-3.1 (EPA ORD)
			Determinative	6020A/6010C (EPA SW-846)	6020A/6010C (EPA SW-846)			IO-3.4/IO-3.5 (EPA ORD)
Chlorpyrifos	2921-88-2	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Crimidine	535-89-7	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	Not of concern
			Determinative	8270D ⁶ (EPA SW-846)	8270D ⁶ (EPA SW-846)	8270D ⁶ (EPA SW-846)	8270D ⁶ (EPA SW-846)	
Cyanide, Total	57-12-5	Spectrophotometry	Sample Prep	ILM05.3 CN (EPA CLP)	Not of concern	ILM05.3 CN (EPA CLP)	335.4 (EPA OW)	6010 (NIOSH)
			Determinative					
Cyanogen chloride	506-77-4	GC-MS	Sample Prep	5035A (EPA SW-846)	3585 (EPA SW-846)	5030C (EPA SW-846)	5030C (EPA SW-846)	TO-15 (EPA ORD)
			Determinative	8260C (EPA SW-846)	8260C (EPA SW-846)	8260C (EPA SW-846)	8260C (EPA SW-846)	
Cyclohexyl sarin (GF)	329-99-7	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
1,2-Dichloroethane (degradation product of HD)	107-06-2	GC-MS	Sample Prep	5035A (EPA SW-846)	3585 (EPA SW-846)	5030C (EPA SW-846)	524.2 (EPA OW)	TO-15 (EPA ORD)
			Determinative	8260C (EPA SW-846)	8260C (EPA SW-846)	8260C (EPA SW-846)		
Dichlorvos	62-73-7	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	525.2 (EPA OW)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)		

Analyte(s)	CAS RN	Determinative Technique	Method Type	Solid Samples	Non-aqueous Liquid/Organic Solid Samples ¹	Aqueous Liquid Samples	Drinking Water Samples	Air Samples
Dicrotophos	141-66-2	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Diesel Range Organics	NA	GC-FID	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	Not of concern
			Determinative	8015C (EPA SW-846)	8015C (EPA SW-846)	8015C (EPA SW-846)	8015C (EPA SW-846)	
Diisopropyl methylphosphonate (DIMP) (degradation product of GB)	1445-75-6	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-10A ⁴ (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Dimethylphosphite	868-85-9	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Dimethylphosphoramidic acid (degradation product of GA)	33876-51-6	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-10A ⁸ (EPA ORD)
			Determinative	8321B ⁷ (EPA SW-846)	8321B ⁷ (EPA SW-846)	8321B ⁷ (EPA SW-846)	8321B ⁷ (EPA SW-846)	
3,5-Dinitroaniline (3,5-DNA)	618-87-1	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	
1,3-Dinitrobenzene (1,3-DNB)	99-65-0	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	
2,4-Dinitrotoluene (2,4-DNT)	121-14-2	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	
2,6-Dinitrotoluene (2,6-DNT)	606-20-2	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	
Diphacinone	82-66-6	HPLC	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	Not of concern
			Determinative	8321B (EPA SW-846)	8321B (EPA SW-846)	8321B (EPA SW-846)	8321B (EPA SW-846)	
1,4-Dithiane (degradation product of HD)	505-29-3	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	Not of concern
			Determinative	8270D ⁹ (EPA SW-846)	8270D ⁹ (EPA SW-846)	8270D ⁹ (EPA SW-846)	8270D ⁹ (EPA SW-846)	

Analyte(s)	CAS RN	Determinative Technique	Method Type	Solid Samples	Non-aqueous Liquid/Organic Solid Samples ¹	Aqueous Liquid Samples	Drinking Water Samples	Air Samples
EA2192 [Diisopropylaminoethyl methylthiophosphonate] (hydrolysis product of VX)	73207-98-4	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-10A ⁸ (EPA ORD)
			Determinative	8321B ⁷ (EPA SW-846)	8321B ⁷ (EPA SW-846)	8321B ⁷ (EPA SW-846)	8321B ⁷ (EPA SW-846)	
Ethyl methylphosphonic acid (EMPA) (degradation product of VX)	1832-53-7	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-10A ⁸ (EPA ORD)
			Determinative	8321B ⁷ (EPA SW-846)	8321B ⁷ (EPA SW-846)	8321B ⁷ (EPA SW-846)	8321B ⁷ (EPA SW-846)	
Ethylchloroarsine (ED)	598-14-1	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-15 (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
N-Ethyldiethanolamine (EDEA) (degradation product of HN-1)	139-87-7	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8321B (EPA SW-846)	8321B (EPA SW-846)	8321B (EPA SW-846)	8321B (EPA SW-846)	
Ethylene oxide	75-21-8	GC-MS	Sample Prep	5035A (EPA SW-846)	3585 (EPA SW-846)	5030C (EPA SW-846)	5030C (EPA SW-846)	TO-15 (EPA ORD)
			Determinative	8260C (EPA SW-846)	8260C (EPA SW-846)	8260C (EPA SW-846)	8260C (EPA SW-846)	
Fenamiphos	22224-92-6	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	525.2 (EPA OW)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)		
Fentanyl	437-38-7	HPLC	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8321B (EPA SW-846)	8321B (EPA SW-846)	8321B (EPA SW-846)	8321B (EPA SW-846)	
Fluoride	16984-48-8	IC	Sample Prep	Not of concern	Not of concern	300.1, Rev 1.0 (EPA OW)	300.1, Rev 1.0 (EPA OW)	Not of concern
			Determinative					
Fluoroacetic acid and fluoroacetate salts (analyze for fluoroacetate ion)	NA	IC	Sample Prep	Analytical Letters, 1994, 27 (14): 2703-2718	Analytical Letters, 1994, 27 (14): 2703-2718	300.1, Rev 1.0 (EPA OW)	300.1, Rev 1.0 (EPA OW)	S301-1 (NIOSH)
			Determinative	300.1, Rev 1.0 (EPA OW)	300.1, Rev 1.0 (EPA OW)			300.1, Rev 1.0 (EPA OW)
Formaldehyde	50-00-0	HPLC	Sample Prep	8315A (EPA SW-846)	Not of concern	8315A (EPA SW-846)	8315A (EPA SW-846)	2016 (NIOSH)
			Determinative					
Gasoline Range Organics	NA	GC-FID	Sample Prep	5035A (EPA SW-846)	3585 (EPA SW-846)	5030C (EPA SW-846)	5030C (EPA SW-846)	Not of concern
			Determinative	8015C (EPA SW-846)	8015C (EPA SW-846)	8015C (EPA SW-846)	8015C (EPA SW-846)	

Analyte(s)	CAS RN	Determinative Technique	Method Type	Solid Samples	Non-aqueous Liquid/Organic Solid Samples ¹	Aqueous Liquid Samples	Drinking Water Samples	Air Samples
Hexahydro-1,3,5-trinitro-1,3,5-triazine (RDX)	121-82-4	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	
Hexamethylenetriperoxidediamine (HMTD)	283-66-9	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	
Hydrogen bromide	10035-10-6	IC	Sample Prep	Not of concern	Not of concern	Not of concern	Not of concern	7903 (NIOSH)
			Determinative					
Hydrogen chloride	7647-01-0	IC	Sample Prep	Not of concern	Not of concern	Not of concern	Not of concern	7903 (NIOSH)
			Determinative					
Hydrogen cyanide	74-90-8	Spectrophotometry	Sample Prep	Not of concern	Not of concern	ILM05.3 CN (EPA CLP)	335.4 (EPA OW)	6010 (NIOSH)
			Determinative					
Hydrogen fluoride	7664-39-3	IC	Sample Prep	Not of concern	Not of concern	Not of concern	Not of concern	7903 ¹⁰ (NIOSH)
			Determinative					
Hydrogen sulfide	7783-06-4	IC	Sample Prep	Not of concern	Not of concern	Not of concern	Not of concern	6013 (NIOSH)
			Determinative					
Isopropyl methylphosphonic acid (IMPA) (degradation product of GB)	1832-54-8	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-10A ⁸ (EPA ORD)
			Determinative	8321B ⁷ (EPA SW-846)	8321B ⁷ (EPA SW-846)	8321B ⁷ (EPA SW-846)	8321B ⁷ (EPA SW-846)	
Kerosene	64742-81-0	GC-FID	Sample Prep	5035A (EPA SW-846)	3585 (EPA SW-846)	5030C (EPA SW-846)	5030C (EPA SW-846)	Not of concern
			Determinative	8015C (EPA SW-846)	8015C (EPA SW-846)	8015C (EPA SW-846)	8015C (EPA SW-846)	
Lewsite 1 (L-1) ¹¹ [2-chlorovinyl]dichloroarsine (analyze for total arsenic)	541-25-3	ICP-MS / ICP-AES	Sample Prep	3050B (EPA SW-846)	3031 (EPA SW-846)	200.8 (EPA OW)	200.8 (EPA OW)	IO-3.1 (EPA ORD)
			Determinative	6020A/6010C (EPA SW-846)	6020A/6010C (EPA SW-846)			IO-3.4/IO-3.5 (EPA ORD)
Lewsite 2 (L-2) [bis(2-chlorovinyl)chloroarsine] (analyze for total arsenic)	40334-69-8	ICP-MS / ICP-AES	Sample Prep	3050B (EPA SW-846)	3031 (EPA SW-846)	200.8 (EPA OW)	200.8 (EPA OW)	IO-3.1 (EPA ORD)
			Determinative	6020A/6010C (EPA SW-846)	6020A/6010C (EPA SW-846)			IO-3.4/IO-3.5 (EPA ORD)
Lewsite 3 (L-3) [tris(2-chlorovinyl)arsine] (analyze for total arsenic)	40334-70-1	ICP-MS / ICP-AES	Sample Prep	3050B (EPA SW-846)	3031 (EPA SW-846)	200.8 (EPA OW)	200.8 (EPA OW)	IO-3.1 (EPA ORD)
			Determinative	6020A/6010C (EPA SW-846)	6020A/6010C (EPA SW-846)			IO-3.4/IO-3.5 (EPA ORD)

Analyte(s)	CAS RN	Determinative Technique	Method Type	Solid Samples	Non-aqueous Liquid/Organic Solid Samples ¹	Aqueous Liquid Samples	Drinking Water Samples	Air Samples
Lewisite oxide (degradation product of Lewisite)	1306-02-1	ICP-MS / ICP-AES	Sample Prep	3050B (EPA SW-846)	3031 (EPA SW-846)	200.8 (EPA OW)	200.8 (EPA OW)	IO-3.1 (EPA ORD)
			Determinative	6020A/6010C (EPA SW-846)	6020A/6010C (EPA SW-846)			IO-3.4/IO-3.5 (EPA ORD)
Mercury, Total	7439-97-6	Thermal Decomposition AA / CVAA / CVAFS	Sample Prep	7473 ¹² (EPA SW-846)	Not of concern	7473 ¹² (EPA SW-846)	245.2 (EPA OW)	IO-5 (EPA ORD)
			Determinative					
Methamidophos	10265-92-6	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-10A ¹³ (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Methomyl	16752-77-5	HPLC	Sample Prep	8318A (EPA SW-846)	8318A (EPA SW-846)	8318A (EPA SW-846)	531.2 (EPA OW)	5601 (NIOSH)
			Determinative					
Methoxyethylmercuric acetate (analyze for total mercury)	151-38-2	Thermal Decomposition AA / CVAA / CVAFS	Sample Prep	7473 ¹² (EPA SW-846)	Not of concern	7473 ¹² (EPA SW-846)	245.2 (EPA OW)	IO-5 (EPA ORD)
			Determinative					
Methyl acrylonitrile	126-98-7	HPLC-UV	Sample Prep	Water extraction	Water extraction	8316 (EPA SW-846)	8316 (EPA SW-846)	PV2004 (OSHA)
			Determinative	8316 (EPA SW-846)	8316 (EPA SW-846)			
Methyl fluoroacetate (analyze for fluoroacetate ion)	453-18-9	IC	Sample Prep	Analytical Letters, 1994, 27 (14): 2703-2718	Analytical Letters, 1994, 27 (14): 2703-2718	300.1, Rev 1.0 (EPA OW)	300.1, Rev 1.0 (EPA OW)	S301-1 (NIOSH)
			Determinative	300.1, Rev 1.0 (EPA OW)	300.1, Rev 1.0 (EPA OW)			300.1, Rev 1.0 (EPA OW)
Methyl hydrazine	60-34-4	GC-MS / Spectrophotometry	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	3510 (NIOSH)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Methyl isocyanate	624-83-9	HPLC	Sample Prep	Not of concern	Not of concern	Not of concern	Not of concern	OSHA 54
			Determinative					
Methyl parathion	298-00-0	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Methyl-2,4,6-trinitrophenylnitramine (Tetryl)	479-45-8	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	

Analyte(s)	CAS RN	Determinative Technique	Method Type	Solid Samples	Non-aqueous Liquid/Organic Solid Samples ¹	Aqueous Liquid Samples	Drinking Water Samples	Air Samples
Methylamine	74-89-5	HPLC	Sample Prep	Not of concern	Not of concern	Not of concern	Not of concern	OSHA 40
			Determinative					
N-Methyldiethanolamine (MDEA) (degradation product of HN-2)	105-59-9	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8321B (EPA SW-846)	8321B (EPA SW-846)	8321B (EPA SW-846)	8321B (EPA SW-846)	
1-Methylethyl ester ethylphosphonofluoridic acid (GE)	1189-87-3	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A ⁴ (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Methylphosphonic acid (MPA) (degradation product of VX, GB, & GD)	993-13-5	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-10A ⁸ (EPA ORD)
			Determinative	8321B ⁷ (EPA SW-846)	8321B ⁷ (EPA SW-846)	8321B ⁷ (EPA SW-846)	8321B ⁷ (EPA SW-846)	
Mevinphos	7786-34-7	GC-MS	Sample Prep	3545A (SW-846)	3580A (SW-846)	3535A (SW-846)	525.2 (OW)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)		
Mustard, nitrogen (HN-1) [bis(2-chloroethyl)ethylamine]	538-07-8	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Mustard, nitrogen (HN-2) [2,2'-dichloro-N-methyldiethylamine N,N-bis(2-chloroethyl)methylamine]	51-75-2	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Mustard, nitrogen (HN-3) [tris(2-chloroethyl)amine]	555-77-1	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Mustard, sulfur / Mustard gas (HD)	505-60-2	GC-MS	Sample Prep	3571 (EPA SW-846)	3571 (EPA SW-846)	3571 (EPA SW-846)	3571 (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D ¹⁴ (EPA SW-846)	8270D ¹⁴ (EPA SW-846)	8270D ¹⁴ (EPA SW-846)	8270D ¹⁴ (EPA SW-846)	
Nicotine sulfate	54-11-5	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	Not of concern
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Nitrobenzene (NB)	98-95-3	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	

Analyte(s)	CAS RN	Determinative Technique	Method Type	Solid Samples	Non-aqueous Liquid/Organic Solid Samples ¹	Aqueous Liquid Samples	Drinking Water Samples	Air Samples
Nitroglycerin (NG)	55-63-0	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	
2-Nitrotoluene (2-NT)	88-72-2	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	
3-Nitrotoluene (3-NT)	99-08-1	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	
4-Nitrotoluene (4-NT)	99-99-0	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	
Octahydro-1,3,5,7-tetranitro-1,3,5,7-tetrazocine (HMX)	2691-41-0	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	
Organophosphate pesticides, NOS	NA	GC-MS / GC-NPD / GC-FPD	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	614 (EPA OW)	507 (EPA OW)	5600 (NIOSH)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)			
Osmium tetroxide (analyze for total osmium)	20816-12-0	ICP-AES / GFAA	Sample Prep	3050B (EPA SW-846)	Not of concern	252.2 (EPA OW)	252.2 (EPA OW)	IO-3.1 (EPA ORD)
			Determinative	6010C (EPA SW-846)				IO-3.4 (EPA ORD)
Oxamyl	23135-22-0	HPLC	Sample Prep	8318A (EPA SW-846)	8318A (EPA SW-846)	8318A (EPA SW-846)	531.2 (EPA OW)	5601 (NIOSH)
			Determinative					
Paraquat	4685-14-7	HPLC	Sample Prep	Not of concern	Not of concern	549.2 (EPA OW)	549.2 (EPA OW)	Not of concern
			Determinative					
Parathion	56-38-2	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Pentaerythritol tetranitrate (PETN)	78-11-5	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	

Analyte(s)	CAS RN	Determinative Technique	Method Type	Solid Samples	Non-aqueous Liquid/Organic Solid Samples ¹	Aqueous Liquid Samples	Drinking Water Samples	Air Samples
Perfluoroisobutylene (PFIB)	382-21-8	GC-MS / GC-NPD	Sample Prep	Not of concern	Not of concern	Not of concern	Not of concern	OSHA 61 ¹⁵
			Determinative					
Phencyclidine	77-10-1	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Phenol	108-95-2	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Phorate	298-02-2	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Phosgene	75-44-5	GC-NPD	Sample Prep	Not of concern	Not of concern	Not of concern	Not of concern	OSHA 61
			Determinative					
Phosphamidon	13171-21-6	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Phosphine	7803-51-2	Spectrophotometry	Sample Prep	Not of concern	Not of concern	Not of concern	Not of concern	6002 (NIOSH)
			Determinative					
Phosphorus trichloride	7719-12-2	Spectrophotometry	Sample Prep	Not of concern	Not of concern	Not of concern	Not of concern	6402 (NIOSH)
			Determinative					
Pinacolyl methyl phosphonic acid (PMPA) (degradation product of GD)	616-52-4	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-10A ⁸ (EPA ORD)
			Determinative	8321B ⁷ (EPA SW-846)	8321B ⁷ (EPA SW-846)	8321B ⁷ (EPA SW-846)	8321B ⁷ (EPA SW-846)	
Propylene oxide	75-56-9	GC-MS / GC-FID	Sample Prep	5035A (EPA SW-846)	3585 (EPA SW-846)	5030C (EPA SW-846)	5030C (EPA SW-846)	1612 (NIOSH)
			Determinative	8260C (EPA SW-846)	8260C (EPA SW-846)	8260C (EPA SW-846)	8260C (EPA SW-846)	
R 33 (VR) [methylphosphonothioic acid, S-[2-(diethylamino)ethyl] O-2-methylpropyl ester]	159939-87-4	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	

Analyte(s)	CAS RN	Determinative Technique	Method Type	Solid Samples	Non-aqueous Liquid/Organic Solid Samples ¹	Aqueous Liquid Samples	Drinking Water Samples	Air Samples
Sarin (GB)	107-44-8	GC-MS	Sample Prep	3571 (EPA SW-846)	3571 (EPA SW-846)	3571 (EPA SW-846)	3571 (EPA SW-846)	TO-10A ⁴ (EPA ORD)
			Determinative	8270D ¹⁴ (EPA SW-846)	8270D ¹⁴ (EPA SW-846)	8270D ¹⁴ (EPA SW-846)	8270D ¹⁴ (EPA SW-846)	
Semivolatile Organic Compounds, NOS	NA	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	525.2 (EPA OW)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)		
Sodium arsenite (analyze for total arsenic)	7784-46-5	ICP-MS / ICP-AES	Sample Prep	3050B (EPA SW-846)	3031 (EPA SW-846)	200.8 (EPA OW)	200.8 (EPA OW)	IO-3.1 (EPA ORD)
			Determinative	6020A/6010C (EPA SW-846)	6020A/6010C (EPA SW-846)			IO-3.4/IO-3.5 (EPA ORD)
Sodium azide (analyze for hydrazoic acid)	26628-22-8	IC	Sample Prep	J. of Forensic Science, 1998, 43(1): 200-202 ¹⁶	3580A ¹⁶ (EPA SW-846)	J. of Forensic Science, 1998, 43(1): 200-202 ¹⁶	J. of Forensic Science, 1998, 43(1): 200-202 ¹⁶	ID-211 (OHSA)
			Determinative	300.1, Rev 1.0 ¹⁷ (EPA OW)	300.1, Rev 1.0 ¹⁷ (EPA OW)	300.1, Rev 1.0 ¹⁷ (EPA OW)	300.1, Rev 1.0 ¹⁷ (EPA OW)	
Soman (GD)	96-64-0	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-10A ⁴ (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Strychnine	57-24-9	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	Not of concern
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Sulfur dioxide	7446-09-5	IC	Sample Prep	Not of concern	Not of concern	Not of concern	Not of concern	6004 (NIOSH)
			Determinative					
Sulfur trioxide	7446-11-9	Titrimetry	Sample Prep	Not of concern	Not of concern	Not of concern	Not of concern	Method 8 (EPA OAQPS)
			Determinative					
Tabun (GA)	77-81-6	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Tear gas (CS) [chlorobenzylidene malonitrile]	2698-41-1	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Tetraethyl pyrophosphate	107-49-3	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	

Analyte(s)	CAS RN	Determinative Technique	Method Type	Solid Samples	Non-aqueous Liquid/Organic Solid Samples ¹	Aqueous Liquid Samples	Drinking Water Samples	Air Samples
Tetramethylenedisulfotetramine	80-12-6	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
Thallium sulfate (analyze for total thallium)	10031-59-1	ICP-MS / ICP-AES	Sample Prep	3050B (EPA SW-846)	3031 (EPA SW-846)	200.8 (EPA OW)	200.8 (EPA OW)	IO-3.1 (EPA ORD)
			Determinative	6020A/6010C (EPA SW-846)	6020A/6010C (EPA SW-846)			IO-3.4/IO-3.5 (EPA ORD)
Thiodiglycol (TDG) (degradation product of HD)	111-48-8	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
1,4-Thioxane (degradation product of HD)	15980-15-1	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	Not of concern
			Determinative	8270D ⁹ (EPA SW-846)	8270D ⁹ (EPA SW-846)	8270D ⁹ (EPA SW-846)	8270D ⁹ (EPA SW-846)	
Titanium tetrachloride (analyze for total Titanium)	7550-45-0	ICP-MS / ICP-AES	Sample Prep	3050B (EPA SW-846)	Not of concern	Not of concern	Not of concern	Not of concern
			Determinative	6020A/6010C (EPA SW-846)				
Triacetone triperoxide (TATP)	17088-37-8	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	
Triethanolamine (TEA) (degradation product of HN-3)	102-71-6	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8321B (EPA SW-846)	8321B (EPA SW-846)	8321B (EPA SW-846)	8321B (EPA SW-846)	
Trimethyl phosphite	121-45-9	GC-MS	Sample Prep	3545A (EPA SW-846)	3580A (EPA SW-846)	3535A (EPA SW-846)	3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D ⁵ (EPA SW-846)	8270D ⁵ (EPA SW-846)	8270D ⁵ (EPA SW-846)	8270D ⁵ (EPA SW-846)	
1,3,5-Trinitrobenzene (1,3,5-TNB)	99-35-4	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	
2,4,6-Trinitrotoluene (2,4,6-TNT)	118-96-7	HPLC	Sample Prep	8330B (EPA SW-846)	8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	3535A/8330B (EPA SW-846)	Not of concern
			Determinative			8330B (EPA SW-846)	8330B (EPA SW-846)	
Vanadium pentoxide (analyze for total vanadium)	1314-62-1	ICP-MS / ICP-AES	Sample Prep	3050B (EPA SW-846)	3031 (EPA SW-846)	200.8 (EPA OW)	200.8 (EPA OW)	IO-3.1 (EPA ORD)
			Determinative	6020A/6010C (EPA SW-846)	6020A/6010C (EPA SW-846)			IO-3.4/IO-3.5 (EPA ORD)

Analyte(s)	CAS RN	Determinative Technique	Method Type	Solid Samples	Non-aqueous Liquid/Organic Solid Samples ¹	Aqueous Liquid Samples	Drinking Water Samples	Air Samples
VE [phosphonothioic acid, ethyl-, S-(2-(diethylamino)ethyl) O-ethyl ester]	21738-25-0	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
VG [phosphonothioic acid, S-(2-(diethylamino)ethyl) O,O-diethyl ester]	78-53-5	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
VM [phosphonothioic acid, methyl-, S-(2-(diethylamino)ethyl) O-ethyl ester]	21770-86-5	GC-MS	Sample Prep	3545A/3541 (EPA SW-846)	3580A (EPA SW-846)	3520C/3535A (EPA SW-846)	3520C/3535A (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	8270D (EPA SW-846)	
VX [O-ethyl-S-(2-diisopropylaminoethyl)methyl-phosphonothiolate]	50782-69-9	GC-MS	Sample Prep	3571 (EPA SW-846)	3571 (EPA SW-846)	3571 (EPA SW-846)	3571 (EPA SW-846)	TO-10A (EPA ORD)
			Determinative	8270D ¹⁴ (EPA SW-846)	8270D ¹⁴ (EPA SW-846)	8270D ¹⁴ (EPA SW-846)	8270D ¹⁴ (EPA SW-846)	
White phosphorus	12185-10-3	GC-NPD / GC-FPD	Sample Prep	7580 (EPA SW-846)	7580 (EPA SW-846)	7580 (EPA SW-846)	7580 (EPA SW-846)	7905 (NIOSH)
			Determinative					

Footnotes

- ¹ An organic solid sample is a solid that completely dissolves in an organic solvent and leaves no solid residue.
- ² If problems occur when using this method, it is recommended that TO-10A be used.
- ³ For this analyte, SW-846 Method 8270D must be modified to include a derivatization step.
- ⁴ If problems occur when using this method, it is recommended that the canister Method TO-15 be used.
- ⁵ If problems occur with analyses, lower the injection temperature.
- ⁶ If problems occur when using this method, it is recommended that SW-846 Method 8321B be used as the determinative method. Sample preparation methods should remain the same.
- ⁷ LC-MS (electrospray) procedures are preferred for these analytes; however, if this technique is not available to the laboratory, GC-MS procedures using derivatization based on SW-846 Method 8270D may be used. Sample preparation methods should remain the same. Both electrospray LC-MS and GC/MS derivatization procedures are currently under development.
- ⁸ For this analyte, Method TO-10A must be modified to include a derivatization step.
- ⁹ If problems occur when using this method, it is recommended that SW-846 Method 8260C and appropriate corresponding sample preparation procedures (i.e., 5035A for solid samples, 3585 for non-aqueous liquid/organic solid samples, and 5030C for aqueous liquid and drinking water samples) be used.
- ¹⁰ If problems occur when using this method, it is recommended that NIOSH Method 7906 be used.
- ¹¹ Laboratory testing is currently under way for speciation of Lewisite 1 using GC-MS techniques.
- ¹² If equipment is not available, use CVAA Methods 7471B (EPA SW-846) for solid samples and 7470A (EPA SW-846) for aqueous liquid samples.
- ¹³ If problems occur when using this method, it is recommended that NIOSH Method 5600 be used.
- ¹⁴ For this analyte, refer to EPA SW-846 Method 8271 for GC-MS conditions.
- ¹⁵ If problems occur when using this method, it is recommended that a method based on the following journal article be used: J. Chrom. A, 1098: (2005) 156-165.
- ¹⁶ Water extraction, filtration, and acidification steps from the Journal of Forensic Science, 1998. 43(1):200-202 should be used for the preparation of solid samples. Filtration and acidification steps from this journal should be used for preparation of aqueous liquid and drinking water samples. The acidification step from the journal should be used with EPA SW-846 Method 3580A for preparation of non-aqueous liquid/organic solid samples.
- ¹⁷ If analyses are problematic, refer to column manufacturer for alternate conditions

Appendix B : Selected Radiochemical Methods

Appendix B: Selected Radiochemical Methods

Analyte Class		Determinative Technique	Drinking Water Samples		Aqueous and Liquid Phase Samples		Soil and Sediment Samples		Surface Wipes		Air Filters	
Gross Alpha		Alpha/Beta spectrometry	900.0 (EPA)		7110 B (SM)		AP1 (ORISE)		FRMAC, Vol 2, pg. 33		FRMAC, Vol 2, pg. 33	
Gross Beta		Alpha/Beta spectrometry	900.0 (EPA)		7110 B (SM)		AP1 (ORISE)		FRMAC, Vol 2, pg. 33		FRMAC, Vol 2, pg. 33	
Gamma		Gamma spectrometry	901.1 (EPA)		Ga-01-R (HASL-300)		Ga-01-R (HASL-300)		Ga-01-R (HASL-300)		Ga-01-R (HASL-300)	
Analyte(s)	CAS RN	Determinative Technique	Drinking Water Samples		Aqueous and Liquid Phase Samples		Soil and Sediment Samples		Surface Wipes		Air Filters	
			Qualitative Determination ¹	Confirmatory	Qualitative Determination ¹	Confirmatory	Qualitative Determination ¹	Confirmatory	Qualitative Determination ¹	Confirmatory	Qualitative Determination ¹	Confirmatory
Americium-241 ²	14596-10-2	Alpha/Gamma spectrometry	D3084 (ASTM)	Am-04-RC (HASL-300)	D3084 (ASTM)	Am-04-RC (HASL-300)	Am-02-RC (HASL-300)	Am-01-RC ³ (HASL-300)	D3084 (ASTM)	Am-04-RC (HASL-300)	D3084 (ASTM)	Am-04-RC (HASL-300)
Californium-252 ²	13981-17-4	Alpha spectrometry	D3084 (ASTM)	Am-04-RC (HASL-300)	D3084 (ASTM)	Am-04-RC (HASL-300)	D3084 (ASTM)	Am-01-RC ³ (HASL-300)	D3084 (ASTM)	Am-04-RC (HASL-300)	D3084 (ASTM)	Am-04-RC (HASL-300)
Cesium-137 ⁴	10045-97-3	Gamma spectrometry	901.1 (EPA)	901.1 (EPA)	7120 (SM)	7120 (SM)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)
Cobalt-60	10198-40-0	Gamma spectrometry	901.1 (EPA)	901.1 (EPA)	7120 (SM)	7120 (SM)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)
Curium-244 ²	13981-15-2	Alpha spectrometry	D3084 (ASTM)	Am-04-RC (HASL-300)	D3084 (ASTM)	Am-04-RC (HASL-300)	D3084 (ASTM)	Am-01-RC ³ (HASL-300)	D3084 (ASTM)	Am-04-RC (HASL-300)	D3084 (ASTM)	Am-04-RC (HASL-300)
Europium-154	15585-10-1	Gamma spectrometry	901.1 (EPA)	901.1 (EPA)	7120 (SM)	7120 (SM)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)
Iridium-192	14694-69-0	Gamma spectrometry	901.1 (EPA)	901.1 (EPA)	7120 (SM)	7120 (SM)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)
Plutonium-238 ²	13981-16-3	Alpha spectrometry	D3084 (ASTM)	EMSL-33 (EPA)	D3084 (ASTM)	EMSL-33 (EPA)	D3084 (ASTM)	EMSL-33 (EPA)	D3084 (ASTM)	EMSL-33 (EPA)	D3084 (ASTM)	EMSL-33 (EPA)
Plutonium-239 ²	15117-48-3	Alpha spectrometry	D3084 (ASTM)	EMSL-33 (EPA)	D3084 (ASTM)	EMSL-33 (EPA)	D3084 (ASTM)	EMSL-33 (EPA)	D3084 (ASTM)	EMSL-33 (EPA)	D3084 (ASTM)	EMSL-33 (EPA)
Polonium-210	13981-52-7	Alpha spectrometry	Po-02-RC (HASL-300)	Po-02-RC (HASL-300)	Po-02-RC (HASL-300)	Po-02-RC (HASL-300)	Po-02-RC (HASL-300)	Po-02-RC (HASL-300)	Method 111 (EPA)	Method 111 (EPA)	Method 111 (EPA)	Method 111 (EPA)

Analyte(s)	CAS RN	Determinative Technique	Drinking Water Samples		Aqueous and Liquid Phase Samples		Soil and Sediment Samples		Surface Wipes		Air Filters	
			Qualitative Determination ¹	Confirmatory	Qualitative Determination ¹	Confirmatory	Qualitative Determination ¹	Confirmatory	Qualitative Determination ¹	Confirmatory	Qualitative Determination ¹	Confirmatory
Radium-226 ⁴	13982-63-3	Alpha Counting	903.0 (EPA)	903.1 (EPA)	7500-Ra B (SM)	7500-Ra C (SM)	D3084 (ASTM)	EMSL-19 (EPA)	D3084 (ASTM)	EMSL-19 (EPA)	D3084 (ASTM)	EMSL-19 (EPA)
Ruthenium-103	13968-53-1	Gamma spectrometry	901.1 (EPA)	901.1 (EPA)	7120 (SM)	7120 (SM)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)
Ruthenium-106 ⁴	13967-48-1	Gamma spectrometry	901.1 (EPA)	901.1 (EPA)	7120 (SM)	7120 (SM)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)
Selenium-75	14265-71-5	Gamma spectrometry	901.1 (EPA)	901.1 (EPA)	7120 (SM)	7120 (SM)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)	Ga-01-R (HASL-300)
Strontium-90 ⁴	10098-97-2	Beta counting by low-background gas flow proportional detector	7500-Sr B (SM)	7500-Sr B (SM)	7500-Sr B (SM)	7500-Sr B (SM)	Sr-03-RC (HASL-300)	Sr-03-RC (HASL-300)	Sr-03-RC (HASL-300)	Sr-03-RC (HASL-300)	Sr-03-RC (HASL-300)	Sr-03-RC (HASL-300)
Uranium-238 ²	7440-61-1	Alpha Counting	908.0 (EPA)	D3972 (ASTM)	7500-U B (SM)	7500-U C (SM)	D3084 (ASTM)	EMSL-33 (EPA)	D3084 (ASTM)	EMSL-33 (EPA)	D3084 (ASTM)	EMSL-33 (EPA)

¹ In those cases where the same method is listed for qualitative determination and confirmatory analysis, qualitative determination can be performed by application of the method over a shorter count time than that used for confirmatory analysis.

² If it is suspected that the sample exists in refractory form (i.e., non-digestable or dissolvable material after normal digestion methods) or if there is a matrix interference problem, use ORISE Method AP11 for qualitative determination or confirmatory analysis.

³ In cases where only small sample volumes (≤ 100 g) are available, use HASL-300 Method Pu-12-RC.

⁴ Methods identified will measure decay product of these isotopes.

Appendix C: Selected Pathogen Methods

Appendix C: Pathogen Methods

Note: If viability determinations are needed (e.g., evaluation of the efficacy of disinfection), a viability-based procedure should be used. Those methods that include viability procedures are indicated as culture, tissue culture, animal infectivity, or embryonation of eggs under the analytical technique column of this appendix, and usually include biochemical and/or immunological differentiation of cultures.

Pathogen(s)	Viability	Analytical Technique	Analytical Method	Solid ¹ (soil, powder)	Particulate ² (swabs, wipes, HEPA)	Liquid/Water (filter, grab)	Drinking Water (filter, grab)	Aerosol (growth media, filter, liquid)
Bacteria								
<i>Bacillus anthracis</i> [Anthrax]	Yes	Culture	LRN This pathogen is included in the HHS/USDA select agent list. Select agent methods are generally provided and controlled by the LRN. In some cases, methods that apply to environmental samples may not be fully developed or validated. If analysis of this agent is required, contact the LRN at (404) 639-2790 for information on the LRN laboratory capable of receiving and processing the sample.					
	No	PCR/Immunoassay						
<i>Brucella</i> spp. (<i>B. abortus</i> , <i>B. melitensis</i> , <i>B. suis</i>)	Yes	Culture	LRN This pathogen is included in the HHS/USDA select agent list. Select agent methods are generally provided and controlled by the LRN. In some cases, methods that apply to environmental samples may not be fully developed or validated. If analysis of this agent is required, contact the LRN at (404) 639-2790 for information on the LRN laboratory capable of receiving and processing the sample.					
	No	PCR/Immunoassay						
<i>Burkholderia mallei</i> [Glanders]	Yes	Culture	LRN This pathogen is included in the HHS/USDA select agent list. Select agent methods are generally provided and controlled by the LRN. In some cases, methods that apply to environmental samples may not be fully developed or validated. If analysis of this agent is required, contact the LRN at (404) 639-2790 for information on the LRN laboratory capable of receiving and processing the sample.					
	No	PCR/Immunoassay						
<i>Burkholderia pseudomallei</i> [Meliodosis]	Yes	Culture	LRN This pathogen is included in the HHS/USDA select agent list. Select agent methods are generally provided and controlled by the LRN. In some cases, methods that apply to environmental samples may not be fully developed or validated. If analysis of this agent is required, contact the LRN at (404) 639-2790 for information on the LRN laboratory capable of receiving and processing the sample.					
	No	PCR/Immunoassay						
<i>Campylobacter jejuni</i>	Yes	Culture	SM 9260 G	Adapted from analytical method	Adapted from analytical method	SM 9260 G	SM 9260 G	Unlikely to be viable
	No	PCR	Molecular and Cellular Probes 20: 269-279	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method	Unlikely to be viable
<i>Chlamydomphila psittaci</i> (formerly known as <i>Chlamydia psittaci</i>)	Yes	Tissue culture	Journal of Clinical Microbiology 38: 1085-1093	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method
	No	PCR/Immunoassay	Journal of Clinical Microbiology 38: 1085-1093	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method
<i>Coxiella burnetii</i> [Q-fever]	Yes	Culture	LRN This pathogen is included in the HHS/USDA select agent list. Select agent methods are generally provided and controlled by the LRN. In some cases, methods that apply to environmental samples may not be fully developed or validated. If analysis of this agent is required, contact the LRN at (404) 639-2790 for information on the LRN laboratory capable of receiving and processing the sample.					
	No	PCR/Immunoassay						
<i>Escherichia coli</i> O157:H7 (<i>E. coli</i>)	Yes	Culture	SM 9260 F	Adapted from analytical method	Adapted from analytical method	SM 9260 F	SM 9260 F	Unlikely to be viable
	No	Immunoassay	SM 9260 F	Adapted from analytical method	Adapted from analytical method	SM 9260 F	SM 9260 F	Unlikely to be viable

Pathogen(s)	Viability	Analytical Technique	Analytical Method	Solid ¹ (soil, powder)	Particulate ² (swabs, wipes, HEPA)	Liquid/Water (filter, grab)	Drinking Water (filter, grab)	Aerosol (growth media, filter, liquid)
<i>Francisella tularensis</i> [Tularemia]	Yes	Culture	This pathogen is included in the HHS/USDA select agent list. Select agent methods are generally provided and controlled by the LRN. In some cases, methods that apply to environmental samples may not be fully developed or validated. If analysis of this agent is required, contact the LRN at (404) 639-2790 for information on the LRN laboratory capable of receiving and processing the sample.					
	No	PCR/Immunoassay						
<i>Leptospira</i> (<i>L. interrogans</i> Serovars Icterohaemorrhagiae, Australis, Balum, Bataviae, Bejro, Pomona)	Yes	Culture	SM 9260 I	SM 9260 I	Adapted from analytical method	SM 9260 I	SM 9260 I	Unlikely to be viable
	No	Immunoassay	SM 9260 I	SM 9260 I	Adapted from analytical method	SM 9260 I	SM 9260 I	Unlikely to be viable
<i>Listeria monocytogenes</i>	Yes	Culture	FDA/Bacteriological Analytical Manual Chapter 10, 2003	FDA/Bacteriological Analytical Manual Chapter 10, 2003	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method
	No	Immunoassay	FDA/Bacteriological Analytical Manual Chapter 10, 2003	FDA/Bacteriological Analytical Manual Chapter 10, 2003	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method
<i>Salmonella</i> spp. (Method not applicable to <i>Salmonella typhi</i>)	Yes	Culture	Method 1682	Method 1682	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method	Unlikely to be viable
	No	Immunoassay	Method 1682	Method 1682	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method	Unlikely to be viable
<i>Salmonella typhi</i> [Typhoid fever]	Yes	Culture	SM 9260 B	Adapted from analytical method	Adapted from analytical method	SM 9260 B	SM 9260 B	Unlikely to be viable
	No	Immunoassay	SM 9260 B	Adapted from analytical method	Adapted from analytical method	SM 9260 B	SM 9260 B	Unlikely to be viable
<i>Shigella</i> spp. [Shigellosis]	Yes	Culture	SM 9260 E	Adapted from analytical method	Adapted from analytical method	SM 9260 E	SM 9260 E	Unlikely to be viable
	No	Immunoassay	SM 9260 E	Adapted from analytical method	Adapted from analytical method	SM 9260 E	SM 9260 E	Unlikely to be viable
<i>Staphylococcus aureus</i>	Yes	Culture	SM 9213	Adapted from analytical method	Adapted from analytical method	SM 9213	SM 9213	Adapted from analytical method
	No	TBD	TBD	TBD	TBD	TBD	TBD	TBD
<i>Vibrio cholerae</i> O1 and O139 [Cholera]	Yes	Culture	SM 9260 H	Adapted from analytical method	Adapted from analytical method	SM 9260 H	SM 9260 H	Unlikely to be viable
	No	Immunoassay	SM 9260 H	Adapted from analytical method	Adapted from analytical method	SM 9260 H	SM 9260 H	Unlikely to be viable
<i>Yersinia pestis</i> [Plague]	Yes	Culture	This pathogen is included in the HHS/USDA select agent list. Select agent methods are generally provided and controlled by the LRN. In some cases, methods that apply to environmental samples may not be fully developed or validated. If analysis of this agent is required, contact the LRN at (404) 639-2790 for information on the LRN laboratory capable of receiving and processing the sample.					
	No	PCR/Immunoassay						
Viruses								
Adenoviruses: Enteric and non-enteric (A-F)	Yes	Tissue culture ³	AEM 71(6): 3131-3136	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴
	No	Real-time PCR	AEM 71(6): 3131-3136	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴

Pathogen(s)	Viability	Analytical Technique	Analytical Method	Solid ¹ (soil, powder)	Particulate ² (swabs, wipes, HEPA)	Liquid/Water (filter, grab)	Drinking Water (filter, grab)	Aerosol (growth media, filter, liquid)
Astroviruses	Yes	Integrated Cell Culture/RT-PCR	Canadian Journal of Microbiology 50: 269-278	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴
	No	Real-time RT-PCR	Canadian Journal of Microbiology 50: 269-278	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴
Caliciviruses: Noroviruses (NoV)	No	Real-time RT-PCR	Journal of Clinical Microbiology 42(10): 4679-4685	Journal of Clinical Microbiology Vol. 42(10): 4679-4685 ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴
Caliciviruses: Sapovirus (SaV)	No	Real-time RT-PCR	Journal of Medical Virology 78(10): 1347-1353	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴
Coronaviruses: SARS-associated human coronavirus (SARS-HCoV)	No	Real-time RT-PCR	Journal of Virological Methods 122: 29-36	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴
Hepatitis E virus (HEV)	No	Real-time RT-PCR	Journal of Virological Methods 131(1): 65-71	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴
Influenza H5N1 virus	No	Real-time RT-PCR	Emerging Infectious Diseases 11(8): 1303-1305	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴
Orthopoxviruses: Variola, Vaccinia, Cowpox, Monkeypox, Camelpox, Ectromelia, and Gerbilpox	Yes	Tissue culture	LRN This pathogen is included in the HHS/USDA select agent list. Select agent methods are generally provided and controlled by the LRN. In some cases, methods that apply to environmental samples may not be fully developed or validated. If analysis of this agent is required, contact the LRN at (404) 639-2790 for information on the LRN laboratory capable of receiving and processing the sample.					
	No	PCR/Immunoassay						
Picornaviruses: Enteroviruses	Yes	Tissue culture	USEPA Manual of Methods for Virology EPA/600/4-84/013, April 2001	Adapted from analytical method ⁴	Adapted from analytical method ⁴	USEPA Manual of Methods for Virology EPA/600/4-84/013, April 2001	USEPA Manual of Methods for Virology EPA/600/4-84/013, April 2001	Unlikely to be viable
	No	RT-PCR	AEM 69(6): 3158-3164	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴	Adapted from analytical method ⁴
Picornaviruses: Hepatitis A virus (HAV)	Yes	Tissue culture	TBD	TBD	TBD	TBD	TBD	TBD
	No	RT-PCR	AEM 69(6): 3158-3164	Adapted from analytical method ⁴	Adapted from analytical method ⁴	AEM Vol. 69 No. 6: 3158-3164 ⁴	AEM Vol. 69 No. 6: 3158-3164 ⁴	Adapted from analytical method ⁴
Reoviruses: Rotavirus (Group A)	Yes	Tissue culture	AEM 69(6): 3158-3164	Adapted from analytical method ⁴	Adapted from analytical method ⁴	AEM Vol. 69 No. 6: 3158-3164 ⁴	AEM Vol. 69 No. 6: 3158-3164 ⁴	Unlikely to be viable
	No	RT-PCR	AEM 69(6): 3158-3164	Adapted from analytical method ⁴	Adapted from analytical method ⁴	AEM Vol. 69 No. 6: 3158-3164 ⁴	AEM Vol. 69 No. 6: 3158-3164 ⁴	Adapted from analytical method ⁴
Protozoa								
Cryptosporidium spp. [Cryptosporidiosis]	Yes	Tissue culture	AEM 65(9): 3936-3941	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method	Unlikely to be found
	No	IMS/FA and/or IFA	Method 1622 and/or Method 1623	Adapted from analytical method	Adapted from analytical method	Method 1623	Method 1622	Unlikely to be found

Pathogen(s)	Viability	Analytical Technique	Analytical Method	Solid ¹ (soil, powder)	Particulate ² (swabs, wipes, HEPA)	Liquid/Water (filter, grab)	Drinking Water (filter, grab)	Aerosol (growth media, filter, liquid)
<i>Entamoeba histolytica</i>	Yes	Culture	Journal of Parasitology 58(2): 306-310	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method	Unlikely to be found
	No	PCR	Journal of Clinical Microbiology 43(11): 5491-5497	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method	Unlikely to be found
<i>Giardia</i> spp.	Yes	Culture	Trans. R. Soc. Trop. Med. Hyg. 77(4): 487-488	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method	Adapted from analytical method	Unlikely to be found
	No	IMS/FA	Method 1623	Adapted from analytical method	Adapted from analytical method	Method 1623	Method 1623	Unlikely to be found
<i>Toxoplasma gondii</i> [Toxoplasmosis]	Yes	Animal infectivity	Emerging Infectious Diseases 12(2): 326-329	Adapted from analytical method	Adapted from analytical method	Emerging Infectious Diseases Vol. 12 No. 2: 326-329	Emerging Infectious Diseases Vol. 12 No. 2: 326-329	Unlikely to be found
	No	PCR	AEM 70(7): 4035-4039	Adapted from analytical method	Adapted from analytical method	AEM Vol. 70 No. 7: 4035-4039	AEM Vol. 70 No. 7: 4035-4039	Unlikely to be found
Helminths								
<i>Baylisascaris procyonis</i>	Yes	Embryonation of eggs	EPA/625/R92/013	EPA/625/R92/013	Adapted from analytical method	EPA/625/R92/013	Adapted from analytical method	Unlikely to be found
	No	Microscopy	EPA/625/R92/013	EPA/625/R92/013	Adapted from analytical method	EPA/625/R92/013	Adapted from analytical method	Unlikely to be found

¹ Solid samples (except those containing viruses) should be (1) prepared for culture according to EPA Method 1680, Fecal Coliforms in Sewage Sludge (Biosolids) by Multiple-Tube Fermentation using Lauryl Tryptose Broth (LTB) and EC Medium, and (2) prepared for PCR according to *Quantification of Bias Related to the Extraction of DNA Directly from Soils*, Frostegard A, Courtois S, Ramisse V, Clerc S, Bernillon D, Le Gall F, Jeannin P, Nesme X, Simonet P, AEM. 2005. 65(12): 5409-5420.

² Particulate samples (except those containing viruses) should be prepared for culture according to *Swab Materials and Bacillus anthracis Spore Recovery from Nonporous Surfaces*, Rose, L., Jensen, B., Peterson, A., Banerjee, S.N., and M.J. Arduino, Emerg. Infec. Diseases. 2004 10(6):1023-1029.

³ Given that adenovirus 40 and 41 can be difficult to grow in culture, additional cell lines such as G293 (Journal of Medical Virology. 1983. 11(3) :215-231) or CaCo-2 (Journal of Medical Virology. 1994. 44(3): 310-315) may be considered when these viruses are suspected to be present.

⁴ Samples should be prepared according to procedures found in USEPA Manual of Methods for Virology EPA/600/4-84/013, April 2001.

Appendix D: Selected Biotoxin Methods

Appendix D: Biotoxin Methods

Analyte(s)	CAS RN	Analysis Type	Analytical Technique	Aerosol (filter/cassette, liquid impinger)	Solid (soil, powder)	Particulate (swabs, wipes, filters)	Liquid/Drinking Water
Protein							
Abrin	1393-62-0	Presumptive	Immunoassay	Adapted from 119th AOAC Annual Meeting & Exposition, 2005, p. 613	Adapted from 119th AOAC Annual Meeting & Exposition, 2005, p. 613	Adapted from 119th AOAC Annual Meeting & Exposition, 2005, p. 613	Adapted from 119th AOAC Annual Meeting & Exposition, 2005, p. 613
		Confirmatory	Competitive binding assay	Adapted from Pharmacology & Toxicology 88(5): 255-260	Adapted from Pharmacology & Toxicology 88(5): 255-260	Adapted from Pharmacology & Toxicology 88(5): 255-260	Adapted from Pharmacology & Toxicology 88(5): 255-260
		Biological Activity	Enzyme activity	Adapted from Analytical Biochemistry 357(2): 200-207	Adapted from Analytical Biochemistry 357(2): 200-207	Adapted from Analytical Biochemistry 357(2): 200-207	Adapted from Analytical Biochemistry 357(2): 200-207
Botulinum neurotoxins (Serotypes A, B, E, F)	NA	Presumptive	Immunoassay	Adapted from Lateral flow immunassay kits	LRN If analysis for this agent is required in solid, particulate, or liquid samples, contact the LRN at (404) 639-2790 for information of the closest LRN laboratory capable of receiving and processing the sample. The terms presumptive and confirmatory as used for LRN methods are described in Section 8.2.1.		
		Confirmatory	Immunoassay (ELISA)	Adapted from FDA Bacteriological Analytical Manual, Chapter 17			
		Biological Activity	Mouse Bioassay	Adapted from FDA Bacteriological Analytical Manual, Chapter 17			
α-Conotoxin	156467-85-5	Presumptive	Immunoassay	Adapted from Biochemical Journal 328: 245-250	Adapted from Biochemical Journal 328: 245-250	Adapted from Biochemical Journal 328: 245-250	Adapted from Biochemical Journal 328: 245-250
		Confirmatory	HPLC-MS	Adapted from Journal of Medicinal Chemistry 47(5): 1234-1241	Adapted from Journal of Medicinal Chemistry 47(5): 1234-1241	Adapted from Journal of Medicinal Chemistry 47(5): 1234-1241	Adapted from Journal of Medicinal Chemistry 47(5): 1234-1241
		Biological Activity	Receptor binding	Adapted from Biochemical Journal 328: 245-250	Adapted from Biochemical Journal 328: 245-250	Adapted from Biochemical Journal 328: 245-250	Adapted from Biochemical Journal 328: 245-250
Ricin	9009-86-3	Presumptive	Immunoassay	Adapted from Lateral flow immunassay kits	LRN If analysis for this agent is required in solid, particulate, or liquid samples, contact the LRN at (404) 639-2790 for information of the closest LRN laboratory capable of receiving and processing the sample. The terms presumptive and confirmatory as used for LRN methods are described in Section 8.2.1.		
		Confirmatory	Immunoassay (ELISA)	Adapted from Journal of Food Protection 68(6): 1294-1301	Adapted from Journal of Food Protection 68(6): 1294-1301	Adapted from Journal of Food Protection 68(6): 1294-1301	Adapted from Journal of Food Protection 68(6): 1294-1301
		Biological Activity	Enzyme activity	Adapted from Analytical Biochemistry 357(2): 200-207	Adapted from Analytical Biochemistry 357(2): 200-207	Adapted from Analytical Biochemistry 357(2): 200-207	Adapted from Analytical Biochemistry 357(2): 200-207
Shiga and Shiga-like Toxins (Stx, Stx-1, Stx-2)	75757-64-1	Presumptive	Optical immunoassay	Adapted from Journal of Clinical Microbiology 45(10): 3377-3380	Adapted from Journal of Clinical Microbiology 45(10): 3377-3380	Adapted from Journal of Clinical Microbiology 45(10): 3377-3380	Adapted from Journal of Clinical Microbiology 45(10): 3377-3380
		Confirmatory	Immunoassay (ELISA)	Adapted from FDA Bacteriological Analytical Manual, Appendix 1	Adapted from FDA Bacteriological Analytical Manual, Appendix 1	Adapted from FDA Bacteriological Analytical Manual, Appendix 1	Adapted from FDA Bacteriological Analytical Manual, Appendix 1
		Biological Activity	Competitive binding assay	Adapted from Pharmacology & Toxicology 88(5): 255-260	Adapted from Pharmacology & Toxicology 88(5): 255-260	Adapted from Pharmacology & Toxicology 88(5): 255-260	Adapted from Pharmacology & Toxicology 88(5): 255-260

Analyte(s)	CAS RN	Analysis Type	Analytical Technique	Aerosol (filter/cassette, liquid impinger)	Solid (soil, powder)	Particulate (swabs, wipes, filters)	Liquid/Drinking Water
Staphylococcal enterotoxins (SEB)	39424-53-8	Presumptive	Immunoassay	Adapted from 993.06 (AOAC)	LRN If analysis for this agent is required in solid, particulate, or liquid samples, contact the LRN at (404) 639-2790 for information of the closest LRN laboratory capable of receiving and processing the sample. The terms presumptive and confirmatory as used for LRN methods are described in Section 8.2.1.		
		Confirmatory	TBD	TBD	TBD	TBD	TBD
		Biological Activity	TBD	TBD	TBD	TBD	TBD
Staphylococcal enterotoxins (SEA, SEC)	37337-57-8, 39424-54-9	Presumptive	Immunoassay	Adapted from 993.06 (AOAC)	Adapted from 993.06 (AOAC)	Adapted from 993.06 (AOAC)	Adapted from 993.06 (AOAC)
		Confirmatory	TBD	TBD	TBD	TBD	TBD
		Biological Activity	TBD	TBD	TBD	TBD	TBD
Small Molecule							
Aflatoxin (Type B1)	1402-68-2	Presumptive	Immunoassay	Adapted from 991.31 (AOAC)	Adapted from 991.31 (AOAC)	Adapted from 991.31 (AOAC)	Adapted from 991.31 (AOAC)
		Confirmatory	HPLC-FL	Adapted from 994.08 (AOAC)	Adapted from 994.08 (AOAC)	Adapted from 994.08 (AOAC)	Adapted from 994.08 (AOAC)
α-Amanitin	23109-05-9	Presumptive	Immunoassay	Adapted from Journal Food Protection 68(6): 1294-1301	Adapted from Journal Food Protection 68(6): 1294-1301	Adapted from Journal Food Protection 68(6): 1294-1301	Adapted from Journal Food Protection 68(6): 1294-1301
		Confirmatory	HPLC amperometric detection	Adapted from Journal of Chromatography 563(2): 299-311	Adapted from Journal of Chromatography 563(2): 299-311	Adapted from Journal of Chromatography 563(2): 299-311	Adapted from Journal of Chromatography 563(2): 299-311
Anatoxin-a	64285-06-9	Presumptive	TBD	TBD	TBD	TBD	TBD
		Confirmatory	HPLC-FL (precolum derivatization)	Adapted from Biomedical Chromatography 10: 46-47	Adapted from Biomedical Chromatography 10: 46-47	Adapted from Biomedical Chromatography 10: 46-47	Adapted from Biomedical Chromatography 10: 46-47
Brevetoxins (B form)	98225-48-0	Presumptive	Immunoassay	Adapted from Environmental Health Perspectives 110(2): 179-185	Adapted from Environmental Health Perspectives 110(2): 179-185	Adapted from Environmental Health Perspectives 110(2): 179-185	Adapted from Environmental Health Perspectives 110(2): 179-185
		Confirmatory	HPLC-MS	Adapted from Toxicon 43(4): 455-465	Adapted from Toxicon 43(4): 455-465	Adapted from Toxicon 43(4): 455-465	Adapted from Toxicon 43(4): 455-465
Cylindrospermopsin	143545-90-8	Presumptive	TBD	TBD	TBD	TBD	TBD
		Confirmatory	HPLC-MS	Adapted from FEMS Microbiology Letters 216(2): 159-164	Adapted from FEMS Microbiology Letters 216(2): 159-164	Adapted from FEMS Microbiology Letters 216(2): 159-164	Adapted from FEMS Microbiology Letters 216(2): 159-164.
Diacetoxyscirpenol (DAS)	2270-40-8	Presumptive	Immunoassay	Adapted from International Journal of Food Microbiology 6(1): 9-17	Adapted from International Journal of Food Microbiology 6(1): 9-17	Adapted from International Journal of Food Microbiology 6(1): 9-17	Adapted from International Journal of Food Microbiology 6(1): 9-17
		Confirmatory	HPLC-MS	Adapted from Rapid Communications in Mass Spectrometry 20(9): 1422-1428	Adapted from Rapid Communications in Mass Spectrometry 20(9): 1422-1428	Adapted from Rapid Communications in Mass Spectrometry 20(9): 1422-1428	Adapted from Rapid Communications in Mass Spectrometry 20(9): 1422-1428

Analyte(s)	CAS RN	Analysis Type	Analytical Technique	Aerosol (filter/cassette, liquid impinger)	Solid (soil, powder)	Particulate (swabs, wipes, filters)	Liquid/Drinking Water
Microcystins (Principal isoforms: LR, YR, RR, LW)	77238-39-2	Presumptive	Immunoassay/ Phosphatase assay	Adapted from Journal of AOAC International 84(4): 1035-1044	Adapted from Journal of AOAC International 84(4): 1035-1044	Adapted from Journal of AOAC International 84(4): 1035-1044	Adapted from Journal of AOAC International 84(4): 1035-1044
		Confirmatory	HPLC-MS	Adapted from Analyst 119(7): 1525-1530	Adapted from Analyst 119(7): 1525-1530	Adapted from Analyst 119(7): 1525-1530	Adapted from Analyst 119(7): 1525-1530
Picrotoxin	124-87-8	Presumptive	Immunoassay	TBD	TBD	TBD	TBD
		Confirmatory	HPLC	Adapted from Journal of Pharmaceutical and Biomedical Analysis 7(3): 369-375	Adapted from Journal of Pharmaceutical and Biomedical Analysis 7(3): 369-375	Adapted from Journal of Pharmaceutical and Biomedical Analysis 7(3): 369-375	Adapted from Journal of Pharmaceutical and Biomedical Analysis 7(3): 369-375
Saxitoxin (STX, NEOSAX, GTX, dcGTX, dcSTX)	35523-89-8	Presumptive	Immunoassay	TBD	TBD	TBD	TBD
		Confirmatory	HPLC-FL (post column derivatization)	Adapted from Journal of AOAC International 78(2): 528-532	Adapted from Journal of AOAC International 78(2): 528-532	Adapted from Journal of AOAC International 78(2): 528-532	Adapted from Journal of AOAC International 78(2): 528-532
T-2 Mycotoxin	21259-20-1	Presumptive	Immunoassay	Adapted from Journal Food Protection 68(6): 1294-1301	Adapted from Journal Food Protection 68(6): 1294-1301	Adapted from Journal Food Protection 68(6): 1294-1301	Adapted from Journal Food Protection 68(6): 1294-1301
		Confirmatory	HPLC-MS	Adapted from Rapid Communications in Mass Spectrometry 20(9): 1422-1428	Adapted from Rapid Communications in Mass Spectrometry 20(9): 1422-1428	Adapted from Rapid Communications in Mass Spectrometry 20(9): 1422-1428	Adapted from Rapid Communications in Mass Spectrometry 20(9): 1422-1428
Tetrodotoxin	9014-39-5	Presumptive	Immunoassay	Adapted from Journal of Clinical Laboratory Analysis 6(2): 65-72	Adapted from Journal of Clinical Laboratory Analysis 6(2): 65-72	Adapted from Journal of Clinical Laboratory Analysis 6(2): 65-72	Adapted from Journal of Clinical Laboratory Analysis 6(2): 65-72
		Confirmatory	HPLC-MS	Adapted from Analytical Biochemistry 290(1): 10-17	Adapted from Analytical Biochemistry 290(1): 10-17	Adapted from Analytical Biochemistry 290(1): 10-17	Adapted from Analytical Biochemistry 290(1): 10-17

