

## **Texas Natural Resource Conservation Commission**

### **Inspection Checklist for Laboratories Operated by Wastewater Permittees**

The purpose of this checklist is to assist TNRCC field investigators who inspect laboratories operated by permittees that provide analyses for self-reporting requirements specified in a wastewater permit. Commercial laboratories that provide analytical services for wastewater permittees are inspected by the Quality Assurance Specialist from the TNRCC Laboratory in Houston.

This checklist is not required to be used at all times when inspecting permittee-operated laboratories. However, it is recommended that the checklist be used as a guide to aid the field investigator in orally reviewing the laboratory's testing procedures and quality control program.

The checklist is set up as follows:

1. Copies of the Appropriate Method References On-hand: This section lists the test procedure references the laboratory should have on hand. A laboratory is not required to have all of those listed, but it must have the ones it is using and referencing at the time of the inspection.

There are seven references listed and each has a "●" with a number inside it next to the reference. After the reference, there is a blank line where the field investigator can identify which references are on hand.

2. General Laboratory Conditions and Laboratory Instruments, Equipment & Reagents: After each item in these sections, a "●" (one or more with a number inside it) is shown. If the field investigator identifies a problem for any of these items, the number inside the bullet will refer the field investigator and the laboratory to the regulations covering that provision.
3. Laboratory Methods: Prior to conducting the inspection, the field investigator must check the permit and identify the self-reporting parameters for which the laboratory is testing.

Each test method/question in this section has a "●" (one or more with a number inside it) followed by a blank line. If a discrepancy is identified, the field investigator can note the problem on the line and the number in the bullet will refer the field investigator and the laboratory to the references covering that method.

# **Texas Natural Resource Conservation Commission**

## **Inspection Checklist for Laboratories Operated by Wastewater Permittees**

**Note:** If laboratory personnel are interested in laboratory training or need technical assistance, they may contact Terry G. Mills in the TNRCC Laboratory at 281-457-5229.

**COMPANY/LABORATORY NAME:** \_\_\_\_\_

—

**PERMIT NO:** \_\_\_\_\_

**INSPECTOR:** \_\_\_\_\_

**DATE:** \_\_\_\_\_

**“\*” INDICATES A RECOMMENDED PROCEDURE OR ACTION**

**“●” with number inside is given as the reference for the specifications listed in each section.**

### **Copies of the Appropriate Method References On-hand:**

Standard Methods For the Examination  
Of Water and Wastewater (Current Edition) \_\_\_\_\_

TAC Chapter 30, Section 319 \_\_\_\_\_  
“General Regulations Incorporated Into  
Permits”, Subchapter A

CFR 29/Part 1910 (current edition) \_\_\_\_\_  
Health & Safety

40 CFR/Part 136 (current edition) \_\_\_\_\_  
Protection of Environment

EPA methods manual 600/4-79-020 \_\_\_\_\_  
“Methods Of Chemical Analysis Of  
Water and Wastes”

EPA QC methods manual 600/4-79-019 \_\_\_\_\_  
“Handbook For Analytical Quality  
Control In Water And Wastewater  
Laboratories”

TNRCC 20/40 Hour Laboratory Manual \_\_\_\_\_

(Copy required if used in lab)

## GENERAL LABORATORY CONDITIONS

Laboratory clean and orderly	_____
All spills cleaned up immediately	_____
All MSDS sheets present and available	_____
Acids and bases separated in storage	_____
Lab properly ventilated and cooled	_____
All lab records kept 3 years	_____
All required licenses/certif. posted	_____
*Voluntary "Certified Laboratory Analyst" Certificates posted	_____
All containers labeled	_____
All chemicals dated rec'd/opened	_____
Methods in use are available in lab	_____
Emergency eye washes/showers present	_____
Fire extinguishes present/accessible	_____

## LABORATORY INSTRUMENTS, EQUIPMENT, & REAGENTS

**"\*" indicates a recommended procedure or technique.**

### Reagents:

American Chemical Society (ACS)/reagent grade or better	_____
--	-------

**Glassware:**

All class "A" and Borosilicate  
(unless otherwise indicated by method)

---

**Water:**

Potable water supply present

---

Distilled or De-ionized water:

Generated on-site

---

Purchased elsewhere

---

Meets Reagent Grade/Type III  
specifications or better

---

**Balances:**

Serviced by technician annually

---

Service certificate on file

---

Calibrated prior to each use

---

Constant wt. chk. prior to each use

---

Daily log & maint. book kept

---

**DO (Dissolved Oxygen) Meter/Probe:**

Temp. compensation auto/manual

---

Meter read-out to 0.1

---

Calibrated prior to use

---

Re-calib. if meter moved

---

Calibrated or standardized against  
Winkler daily

---

Winkler titrant normality chk weekly

---

Daily log & maint. book kept

---

**pH Meter:**

Temp. compensation auto/manual	_____
Meter read-out to 0.1	_____
Calib. with 2 buffers prior to use	_____
Do buffers bracket pH of samples	_____
Re-calib. if meter moved	_____
*Buffer/samples near same temperature (Follow manufacturer's instructions)	_____
Daily log & maint. book kept	_____

**BOD Incubators:**

Interior clean/no spills	_____
No light in incubator	_____
Temp. at $20^{\circ}\text{C} \pm 1^{\circ}\text{C}$	_____
Temp. checked daily when in use	_____
Temp. adjustments noted	_____
Temp. measured in $\text{H}_2\text{O}$ filled btl	_____
Daily log & maint. book kept	_____

**Fecal Coliform Incubators:**

Water/interior clean/no spills	_____
Temp. checked daily when in use	_____
Temp. at $44.5^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$	_____
Temp. adjustments noted	_____
Thermometer submerged in $\text{H}_2\text{O}$	_____

Water level maint. daily \_\_\_\_\_

Daily log & maint. book kept \_\_\_\_\_

**Ovens/Furnaces:**

Interior clean/no spills \_\_\_\_\_

Temp. checked daily when in use \_\_\_\_\_

TSS oven at 103°C to 105°C \_\_\_\_\_

VSS furnace at 550°C  $\pm$  50°C \_\_\_\_\_

Temp. adjustments noted \_\_\_\_\_

Thermometer in oven when in use \_\_\_\_\_

Daily log & maint. book kept \_\_\_\_\_

**Refrigerators:**

Interior clean/no spills \_\_\_\_\_

Temp. checked daily when in use \_\_\_\_\_

Temp. at 4°C (\* $\pm$  2°C) \_\_\_\_\_

Temp. adjustments noted \_\_\_\_\_

Temp. measured in H<sub>2</sub>O filled btl \_\_\_\_\_

Daily log & maint. book kept \_\_\_\_\_

**Thermometers:**

Calib. checked annually \_\_\_\_\_

All thermometers "NITS" traceable \_\_\_\_\_

Broken units (Hg) disp. of properly \_\_\_\_\_

Daily log & maint. book kept \_\_\_\_\_

## LABORATORY METHODS

**\*\*\* indicates a recommended procedure or technique.**

**\*\*NOTE: THE TERM "BATCH" MEANS 10 SAMPLES OR LESS ANALYZED\*\*  
(EXCEPTION TO ABOVE IS FECAL COLIFORM WHERE A BATCH IS 20  
SAMPLES OR LESS)**

### **Chlorine Residual:**

Name/Number of method used \_\_\_\_\_

FAS normality checked weekly \_\_\_\_\_

Cl<sub>2</sub> Equiv. Std. each use \_\_\_\_\_

Sample volume adequate (≥ 500 mL) \_\_\_\_\_

Each measurement "MnOx" corrected \_\_\_\_\_

Each measurement "Blank" corrected \_\_\_\_\_

Daily log & maint. book kept \_\_\_\_\_

### **Laboratory BOD Chlorine Neutralization:**

Sodium Sulfite used \_\_\_\_\_

Sodium Sulfite prepared fresh daily \_\_\_\_\_

Sample volume adequate (≥ 2 liter) \_\_\_\_\_

Daily log & maint. book kept \_\_\_\_\_

### **BOD<sub>5</sub> (Biochemical Oxygen Demand) at 20°C ± 1°C**

Name/Number of method used \_\_\_\_\_

Samples at **20°C ± 1°C** when analyzed \_\_\_\_\_

Dilution water at **20°C ± 1°C** when used \_\_\_\_\_

## Controls

### Dilution Water Blanks

Dilution water blank with each day  
(\*2 recommended) \_\_\_\_\_

Dilution water @  $20^{\circ}\text{C} \pm 1^{\circ}\text{C}$  prior to use \_\_\_\_\_

Depletion  $\geq 0.2$  mg/L \_\_\_\_\_

Nutrients:

In-house nutrients used \_\_\_\_\_

Commercially prepared used \_\_\_\_\_

If “yes”, is certification on file \_\_\_\_\_

Reagents dated Rec'd/Opened \_\_\_\_\_

### Seed Controls

Seed controls with each day \_\_\_\_\_

Prepared fresh daily \_\_\_\_\_

1 control used (required)  
(\*3 Recommended at 5%, 10%, 15%) \_\_\_\_\_

40% to 70% DO depletion in control \_\_\_\_\_

Seed Corr formula  $[(B1 - B2) \times F]$  \_\_\_\_\_

Seed Corr  $\leq 0.6$  mg/L and  $\geq 1.0$  mg/L  
(\* .5 to 1.1 acceptable on occasion) \_\_\_\_\_

If commercial seed used: \_\_\_\_\_

Certification on file \_\_\_\_\_

Used immed. after aeration \_\_\_\_\_

Reagents dated Rec'd/Opened \_\_\_\_\_

Prepared fresh daily \_\_\_\_\_



If in-house (raw/influent) seed used:	_____
BOD tested on in-house seed	_____
Toxic effects observed	_____
pH adjusted 6.5 to 7.5	_____
Used immed. after aeration	_____
Prepared fresh daily	_____

**Glucose/Glutamic Acid Standards (GGA's)**

GGA's set up each day	_____
GGA prepared fresh daily	_____
Oven dried/cooled desic. prior to wt.	_____
GGA set up as duplicates of 2% conc.	_____
GGA BOD = 198 mg/L $\pm$ 30.5 mg/L	_____

**BOD SAMPLES**

48 hr. holding time observed	_____
Sx @ 20°C $\pm$ 1°C prior to analysis	_____
Sample volume adequate ( $\nless$ 2 liter)	_____
At least 3 <u>different</u> conc. used	_____
All particle sizes < 3 mm (homogenized)	_____
If Sx homogenized, what device is used	_____
if particle sizes > 3 mm	_____
BOD <u>SAMPLE</u> pipet tips $\nless$ 3 mm bore	_____
Was sample seeded if:	
Sample not analyzed within 2 hrs.	_____
Sample temperature changed	_____
pH adjusted to 6.5 - 7.5	_____
Sample disinfected	_____

BOD btl's have ground glass stoppers	_____
Distilled water seal on bottles	_____
(* <i>Dilution</i> water seal recommended)	
*BOD verified with BOD/VSS ratios	_____
True duplicate test performed	_____
1 dup per 10 or per batch	_____
Mean, WL, & CL charted	_____
Reportable BOD calculated	_____
5 day DO $\geq$ 1.0 mg/L	_____
Corrected depl. $\geq$ 2.0 mg/L	_____
Statistical analysis used	_____
Toxic effects indicated	_____
Daily log & maint. book kept	_____

### **Carbonaceous BOD<sub>5</sub> @ 20°C $\pm$ 1°C**

TCMP inhibitor used	_____
Sample volume adequate ( $\nless$ 2 liter)	_____
All controls with inhibitor added	_____
TCMP added to all sample dilutions	_____
All above rules for BOD applied	_____
Daily log & maint. book kept	_____

### **Fecal Coliforms**

Name/Number of method used	_____
Sx Container pre-sterilized with Sodium Thio. in container	_____
Sample volume adequate ( $\nless$ 125 mL)	_____

6 hour holding time observed	_____
Sterilized grad cyl/pipets used	_____
Culture media prepared fresh	_____
Commercially prepared media used (If "Yes", documentation on file)	_____ _____
Commercial sterile culture plates used (If "Yes", documentation on file)	_____ _____
Commercial sterile filters used (If "Yes", documentation on file)	_____ _____
0.1 mL dilutions done in serial dilu's	_____
Sterile buffer wash water used	_____
Apparatus sterilized between each use 3 different dilutions used	_____ _____
All plates/dilutions counted	_____
Blank analyzed daily	_____
If more than two analysts performing the test, are both analysts counting plates at least once per month	_____
Is Presumptive & Confirmed QC using MPN being done once per month on <b><u>positive &amp; negative</u></b> cultures	_____
True duplicate test performed 1 dup per 5% of tests Mean, WL, & CL charted	_____ _____ _____
Daily log & maint. book kept	_____

## TSS/VSS:

Name/Number of method used \_\_\_\_\_

7 Day holding time observed \_\_\_\_\_

Blank analyzed \_\_\_\_\_  
(\*500 mL recommended)

\*Standard Used (Celite or other) \_\_\_\_\_

Sample volume adequate (≠ 500 mL) \_\_\_\_\_

Filter used with wrinkled side up \_\_\_\_\_  
("N/A" if no wrinkled side present)

Filter treated prior to use \_\_\_\_\_

All dried/cooled until at constant wt. \_\_\_\_\_  
± 0.5 mg 1st wt. to 2nd wt. \_\_\_\_\_

Volume of sample filtered achieves: \_\_\_\_\_  
≥ 1 mg wt increase over tare wt \_\_\_\_\_  
< 200 mg increase over tare wt \_\_\_\_\_

Grad. cylinder rinsed properly \_\_\_\_\_  
cyl. 3 times-each through filter \_\_\_\_\_  
funnel walls-3 rinses \_\_\_\_\_  
filter-3 rinses \_\_\_\_\_

True duplicate test performed \_\_\_\_\_  
1 dup per batch \_\_\_\_\_  
Mean, WL, & CL charted \_\_\_\_\_  
All QC > CL re-analyzed \_\_\_\_\_

If VSS analyzed: \_\_\_\_\_  
Filter ashed at 500°C ± 50°C prior to use \_\_\_\_\_

All VSS values ≤ TSS values \_\_\_\_\_

True duplicate test performed \_\_\_\_\_  
1 dup per batch \_\_\_\_\_  
Mean, WL, & CL charted \_\_\_\_\_

All QC > CL re-analyzed

\_\_\_\_\_

Daily log & maint. book kept

\_\_\_\_\_

### **pH @ 25°C (Laboratory)**

Name/Number of method used

\_\_\_\_\_

BOD sample pH measure @ 20°C

\_\_\_\_\_

pH meter approved for self reporting  
(See "Laboratory Instruments, Equipment,  
& Reagents" above for spec's)

\_\_\_\_\_

All pH values within calib. range

\_\_\_\_\_

Sample volume adequate ( $\neq$  50 mL)

\_\_\_\_\_

True duplicate test performed

\_\_\_\_\_

1 dup per batch

\_\_\_\_\_

Mean, WL, & CL charted

\_\_\_\_\_

All QC > CL re-analyzed

\_\_\_\_\_

Daily log & maint. book kept

\_\_\_\_\_

### **DO (Winkler Titration)**

Name/Number of method used

\_\_\_\_\_

Sample collected in 300 mL BOD bottle

\_\_\_\_\_

If volume titrated is 200 mL:

\_\_\_\_\_

Is volume adj to 201.3 mL:

\_\_\_\_\_

Titrant Normality = 0.025N

\_\_\_\_\_

If volume titrated 300 mL:

\_\_\_\_\_

Titrant Normality = 0.0375N

\_\_\_\_\_

If sodium thio titrant purchased:

\_\_\_\_\_

Normality checked weekly:

\_\_\_\_\_

Sodium Thio titrant prepared from stock:

\_\_\_\_\_

Thio. Normality check weekly: \_\_\_\_\_

Daily log & maint. book kept \_\_\_\_\_

### **Ammonia Nitrogen (NH<sub>3</sub>-N)**

Name/Number of method used \_\_\_\_\_

Holding times observed as follows:

Non-preserved analyze immediately  
(\*\*\*\***WITHIN 15 MINUTES**\*\*\*\*) \_\_\_\_\_

Preserved analyze within 28 days \_\_\_\_\_

Sample volume adequate (≥ 1 Liter) \_\_\_\_\_

Sample distilled prior to analysis \_\_\_\_\_

Chlorine neutralized before distill. \_\_\_\_\_

pH neutralized to pH 9.5 before distill. \_\_\_\_\_

Distillation apparatus all glass \_\_\_\_\_

Ammonia free distilled water used \_\_\_\_\_

If specific ion probe used:

Sample not distilled \_\_\_\_\_

distilled/non-distilled com-  
parison in files \_\_\_\_\_

Comparison performed annually  
(or if meter/probe repaired or replaced) \_\_\_\_\_

Blank and standards:

Analyzed with each batch \_\_\_\_\_

Distilled with samples \_\_\_\_\_

True dup/spk tests performed \_\_\_\_\_

1 dup/spk per batch \_\_\_\_\_

Mean, WL, & CL charted \_\_\_\_\_

All QC > CL re-analyzed \_\_\_\_\_

Daily log & maint. book kept \_\_\_\_\_

## TKN (Total Nitrogen)

Name/Number of method used \_\_\_\_\_

Digestion apparatus all glass \_\_\_\_\_

TKN values **always**  $\geq$  Ammonia \_\_\_\_\_

Proc. for "Ammonia Nitrogen (NH<sub>3</sub>-N)"  
followed after digestion \_\_\_\_\_

## O&G (Oil and Grease, Total Recoverable)

Name/Number of method used \_\_\_\_\_

Sample volume adequate ( $\leq$  1 Liter) \_\_\_\_\_

Sample preserved with H<sub>2</sub>SO<sub>4</sub> \_\_\_\_\_

28 day holding time observed \_\_\_\_\_

Blank analyzed each day O&G tested \_\_\_\_\_

\*Blank value  $\geq$  half of MDL \_\_\_\_\_

Field Dups **OR** Stds each batch \_\_\_\_\_

Type of standard material used \_\_\_\_\_

Widemouth, glass jar with TF liner \_\_\_\_\_

Jar/liner rinsed with freon before use \_\_\_\_\_

All are grab samples-NO composites \_\_\_\_\_

Daily log & maint. book kept \_\_\_\_\_

## COD (Chemical Oxygen Demand)

Name/Number of method used \_\_\_\_\_

Open Reflux: \_\_\_\_\_

Closed Reflux: \_\_\_\_\_

Sample preserved with H<sub>2</sub>SO<sub>4</sub> \_\_\_\_\_

Sample volume adequate (≧ 1 Liter) \_\_\_\_\_

28 day holding time observed \_\_\_\_\_

Blank and standard each batch \_\_\_\_\_

Dup and Spk each batch \_\_\_\_\_

Cl<sup>-</sup> < 2000 mg/L per sample \_\_\_\_\_

Cl<sup>-</sup> levels not measured \_\_\_\_\_

Method used to measure Cl<sup>-</sup> \_\_\_\_\_

\*COD/TOC ratio used to verify COD \_\_\_\_\_

Daily log & maint. book kept \_\_\_\_\_

### **TOC (Total Organic Carbon)**

Name/Number of method used \_\_\_\_\_

Sample preserved with H<sub>2</sub>SO<sub>4</sub> \_\_\_\_\_

Sample volume adequate (≧ 250 mL) \_\_\_\_\_

28 day holding time observed \_\_\_\_\_

Blank and standard each batch \_\_\_\_\_

Dup and Spk each batch \_\_\_\_\_

Daily log & maint. book kept \_\_\_\_\_

### **Phenol (Total Recoverable)**

Name/Number of method used \_\_\_\_\_ \*

(\* use 14th Edition **Only**)

Sample preserved with H<sub>2</sub>SO<sub>4</sub> \_\_\_\_\_

Widemouth, glass jar with TF liner \_\_\_\_\_



Use brown glass or store in dark \_\_\_\_\_

Jar/liner rinsed with freon before use \_\_\_\_\_

Sample volume adequate (≠ 1 Liter) \_\_\_\_\_

28 day holding time observed \_\_\_\_\_

If Cl<sub>2</sub> present, FAS added \_\_\_\_\_

Sample/blk/std all distilled \_\_\_\_\_

Blank and standard each batch \_\_\_\_\_

Dup each batch \_\_\_\_\_

\*Spk with batch \_\_\_\_\_

Chloroform extracts filtered \_\_\_\_\_

All grab samples-NO composites \_\_\_\_\_

Daily log & maint. book kept \_\_\_\_\_

## AA (METALS ANALYSIS)

Name/Number of method used \_\_\_\_\_

If SW846 used, that Method # \_\_\_\_\_

Sample volume adequate (≠ 1 Liter) \_\_\_\_\_

Sample preserved with HNO<sub>3</sub> \_\_\_\_\_

    Is plastic container used \_\_\_\_\_

    Is **clear**, glass jar used \_\_\_\_\_

    If “No”, what type jar \_\_\_\_\_

6 mth holding time observed (excpt. Hg) \_\_\_\_\_

Hg holding time 28 days (sol/liq) \_\_\_\_\_

Cr<sup>6</sup> not pres/run within 24 hrs \_\_\_\_\_

Blank and standard each batch \_\_\_\_\_

Sample/blk/std all digested \_\_\_\_\_

Dup and spk each batch \_\_\_\_\_

Standard curve/calib. check daily \_\_\_\_\_

Daily log & maint. book kept \_\_\_\_\_