

Micro-liquid Chromatography- Electrospray/Ion Trap Mass Spectrometry for Detection of - Dibutyltin - in Natural Waters



**T.L. Jones-Lepp and K.E. Varner
U.S. Environmental Protection Agency
Office of Research and Development
National Exposure Research Laboratory
Las Vegas, NV 89119**





U.S. EPA Notice

The U.S. Environmental Protection Agency (EPA), through its Office of Research and Development (ORD), funded this research and approved the materials that formed the basis for this presentation. The actual presentation has not been peer reviewed by EPA.



Introduction

Pharmaceutical Sources of Dibutyltin (DBT)

- ❑ Feed additive for chickens and turkeys - used as an antihelminthic.

Other Possible Environmental Sources of Dibutyltin

- ❑ Widely used as a plastic stabilizer in PVC products: PVC pipe (i.e., 1 lb of DBT for every 100 lbs of PVC), packaging materials, baking parchments.
- ❑ Breakdown product of Tributyltin.



Adverse Health Effects

DBT

- Significant neurotoxicity exhibited in immature brain cell cultures at 30 ppb.¹
- 100% Suppression of natural killer lymphocytes function at approx. 3 ppm - 30% suppression at 150 ppb.²
- Potent teratogen if exposure is during period of organogenesis.³



Experimental⁴

Water Extractions

- water samples (1-2 L) are pH adjusted to approximately pH 2.5 with HCl (12N)
- The pH-adjusted waters are poured through prepped 48-mm nu*phase SPE C18 discs [using the Accuprep 7000TM manifold (CPI International)]
- Dried discs are then extracted with four 10-mL volumes of 99% methanol/1% acetic acid
- The extracts are then evaporated to 0.5 mL using nitrogen in a water bath 30 °C (TurboVap® II) (Zymark Corporation, Hopkinton, MA, USA).



Experimental - cont.

Analysis

- chromatography - micro-LC columns packed in-house: 160- μ m i.d. x 30-cm fused-silica columns (Polymicro Technologies, Phoenix, AZ). Packed with 10 to 12 cm of 5- μ m ODS-Hypersil (Shandon, Astmoor, England)
- mobile phase: isocratic mobile phase:
80% methanol : 14% water : 6% acetic acid :
0.1% tropolone (v/v/v/m)
- flow rates: 4 to 6 : L/min




Experimental - cont.

Analysis - cont.

- Mass Spectrometry - Finnigan LCQ™, configured with an electrospray (ES) ion source
- positive ion mode
- ES needle: 4.5 kV to 5.2 kV
- 150 to 430 amu (full-scan mode) in 3 microscans with an ion injection of 200 ms.

Results and Discussion

Spectrum of Dibutyltin



Results & Discussion

Amount of DBT found in Natural Waters

State	DBT avg. : g/L
North Dakota	ND ^a
Texas	ND
Maine	ND
Alaska	ND
North Carolina Site #1	1.3
North Carolina Site #2	2.0
North Carolina Site #3	2.6

^aND = not detected, below LOD (LOD = 970 pg – on-column)

Map of North Carolina


Sampling Site #1



★ → Sampling site

Map of North Carolina

Sampling Site #2



Map of North Carolina

Sampling Site #3

★ → Sampling site





Conclusions

- **LC/MS provides a low solvent usage (green), fast, and definitive method for determining dibutyltin in natural waters.**
- **Trace levels (ppb) of dibutyltin are present in some natural waters from various regions in the United States.**



References

1. Eskes, C.; Honegger, P.; Jones-Lepp, T.; Varner, K.; Matthieu, J.-M.; and Monnet-Tschudi, F. “Neurotoxicity of Dibutyltin in Aggregating Brain Cell Cultures,” *Toxicol. In Vitro* **13** (1999) 555-560.
2. Whalen, M.M.; Loganathan, B.G.; Kannan, K. “Immunotoxicity of Environmentally Relevant Concentrations of Butyltins on Human Natural Killer Cells *in Vitro*,” *Environ. Res. Section A* **81** (1999) 108-116.
3. Ema, M.; Kurosaka, R.; Amano, H.; Ogawa, Y. “Comparative Developmental Toxicity of Di-, Tri- and Tetrabutyltin Compounds after Administration during Late Organogenesis in Rats,” *J. Appl. Toxicol.* **16** (1) (1996) 71-76.
4. Jones-Lepp, T.L.; Varner, K.E.; McDaniel, M.; Riddick, L. “Determination of Organotins in Water by Micro Liquid Chromatography–Electrospray/Ion Trap Mass Spectrometry,” *Appl. Organometal. Chem.* **13** (1999) 881-889.