



William C. Brumley and William H. Matchett
U.S. Environmental Protection Agency,
National Exposure Research Laboratory
Characterization Research Division
P.O. Box 93478, Las Vegas, NV 89193-3478

CE experiments were carried out using a Beckman P/ACE 2100 or 5000 instrument with indirect UV detection or laser-induced fluorescence detection (Ar ion laser, 488 nm excitation, 520/20 nm emission).

Solid phase extraction (SPE) studies used SCX cartridges, C18 extraction disks with ion-pairing, or cation exchange extraction disks. Elution was 4% ammonia in methanol.

Amines were derivatized using fluorescein isothiocyanate at 55° C with sodium bicarbonate buffer, usually in aqueous solution. An alternate solvent was acetonitrile. CE separations using indirect detection under free zone conditions used imidazole at 5 mM at pH 5.0 or as otherwise indicated in the tables.

Separations under LIF detection were free zone at pH 7.0 of 5:1 buffer (50 mM phosphate:methanol). Separations under MEKC were at pH 7.0 using urea, SDS, and methanol as indicated on figure.

Aliphatic amines are toxic substances and irritants to mucous membranes that are among the common chemicals of commerce. They are used as corrosion inhibitors in steam boilers and as starting materials in the manufacture of pharmaceuticals, insecticides, herbicides, fungicides, polymers, surfactants, and rubber accelerators. The related alkanolamines function as solvents and starting materials for surfactants, but they appear to be less toxic than the aliphatic amines.

The many commercial uses and natural occurrence of aliphatic amines (here we refer primarily to C_1 to C_4 alkyl substituted primary, secondary, and tertiary amines) suggest that ultimately they will appear in the environment as pollutants. Thus, they are target analytes of U.S. EPA Method 8260 (also Method 624) where they are classified as volatiles. The U.S. EPA, EMSL-LV, maintains a continuing interest in analytical methods for amines because of their wide occurrence. In addition, there is need for determinative methods for amines as a result of the listing activities for various hazardous wastes under RCRA when amines are suspected to be present.

		percent recovery					
		propy- lamine		dipropy- lamine		tripropy- lamine	
sample ¹		Exp.1	Exp.2	Exp.1	Exp.2	Exp.1	Exp.2
1		70	70	77	77	61	61
2		86	86	88	88	65	65
3		92	92	103	100	70	70
Average		82		89		65	


Samples 1, 2, and 3 consisted of 300 μmol of dodecylbenzenesulfonate (Na^+); 37 μmol propylamine; 33 μmol dipropylamine; 31 μmol tripropylamine; 5 mL methanol in a final volume of 1 L. Each of the three samples was passed through a separate disk prepared as suggested by the manufacturer (Empore Bakerbond Octadecyl C₁₈). The 100% level of recovery was established by an appropriate portion of the original solutions used in preparation of the samples.

disk ¹	percent recovery ²		
	propylamine	dipropylamine	tripropylamine
1	105	107	108
2	98	112	107
3	96	96	86
4	0	0	0

¹ Empore Bakerbond Octadecyl (C₁₈) extraction disks were prepared following the manufacturer's instructions.

² The test solution contained per L, 45 μ mol of total amines; 150 μ mol of dodecylbenzene-sulfonate, sodium salt; and 5.0 mL of methanol. The filtrate from disk 3 was passed through disk 4 to check for amines remaining in solution after the first passage.

Figure 1. Electropherogram of propylamines using indirect detection.



Compound	Migration Time Seconds	Mobility $10^4 \text{ cm}^2 \text{ V}^{-1} \text{ s}^{-1}$	Response Factor Per mM	10^4 Theoretical Plates
propylamine	119	3.76	0.385	5.48
dipropylamine	135	2.83	0.421	3.13
tripropylamine	144	2.39	0.395	5.55
tetrabutylammonium bromide	155	1.81	-	0.75

compound ²	sample 1	sample 2	sample 3	average
butylamine	49	67	60	59
dibutylamine	50	66	76	64
tributylamine	42	58	46	49
dimethylamine	12	15	23	17
diethylamine	10	14	18	14
triethylamine	12	16	20	16
propylamine	79	85	83	82
dipropylamine	72	84	78	78

¹ Cation exchange disks (47-mm, sulfonic acid bonded to poly (styrene-divinylbenzene) copolymer, hydrogen form) were a gift of 3M, new products department, St. Paul, MN 55144. Disks were prepared by washing with acetone, methanol, water, very dilute H₂SO₄ (2 drops/100 mL), and finally with H₂O until the pH reached 5.5. Disks were eluted with four 5.0-mL portions of 4% NH₄OH in methanol (v/v); following elution, the disks were regenerated as indicated in the washing procedure.

² Samples consisting of 1L were passed through the disk at a flow rate of approximately 200 mL/min.

Compound	Lowest Detected Concentration mMolar	Mass Detection Limit 10^{-4} nMoles
Propylamine	0.033	1.65
Dipropylamine	0.020	1.0
Tripropylamine	0.014	0.7
Ethanolamine	0.039	1.9
Diethanolamine	0.026	1.3
Triethanolamine	0.018	0.9

Relative Response Arbitrary Units

Concentration (mM)

Dipropylamine Concentration (mM)

Relative Response Arbitrary Units

Concentration (mM)

Figure 4. LIF detection of amines as fluorescein isothiocyanate derivatives; micellar conditions.

The chemical structure of fluorescein isothiocyanate (FITC) is shown in the inset. It consists of a fluorescein core, which is a tricyclic system with a central oxygen atom, fused benzene and pyrone rings, a hydroxyl group, and a carboxylic acid group. The 2-position of the pyrone ring is substituted with an isothiocyanate group (-N=C=S).

The chromatogram displays the following retention times (min): 4.50, 8.82, 9.77, 10.86, 11.80, 12.09, 12.68, 13.35, 13.89, 15.99, and 16.77. The peak at 16.77 min is the most intense, reaching an RFI of approximately 20.0000.

Separation of aliphatic and other amines with LIF detection; buffer and conditions: phosphate (0.039 M), SDS (0.070 M), urea (2.0 M), 23 methanol, pH 7.0, 47 cm X 0.05 mm ID LIF detection (488 nm excitation and 520/20 nm emission, 40 cm to detector). Peak 1 isothiocyanate derivatives): MT =9.77, 2-(2-aminoethyl)pyridine; MT =11.80, 5-diaminopentane; MT =12.68, butylamine; MT =13.35, diethylamine; MT =13.89, propylamine.

Figure 5. LIF detection of amines under free zone separation; 1,5-diaminopentane is well separated.

The aliphatic amines may be determined in aqueous matrices by indirect methods provided the background interferences are minimal. SPE may be used to afford 100 to 1000 fold concentrations of selected analytes. Electrophoretic mobilities were correlated with the Stokes radius of each analyte ion. Responses were a linear function of concentration. Separations of closely migrating ions was effected by use of optimal background electrolyte, additives, or pH adjustment to the pK_b of the amines.

IF detection offers improved sensitivity and robustness with respect to inorganic interferences. Limitations include derivatization of matrix coextractives and by-products of derivatization for ng and sub-ng amounts of amines. Two orthogonal chromatographies are being investigated to address these limitations. GPC provides some removal of many matrix coextractives prior to derivatization. Peak selection from reverse-phase HPLC of the fluorescein isocyanate derivatives provides a cleanup step prior to final separation/determination by CZE or MEKC. Derivatization of tertiary amines is under investigation by other reagents.