

Procedure for Importing DSSTox SDF file into ACD ChemFolder (ver 6.0 for Windows) *(16 April 03)*

1. Open ACD ChemFolder v 6.0 application.
2. **MainMenu → Database → Update Mode**
3. **Database → New →** type name and location of new *.CFD file being created, click **Save**
4. **Enter New Password for Database Update → Cancel**
5. **Database → Import →** *.SDF browse, find, and select SDF file for import
6. Pop up menu for **Import Options** will appear.
 - a. To check to see if all structures are imported and considered valid, click **OK**. Continue procedure and make sure that the number of chemical records imported into ChemFolder equals the known count for valid structures.
 - b. To ensure that all structures are imported, regardless of whether they are considered valid in ChemFolder, or to check to see which structures were considered invalid on previous import, enter a file name in the **Log File** box (browse to select convenient location for Log File), click **Other Options button**, click **Import Invalid Structures** box (box should have check mark), and click **OK**.
7. If only Structure appears in form view, change default view: **MainMenu → View → Screen Forms → Default with User Data**
8. In ACD ChemFolder, a record lists only fields with non-blank entries. Hence, to see all fields, one must scan through database to find examples of records with non-blank entries.
9. **View → Tile** shows all structures in the database. Right click on screen, double click on desired field names to display these field contents in Tiling View along with structures.
10. **View → Table** shows an Access-like table view of all records.