EPA/OPP MICROBIOLOGY LABORATORY ESC, Ft. Meade, MD

Standard Operating Procedure for

Media and Reagents: Examination and Expiration Time

SOP Number: QC-10-02

Date Revised: 08-28-02

Prepared By:		Date:	/	/
	Print Name:			
Reviewed By	:	Date: _	/	_/
	Print Name: Technical Staff			
		Date:	/	/
	Print Name: QA Officer			
		Date:	_/	/
	Print Name: Laboratory Director			
Date Issued:	//			
Withdrawn B	y:	Date:	_//	/
Controlled Co	ppy No.:			

TABLE OF CONTENTS

	<u>Contents</u>	Page Number
1.0	SCOPE AND APPLICATION	2
2.0	DEFINITIONS	2
3.0	HEALTH AND SAFETY	2
4.0	CAUTIONS	2
5.0	INTERFERENCES	2
6.0	PERSONNEL QUALIFICATIONS	2
7.0	SPECIAL APPARATUS AND MATERIALS	2
8.0	INSTRUMENT OR METHOD CALIBRATION	2
9.0	SAMPLE HANDLING AND STORAGE	2
10.0	PROCEDURE AND ANALYSIS	3
11.0	DATA ANALYSIS/CALCULATIONS	5
12.0	DATA MANAGEMENT/RECORDS MANAGEMENT	5
13.0	QUALITY CONTROL	5
14.0	NONCONFORMANCE AND CORRECTIVE ACTION	5
15.0	REFERENCES	5
16.0	FORMS AND DATA SHEETS	5

1.0 SCOPE AND APPLICATION:

1.1 This protocol establishes expiration dates for media and reagents used in conjunction with efficacy testing of disinfectants, and describes the procedure for visual assessment of autoclaved stored media.

2.0 DEFINITIONS:

- 2.1 DI = de-ionized water
- 2.2 Target Volume = required volume of media
- 3.0 HEALTH AND SAFETY: Not applicable
- 4.0 CAUTIONS:
 - 4.1 Products will not be used past the manufacturer's recommended expiration date.

5.0 <u>INTERFERENCES</u>:

5.1 If storage conditions for the media have been compromised in any way, then the media will be discarded.

6.0 PERSONNEL QUALIFICATIONS:

- 6.1 Personnel are required to be knowledgeable of the recommended expiration times for the various media used, and the procedures for visual assessment of media. Documentation of training and familiarization with this SOP can be found in the training file for each employee.
- 7.0 SPECIAL APPARATUS AND MATERIALS: None
- 8.0 <u>INSTRUMENT OR METHOD CALIBRATION</u>: Not applicable
- 9.0 <u>SAMPLE HANDLING AND STORAGE</u>:
 - 9.1 Refer to section 10.0 for storage conditions.

10.0 PROCEDURE AND ANALYSIS:

10.1 <u>Liquid Media/Reagents</u>:

- Liquid Volume Check Prior to and After Autoclaving. To determine if liquid media is lost during autoclaving, mark a random sampling (approximately 5%) of tubes "V" (for volume check) prior to autoclaving. The volume of the media dispensed and autoclaved should be checked after autoclaving. Record the target volume, and the volume after autoclaving on the Liquid Media/Reagent Volume Check Before and After Autoclaving form. Calculate the percent difference. The percent difference should be within ±5% of the target volume for liquid media not used for enumeration. For liquid used for enumeration, the percent difference should be within ±2% of the target volume. If the percent difference exceeds the established difference, the batch will be discarded.
- 10.1.2 <u>Liquid Media/Reagents in Tubes with Morton Closures</u>. Liquid tubed media and reagents with Morton Closures can be stored at room temperature for up to 2 months in plastic bags. The media will be inspected for dehydration according to the instructions stated in 10.3.
- 10.1.3 <u>Liquid Media in Capped Tubes or Bottles</u>. Liquid media in tightly capped bottles or tubes may be stored for up to 3 months. The media in tubes will be inspected for dehydration according to the instructions stated in 10.3.
- 10.1.4 Reagents that are not very labile or are not harmed by long storage periods may be used as long as it can be demonstrated that they perform as they should. Reagents for which a recommended shelf life is known will be used within that time. Commercially prepared media and reagents will not be used beyond the manufacturer's expiration date.

10.2 Solid Media:

10.2.1 Solid Media in Plates. All media poured into Petri dishes will

SOP No. QC-10-02 Date Revised 08-28-02 Page 4 of 8

be stored in plastic bags at 2-5°C. The plastic bags will be punctured in a few areas to allow air exchange and to reduce excessive condensation. Both selective and non-selective media may be stored for up to 3 months. After 3 months the media must be discarded.

10.2.2 <u>Solid Media in Tubes</u>. Tubed agar media with screw caps or other sealed containers may be stored at 2-5°C for up to 3 months. Both selective and non-selective media may be stored for this period of time. After 3 months the media must be discarded.

10.3 <u>Inspection of Stored Media</u>:

- 10.3.1 <u>Solid Media in Tubes and Plates</u>. Approximately 2 to 3 days prior to use, solid media will be inspected. If any discernible drying, cracking or contamination is present, the media will be discarded.
- Liquid Media. Approximately 2 to 3 days prior to use, liquid media in tubes with Morton Closures stored longer than one month or two months for screw caps will be inspected. The fluid level will be checked in approximately 5% per batch and/or per rack of media to be used in testing to determine if more than 5% of the fluid has evaporated. If less than 5% of the solution has evaporated, the media can be used. If more than 5% has evaporated, the media will be discarded.

10.4 Performance and Sterility of all Purchased or Prepared Media/Reagents:

10.4.1 All prepared and purchased media and reagents will undergo sterility and performance verification according to guidelines stipulated in SOP QC-11, Performance and Sterility of Media and Reagents.

11.0 DATA ANALYSIS/CALCULATIONS:

11.1 Calculate the percent difference due to evaporation or loss during autoclaving. The percent difference is determine using the formula:

Percent Difference = 100% - [(Measured Volume / Target Value) X 100]

12.0 DATA MANAGEMENT/RECORDS MANAGEMENT:

12.1 Data will be recorded promptly, legibly, and in indelible ink on the appropriate forms. Completed forms are archived in notebooks kept in locked file cabinets in the file room D217. Only authorized personnel have access to the locked files. Archived data is subject to OPP's official retention schedule contained in SOP ADM-03, Records and Archives.

13.0 QUALITY CONTROL:

13.1 The OPP Microbiology Laboratory conforms to 40CFR Part 160, Good Laboratory Practices. Appropriate quality control measures are integrated into each SOP.

For quality control purposes, the required information is documented on the appropriate record form(s) (see 16.0).

14.0 NONCONFORMANCE AND CORRECTIVE ACTION:

- 14.1 Outdated media or that which has measurably lost more than 5% of the target volume is not to be used.
- 15.0 REFERENCES: None

16.0 FORMS AND DATA SHEETS:

- 16.1 Check of Stored Solid Media
- 16.2 Liquid Media/Reagents Volume Check Before and After Autoclaving
- 16.3 Volume Check for Stored Liquid Media and Reagents

SOP No. QC-10-02 Date Revised 08-28-02 Page 6 of 8

Check of Stored Solid Media OPP Microbiology Laboratory

Date	Initials	Media	Media Preparation No.	Observations (signs of contamination and/or dehydration)	

^{*}Save or Discard

SOP No. QC-10-02 Date Revised 08-28-02 Page 7 of 8

Liquid Media/Reagents Volume Check Before and After Autoclaving OPP Microbiology Laboratory

Date	Initials	Media/Reagent	Media Preparation No.	Target Volume (before autoclaving)	Volume after autoclaving	

^{*}Save or Discard

SOP No. QC-10-02 Date Revised 08-28-02 Page 8 of 8

Volume Check for Stored Liquid Media and Reagents OPP Microbiology Laboratory

Date	Initials	Media/Reagent	Media Preparation					
Date	initials	would, reagent	No.	No.	Measured Volume	Target Volume	% Dehydration	

^{*}Save or Discard