

Laboratory Evaluation of Amalgam Separators

ADA Laboratories

Introduction

- **POTW effluent limits for Hg**
- **Binational Toxics Strategy**
- **EPA Virtual Elimination**
- **Discharge into the environment**
- **Dental office wastewater**

Introduction

- **Environmental impact ?**
- **Hg Amalgam**
- **Identifiable source**
- **Source reduction**
- **Amalgam separators**

Aim

- **Evaluate amalgam removal efficiency**
- **Determine mercury concentration in effluent**

Amalgam Separators Tested

A1000

Amalgam Collector

Asdex

BullfroHg

Durr 7800

ECO II

Hg 5

Hg 10

MRU

MSS 2000

Rasch 890-4000

RME 2000

Removal Technologies

- **Sedimentation**
- **Filtration**
- **Chemical Removal**
- **Centrifugation**
- **Combinations**

ISO Protocol: Amalgam Removal Efficiency

ISO Standard 11143

10 g amalgam with size distribution:


$60\% \leq 3.15 \text{ mm} > 0.5 \text{ mm}$

$10\% \leq 0.5 \text{ mm} > 0.1 \text{ mm}$

$30\% \leq 0.1 \text{ mm}^*$

* meets distribution curve

ISO Protocol: Amalgam Removal Efficiency



ISO Protocol: Amalgam Removal Efficiency

- 1L filtered water
- 1g sodium pyrophosphate
- Slurry + filtered water to achieve maximum flow rate
- Collect effluent

ISO Protocol: Amalgam Removal Efficiency

- Collect amalgam in effluent using pre-weighed filters (12 μm , 3 μm , 1.2 μm)
- Dry filters with collected amalgam to constant weight

Removal Efficiency (E)

$$E (\%) = 100 (W_I - W_F)/W_I$$

W_I = Wt. of amalgam sample

W_F = Wt. of amalgam in filters

Removal Efficiency (E)

- Amalgam separators tested
“Empty” & “Full”
except for Hg 10
- Full = 70% glass beads (1mm)
+ 30% amalgam ($\leq 0.3\text{mm}$)

Removal Efficiency (E)


- Tests done in triplicate
- Lower mean value for “Empty” or “Full” is the removal efficiency

Statistical Analysis for Removal Efficiency

- **Among amalgam separators:
ANOVA & Multiple
Comparison (Tukey)**
- **Each amalgam separator
“Empty” & “Full”: t-test**

Total Mercury Concentration in Effluent

- Mercury from amalgam in pre-weighed filters ($\geq 1.2 \mu\text{m}$) by calculation from weight/volume
- Mercury in effluent ($< 1.2 \mu\text{m}$) by EPA 245.1: acid digestion & AA
- Total = Sum of the above



Results: Removal Efficiency

ANOVA and Tukey Multiple Comparison showed differences in the amalgam removal efficiencies among the amalgam separators tested using the ISO 11143 protocol

Hg 10	99.99%
MRU	99.95%
Rasch 890-4000	99.90%
Amalgam Collector	99.89%
RME 2000	99.66%
Hg 5	99.28%
Asdex	99.10%
MSS 2000	98.94%
BullfroHg	98.88%
Durr 7800	97.66%
ECO II	97.51%
A1000	96.09%

Removal Efficiency

Except for MSS 2000, for each amalgam separator there are no differences (t-test) between the “Empty” and “Full” amalgam removal efficiencies as measured using the ISO protocol

Total Mercury Concentration in Effluent

**Total mercury concentration in
the effluent from laboratory
testing varied widely for each
amalgam separator**

Total Mercury Concentration in Test Effluent (ppb)

A1000	30,200 – 34,899
Amalgam Collector	689 – 3,349
Asdex	9,474 – 17,953
BullfroHg	3,498 – 16,269
Durr 7800	967 – 4,045
ECO II	16,307 – 39,770

Total Mercury Concentration in Test Effluent (ppb)

Hg 5	6,429 – 15,150
Hg 10	20 -- 100
MRU	196 – 697
MSS 2000	732 – 6,168
Rasch 890-4000	601 – 1,757
RME	767 – 4,449

Total Mercury Concentration in 55 L of Effluent (ppb)

A1000	3,283 - 4,671
Amalgam Collector	115 - 351
Asdex	784 - 8,429
BullfroHg	696 - 2,913
Durr 7800	887 - 3,615
ECO II	1,400 - 2,076

Total Mercury Concentration in 55 L of Effluent (ppb)

Hg 5	557 - 6,680
Hg 10	9 - 13
MRU	93 - 169
MSS 2000	156 - 1,168
Rasch 890-4000	67 - 157
RME	196 - 629

Conclusions

- **All 12 amalgam separators tested exceeded the ISO 11143 requirement of 95% amalgam removal efficiency**

Conclusions

- Amalgam effluent from separators in the laboratory test, when chemically digested, had total mercury levels in the ppb range.

Publications

**Laboratory evaluation of amalgam
separators.**

JADA 2002; 133:577-584.

Amalgam in dental office wastewater.

JADA 2002; 133:585-589.