

# **Rapid, Quantitative Biological Indicator System with *Bacillus thuringiensis* Al Hakam Spores**

---

**May 9, 2018**

**Yoojeong Kim<sup>1</sup>, John Lovaasen<sup>1</sup>, Kristopher Daly<sup>2</sup>, Anant Singh<sup>1</sup>**

<sup>1</sup> Triton Systems, Inc., 200 Turnpike Road, Chelmsford, MA 01824

<sup>2</sup> KBioSim, Worcester, MA 01606



# Current Practice of Decontamination Assurance




---

- Surface sampling after decontamination - swabs
- Spore strips or coupons with a known population (biological indicators) placed before decontamination and retrieved afterward
- They are extracted, serially diluted, plated, and enumerated.
- It requires considerable labor.
- Results typically cannot be obtained before 24 – 48 hours, up to 7 days.



# Triton's BIs and Detector

- Self-contained BIs and an incubator/detector system
- Quantitative
- *Bacillus thuringiensis* (Bt) Al Hakam
  - A close relative of *B. anthracis*, with the same growth and spore production properties
  - Without pathogenicity: do not have pXO1 or pXO2 plasmids found in *B. anthracis* that encodes for virulence genes and the anti-phagocytic capsule





# Principle of Assay to Detect Spore Viability

---

- Evaluates the ability of the spores to germinate and carry out protein synthesis as a measure of the viability of the spores
- Fluorogenic Substrate A
  - Based on the enzyme activity packaged in dormant spores of Bt AI Hakam
  - The enzyme is either not active or not accessible to Substrate A in dormant spores.
  - When the spores germinate, Substrate A is taken up by the spores and hydrolyzed into a highly fluorescent compound by the enzyme.
  - The fluorescence yield is further increased by promoting spore outgrow and vegetative growth.


# Quantitative Assay

- 10 different batches were tested.



# Principle of Assay to Detect Spore Viability


- Evaluates the ability of the spores to germinate and carry out protein synthesis as a measure of the viability of the spores
- Fluorogenic Substrate A
  - Based on the enzyme activity packaged in dormant spores of Bt Al Hakam
  - The enzyme is either not active or not accessible to Substrate A in dormant spores.
  - When the spores germinate, Substrate A is taken up by the spores and hydrolyzed into a highly fluorescent compound by the enzyme.
  - The fluorescence yield is further increased by promoting spore outgrow and vegetative growth.
  - **Substrate A is heat and moisture sensitive and cannot be autoclaved to sterilize after glass ampoule production.**




# Principle of Assay to Detect Spore Viability

- Fluorogenic Substrate B
  - Based on the enzyme activity synthesized during germination, outgrowth, and vegetative growth
  - Heat and moisture stable – does not require two compartment ampoules

## Fluorogenic Substrate B




## Fluorogenic Substrate C




# Comparison between Substrates A and B

- Both give linear relationships.
- Fluorescence generation is slower with Substrate B.




# Growth Medium Optimization for Substrate B


## Growth Medium 1




## Growth Medium 2




## Growth Medium 3



## Growth Medium 4




# Correlation for Growth Medium 1



# Production of Glass Ampoules


- 7 mm OD × 40 mm H
- Fill volume: 0.6 mL
- 5000 glass ampoules were made and autoclaved to sterilize.
- 6000 more plastic vials and caps were produced.
- Spore strip: inoculated with  $10^7$  CFU Bt Al Hakam spores



# Triton's Detector

- Simultaneous incubation (37°C) and fluorescence detection of 12 BIs
- Touch screen user interface

Well Status Indicator Lights




Touch Screen




# Triton's BIs and Detector

- 10 – 12 hours to detect a single spore → 12 hour run time
- Data processing to display estimated spore population is being implemented.

**Killed  
spores**



**10<sup>7</sup>  
spores**




**10<sup>0</sup>  
spores**




# Assay Results with Glass Ampoules

- Glass ampoules cause more scatter in the data.



# Spores with Organic Debris



- Spore strip: Bt Al Hakam spores + humic acid
- Buhr et al.<sup>1</sup> – 5 g/L humic acid for  $1-2 \times 10^8$  spores/mL
- $10^7$  spores with 250  $\mu$ g humic acid



<sup>1</sup> *J Appl Microbiol* 2015. 119:1263–1277

# Spores with Organic Debris


- Spore strip: Bt Al Hakam spores + humic acid
- Buhr et al.<sup>1</sup> – 5 g/L humic acid for  $1-2 \times 10^8$  spores/mL
- $10^7$  spores with 250  $\mu$ g humic acid



<sup>1</sup> *J Appl Microbiol* 2015. 119:1263–1277

# Spores with Organic Debris

- Spore strip: Bt Al Hakam spores + humic acid
- Buhr et al.<sup>1</sup> – 5 g/L humic acid for  $1-2 \times 10^8$  spores/mL
- $10^7$  spores with 250  $\mu$ g humic acid




<sup>1</sup> *J Appl Microbiol* 2015. 119:1263–1277

# Wet Heat Exposures of BIs


- Tested whether the presence of the ampoule would hinder the action of a sterilant/decontaminant in the BI



121°C, 20 min





100°C, 10 min



# Hot, Humid Air Exposures of BIs

- Tested whether the presence of the ampoule would hinder the action of a sterilant/decontaminant in the Bi
- 75 °C, 80% RH



# Summary

---

- Quantitative biological indicator system with Bt AI Hakam spores
- 10 – 12 hours to detect a single spore and < 3 hours for  $10^7$  spores
- Allows determination of the decontamination kinetics for modeling to plan decontamination schemes for emergency responses or facilitate developing new decontamination systems for biological agents
  
- Inclusion of data processing in the detector
- Conversion of the current prototype detector into a manufacturable form with standards compliance



# Acknowledgments

- DTRA
  - Dr. Charles Bass
  - Dr. Glenn Lawson
  - Mark Morgan
  - Kelly Crigger
- Dr. Tony Buhr at NSWCCD
- ARO
  - Dr. Jennifer Becker
  - Dr. Dawanne Poree
  - Wendy Mills
- Funding (JSTO-CBD)
  - DoD CBD SBIR Phase I Contract No. W911NF-11-C-0096
  - DoD CBD SBIR Phase II Contract No. W911NF-12-C-0048
  - DoD CBD SBIR Phase II Contract No. W911NF-16-C-0074
- Triton
  - John Lovaasen
- KBioSim
  - Kristopher Daly
- UConn Health Center
  - Dr. Peter Setlow
  - Dr. Barbara Setlow

The effort depicted here has been funded by the Joint Science and Technology Office for Chemical and Biological Defense (JSTO-CBD) through CBD SBIR Phase I and II programs under U.S. Army Research Office contracts. The content of the information does not necessarily reflect the position or the policy of the Government, and no official endorsement should be inferred.

