



Determination of Ethylene Oxide at Ultra Trace Concentrations in Ambient Air Using EPA Method TO-15A: Optimization of VOC Preconcentrator and GC-MS Analytical Method Parameters

*Air & Waste Management Association
Air Quality Measurement Methods and Technology Conference
November 14 – 16, 2023
Durham, NC*

Tamira Cousett, EPA Office of Research and Development (ORD)

Karen Oliver, EPA ORD

Carlton Witherspoon, Jacobs Technology

Andrew Whitehill, EPA ORD

Disclaimer


The views expressed in this presentation are those of the authors and do not necessarily represent the views or policies of the U.S. Environmental Protection Agency (EPA). Any mention of trade names, products, or services does not imply an endorsement by the U.S. Government or EPA. EPA does not endorse any commercial products, services, or enterprises.

Ethylene Oxide – An EPA Priority for Air

- EPA updated IRIS assessment for ethylene oxide (EtO) in December 2016
 - Cancer risk value is now ~60x higher than previous estimate
 - Based on lifetime inhalation unit risk estimate, 100-in-million cancer risk is 11 pptv
- EtO is a hazardous air pollutant of concern for communities near sources, including communities with environmental justice concerns
- ORD and Office of Air and Radiation developing methods and guidance for monitoring
- Our team evaluating measurement using EPA Method TO-15A

EtO Analytical Challenges

- IRIS value necessitates ultra-low method detection limit (MDL) of ~10 pptv
- National Air Toxics Trends Stations (NATTS) Technical Assistance Document (TAD):
 - EtO vendor gas standards instability
 - Some cleaned canisters may exhibit growth or contamination
 - Difficulty quantitating due to coelutions and instrument sensitivity





<https://www.epa.gov/system/files/documents/2022-08/NATTS-TAD-Revision-4-Final-July-2022-508.pdf>

Research Plan

- Evaluating state-of-the-art volatile organic compound (VOC) preconcentrator GC-MS systems
- Evaluating silico-ceramic lined canisters
- Optimizing settings to achieve robust measurements of EtO at low-pptv levels – cleanliness target for TO-15 of 200 pptv reduced to 20 pptv for TO-15A (our target EtO MDL is even lower)

Canister Cleaning

- Precleaning step with four fill/evacuation cycles using ultra-zero air and rough pump
- Xontech/Xonteck Model 960 systems
 - 4 hours total cycle time currently, may be modified
 - Research-grade nitrogen for fill gas
 - Fill gas humidified to 95% RH
 - Canisters heated to ~85 °F
 - Final high-vacuum finish (8 microns Hg)



Standards Preparation


- **Entech 4700 Static Diluter**

- Lower-concentration cylinder standards (<1 ppm) minimize serial dilutions
- 100 ppb TO-14/50 ppb EtO
- 100 ppb 4-component internal standard



- **Exposure Chamber**

- Canister-free approach to calibration



A Tale of Two Systems – “System 1”

Preconcentrator: Entech 7200A

GC: Agilent 8890

MS: 5977B Turbo EI

Source: Inert Plus Extractor



System 1 – Background Information

- Chosen to meet the low-pptv MDL and cleanliness requirements of Method TO-15A—and specifically as applicable for EtO
- Sensitive, inert, reasonably priced system
 - Preconcentrator introduced in 2020 offers improvements in both water management and system cleanliness compared to legacy 7200
 - Inert Plus MSD with extractor source offers better sensitivity and inertness than our aging Agilent systems

System 1 - Initial Method

Analytical Method

Parameter	Specifications
Column	DB1 30-m x 0.25-mm ID x 0.50- μ m column
Source/quad	350°C/200°C
Draw-out plate	3 mm
Tune	Etune (10x sensitivity)
Acquisition	Selected ion monitoring (SIM), low resolution
Gain	10 (tested multiple settings)
Column flow	0.6 mL/min


8890 GC Oven Program

	Rate (°C/min)	Value (°C)	Hold (min)
Initial	–	35	3
Ramp 1	20	80	0
Ramp 2	25	160	0
Ramp 3	40	250	1
Ramp 4	50	180	4
Post bake		35	0

Entech 7200A Preconcentrator settings specified in Entech Application Note [A-3742-01](#)

System 1 – Initial Method (cont.)


- Literature reviews/technical discussions
- Not sufficient
 - 80 pptv MDL* for EtO
 - Separation of coeluting compounds not achieved
 - GC oven program 16.35 min
- Further evaluated 12 analytical methods to reach our intermediate method



*Tentative MDL: 40 CFR, Part 136, Appendix B, 1984 – “best-case” using 7 replicates


System 1 – Method Refinement Highlights

10 ppbv TO14/EtO Blend




3/31/22	Rate (°C/min)	Value (°C)	Hold (min)
Initial	–	-30	5
Ramp 1	5	130	0
Ramp 2	20	220	3

5 ppbv TO14/EtO Blend



4/12/22	Rate (°C/min)	Value (°C)	Hold (min)
Initial	–	-30	2
Ramp 1	3	69	0
Ramp 2	4	141	0
Ramp 3	40	240	3.52

5 ppbv TO14/EtO Blend



4/12/22	Rate (°C/min)	Value (°C)	Hold (min)
Initial	–	-50	2
Ramp 1	3	69	0
Ramp 2	4	141	0
Ramp 3	40	240	3.52

System 1 – Intermediate Method

Analytical Method

Parameter	Specifications
Column	Rxi-624Sil MS 60-m x 0.25-mm ID x 1.4- μ m GC column
Source/quad	300°C/180°C
Draw-out plate	3 mm
Tune	Etune (10x sensitivity)
Acquisition	SIM, low resolution
Gain	10/25; time segmented gain
Column flow	1.5 mL/min

8890 GC Oven Program

	Rate (°C/min)	Value (°C)	Hold (min)
Initial	–	-20	2
Ramp 1	3	69	0
Ramp 2	4	141	0
Ramp 3	40	240	3.52
Post bake		180	6

Modified Entech 7200A Preconcentrator settings specified in Entech Application Note [A-3742-01](#)


System 1 – Intermediate Method (cont.)

- **Refinements to reach intermediate method:**

- Optimized Entech/GC-MS settings

- **Achievements:**

- Desired chromatographic separation
- Identified potential interferents
- MDL* range of 12–20 pptv



*Tentative MDL: 40 CFR, Part 136, Appendix B, 1984 – “best-case” using 7 replicates

System 1 – Intermediate Method (cont.)


- **Caveat:**

- Oven program 55+ min

- **Continued method refinements:**

- Checked coelutions: TO-15 subset, TO-14/EtO blend, and Photochemical Assessment Monitoring Stations (PAMS) VOCs
- Further adjust preconcentrator/MS settings/GC oven programs

Co-elutions with PAMS, TO14, EtO, and TO15 subset



System 1 – Optimized Method

Analytical Method

Parameter	Specifications
Column	Rxi-624Sil MS 60-m x 0.25-mm ID x 1.4- μ m GC column
Source/quad	230°C/150°C
Draw-out plate	6 mm
Tune	Atune/Etune evaluations
Acquisition	SIM, low resolution
Gain	Testing overall and segmented
Column flow	1.5 mL/min

8890 GC Oven Program


	Rate (°C/min)	Value (°C)	Hold (min)
Initial	–	-20	3
Ramp 1	1.5	1	0
Ramp 2	18	150	0
Ramp 3	12	240	3
Post bake		180	18

Modified Entech 7200A Preconcentrator settings specified in Entech Application Note [A-3742-01](#) to hybrid mode

System 1 – Optimized Method (cont.)

Method refinements to reach optimized method:

- Varied temperatures, ramps, hold times:
 - Dichlorotetrafluoroethane (F-114) and chloromethane
 - Styrene and *o*-xylene
 - Bromomethane and EtO
 - All compounds separated in the quantitation software
- Optimized oven program to 35.1 min



System 1 – Optimized Method (cont.)

- Implemented standard source and quad temperatures (230°C/150°C)
- 6 mm draw-out plate and new column (same phase) installed
- Evaluated Atune vs Etune
- Tentative MDLs* of 5–11 and 7–18 pptv were calculated for Atune and Etune, respectively
- Analyzed ambient air samples – data reduction in progress

*Tentative MDL: 40 CFR, Part 136, Appendix B, 1984 – “best-case” using 7 replicates

System 1 – Tentative MDLs*

MDL Studies - Entech 7200A/8890A GC - 5977B MS						
Study	Tune	Source/Quad, °C	Lens (mm)	Gain	EtO MDL (pptv)	Timeline
03/07/22	Etune	350/200	3	10	79	initial
05/04/22	Etune	350/200	3	10	24	initial
05/25/22	Etune	350/200	3	10	9	initial
08/08/22	Etune	350/200	3	10	17	intermediate
10/25/22	Etune	350/200	3	10	20	intermediate
03/07/23	Etune	300/180	3	10	12	intermediate
03/30/23	Etune	300/180	3	10	19	intermediate
05/23/23	Etune	300/180	6	25	7	optimized
06/05/23	Etune	300/180	6	25	10	optimized
06/08/23	Etune	300/180	6	25	8	optimized
06/15/23	Etune	300/180	6	25	11	optimized
08/17/23	Atune	230/150	6	default	10	optimized
08/23/23	Etune	230/150	6	default	18	optimized
08/30/23	Etune	230/150	6	25	8	optimized
09/13/23	Atune	230/150	6	default	10	optimized
10/17/23	Atune	230/150	6	default	5	optimized
10/30/23	Atune	230/150	6	25	11	optimized

*Tentative MDL: 40 CFR, Part 136, Appendix B, 1984 – “best-case” using 7 replicates

System 1 – NATTS PT Results

- Solid NATTS CY2023 Q2 results – met TO-15A $\pm 30\%$ criteria – our goal
- Analytical parameters:
 - Etune
 - SIM mode
 - 230°C/150°C
 - Default gain factor of 1
 - MDL* for EtO: 18 pptv

*Tentative MDL: 40 CFR, Part 136, Appendix B, 1984 – “best-case” using 7 replicates

Percent Difference against the NATTS Lab Average, per laboratory (%)		
Analyte	Target Conc. (ppbv)	AMCD AAB
1,1,2,2-Tetrachloroeth.	0.034	-2.9%
1,2-Dibromoethane	0.29	-3.1%
1,2-Dichloroethane	0.22	-3.6%
1,2-Dichloropropane	0.20	0.5%
1,3-Butadiene	0.16	-25.6%
1,3-Dichloropropene -	0.16	7.5%
1,3-Dichloropropene -	0.15	8.0%
Acrolein	0.16	NR
Benzene	0.17	-6.5%
Carbon Tetrachloride	0.056	-26.8%
Chloroform	0.26	-6.2%
Dichloromethane	0.23	-9.6%
Ethylene Oxide	0.11	16.4%
Tetrachloroethylene	0.12	7.5%
Trichloroethylene	0.14	-0.7%
Vinyl Chloride	0.15	-8.7%
overall average		-3.6%

System 1 - Challenges

- Residual water may be impacting MS internal components
- Instrument stability issues potentially due to temperature stress on source/quad
- Etune vs. Atune operational issues
- Poor Etunes are a first indicator of system beginning to fail; need to check voltages; source cleanings not necessarily sufficient
- Considerable downtime over last year

System 1 - Conclusions

- Excellent collaborations on optimization of instruments and methods
- Thoroughly evaluated/optimized method parameters, completed SOP
- Etune with gain plus SIM mode provides optimal sensitivity
- Settled on standard 230°C/150°C source/quad temperatures
- 6 mm draw-out plate is viable option vs 3 mm draw-out plate
- Tentative MDLs in 5-18 pptv range (combined Atune and Etune results)
- Solid audit results for EtO

A Tale of Two Systems – “System 2”

Preconcentrator: Entech 7200
(refurbished)

GC: Agilent 7890A

MS: 5975C with diffusion pump

Source: Inert with triple-axis
detector



System 2 – Background Information

- Evaluating legacy system viability for determination of ultra-low EtO concentrations
- Refurbishment of 9-year-old Entech 7200 completed
- Diffusion pump rather than turbopump system for MS

System 2 – Current Method

Analytical Method


Parameter	Specifications
Column	Rxi-1ms 60-m x 0.32-mm x 1- μ m column
Source/quad	230°C/150°C
Draw-out plate	Standard
Tune	HiSense/Atune
Acquisition	SIM, low resolution
Gain	Various settings tested
Oven program	Ambient starting temp

7890A GC Oven Program

	Rate (°C/min)	Value (°C)	Hold (min)
Initial	–	35	5
Ramp 1	5	130	0
Ramp 2	20	220	3
Post bake		220	7


System 2 – Current Method (cont.)

- Legacy ambient GC oven program and MS method in place, GC oven program allowed solid separations
- Evaluated Trace Ion Detection (TID) and HiSense modes
- Challenges with peak integration for EtO below 125 pptv
- Length of oven program 31.5 min




125 pptv EtO

125 pptv Perc



System 2 - Conclusions

- EtO viable on Rxi-1MS column; <0.25 ppbv difficult to integrate; may try 624 phase column
- HiSense tune improved sensitivity but not sufficiently
- EtO tentative MDLs*:
 - Atune – 96 pptv
 - Atune TID – 60 pptv
 - HiSense – 27 pptv



*Tentative MDL: 40 CFR, Part 136, Appendix B, 1984 – “best-case” using 7 replicates

“System 3”

Preconcentrator: Markes CIA Advantage/Kori-xr/Unity-xr

Agilent GC: 7890B

MS/source: Bench-TOF Select
time-of-flight MS (TOF MS)



System 3 – Background Information

- GC-TOF MS system – research-grade instrument fully operational in lab
- Different preconcentration/water management system and MS
- As time allows, implement optimized method (modified as appropriate) using TO-15A sorbent traps
- Technical Note “Advances in ultra-trace air toxics analysis: Cryogen-free canister pre-concentration and GC–MS for US EPA Method” available (<https://markes.com/content-hub/application-notes/application-note-169>)

Next Steps

- Installation and evaluation of 5977C diffusion pump MS
- Continue analysis of ambient air samples
- Continue sample integrity studies
- Move to Method Update Rule MDL procedure
- Purchase new canisters
- Evaluate and optimize canister cleaning processes
- GC-TOF MS system experiments and methods development
- Conduct interlaboratory comparison

Acknowledgments

For valuable technical discussions regarding optimization of analytical methodology...

- Jeff Schroeder, Entech Instruments
- Yvonne Brown, Agilent Technologies, Inc.
- Ericka Hachmeister, Hannah Calder, and Nathan Shafer, Markes International
- Jason Hoisington and Jason Herrington, Restek Corporation

For technical and content editing...

- Stacy Henkle, Jacobs Technology



References

Agilent Technologies. 2006. The 5975C Series MSDs: Normalized Instrument Tuning: Technical Overview. Agilent Technologies, Inc., 5301 Stevens Creek Blvd, Santa Clara, CA 95051. <https://www.agilent.com/Library/technicaloverviews/Public/5989-6050EN.pdf>

Agilent Technologies. 2007. The 5975C Series MSDs: Method Optimization and Trace Ion Detection Technical Overview. Agilent Technologies, Inc., 5301 Stevens Creek Blvd, Santa Clara, CA. <https://www.agilent.com/cs/library/technicaloverviews/public/5989-6425EN.pdf>

Agilent Technologies. 2013. A Quick Start Guide to Optimizing Detector Gain for GC/MSD. Agilent Technologies, Inc., 5301 Stevens Creek Blvd, Santa Clara, CA. <https://www.agilent.com/cs/library/technicaloverviews/public/5991-2105EN.pdf>

Entech Instruments. 2018. Application Note A-3742-01, Fast US EPA TO15 Analysis Providing Higher Productivity Using the New Entech 7200A and Agilent 7890B/5977 GCMS. Entech Instruments, 2207 Agate Court, Simi Valley, CA 93065.

Hoisington, J. and J.S. Herrington. 2021. Rapid determination of ethylene oxide and 75 VOCs in ambient air with canister sampling and associated growth issues. Separations 8, 35. <https://doi.org/10.3390/separations8030035>.

U.S. EPA. 2019. Method TO-15A: Determination of Volatile Organic Compounds (VOCs) in Air Collected in Specially Prepared Canisters and Analyzed by Gas Chromatography–Mass Spectrometry (GC-MS). Research Triangle Park, NC: Office of Research and Development, National Exposure Research Laboratory, and Office of Air Quality Planning and Standards, Air Quality Assessment Division. https://www.epa.gov/sites/production/files/2019-12/documents/to-15a_vocs.pdf (last accessed August 30, 2023).

U.S. EPA. 2023. J-AMCD-AAB-SOP-5188-0. Standard Operating Procedure for Analysis of Volatile Organic Compounds in Whole Air Samples Using the Entech 7200A Preconcentrator and Agilent 8890/5977B Gas Chromatograph–Mass Selective Detector. U.S. Environmental Protection Agency, Center for Environmental Measurement and Modeling.

U.S. EPA. 2023. J-AMCD-0033182-QP-1-2. Quality Assurance Project Plan: Guidance Document for Determination of Ethylene Oxide in Ambient Air Using EPA Method TO-15A. U.S. Environmental Protection Agency, Center for Environmental Measurement and Modeling.

<https://www.epa.gov/hazardous-air-pollutants-ethylene-oxide>

<https://www.epa.gov/sites/default/files/2021-03/documents/eto-method-for-natts-labs-2019.pdf>

<https://www.agilent.com/cs/library/eseminars/public/mass-spectrometer-optimization-agilent-webinar-janyary232020.pdf>

<https://www.epa.gov/hazardous-air-pollutants-ethylene-oxide>

https://www.agilent.com/cs/library/eseminars/public/AFO_Simple_Strategies_for_GCMS_04112018.pdf