

Manufacturer: Intini Marketing Inc.  
Model: Medi-Lock Slider Pack  
ASTM Type: IXF

Description:

This package consists of two pieces, an outer plastic shell and a permanently attached inner plastic dispensing slider (Figure 1A). The inner slider tray has a rectangular dispensing compartment that holds one tablet when opened (Figure 1B).



Figure 1A



Figure 1B

The package is opened by pushing the button on the top of the package (black arrow, Figure 1A), while simultaneously pushing the slider tray in (blue arrow Figure 1A) to reveal the dispenser compartment (Figure 1B). The slider tray is spring loaded and will return to the closed position when the slider tray is no longer pressed. (Note the tray is fixed in the open position in Figure 1B for illustrative purposes).

The manufacturer's instruction on the button of the package is below (Figure 2).

PUSH



Figure 2

The manufacturer's instruction on the slider tray is below (Figure 3).

Push



Figure 3